

SOFTWARE

Educational Training Guide



NMI ParaVision 360 + pmod 4.5 Workflow Basics

SOFTWARE

Educational Training Guide

PET/CT & PET/MR: Animal Transport System Workflow

PET/MR & PET/CT Animal Transport System (ATS)* Workflows



- 1. Camera.** After placing the animal in the cradle set the ATS to camera.
- 2. Position.** Move the crosshair to the center of your target object, noting the boundaries of the FOV(s) displayed.
- 3. Set the position for transport.** **Tip: Set the position for each subject and study.**
- 4. Select the Modality.** Select the modality button to drive the transport to PET, CT or MR.
- 5. ATS Reference** box may be adjusted to move transport during study.

**For systems equipped with the ATS only.*

SOFTWARE

Educational Training Guide

PET/CT & PET/MR: ParaVision 360 Study Workflows


PET/CT & PET/MR: ParaVision 360 Study Workflows

Key Menus

- 1. Study Registration.** The process for entering study registration details (intended for a single animal at a single time point).
- 2. Protocols & Multimodal Scan Programs.** Multimodal Scan Programs are predefined (or customized) PET + CT (or MR) scans programs with predefined scan protocol order, processing (stitching, attenuation correction, etc.).
- 3. Palette.** Most centralized location for Scan Programs, Image Fusion, and Dataset Filtering. **Tip: Return to Palette for access to functions in central location.**
- 4. Dataset Browser.** Location for accessing/searching data with search menus, DICOM export, and other study management tools.

PET/CT & PET/MR: ParaVision 360 Study Workflows

1. Study Registration



ANIMAL

2 **Animal Name** Animal1

Animal ID Animal1

Specimen Rodent

Gender UNKNOWN

Date of Birth / Transfer December 2, 2023

Age 90 Days

Additional Info

STUDY

4 **Study Name** XenograftTumor

Operator nmrsu

Location Rat/Body

Position HEAD_PRONE

Animal Weight 250 g

Animal Cradle Use Animal Transport System ☒

Additional Info Day 1

5

Modality

1 ☒ PET ☐ CT

Isotope F-18

Compound FDG

Administration Route Unknown Route

6 **Initial Activity** 200 µCi

Calibration Time March 1, 2024 08:12:58 Now

Injection Time March 1, 2024 08:12:58 Now

Residual Activity 0 µCi


Residual Time March 1, 2024 08:12:58 Now

7 **Create** **Exam** **Close**

- Modality.** Check the PET and CT or MR modality boxes. **Tip: required to be set before setting other parameters.**
- Animal.** Enter the Animal Name & Animal ID.
- Study Name.** Enter a Study Name.
- Select the Location.** Select the Location defines the scan program and protocols for your acquisitions.
- Additional Info.** Enter additional information (e.g. study day in a series). **Tip: This variable is useful when sorting data in a PMOD.**
- Isotope/Compound/Activity.** Enter the activity and calibration time.
- Create, Exam.** Select Create and Exam to initiate a study.

PET/CT & PET/MR: ParaVision 360 Study Workflows




2. Protocols & Multimodal Scan Programs



1. **Instruction Card.** Scan Programs & Protocols can be drag and dropped to the instruction queue for acquisitions.
2. **Multimodal Scan Programs** are found under Palette > Explorer. These are interleaved PET & CT or MR scan protocols with pre-optimized with predefined AC processing.
3. **Protocols.** Individual PET, CT or MR scan protocols that can be added to the instructions. **Tip: Always perform the PET Daily QA workflow at the beginning of the day. Include the CT protocol that will be used to allow for CT warmup & flat-fielding.**
4. **Continue.** Select Continue to initiate scans.

PET/CT & PET/MR: ParaVision 360 Study Workflows


2. Custom Multimodal Scan Programs



Caution! Please do not Overwrite Existing Scan Programs

PET/CT & PET/MR: ParaVision 360 Study Workflows


3. Palette



1. **Datasets.** Datasets for the active study are found under Palette > Explorer. Thumbnails for images appear after reconstruction.
2. **Right Click** on image thumbnails for common processing tools such as: 3. **Image Fusion**, 4. **Create Image Series**, and 5. **Export to DICOM**.

PET/CT & PET/MR: ParaVision 360 Study Workflows

4. Database Browser




1. **Database Browser.** Use database browser to filter data by Subject and Study.
2. **Export to DICOM.** Select Export to DICOM or PMOD server.
3. Select **Add Study** to initiate to add a new study registration in preparation for your next study in a longitudinal time series. Modify the Study Name and Additional Info.

SOFTWARE

Educational Training Guide

pmod v4.5 Multi Animal Image Cropping – Preprocessing in View (Skip this Step if Single Animal)

PET/MR & PET/CT: pmod Multi Animal Image Cropping – Preprocessing in View



1. Open the **View** module.

2. From "File" menu, select option: Load image data > Database.

3. Select the **Subject Name** from the menu (in this example 4 Mouse Demo).

4. Highlight a PET and Anatomical images in the **Series** tab.

5. Select **Add Selected series**.

6. Select **Open**.


Subject Name	Subject ID	Modification Date	Sex	Date of Birth
4 Mouse Demo	4 Mouse Demo	2025-08-08 14:05:52.992	O	2024.07.30

Subject Name	Study Date	Study Time	Series Date	Series Time	Study Description	Series Description	Modification	Last Use	Mod
4 Mouse Demo	2024.09.17	08:16:24	2025.07.31	06:40:46	Static_20min 0.25 36it	Static_20min 0.25 36it	2025-08-08 14:00	2025-08-11 11: PT	PT
4 Mouse Demo	2024.09.17	08:16:24	2025.07.31	06:40:46	AnatomicalRef_4mice	AnatomicalRef_4mice	2025-08-08 14:00	2025-08-11 11: CT	CT

Subject Name	Study Date	Study Time	Series Date	Series Time	Study Description	Series Description	Modification	Last Use	Mod	nt	nz	nx	ny	Organ	At
4 Mouse Demo	2024.09.17	08:16:24	2025.07.31	06:40:46	Static_20min 0.25 36it	Static_20min 0.25 36it	2025-08-08 14:00	2025-08-11 11: PT	PT	1	728	384	384	WHOLEBOD	0.
4 Mouse Demo	2024.09.17	08:16:24	2025.07.31	06:40:46	AnatomicalRef_4mice	AnatomicalRef_4mice	2025-08-08 14:00	2025-08-11 11: CT	CT	1	728	384	384	WHOLEBOD	0.


1. Open the **View** module.
2. From "File" menu, select option: Load image data > Database.
3. Select the **Subject Name** from the menu (in this example 4 Mouse Demo).
4. Highlight a PET and Anatomical images in the **Series** tab.
5. Select **Add Selected series**.
6. Select **Open**.

PET/MR & PET/CT: pmod Multi Animal Image Cropping – Preprocessing in View



1. Split Beds in the Select “Tools” menu. Select 4 beds or 3 beds from the pull down and select Yes. 4 (or 3) cropped images will be created.
2. Select the Anatomical Image and select Split Beds in the “Tools” menu again to split the anatomical images.

PET/MR & PET/CT: pmod Multi Animal Image Cropping – Preprocessing in View




SOFTWARE

Educational Training Guide

pmod v4.5 Multimodal Image Fusion & Display Workflows

PET/MR & PET/CT: pmod Multimodal Image Fusion & Display Workflow Basics


- 1. Multimodal Image Fusion.** Workflow for fusion of PET data to CT or MR data.
- 2. Multimodal Image Masking/Cropping.** Process for subtraction of hardware components in CT image.
- 3. Multimodal Image Display.** Workflow for Image Capture of Linear and/or MIP image display.



Tip: Begin by referencing the “PMOD File Management & Starting Reference For Bruker PET Data” to configure default application menus for simple workflows.

PET/MR & PET/CT: pmod Multimodal Image Fusion & Display

1. Multimodal Image Fusion



The screenshot displays the Bruker pmod software interface with the following components and steps:

- Step 1:** The **Fusion** module is selected in the left sidebar.
- Step 2:** The **File** menu is open, showing the option **Load image data > Database**.
- Step 3:** The **Subject Name** field is selected in the top form.
- Step 4:** A PET image is highlighted in the **Series** tab table.
- Step 5:** The **Add Selected series** button is selected.
- Step 6:** The **Open** button is selected.

Series [2] Table:

Subject Name	Study Description	Mod	Series Description	Study Date	Study Time	Series Date
Animal1	Day 1		AnatomicalReferenceStandardContrastMultiPositionAC	2024.07.17	09:47:23	
Animal1	Day 1		ImReg(int).Src E5:P20, ref. E6:P6	2024.07.17	09:47:23	2025.03.06


Selected for loading [1] Table:

Subject Name	Study Description	Mod	Study Date	Study Time	Series Date	Series Time	Series Description	Modification	Last Use
Animal1	Day 1	PT	2024.07.17	09:47:23	2025.03.06	08:32:49	ImReg(int).Src E5:P20, re 2025-03-25 12:14: 2025-03-25		

1. Open the **Fusion** module.
2. From "File" menu, select option: Load image data > Database
3. Select the **Subject Name** from the menu (in this Example Animal 1).
4. Highlight a PET image in the **Series** tab.
5. Select **Add Selected series**.
6. Select Open.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

2. Multimodal Image Fusion



1. From "File" menu, select option: Load reference data > Database

2. Select the **Subject Name** from the menu (in this Example Animal 1).

3. Highlight the CT (or MR) image in the **Series** tab.


4. Select Open.

Subject Name	Study Description	Mod	Series Description	Study Date	Study Time
Animal1	Day 1	CT	AnatomicalReferenceStandardContrastMultiPositionAC	2024.07.17	09:47:23
Animal1	Day 1	PT	ImReg(int)Src E5.P20, ref. E6.P6	2024.07.17	09:47:23

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


2. Multimodal Image Masking

1. In workflow controls select **VOI ANALYSIS** sub-page.



PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


2. Multimodal Image Masking



1. Place the orthogonal crosshairs in the image center. To create a VOI based on the Bruker MMPF cradle, first select **Create Regular VOI > Organs (Predefined VOI)**.
2. Select **Rodent Beds > Bruker MMPF Mouse**.
3. Select the **Operation on Entire VOI** button.
4. Drag the VOI at the crosshair to align with the cradle.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


2. Multimodal Image Masking



1. Select **Tab B** to set the current active working image data to the CT (or 3D MR) data.
2. Select the “**Masking...**” tab.
3. Select the “**Mask voxels inside selected VOI(s)**” button.
4. For Hounsfield calibrated CT enter “-1000” in the dialogue, and select yes.
5. Finally select **Crop**. Then, change to **Tab A** and select **Crop**.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


3. Multimodal Image Display



1. Select the **Image Display Layout** tab.
2. Select the **Linear arrangement** display selection, common for display in figures.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


3. Multimodal Image Display



PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

3. Multimodal Image Display

1. Set the first image to PET and second image to MR or CT using the pulldown menus A and B.



SOFTWARE

Educational Training Guide



PET/MR & PET/CT Software Workflows: pmod 4.5 VOI Basic Workflows

PET/MR & PET/CT: pmod 4.5 VOIs Basics

VOIs by Iso-Contouring and Contour Interpretation

1. VOIs & Iso-Contouring by Region Growing

2. VOIs & Freehand with Contour Interpretation




Study data (courtesy Virginia Tech) was acquired in a BioSpec 9.4 MRI with a PET Insert Si103 and analyzed using PMOD v4.4

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

1. pmod PFUS VOIs & Iso-Contouring by Region Growing


1. Load multi-modal data and view at the VOI ANALYSIS sub-page menu.



Tip: The Iso-Contouring By Region Growing VOI method is suitable with targets with well defined anatomical or PET signal boundaries

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

1. pmod PFUS VOIs & Iso-Contouring by Region Growing




1. Select **Tab B** to set the current active working image to the CT or MR reference.
2. Set the **Fusion Balance** slider all the way to the right, to visualize only Tab B (fusion slider value 1.0).
3. Select **Create regular VOI**. Select and choose SPHERE (analytic object).
4. Select the **Operation on Entire VOI** button.
5. Adjust the location and boundaries of the sphere to outside the tumor margins.
6. Select **Save VOI**. Tip: If additional studies will be analyzed, it may be useful to save the sphere VOI to the database for recall.

Tip: The Iso-Contouring By Region Growing VOI method is suitable with targets with well defined anatomical or PET signal boundaries


PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

1. pmod PFUS VOIs & Iso-Contouring by Region Growing



PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

1. pmod PFUS VOIs & Iso-Contouring by Region Growing




1. **Tab A.** Select Tab A to set the active image to the PET input data.

2. **VOI Statistics.** Select the VOI Statistics button to view VOI statistic. Tip: Choose “Selected Statistics” or “Default Statistics” when prompted to view your predefined subset of statistics or an interactive menu of statistics respectively.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


2. pmod PFUS VOI & Freehand with Contour Interpretation



1. Select **Tab B** to set the active working image to CT or MR reference data.
2. Select the **Draw Polygon with Dense Vertices** button.
3. Tracer the ROI boundaries on a selected slice.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

2. pmod PFUS VOI & Freehand with Contour Interpretation




1. **VOI.** Trace the ROI boundaries on additional non-contiguous slices. Slices must be drawn in a single axis only (in this example: Z axis – indicated by red arrow).

Tip: The Freehand with Region Growing VOI method is useful for targets with faint boundaries.

PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display


2. pmod PFUS VOI & Freehand with Contour Interpretation



1. **Masking Tab.** Select the “Masking...” tab.
2. **Contour Interpolation.** Select the **Contour Interpolation** button to fill the non-contiguous spaces.


PET/MR & PET/CT: pmod Multimodal Imaging Fusion & Display

2. pmod PFUS VOI & Freehand with Contour Interpretation



1. **Tab A.** Select Tab A to set the active image to the PET input data.

2. **VOI Statistics.** Select the VOI Statistics button to view VOI statistic.



Innovation with Integrity