

# Top-Down Mapping of Protein Modifications by Multimodal MS<sup>n</sup> on the New timsOmni Platform

Francis Berthias<sup>1</sup>, Nurgül Bilgin<sup>2</sup>, Athanasios Smyrnakis<sup>3</sup>, Mariangela Kosmopoulou<sup>3</sup>, Christian Albers<sup>4</sup>, Jasmin Mecinović<sup>2</sup>, Dimitris Papanastasiou<sup>3</sup>, and Ole N. Jensen<sup>1</sup>

<sup>1</sup> Department of Biochemistry and Molecular Biology; <sup>2</sup> Department of Physics, Chemistry and Pharmacy, University of Southern Denmark, Odense, Denmark  
<sup>3</sup> Fasmatech Science & Technology SA, NCSR Demokritos, Athens-Chalandri, Greece, <sup>4</sup> Bruker Daltonics GmbH & Co. KG, Bremen, Germany

## Aims

- Evaluate multimodal Omnitrap™ (MS<sup>2</sup>/MS<sup>3</sup>) workflows for rapid top-down proteoform identification and high-coverage amino acid sequencing of intact proteins.
- Precisely localize Post-translational modification (PTM) sites on intact histone proteins, capturing co-occurring modifications and resolving positional isomers.

## Introduction

**Biological context:** Histone proteins carry dense array of PTMs that orchestrate chromatin structure and gene expression. Their spatial patterns, not merely their presence, drive functional outcomes and are frequently disrupted in disease.


- Analytical challenge:** determining the presence and exact location of PTMs on intact proteins remains challenging, especially when isomers (same m/z) are involved[1].
- Complementary fragmentation:** Combining electron- and collision-based dissociation in one device broadens sequence coverage and enables residue-level discrimination of PTM positional isomers[2].
- This study:** We deploy eXd MS<sup>2</sup>/MS<sup>3</sup> on the new timsOmni platform to localise multiple and co-occurring lysine acetylations on histones H3/H4.

## timsOmni Platform


timsOmni™ integrates TIMS-Q-TOF with a Omnitrap™ platform for advanced top-down proteomics[3].

**Omnitrap™ eXd:** The Omnitrap™ segmented linear ion trap integrate multimodal fragmentation workflow delivering the rich spectral detail required for precise PTM mapping on intact proteins.


## Multimodal MS<sup>n</sup> eXd workflows:



## Proteoform Mapping of Multiply Acetylated Histone H4 (11.2 kDa)




## In vitro acetylation reaction with p300




Mass spectrum reveals H4 proteoforms bearing 1 to 7 acetylations following p300 treatment (1h)

## eXd MS<sup>2</sup> and MS<sup>3</sup> of H4(Kac)6 (13+)




## Proteoform identification

### H4K5acK8acK12acK16acK20acK44ac - (13+)



Fragment ion types	Sequence coverage
ECD	89.1 %
EID	94.1 %
ECD-rCID	89.1 %
rCID	51.5 %
Combined	100 %

### H4K5acK8acK12acK16acK20acK31ac - (13+)



Fragment ion types	Sequence coverage
ECD	87.1 %
EID	91.1 %
ECD-rCID	89.1 %
rCID	49.5 %
Combined	99 %

### H4K5acK8acK12acK16acK31acK44ac - (13+)




Fragment ion types	Sequence coverage
ECD	86.1 %
EID	91.1 %
ECD-rCID	85.1 %
rCID	52.5 %
Combined	99 %

(Error of 5-6 ppm and intensity deviation of 0.15).


- ~100% amino acid sequence coverage achieved using eXd workflow, enabling residue-level PTM mapping of complex proteoforms.
- eXd MS<sup>2</sup>/3 of H4(Kac)<sub>6</sub>-(13<sup>+</sup>) localized 7 acetylation sites: K5–K20, K31, and K44. Results were consistent at z = 16<sup>+</sup>.
- For lower acetylation states (Kac1–Kac4), acetylation preferentially targeted N-tail lysine residues (K5–K16).

## Site-Specific Monoacetylation of H3.1 at K14 by GCN5 and PCAF


### H3.1 + GCN5 (10 min)



### H3.1 + PCAF (30 min)



Enzymes GCN5 and PCAF are highly specific and K14 is identified as acetylation site for H3.1(Kac)<sub>1</sub>



- H3.1K14ac is confidently identified as the dominant monoacetylated proteoform, even with minimal MS/MS averaging.
- Sequence coverage increase logarithmically with the number of MS/MS spectra.

## Conclusion

- Multimodal MS<sup>n</sup> (eXd) on the timsOmni™ platform enables residue-level mapping of co-occurring lysine acetylations on intact histones H3.1 and H4.
- Nearly ~100 % sequence coverage is achieved.
- For heavily acetylated H4, we confidently localized up to 7 acetylation sites using combined MS<sup>2</sup>/MS<sup>3</sup> strategies, with consistent results across charge states.
- Monoacetylated H3 proteoforms from GCN5 and PCAF treatments were also characterized, confirming enzyme specificity (e.g. H3K14ac).

## Funding

PLATO: Danish National Mass Spectrometry Platform for Proteomics and Biomolecular Imaging is supported by a generous grant from the Danish Agency of Higher Education and Science (grant no. 5229-00012B, [www.sdu.dk/PLATO](http://www.sdu.dk/PLATO)).  
+ Jasmin funding

## Conflict of Interest Disclosure

C. A., A. S., M. K., D. P. are employees of Bruker Daltonics, which develops and manufactures timsOmni.