



Single atom chemical identification of TMD defects in ambient conditions

by E. J. Dunn, A. J. Robson, R. J. Young, S. P. Jarvis

Nanotechnology 37, 13 (2026): 135705

OPEN ACCESS ARTICLE

Hello {{lead.First Name:default=researcher}},

For this month's [AFM Materials Journal Club](#), I would like to share with you a paper titled, "[Single atom chemical identification of TMD defects in ambient conditions](#)."

Observation of atomic point defects with AFM remains challenging, particularly under ambient conditions. In this paper, researchers from Lancaster University and University of Sheffield demonstrate that [conductive atomic force microscopy \(C-AFM\)](#) can achieve single-atom resolution imaging of defects in 2D materials while providing the sensitivity to the local electronic structure needed to infer defect identity.

The authors use a Bruker [MultiMode AFM](#) to investigate mechanically exfoliated MoS₂, WSe₂, and WS₂ monolayers on conductive HOPG substrates under ambient conditions. The atomic resolution C-AFM reveals two dominant defect classes:

1. point defects localized to single atomic sites, and
2. less common diffuse defects extending several nanometers and affecting multiple lattice sites.

Defects were observed and densities quantified in all three materials, but characterization was challenging in WS₂ due to strong moiré patterns arising from the underlying HOPG substrate.

The chemical sensitivity of the method was then investigated by intentionally doping WSe₂ using nitrogen plasma. C-AFM distinguishes oxygen substitutions (OSe) from nitrogen substitutions (NSe) based on subtle differences in current contrast between defect sites.


Finally, this work provides insight into the mechanism of ambient atomic resolution C-AFM, attributing it to stable nanoscale asperities at the diamond coated probe apex rather than insulating tunneling barriers.

I hope you find this article interesting. For questions, please feel free to [reach out to me](#) or anyone in our [AFM applications team](#).

**Bede Pittenger, Ph.D.**


Senior Staff Applications Scientist
Bruker Nano Surfaces and Metrology
[\[email protected\]](#)

Bruker Technology Featured in this Article

**Conductive AFM (C-AFM)**

High-resolution current and topography mapping.

[Learn More](#)

**MultiMode 8-HR AFM**

The benchmark advanced-research small-sample AFM for highest resolution and extendable versatility.

[Learn More](#)

Engage with the AFM Materials Community

**AFM Summer Savings Event**

Take advantage of special savings on select Bruker AFM and BioAFM products and accessories through September. Save up to 50% on a wide range of components, accessories, and upgrades designed to support your next research project.

[AFM Summer Savings](#)

**AFM Materials Journal Club Archive**

Browse the complete archive of articles and write-ups featured in our AFM Materials Journal Club series.

[Browse Archive](#)

Contact a Bruker AFM Expert



Get in touch with an expert to discuss your challenges and schedule a personalized virtual demo.

[Contact Expert](#)



Innovation in AFM: A Timeline of Scientific and Technical Progress

See how AFM developed from its origins to today – and how Bruker innovation helped drive major advances along the way.

[READ MORE](#)

Bruker Nano Surfaces and Metrology

Atomic Force Microscopy

112 Robin Hill Road

Santa Barbara, CA 93117 USA

+1 (805) 967-2700

[\[email protected\]](#)



© 2026 Bruker Corporation