Introduction

Trapped ion mobility spectrometry coupled with quadrupole time-of-flight mass spectrometry (timsTOF Pro) offers a unique dimension of characterization and separation in complex samples.[1] The timsTOF Pro instrument also enables the "Parallel Accumulation Serial Fragmentation" (PASEF) method that allows continuous ion accumulation, sorting and eluting, therefore reaching a duty cycle of nearly 100%.[1] Notably, the timsTOF Pro powered by PASEF offers sequencing speed of more than 100 Hz.[1] In this work, we tested the effects of PASEF in shotgun proteomic analysis of a HeLa digest sample using the timsTOF Pro instrument. Proteomic data from the same sample were also collected using a Thermo Orbitrap Fusion mass spectrometer. The results from these two instruments were compared.

[1] https://www.bruker.com

Methods

A tryptic peptide mixture derived from HeLa cells (Pierce) mixed iRT (Biognosys) (100 ng or 200 ng) was loaded onto a C18 column (25 cm X 75 μm, 1.6 μm, lonOpticks, Australia) using a nanoElute UHPLC that was coupled to the timsTOF Pro mass spectrometer (Bruker), or a self packed column (5 cm x 150 μm, 1.9 μm) (abbreviaed as selfpack hereafter) using a NanoLC coupled to the Orbitrap Fusion mass spectrometer (Thermo Fisher Scientific). Chromatographic separation was carried out using a linear gradient of 2-37% of acetonitrile with 0.1% formic acid (FA) in water with 0.1% FA at a flow rate of 400 nL/min over 90 min (lonopticks) or at a flow rate of 800 nL/min over 75 min (Selfpack). The acquired data were analyzed for peptide and protein identifications using the Peaks Studio 8.5 (Bioinformatics Solutions).

Data analysis scheme:



Data highlighted were listed in Tables 1 and 2.

Performance of a timsTOF Pro mass spectrometer in shotgun proteomics of a tryptic HeLa digest mixture

Guoting Qin,^{1,2} JongMin Choi,³ Hong Lin,⁴ Sung Jung,³ Chengzhi Cai¹

¹Department of Chemistry, Mass Spectrometry Laboratory, ²College of Optometry, University of Houston, Houston, TX 77204 USA ³Proteomic Core, Baylor College of Medicine, Houston, TX 77030 ⁴Department of Computer and Mathematical Sciences, University of Houston - Downtown, Houston TX 77002





Figure 2. The intensity of iRT peptides determined by A) timsTOF Pro using a commercial column from lonopticks, B) timsTOF Pro using a self packed column, and C) Fusion using a self packed column. One hundred nanograms of HeLa-iRT sample were used for all experiments. Data were expressed as mean ± SD (n=3). Raw data were processed by Peaks Studio 8.5. Detailed column information was in Methods section.

Table 1. Comparison between timsTOF Pro and Fusion using selfpack columns.

| | HeLa + iRT (ng) | Matched peptide [2] | Peptide within defined RT & m/z window [3] | Peptide only found by timsTOF Pro or Fusion [4] | Peptide in common [4] | Percentage [5] |
|-------------|--------------------|------------------------|-----------------------------------------------|----------------------------------------------------|--------------------------|----------------|
| timsTOF Pro | 100 | 21377 ±363 | 2129±67 | 8407 | 155 ± 11 | 1.8 ± 0.1 |
| | 200 | 23166 ± 629 | 2512 ± 37 | 8678 | 195 ± 10 | 2.2 ± 0.1 |
| Fusion | 100 | 14434 ± 418 | 250 ± 33 | 3202 | 5 ± 2 | 0.2 ± 0.1 |
| | 200 | 17975 ± 345 | 228 ± 31 | 5451 | 7 ± 1 | 0.1 ± 0.0 |

Table 2. Comparison between a commercial column from lonopticks and a selfpack column using the timsTOF Pro.

| | HeLa + iRT (ng) | Matched peptide [2] | Peptide within defined RT & m/z window [3] | Peptide only found by Ionopticks or Selfpack [4] | Peptide in common [4] | Percentage [5] |
|------------|--------------------|------------------------|-----------------------------------------------|-----------------------------------------------------|--------------------------|----------------|
| lonopticks | 100 | 25849 ± 211 | 2275 ± 68 | 8407 | 144 ± 15 | 1.9 ± 0.2 |
| | 200 | 31776 ±684 | 2957 ± 72 | 11746 | 230 ± 6 | 2.0 ± 0.1 |
| Selfpack | 100 | 21377 ± 363 | 2129 ± 67 | 4137 | 128 ± 55 | 1.5 ± 0.7 |
| | 200 | 23166 ± 629 | 2512 ± 37 | 4515 | 110 ± 10 | 2.4 ± 0.2 |
| | | | | | | |

Results

Figure 1. Correlation between measured retention time (min) and iRT hydrophobicity using the timsTOF Pro and Fusion systems. Data were collected using selfpack columns over a 75 min gradient. (0-71 min, 4-25% B; 71-73.5 min, 25-90%B; 73.5-75 min, 90-2%B). Data were expressed as mean ± SD (n=3). Raw data were processed by Peaks Studio 8.5.

[2] Matched peptide list was obtained from Peaks Studio 8.5. [3] Peptide search within the defined RT and m/z windows was achieved using a custom program written in R. The defined RT and m/z windows were 1.2 s and 0.015, respectively.

[4] Peptides that were uniquely found in each scheme was searched against the matched peptide lists using a custom program written in R. Peptides in common were found by the same program. [5] Percentage = 100*Peptide in common/peptide only found by timsTOF or Fusion (Table 1) or peptide only found by lonopticks or Selfpack (Table 2).

gun proteomics.

This work was supported by National Science Foudation Grant (DMR-1508722) and the Technology Gap Fund (University of Houston).



Summary

The timsTOF Pro identified 30-50% more peptide sequences compared to Fusion.

Both timsTOF Pro and Fusion exhibited good LC performance.

The fast scanning rate of the timsTOF Pro may be the major contributor for its superb performance in peptide identification.

Ion mobility of the timsTOF Pro may contribute to its superb performance in peptide identification, however, it's unlikely to be the major contributor.

Conclusion

The PASEF method utilized in the timsTOF Pro represents a robust approach for shot-

Acknowledgment

