

Mouse fatty liver proteomes in dda-PASEF and dia-PASEF from low input using the timsTOF SCP

Elizabeth Gordon¹, Matthew Willetts¹, Diego Assis¹, Liwen Zhang², Yanqiao Zhang³

¹Bruker Daltonics Inc., 40 Manning Road., Billerica, MA 01821, USA

²Proteomic Shared Resources, CCC Mass Spectrometry and Proteomics Facility, 246 Biomedical Research Tower, Ohio State University

³Northeast Ohio Medical University (NEOMED) 4209 St, OH-44, Rootstown, OH 44272

Introduction

Liver diseases are a global health issue and especially pronounced with a western diet, including a continually increasing disease burden from non-alcoholic fatty liver disease. There is both a need for biomarkers and ability to test potential drugs for this condition and thus a need for translational mouse studies. The mice in this study were fed a Western diet for 16 weeks and their livers were collected and flash frozen. The liver samples were from mice over-expressing the protein of interest and on the high fat diet or a control group of matched wild type were used for proteomic analysis with liquid chromatography and tandem mass spectrometry analysis.

Methods

Mouse livers were harvested and frozen with liquid nitrogen. Frozen cells were thawed, reduced, and alkylated. 10ng of the protein was digested with trypsin and cleaned up with STRAP purification. Chromatography was performed on a nanoElute (Bruker Daltonics) using an Aurora nano column (25 cm x 75 µm ID, C18 - IonOpticks, Australia) at 400 nl/min with a 90 min gradient. LC-TIMS MS/MS data were acquired on a timsTOF SCP instrument operated in dda-PASEF and dia-PASEF modes. DDA data were analyzed using PEAKS OnLine (Bioinformatics Solutions), and dia-PASEF and raw data were analyzed with Spectronaut 17 (Biognosys) at a 1% FDR.



Fig. 1. A) Bruker NanoElute and timsTOF SCP B) dia-PASEF windows

Results




Fig. 2 Total Ion Chromatogram of the MS2 of replicate runs

| Sample Name | # Precursor | #PSM | #Peptides | # Protein Groups | # PSM / # Precursor |
|-----------------|-------------|---------|-----------|------------------|---------------------|
| All (8 samples) | 1,871,262 | 718,290 | 92,442 | 5,704 | 38% |
| WT 1 | 218,285 | 84,319 | 42,505 | 3,779 | 39% |
| WT 2 | 221,475 | 84,998 | 42,856 | 3,804 | 38% |
| WT 3 | 230,860 | 86,599 | 43,620 | 3,859 | 38% |
| WT 4 | 185,550 | 69,441 | 36,430 | 3,189 | 37% |
| Average | 214,043 | 81,339 | 41,353 | 3,658 | 38.00% |
| FL 1 | 241,723 | 98,997 | 49,228 | 4,113 | 41% |
| FL 2 | 265,269 | 94,947 | 48,143 | 4,148 | 36% |
| FL 3 | 239,567 | 98,859 | 49,118 | 4,047 | 41% |
| FL 4 | 268,533 | 100,130 | 50,561 | 4,271 | 37% |
| Average | 253,773 | 98,233 | 49,263 | 4,145 | 38.75% |

Table 2: PEAKS Results from 10ng of sample on column




Fig. 3 dia-PASEF A) Ion cloud coverage B) base peak chromatogram in dark blue over the blank run in green.

| Parameter | Value | Parameter | Value |
|------------------------|----------------------|-----------------------------|----------------|
| Precursors | 75,303 of 75,313 | Average XIC Width | 1.47 min |
| Modified Peptides | 59,014 of 59,017 | Median Absolute Delta iRT | 0.19 |
| Peptides | 55,113 of 55,115 | Median Absolute Delta RT | 0.08 |
| Proteotypic Peptides | 11,766 | MS1 Average Delta (ppm) | 2.32 |
| Protein Groups | 5,996 (5,996) | MS2 Average Delta (ppm) | 2.04 |
| Proteins | 12,586 (12,543) | MS1 Average Tolerance (ppm) | 6.5 |
| Single Hits | 1,042 (1,159) | MS2 Average Tolerance (ppm) | 6.5 |
| Fragments | 450,309 of 1,701,688 | Average IM Width | 0.05 |
| Library Source | Spectronaut | | |
| Software Version | 17.4.230317.55965 | | |
| Merged | False | | |
| Search Engine Platform | Pulsar | | |
| Digest Type | Specific | Peak Capacity | 544.72 (NaN) |
| Missed Cleavage | 3 | Median FWHM | 0.092026 (NaN) |
| Min Peptide Length | 7 | Median Peak Width | 0.155975 (NaN) |
| Max Peptide Length | 52 | Data Points per Peak (MS1) | 10 (NaN) |
| Toggle N-terminal M | True | Data Points per Peak (MS2) | 10 (NaN) |
| Digest Rule | Trypsin/P | Gradient Length | 85 |
| Protein Databases | | | |

Table 3: Spectral library and characteristics




Fig. 4. dia-PASEF A) protein groups per condition B) precursors per condition




Fig. 5. dia-PASEF volcano plot highlighting genes that are known in fatty liver early disease




Fig. 6. DIA heat map with independent clustering separates high fat diet wild type control from wild type control

Summary

dda-PASEF Results

Control samples had good technical reproducibility as well as consistency of detected protein groups between biological replicates.

Wild Type directly loaded at 10ng on column gave an average of 3,365 protein groups and 41k peptides. The fatty liver that was 16 weeks on a western diet 10ng averaged 4,415 protein groups at an FDR of 1% and 49k peptides.

dia-PASEF Library

From a small sample size [n=4] biological with 2 technical replicates a spectral library representing 5,996 Protein groups and 55k peptides and 75 k precursors was generated.

dia-PASEF Results

Results from just 1ng on column averaged 4,500 protein groups in the wild type samples, with 40k peptides, 45k sequences, and over 5,500 protein groups in the fatty liver with western diet, with 50k peptides, and 60k sequences.

Conclusion

TIMS-PASEF in DDA mode was able to identify over five thousand protein groups from 10ng of mouse liver and in DIA over four thousand protein groups were identified from just 1ng of protein. This methodology allows deep proteome coverage from very small amounts of tissue and enables isolation of proteomes from specific regions of tissue vs whole tissue.

timsTOF SCP Technology