

**Optimal Sample Preparation for Large-Area Transmission Kikuchi Diffraction Analysis** 

Dr. Laurie Palasse Bruker Nano Analytics, Berlin (DE)

Dr. Pawel NOWAKOWSKI E. A. Fischione Instruments, Inc., Pittsburgh (US)

#### **Presenters**

Dr. Laurie Palasse

*Global Application Manager EMA* Bruker Nano Analytics, Berlin (DE)

- Dr. Pawel NOWAKOWSKI
- Sr. Applications Scientist
- E. A. Fischione Instruments, Inc., Pittsburgh (US)









01 What is on-axis TKD? A review vs. EBSD vs. conventional (off-axis) TKD in SEM

# 02 On-Axis TKD – Capabilities and Benefits

Optimal Sample Preparation for Large-Area TKD

# Q&A / Discussion



# 01 What is on-axis TKD?

Innovation with Integrity

### **Spatial resolution** EBSD vs. conventional TKD vs. On-Axis TKD

- **Physical** Spatial Resolution (PSR) is set by the interaction volume
- Effective Spatial Resolution (ESR) is given by the SW ability to correctly index patterns produced by multiple crystals

Effective spatial resolution values:

- EBSD at 20kV: down to ~50nm (low probe current)
- Low-kV EBSD: down to ~20nm
- Conventional TKD : down to ~8nm
- On-Axis TKD: down to ~2nm



5





### **Experimental setup** EBSD vs. conventional TKD vs. On-Axis TKD



Conventional (Off-Axis) TKD





- EBSD and standard TKD use the same hardware & software
- Non-optimum sample-detector geometry for TKD ⇔ weak signal



- Additional hardware OPTIMUS 2
- Provides optimum sample-detector geometry for TKD ⇔ strong signal

"Orientation mapping by transmission-SEM with an on-axis detector" J.J. Fundenberger et all, Ultramicroscopy, 161, 17–22, 2016. "A systematic comparison of on-axis and off-axis transmission Kikuchi diffraction" F. Niessen et all, Ultramicroscopy, 186, 158-170, 2018.

#### **On-Axis TKD in SEM** DF imaging with OPTIMUS



Off-axis TKD



- 2.6nA and 10ms/point
- 43:05 min





vs. 1.75nA and 3ms/point

vs. 14:04min

Lower probe currents and faster data acquisition – less prone to beam instability



We all want to acquire orientation maps with:

- 1. Best spatial resolution to resolve even the finest crystals/features
- 2. Fastest speed possible to use lab's resources efficiently
- 3. Highest indexing rate, i.e. reliable data to help us get a realistic understanding of sample's properties

All the above is possible but it "costs" signal:

- 1. High spatial resolution requires low probe diameter (low current) => to keep the pattern quality constant we need signal
- 2. High speed  $\Leftrightarrow$  short exposure times => to keep the pattern quality constant we need signal
- 3. High indexing rates require high quality patterns which can only be obtained by having a **strong signal** (SNR)

### **On-Axis TKD** Best spatial resolution, speed and data quality

- Important parameters:
- EHT: 30kV
- Probe current: 2nA
- Step size: 3nm
- Mapping speed: 320fps (3ms/point)
- Zero solutions: 11.5%
- Annealing twin: <6 nm wide
- No data cleaning!

20nm Au film on 5nm Si<sub>3</sub>N<sub>4</sub> membrane







## **On-Axis TKD** Best spatial resolution, speed and data quality

- Important parameters:
- EHT: 30kV
- Probe current: 2nA
- Step size: 3nm
- Mapping speed: 320fps (3ms/point)
- Zero solutions: 11.5%
- Annealing twin: <6 nm wide
- No data cleaning!

#### 20nm Au film on 5nm $Si_3N_4$ membrane





### **On-Axis TKD** Best spatial resolution, speed and data quality

- Important parameters:
- EHT: 30kV
- Probe current: 2nA
- Step size: 3nm
- Mapping speed: 320fps (3ms/point)
- Zero solutions: 11.5%
- Annealing twin: <6 nm wide
- No data cleaning!

BRUKER

20nm Au film on 5nm Si<sub>3</sub>N<sub>4</sub> membrane





# 02

# on-axis TKD – capabilities and benefits

Innovation with Integrity



#### **OPTIMUS 2 detector head for On-Axis TKD** Features & benefits

© 2023 Bruker



- **OPTIMUS-Vue screen** BF-like imaging in mapping position (more details later)
- e- beam friendly materials *improved spatial resolution and minimized drift*
- Additional layer in screen structure *minimized beam interference*
- Redesigned screen frame easier/safer to use



### **High resolution DF imaging with OPTIMUS 2** 15 nm Au film deposited on 10 nm Si<sub>3</sub>N<sub>4</sub> membrane



High resolution DF image acquired at 30kV with 0.8nA probe current and a **pixel size of 1 nm** 

Image is courtesy of Hong Zhang from Eurofins in Santa Clara CA, USA



#### DF & BF imaging with OPTIMUS 2 Ruthenium film Appl



Applications & benefits:

- Qualitative characterization of microstructures
- Finding area/features of interest
- Refining beam focus and astigmatism
- Ideal for drift correction
- Essential for three new SW features:
  - ESPRIT FIL-TKD
  - ESPRIT TRM
  - ESPRIT MaxYield







# DF & BF imaging with OPTIMUS 2 color-coded dark field ARGUS<sup>™</sup> imaging







#### ESPRIT MaxYield Productivity boost



Argus MAG: 73.2kx HV: 30 kV WD: 4.3 mm Px: 5 nm

ARGUS image of nanoparticles on C-lacey support film



MaxYield

Binarized image to be used as a mask

ESPRIT MaxYield => key details / benefits:

- Acquire and binarize ARGUS/SEM images
- Use such images as masks to map sparse samples, e.g. nanoparticles, nanorods, nanotubes, etc.
- Acquire data only from the area of interest:
  - Productivity boost
  - Reduced drift induced artifacts

#### ESPRIT MaxYield Productivity boost



 Use masks to map sparse samples, e.g. nanoparticles, nanorods, nanotubes, etc.

MAG: 25.0kx HV: 30 kV WD: 5.0 mm Px: 5 nm

Acquire data only from the area of interest



#### OPTIMAL SAMPLE PREPARATION FOR LARGE-AREA TRANSMISSION KIKUCHI DIFFRACTI ESPRIT TRM Time Resolved Measurements

- Automatic and repetitive acquisition of time resolved images & maps
  - close to real time visualization of samples during in-situ experiments
- Applicable to EDS/EBSD/TKD on SEM and EDS on TEM
- Works on same location or multiple user defined locations

See: "Elevated temperature transmission Kikuchi diffraction in the SEM" Fanta et all, Materials Characterization, Vol ´.139, May 2018, Pages 452-462



## Integration of multiple instruments OPTIMUS 2 and PicoIndenter PI 89 with Push-to-Pull device



Sample and detectors setup for Time Resolved Measurements during in-situ tensile testing of nanomaterials



Time Resolved Measurements during in-situ experiments:

- Imaging at **125,000 pixels/second** during in-situ tensile testing with Push-to-Pull device on PI 89
- Time resolved mapping at up to 600 points/second quantify every microstructural change in the sample
- Applicable for in-situ heating and electrical biasing experiments

### ESPRIT FIL-TKD Enabling nano-scale TKD mapping







- Immersion mode ⇔ Strong Magnetic Field inside SEM chamber
  - Electrons are constrained within a narrow space around the optical axis

of SEM – limits their spread laterally

- Electron trajectories are affected distorted patterns
- OPTIMUS 2 and ESPRIT FIL-TKD (EP000003835768B1)





FIL-TKD

#### Goal: Achieve the best spatial resolution possible!

### **ESPRIT FIL-TKD** Enabling nano-scale TKD mapping

- Available for certain SEMs using Full Immersion Lens technology
- Fully integrated in the pattern analysis process of ESPRIT 2
- Works for all applications except residual strain analysis
- Results show here were acquired using a +10 years old NovaNano SEM
- **DISCLAIMER**: Resolution difference between the two functioning modes is very likely not as dramatic on latest gen. e- columns

Data courtesy of Alice Da Silva Fanta from DTU Nanolab in Copenhagen, Denmark

No magnetic field  $\Leftrightarrow$  analytical mode

22









# **Application examples**

Innovation with Integrity

## Integration of multiple instruments OPTIMUS 2 and FlatQUAD EDS



Sample and detectors setup for simultaneous on-axis TKD and EDS measurements on nanomaterials



Combined on-axis TKD and EDS => benefits

- Fastest simultaneous TKD&EDS (enabled by ultra-high solid angle of FlatQuad 1sr)
- X-Rella mask for blocking most stray X-Rays from reaching the FlatQUAD
- Applications:
  - Phase identification & distribution (same as for bulk samples)
  - Characterization of complex samples with crystalline, semicrystalline and amorphous regions
  - Discrimination of phases creating similar Kikuchi patterns using EDS information



#### Cu(InGa)Se<sub>2</sub> solar cell FIB lamella



Sample courtesy: M. Raghuwanshi, RWTH Aachen University, Aachen, Germany



#### Cu(InGa)Se<sub>2</sub> solar cell FIB lamella

Pattern quality map









#### Cu(InGa)Se<sub>2</sub> solar cell FIB lamella

#### IPF-Z map





#### Mapping parameters:

- EHT: 30kV
- Probe current: 1.6 nA
- Step size: 4nm
- Acq. speed: 330 fps
- Map size: 1.93 million pixels
- Zero sol.: 19.3% (whole map)
- Zero sol.: 1.18% on CIGS

Indexing rate in CIGS: 98.82%

### **On-Axis TKD in SEM** Orientation mapping of Mo thin film – plane view

- Map size: +1Mpixels
- Acquisition speed: **626 fps**
- Measurement time: 28:52 minutes
- Pixel size: 5 nm
- Field of view: 6x4.4 μm<sup>2</sup>
- Pattern Quality Map







### **On-Axis TKD in SEM** Orientation mapping of Mo thin film – plane view

- Map size: +1Mpixels
- Acquisition speed: **626 fps**
- Measurement time: 28:52 minutes
- Pixel size: 5 nm
- Field of view: 6x4.4 μm<sup>2</sup>
- IPF map (orientation distribution)
- Map shows raw indexing results
- High indexing rate no misindexing





Molybdenum

[111]

[001]

[101]



### **On-Axis TKD in SEM** Orientation mapping of Mo thin film – plane view

- Map size: +1Mpixels
- Acquisition speed: **626 fps**
- Measurement time: 28:52 minutes
- Pixel size: 5 nm
- Field of view: 6x4.4 μm<sup>2</sup>
- 9284 grains! Average size 70 nm
- Map shows raw indexing results



#### Summary



- Best spatial resolution period Don't miss any detail in your samples
  - Down to 1.5nm effective spatial resolution\*
  - *New* ESPRIT FIL-TKD for orientation mapping using Full Immersion Lens mode on certain e- columns
  - *New* BF-like imaging in mapping position for ideal beam focus optimization
  - New Use of high-performance materials for minimized interference with electron beam
- Understand your samples better using the best data quality & integrity
  - Maximized SNR during imaging & mapping due to unique on-axis geometry
  - Ideal also for low kV and low probe current applications
  - *New* Improved drift correction due to BF-like imaging in mapping position
  - *New* Sample holders made of high-performance alloys to minimize thermal expansion induced drift

\* Dependent on SEM type, vacuum quality & room environment, e.g. vibrations, acoustics, etc.



# 03

# **Optimal Sample Preparation for Large-Area TKD**

Innovation with Integrity



# 04 Questions & Answers

please type your questions in the Q&A window and press SEND

Innovation with Integrity



## on-axis TKD in SEM Effect of plasma cleaning on data quality









20nm Au thin film on  $Si_3N_4$  membrane

- Acquisition at 4ms
- Step size 1.6nm
- e<sup>-</sup>Flash<sup>FS</sup> with OPTIMUS TKD

Beam parameters:
 EHT: 30kV / probe current: 1.75nA

Indexing rate:46.5% vs. 87.6%

[001]

[101]

[111]

Raw data – no cleaning no processing!

Sample is courtesy of Alice Da Silva Fanta from CEN, DTU, Denmark

## **Off-axis TKD vs on-axis TKD**



#### Stainless steel TEM foil (10 nm step, 6nA)



Pattern quality+IPFYMap MAG: 15,0kx HV: 29 kV WD: 4,8 mm

Px: 10.0 nm

#### < OFF-AXIS TKD >

#### < ON-AXIS TKD >

Sharp grain/twin boundaries Almost no misindexing Higher indexing rate Faster speed>>less beam instabilities

better data quality in less time!

Au thin film (2nm step, 1,6nA)





### **Productivity boost** TKD sample holders

TKD holder (EP000002824448B1) in various configurations:

- Designed for increased productivity & data integrity
- Various concepts fitting most SEM models and airlocks
- Flexible design to allow quick upgrades based on customer requirements





Adapter for Zeiss, TFS\* and TESCAN SEMs

\* Only for models with standard DC stages



Adapter for TFS SEMs with 6in Piezo stages

X-Rella mask mounted on TKD holder for simultaneous EDS & TKD mapping



Jeol adapter



Innovation with Integrity

Optimal sample preparation for large-area transmission Kikuchi diffraction analysis



## **Introduction** FIB sample preparation



FIB milling is a versatile, powerful technique for electron-transparent sample preparation

#### **Disadvantages**

- Requires significant polishing of the sample surface to remove:
  - > Amorphization
  - > Ga implantation
  - > Structural damage
- Control of sample thickness can be complex and time consuming



## **Introduction** FIB sample preparation disadvantages



Excessive FIB sample thinning can lead to lamella banding



Non-uniform sample thickness. Sample very thin on front side and thick on back side.

#### Neither situation is ideal for large-area TKD acquisition



excellence...magnified

## Introduction FIB-induced structural damage



FIB-induced amorphization and ion implantation.\*

\*Bonifacio, C., ISTFA Proceedings, ASM (2018)



FIB milling introduced enough atom displacement and vacancies to cause dynamic strain-induced phase transformation of austenite to martensite. Ga ions can penetrate austenite steel up to 300 nm at 30 keV, depending on grain crystal orientation.<sup>†</sup>

<sup>†</sup>R. Prasath Babu, et al. *Acta Materialia* 120 (2016)



excellence...magnified

## Introduction FIB-induced structural damage



Sample damage after 1 keV Ga FIB





In the area affected by the Ga ion beam, the EBSD pattern is of very poor quality (Kikuchi bands shift and overlap)

P. Nowakowski, et al. Microsc. Microanal 28 (2022)

To accurately characterize a material, a sample preparation technique that preserves the microstructure of the material in its native state is necessary.



## Model 1040 NanoMill<sup>®</sup> TEM specimen preparation system Post-FIB sample preparation system



Low-energy, inert, concentrated ion beam avoids sample damage

- Inert gas (Ar)
- Ultra-low ion energies (~50 eV)
- Small spot size (1 µm)
- Targeted ion beam positioning
- Imaging of ion-induced secondary electrons

Application of low ion beam energies and currents allow improvement of FIB-prepared specimens for electron microscopy applications



## NanoMill system







## NanoMill system Milling geometries



- Ion milling at low energy (50 eV to 2 keV) and low beam current
- Thinning of FIB lamella from the front side



## NanoMill system Milling geometries



- Thinning of FIB lamella from the back side
- Milling angles up to 30°



## NanoMill system



- Specimen visualization using induced secondary electrons
- Region of interest is precisely targeted by the user
- Micron-sized beam moves in a raster pattern



# NanoMill system



- Removal of implanted Ga ions
- Removal of FIB-induced amorphization
- Controllable sample thinning



## NanoMill system Large-area TKD sample preparation workflow

#### **FIB** sample preparation



Lamella after standard lift-out

Wedge-shaped FIB lamella milled at high angle from 7 to 2°

0.5 um



NanoMill system back side thinning at high angle from 30 to 15°

P. Nowakowski, et al. Microsc. Microanal 29 (2023)



excellence...magnified

## NanoMill system Large-area TKD sample preparation workflow



- Ion milling at:
  - > low energy (900 eV 300 eV) and
  - low beam current (100 pA)
- Lamella is thinned from the back side
- Region of interest is precisely targeted; user-defined milling box is 10 x 3 µm
- Specimen is tilted up 30°



# **Highly deformed Ni alloy** (cold rolled; total sheet thickness reduced by 95%) FIB tool vs. NanoMill system

# After 5 kV FIB milling standard preparation

STRUMENTS



#### After 500 eV high-angle back side NanoMilling



Large area can be easily and accurately measured by TKD

Indexing rate: 97%

#### Acquisition conditions

- 30 kV acceleration voltage
- 6.4 nA beam current
- 5 nm step size



excellence...magnified

## Highly deformed Ni alloy Entire lamella is milled uniformly

#### TKD IPF map



# Acquisition conditions

- 30 kV acceleration voltage
- 6.4 nA beam current
- 10 nm step size
- 97 % indexing rate

Nickel (New)

P. Nowakowski, et al. Microsc. Microanal 27 (2021)



# Highly deformed Ni alloy

Dark field electron contrast image 400 nm MAG: 120kx Px: 1 nm

#### TKD KAM map



Deformation defects can be observed and correlated with kernel average misorientation TKD measurements

P. Nowakowski, et al. Microsc. Microanal 27 (2021)



### Front side standard FIB 5 keV preparation

#### Sample thickness ~60 nm



- The analysis challenge is the exceptionally fine grain size – significantly less than 50 nm.
- Only grains > 100 nm are indexed.



#### **Acquisition conditions**

- 30 kV acceleration voltage
- 1.6 nA beam current
- 2 nm step size



Optimum sample thickness for TKD measurements





Structural damage by Ga ion implantation

Sample milled by FIB at 30 and 5 keV





INSTRUMENTS

Ga ion implantation / damage is removed after NanoMilling

#### Sample back side milled by NanoMill system at 500 eV



Thickness is not the only concern. Thin lamella must also be free from structural damage for nanostructural material characterization.

Optimal sample preparation for TKD analyses

#### Sample NanoMilled at 500 eV; final thickness ~40 nm



Very fine Ni structure is revealed after NanoMilling Acquisition conditions30 kV acceleration voltage

- 1.6 nA beam current
- 2 nm step size



Optimal sample preparation for TKD analyses

Sample NanoMilled at 500 eV; final thickness ~40 nm





#### Sample NanoMilled at 500 eV; final thickness ~40 nm

The NanoMill system enables large area, electron-transparent sample preparation for accurate TKD analyses



- Average grain size: 30 nm
- 10% Σ3 <111>60° twin boundary

#### **Acquisition conditions**

- 30 kV acceleration voltage
- 1.6 nA beam current
- 2 nm step size



#### Sample NanoMilled at 500 eV; final thickness ~40 nm

The NanoMill system enables large area, electron-transparent sample preparation for accurate TKD analyses



- Average grain size
  30 nm calculated
  over 5828 grains
- 70% Σ3 <111>60° twin boundary

P. Nowakowski, et al. *Microsc. Microanal 29 (2023)* 



## Conclusions

- A new approach for large area electron-transparent samples for TKD measurement was demonstrated
- Proposed workflow allows:
  - Significantly reduced FIB milling time
  - Removal of FIB-induced structural damage
  - > Control of sample thickness
  - Improved overall TKD signal and indexing rate

Accurate sample preparation makes TKD analyses easier and provides improved results



# E.A. Fischione Instruments, Inc.

**9003 Corporate Circle, Export, PA 1563 USA** Tel **+1 724.325.5444** Fax **+1 724.325.5443** www.fischione.com