

Pesticides Target Screening with GC-APCI coupled to high-resolution Q-TOF-MS

ASMS 2015, WP 062

Thomas Arthen-Engeland, Dorith Claes, Petra Decker, Oliver Raether, Carsten Baessmann

Bruker Daltonik GmbH, Bremen, Germany

Introduction

Pesticides are frequently found as contaminations in food and environmental matrices. They are analyzed by GC-MS or LC-MS. Both are complementary, "orthogonal" methods, comprising in total >1100 known pesticides. While they overlap in scope, each method alone covers exclusively a certain range of pesticides: GC-MS is more common for semi-volatile compounds, LC-MS is favorable for polar and thermo-labile pesticides.

Full-scan accurate mass screening with atmospheric pressure inlet LC-MS became very common in the last years. It is able to cover hundreds of target compounds in a single run and additionally enables the identification of unknowns and retrospective analysis. Target compounds are identified by their retention time, mass accuracy and isotope pattern. Reliability of identification is improved by using diagnostic ions generated by broad-band CID alternating with full scan data acquisition. Diagnostic ions are valuable in complex matrices as they support the differentiation of target analytes signals from the matrix background.

Here we describe for the first time the application of bbCID data acquisition for pesticides target screening by coupling a GC to an atmospheric pressure chemical ionization source (GC-APCI) and a high-res QTOF-MS.




Fig. 1 left: Base Peak (grey BPC) and Extracted Ion Chromatograms (EIC) of 48 pesticides in peach QueCHERS extract (50 ng/mL each): BPC shows a complex peak pattern of non-target matrix substances. Right: Excerpt of BPC (grey) and EICs at 13.2 min, EIC peaks with symmetric and sharp peak shapes with typical FWHM of 1.2 to 1.4 s. Peak shapes are well described when using scan rates faster > 8 Hz.

Methods

For this GC-APCI-MS screening study a mix of 60 representative pesticides were chosen with regard to their relevance to routine screening. The mix was diluted in dichloromethane to appropriate concentrations for the generation of calibration curves between 0.01 ng/mL (0.01 µg/µl) and 1000 ng/mL. Samples were prepared by spiking 10 µg and 50 µg of the pesticide mix into QueCHERS extracts of orange, peach and tomato. For all analyses 2 µl of each sample were injected into the GC.

GC-MS analysis was performed using a 450-GC with PAL Combi-xt Autoinjector coupled with a GC-APCI II source to an impact II Q-TOF mass spectrometer (all Bruker Daltonics). The GC was operated with a 30 m Rxi-5ms capillary column (0.25 mm ID, 0.25 µm film thickness), operated at 1.2 ml/min constant helium flow and a GC oven temperature program at 70°C (1 min) - 25°C/min - 180°C - 15°C/min - 300°C (8.1 min). Pulsed splitless injection was at 280°C (40 psi for 0.25 min, 1 min splitless time).

MS Data were acquired from 50 - 1000 m/z in alternating full-scan and broad-band CID acquisition mode at 8 spectra per second operated in the positive ionization mode. All spectra were calibrated using PFTBA as external calibration gas injected automatically into the APCI source at the beginning of each MS run. Data were evaluated using DataAnalysis and TASQ software (Bruker Daltonics) for target analysis and quantification (TASQ settings: ΔRT 0.7 min, EIC width ±3 mDa of M⁺, [M+H]⁺, [M+H+2]⁺ and fragment ions).




Fig. 2 Representative calibration curve of Indoxacarb: LOQ 5 µg/µl, linearity 5-1000 µg/µl, R² = 0.9979

Results

The purpose of this study was to evaluate the applicability of GC-APCI-QTOF-MS for pesticides target screening following the concepts known from LC-MS: mass spectra were acquired with full and alternating bbCID scan. (1) In a first step a fast GC program for the set of 60 pesticides was optimized. (2) 48 of the 60 pesticides were assigned by retention time and nominal masses of the ions (which is an acceptable result because of the varieties of pesticide properties as polarity, volatility, etc.).

Typical BPC plus EICs of 50 ng pesticides in 1 mL peach matrix sample are shown in fig.1. Most of the GC peaks are very sharp, e.g. 79% show Full Width at Half Peak Maximum (FWHM) below 2 s and 95% below 3 s FWHM. Therefore fast scan rates are very important to generate enough MS data points and describe peaks properly and reproducibly. Although some GC peaks are overlapping, they are separated by mass. An excerpt of EICs at RT 13.2 min is shown in fig. 1 right where MS data points are drawn for the [M+H]⁺ ion (322.1442 Da) of Pyriproxyfen.

(3) In the 3rd step master list files for TASQ import were created for all target pesticides containing retention time, nominal mass of quantification and diagnostic ions, type of ion, concentration of calibration standards. A summary of TASQ quantification results is shown in the table. Typical LODs are between 0.1 and 50 ng/mL for GC peaks with low response. Linear range of quantification typically was 2 orders of magnitude or higher between 1 - 1000 ng/mL, while some analytes with very good GC response show 0.5 - 500 ng/mL. Good R² linearity values with >0.99 are significant for most of the 48 pesticides.

Analyte	RT [min]	LOD [ng/mL]	Working range	R ² Linearity	Quant m/z [Da]	Ion type	Δm/z [ppm]	Δm/z [mDa]
Pyriproxyfen	13.20	0.1	0.5-500	0.9985	322.1438	MH	0.61	0.20
Quinoxiphen	12.08	0.5	0.5-500	0.9961	308.0040	MH	0.81	0.25
Azoxystrobin	16.47	0.5	1-500	0.9987	404.1241	MH	0.63	0.25
EPN	12.72	0.5	1-500	0.9979	324.0454	MH	0.55	0.18
Myclobutanil	11.18	0.5	1-500	0.9988	289.1215	MH	0.62	0.18
Trifloxystrobin	11.88	0.5	1-500	0.9994	408.1370	MH	0.55	0.22
Chlorpyrifos	9.89	0.5	1-1000	0.9938	351.9306	MH+2	0.61	0.22
Pendimethalin	10.30	0.5	1-1000	0.9974	212.0666	Fragment	0.41	0.09
Metaxyl	9.48	0.5	5-500	0.9970	280.1543	MH	0.54	0.15
Pririmicarb	9.00	0.5	5-500	0.9967	239.1503	MH	0.53	0.13
Triazophos	11.80	0.5	5-500	0.9985	314.0723	MH	0.69	0.22
Boscalid	14.68	1	5-500	0.9991	343.0399	MH	1.13	0.39
Cyhalothrin.lambda	13.24	1	5-500	0.9962	225.0289	Fragment	0.61	0.14
Cyprodinil	10.32	1	5-500	0.9965	226.1339	MH	0.97	0.22
Diazinon	8.68	1	5-500	0.9995	305.1083	MH	0.44	0.14
Fludioxonil	11.01	1	5-500	0.9938	248.0392	M	0.48	0.12
Propargite	22.28	1	5-500	0.9971	350.1546	M	0.54	0.22
Chlorpyrifos-methyl	9.34	1	5-1000	0.9936	323.8993	MH+2	0.52	0.17
Dimethomorph Peak 1	16.67	1	5-1000	0.9949	388.1310	MH	1.16	0.45
Dimethomorph Peak 2	27.09	1	5-1000	0.9974	388.1310	MH	1.25	0.48
Ethion	11.62	1	5-1000	0.9975	384.9949	MH	0.52	0.20
Indoxacarb	16.02	1	5-1000	0.9979	528.0780	MH	0.67	0.36
Penconazole	10.40	1	5-1000	0.9971	284.0716	MH	0.42	0.12
Phosmet	12.71	1	5-1000	0.9973	169.0393	Fragment	0.69	0.11
Profenophos	11.09	1	5-1000	0.9972	374.9402	MH+2	0.96	0.36
Tebuconazole	12.31	1	5-1000	0.9968	308.1524	MH	0.44	0.13
Tolylfluand	10.44	5	5-1000	0.9988	237.9555	Fragment	0.51	0.12
Kresoxim-methyl	11.18	5	5-500	0.9885	205.0812	Fragment	1.86	0.38
Etifopros	14.83	5	10-500	0.9993	359.2006	Fragment	0.44	0.16
Carbafuran	8.40	5	10-1000	0.9915	222.1125	MH	0.79	0.17
Difenoconazole Peak 1	15.56	5	10-1000	0.9958	406.0720	MH	0.57	0.23
Difenoconazole Peak 2	16.03	5	10-1000	0.9932	406.0720	MH	0.57	0.23
Fenhexamid	12.18	5	10-1000	0.9982	302.0709	MH	0.79	0.24
Prochloraz	14.06	5	10-1000	0.9959	376.0381	MH	0.41	0.16
Pyrimethanil	8.81	1	50-500	0.9972	200.1182	MH	0.69	0.14
Carbendazim	10.30	10	50-500	0.9883	191.0689	M	1.18	0.23
Chlorpropham	7.91	10	50-500	0.9955	215.0551	M	0.92	0.20
Bifenthrin	12.63	1	50-1000	0.9906	181.3012	Fragment	0.68	0.12
Azinphos-methyl	13.22	5	50-1000	0.9985	132.0444	Fragment	1.01	0.13
Dimethoate	8.39	5	50-1000	0.9973	230.0669	MH	0.62	0.14
Carbaryl	9.51	10	50-1000	0.9948	145.0648	Fragment	0.82	0.12
Cypermethrin I	14.54	10	50-1000	0.9971	416.0815	MH	1.07	0.44
Cypermethrin II	14.62	10	50-1000	0.9940	416.0815	MH	0.77	0.32
Cypermethrin III	14.69	10	50-1000	0.9927	416.0815	MH	0.73	0.30
Imazalil	11.07	10	50-1000	0.9958	297.0556	MH	0.54	0.16
Linuron	9.81	10	50-1000	0.9935	249.0192	MH	0.53	0.13
Thiadicoprid	15.70	10	50-1000	0.9969	253.0309	MH	0.33	0.08
Triadimenol I	10.55	10	50-1000	0.9965	296.1160	MH	0.76	0.21
Triadimenol II	10.65	10	50-1000	0.9982	296.1160	MH	0.67	0.20
Monocrotophos	8.03	50	50-1000	0.9957	159.0260	Fragment	0.33	0.06
Qmethoate	7.59	50	50-1000	0.9675	214.0297	MH	0.52	0.11
Thiabendazole	10.80	50	50-1000	0.9928	202.0433	MH	0.64	0.13
Acetamiprid	12.67	<50	50-1000	0.9977	126.0105	Fragment	0.50	0.06

Table 1: Analytical results for 48 Pesticides sorted according to LOQ

OR PICTURE OF TASQ ??




Fig.3. Comparison of mass accuracy over all EICs (±3 mDa): a) only for calibrant runs, b) for calibrant + matrix sample runs. For both groups nearly identical results are observed below 3 ppm (98% and 96%), even if matrix sample runs are included. For latter group a slightly increased % of mass deviations <8.5 ppm are observed (4% over 1% for calibrant runs).

Table 1 also shows mass accuracy data for each calibrant averaged over all concentrations: low mass accuracy for all calibrants is below 0.4 mDa (w/o Dimethomorph), 76% of all calibrants are below 1 ppm. Average mass accuracy over all calibration runs is 0.78 ppm and 0.88 ppm when the 9 matrix sample runs were included. Six pesticides show accuracies at ca. 1 ppm related to low response or noisy peaks, esp. for low concentration calibrant samples. Carbendiazim, Kresoxim and Boscalid are co-eluting with other pesticides. Following fig. 3 we observe only minor matrix effects with GC-APCI pesticides target screening, even if samples include substantial amounts of matrix compounds. In future we will extend number of pesticides for GC-APCI target screening.

Conclusions

- GC-APCI-QTOF-MS was successfully applied for pesticides target screening acquiring fast (8 Hz) alternating full scan and bbCID MS data
- Linearities are R² >0.99 between 1-1000 ng/mL (LOD down to <1 ng/mL)
- All MS acquired with automated external PFTBA mass calibration during each GC/MS run
- Average mass accuracies of analytes were <0.8 ppm and for each of the 48 analytes in all 36 analyses typically <1 ppm or < 0.4 mDa

GC-APCI hr-QTOF-MS