

## X-RAY MICROCOMPUTED TOMOGRAPHY

# X4 POSEIDON – Mouse bone osteocytes

### Innovation with Integrity

Submicron imaging of mineralised bone opens the field of ultra-structural analysis of bone including the distribution of osteocyte lacunae (spaces) and blood vessel canals. Both of these are highly informative of the metabolic state and signalling responses of bone in a rodent research model.

In the image above, submicron pixel scanning of the midshaft mouse femur diaphysis clearly resolves both blood vessel canals and osteocyte lacunae. 3D morphometric analysis of both extends the knowledge of the status of the bone, illustrating physiological and disease processes.

### References

Javaheri et al. (2015) *Bone* 81: 277-291.


Ahmed et al. (2023) *Annals of Anatomy-Anatomischer Anzeiger* 250: 152142.

Azari et al. (2024) *Calcified Tissue International* 115(1): 78-84

### Scan parameters


- Detector: 16MP sCMOS
- Voxel size: 550 nm (1x1 binning mode)
- Source: Transmission type
- Source power: 50 keV, 6 W
- X-ray filter: 0.25 mm Al
- Rotation step: 0.12° over 360°
- Phase retrieval (Paganin) reconstruction,  $\beta/\delta$  ratio 100






**Figure 1:** Mouse femur diaphyseal cortical bone, imaged at 550 nm voxel size, with absorption reconstruction. Osteocyte lacunae and blood vessel canals are clearly resolved.

The high resolution and low noise active pixel performance of the sCMOS camera in the X4 POSEIDON, as well as its large 16Mp format, allows ultrastructural analysis of a whole mouse femur with no cutting required. Osteocyte lacunae and blood vessel canals are clearly resolved (figure 1, 2). Phase retrieval (Paganin) allows spatial mapping of mineralization density (figure 3) revealing the history of remodeling events, possibly indicating the relative age of bone regions based on nanodensitometry.



**Figure 2:** Mouse femur diaphyseal cortical bone, imaged at 550 nm voxel size, allowed segmentation of osteocyte lacunae (green) and blood vessel canals (red)



**Figure 3:** Phase retrieval processing of the cortical bone scan reveals mineral density inhomogeneities

Bruker BioSpin is continually improving its products and reserves the right to change specifications without notice. BS-100034 01/2025 Bruker BioSpin.

**Bruker BioSpin**  
info@bruker.com

**Worldwide Offices**  
bruker.com/

**X4 POSEIDON**  
bruker.com/

bruker.com

