

Profiling and quantitation of bile acids in human biofluids by LC-TIMS-MS

Results

Aiko Barsch¹, Surendar Tadi², Xuejun Peng³, Narayanaganes Balasubramanian², Elena Checkmeneva⁴, Katie E. Chappell⁴, Niels Goedecke¹, Matthew R. Lewis¹

¹Bruker Daltonics GmbH & Co. KG, Bremen, Germany; ²Bruker Scientific, LLC, Billerica, MA; ³Bruker Scientific, LLC, San Jose, CA; ⁴National Phenome Centre & Section of Analytical Chemistry, Division of Systems Medicine, Department of Metabolism, Digestion & Reproduction, Imperial College London, London, United Kingdom

Introduction

The physiological relevance of bile acids (BA) resulted in extensive research and the recent discovery of novel BA conjugates shown to be related to human health and the microbiome [1]. BAs are derived from cholesterol and their chemical complexity arises from the occurrence of hydroxyl or carbonyl groups in different stereochemistries. This renders bile acids analysis challenging as several isobaric and isomeric BAs can't be separated by RP-LC.

We optimized a RP-LC-TIMS-MS method providing trapped ion mobility separation (TIMS) capabilities orthogonal to RP-LC. This allows separating co-eluting unconjugated BAs, among others. This is crucial as they lack characteristic MS/MS fragments, rendering their annotation using common, low-specificity SIM LC-QQQs methods ambiguous.

The established method provides quantitation and profiling capabilities for BAs extracted from complex human biofluids.

Methods

Bile acid standards were obtained from Cambridge Isotope Labs (Tewksbury, USA), Steraloids (Newport, USA) and Medical Isotopes (Pelham, USA). NIST SRM 1950 human reference plasma was obtained from Sigma-Aldrich (Taufkirchen, Germany). C8-RP chromatographic conditions and plasma extraction were conducted following a protocol from Sarafian et al. 2015 [2]. LC-trapped ion mobility (TIMS) MS analyses were performed using an Elute UHPLC connected to a timsTOF Pro 2 system equipped with a VIP-HESI source operated in negative ion mode (Bruker Daltonics, Bremen, Germany). The resulting four-dimensional data (*m/z*, RT, mobility, and MS/MS) was processed using TASQ[®] 2023 and MetaboScape[®] 2023 (Bruker Daltonics, Bremen, Germany) for targeted quantitation and non-targeted profiling, respectively. In both software workflow solutions, raw data was automatically recalibrated for mass and mobility.

1) LC-TIMS-MS method for profiling of bile acids in biological samples




Fig. 1 LC-TIMS-MS analysis of >70 bile acid standards

In total 71 bile acids can be separated by the LC-TIMS-MS method and hence confidently annotated utilizing the established assay. The optimized LC gradient and HRMS used in our LC-TIMS-MS analysis allowed many bile acid species to be separated by either LC or MS alone (see Figure 1 A). Even in this simplified system of 71 synthetic reference materials, co-elution of several isobaric bile acid pairs was observed. Lithocholic acid and Allolithocholic acid (see Figure 1 B), for example, were further separable using TIMS of their [M+Acetate]⁻ ion species as highlighted in Figure 1 C. This clearly demonstrates the additional benefit of TIMS for increasing the specificity of measurement and enhancing the confidence in annotation for these bile acids which can not be achieved by common LC-QQQ methods.

2) Quantitation and comparison to reference values




Fig. 2 Taurocholic acid quantitation

Figure 2 shows the calibration curve for Taurocholic acid with a linear dynamic range of 3.7 orders of magnitude (1nM – 5000 nM) an *R*² value of 0.9987 and residuals below 20%. This demonstrates the quantitative capabilities of the established assay. We extracted bile acids from human reference plasma (SRM 1950) and compared the determined concentrations to the MEDM Locations established by Bowden *et al.* [3]. Table 1 lists the bile acid that were detected within the dynamic range of the LC-TIMS-MS method. All determined concentrations are within the range of the published standard uncertainty locations.

Table 1: Bile acid concentrations determined by LC-TIMS-MS are within the reference concentration range

Name	MEDM Location +/- Standard uncertainty [nmol/ml]; Bowden L, <i>et al.</i> 2017 [3]	LC-TIMS-MS determined concentration [nmol/ml] n=2
Chenodeoxycholic acid	0.30 +/- 0.11	0.28
Cholic acid	0.12 +/- 0.034	0.14
Deoxycholic acid	0.35 +/- 0.083	0.3
Glycochenodeoxycholic acid	1.1 +/- 0.18	1.03
Glycodeoxycholic acid	0.43 +/- 0.069	0.38
Glycoursodeoxycholic acid	0.15 +/- 0.024	0.17
Glycocholic acid	0.24 +/- 0.069	0.18
Taurocholic acid	0.026 +/- 0.0056	0.022
Ursodeoxycholic acid	0.11 +/- 0.024	0.11

3) Intertechnology and interlaboratory CCS value comparability

We measured bile acid reference mixtures using two LC-TIMS-MS setups. The same MS method was applied, but different LC gradients and column chemistries were used (a 15 min gradient in Bremen, Germany and a 12.5 min gradient in Billerica, US). Table 2 shows the CCS values deviated on average by only 0.1% permitting confident interlaboratory annotations.

The average deviation of 0.3% in CCS values determined by TIMS vs. Drift Tube (DT) technology underlines the ability to match CCS values acquired by TIMS to reference values in public repositories like the Unified CCS Compendium [4].

Table 2–TIMSCCS values are comparable between laboratories and to reference DTCCS values

Name	CCS (Å ²) [M-H] ⁻ Reference			% CCS deviation	
	Poland C., <i>et al.</i> 2020 [5]	TIMS Bremen [n=6]	TIMS Billerica [n=6]	TIMS Bremen vs. Poland <i>et al.</i> [5]	TIMS Bremen vs. Billerica
Glycolithocholic acid	199.5 +/- 0.2	198.8 +/- 0.1	199.1 +/- 0.2	0.3	-0.1
Glycodeoxycholic acid	199.9 +/- 0.2	199.2 +/- 0.1	199.4 +/- 0.1	0.4	-0.1
Glycoursodeoxycholic acid	201.1 +/- 0.1	200.5 +/- 0.1	200.6 +/- 0.1	0.3	0.0
Glycocholic acid	202.2 +/- 0.1	201.6 +/- 0.0	201.8 +/- 0.1	0.3	-0.1
Taurolithocholic acid	206.4 +/- 0.1	206.0 +/- 0.1	206.2 +/- 0.2	0.2	-0.1
Taurochenodeoxycholic acid	207.2 +/- 0.2	206.7 +/- 0.1	206.9 +/- 0.1	0.2	-0.1
Tauroursodeoxycholic acid	207.6 +/- 0.1	207.2 +/- 0.1	207.4 +/- 0.2	0.2	-0.1
Taurocholic acid	207.6 +/- 0.2	207.0 +/- 0.1	207.3 +/- 0.1	0.3	-0.2
Average				0.3	0.1

4) Bile acid profiling in SRM 1950

Index	m/z	mob.	Δm/z [mDa]	RT [min]	ΔRT	mSigma	Mod. VLD	CCS (Å ²)	ACC (%)	Name	Molecular For.	Annotations	AQ
1	389.2557	-0.164	10.27	0.02	24.5	0.961	199.2	0.1	5-beta-Cholanic Acid 12-alpha-ol-3-one	C ₂₇ H ₄₆ O ₆	1950	100	
2	391.2503	-0.353	9.35	0.03	2.8	0.972	201.5	-0.0	5-beta-Cholanic Acid 3-beta, 12-alpha-diol	C ₂₇ H ₄₈ O ₇	1950	100	
3	448.30038	-0.407	5.22	0.04	8.2	0.971	200.4	-0.0	Glycoursodeoxycholic Acid	C ₂₇ H ₄₆ NO ₇	1950	100	
4	464.30126	-0.502	5.75	0.06	10.8	0.977	201.5	-0.0	Glycocholic Acid	C ₂₇ H ₄₆ NO ₇	1950	100	
5	498.29943	0.526	6.97	0.03	10.3	0.999	205.6	1.0	Taurodeoxycholic Acid	C ₂₇ H ₄₆ NO ₈	1950	100	
6	391.28482	-0.559	10.66	0.00	8.0	0.972	201.5	0.1	Deoxycholic Acid	C ₂₇ H ₄₆ O ₆	1950	100	
7	391.28482	-0.562	7.70	0.02	3.3	1.000	207.2	-0.8	Murocholic Acid	C ₂₇ H ₄₆ O ₆	1950	100	
8	391.28477	-0.612	10.55	0.01	17.9	1.001	207.6	-0.2	Chenodeoxycholic Acid	C ₂₇ H ₄₆ O ₆	1950	100	
9	448.30042	-0.244	8.09	0.01	16.4	0.964	199.1	-0.1	Glycodeoxycholic acid	C ₂₇ H ₄₆ NO ₇	1950	100	
10	432.31161	-0.323	9.87	0.04	6.5	0.961	198.5	-0.2	Glyolithocholic acid	C ₂₇ H ₄₆ NO ₇	1950	100	
11	498.29993	0.399	6.55	0.02	8.7	1.004	206.5	-0.1	Taurchenodeoxycholic Acid	C ₂₇ H ₄₆ NO ₈	1950	100	
12	407.27967	-0.397	9.79	0.03	14.8	0.979	202.7	-0.2	Cholic Acid	C ₂₇ H ₄₆ O ₆	1950	100	
13	514.28392	-0.480	4.82	0.01	16.7	1.007	207.1	0.1	Taurchenolic Acid	C ₂₇ H ₄₆ NO ₈	1950	100	
14	528.26340	-0.257	2.83	0.05	33.5	1.049	215.7	0.0	Glycoursodeoxycholic Acid-3-Sulfate	C ₂₇ H ₄₆ NO ₇ S	1950	100	
15	578.24657	0.270	3.61	0.02	45.9	1.081	217.5	0.1	Taurchenodeoxycholic Acid-3-Sulfate	C ₂₇ H ₄₆ NO ₇ S	1950	100	
16	562.23165	0.471	4.87	0.03	42.7	1.055	216.5	-0.2	Taurolithocholic Acid-3-Sulfate	C ₂₇ H ₄₆ NO ₇ S	1950	100	
17	391.28473	-0.652	9.21	0.17	12.3	1.001	207.5	0.5	Hyodeoxycholic Acid	C ₂₇ H ₄₆ O ₆	1950	100	
18	391.28426	-1.125	8.55	-0.06	3.2	0.998	206.8	-0.2	Ursodeoxycholic acid	C ₂₇ H ₄₆ O ₆	1950	100	
19	448.30094	-0.905	7.64	-0.01	5.8	0.968	199.9	-0.1	Glychenodeoxycholic acid	C ₂₇ H ₄₆ NO ₇	1950	100	
20	528.26319	-0.911	5.22	0.14	34.4	1.041	214.2	-0.0	Glycodeoxycholic Acid-3-Sulfate	C ₂₇ H ₄₆ NO ₇ S	1950	100	
21	528.26332	-0.341	4.94	0.13	25.8	1.041	213.8	0.0	Glycochenodeoxycholic Acid-3-Sulfate	C ₂₇ H ₄₆ NO ₇ S	1950	100	
22	464.30078	-0.977	4.89	0.06	41.7	0.977	201.4	0.1	Glycocholic acid	C ₂₇ H ₄₆ NO ₇	1950	100	
23	467.30254	1.110	8.22	0.02	45.1	0.996	205.3	0.2	Hyocholic acid	C ₂₇ H ₄₆ O ₆	1950	100	
24	432.34883	1.946	9.02	0.04	38.4	1.040	214.5	0.0	Lithocholic Acid 3-Sulfate	C ₂₇ H ₄₄ O ₆ S	1950	100	
25	498.28848	-1.001	4.33	-0.02	61.4	1.005	206.8	-0.1	Tauroursodeoxycholic acid	C ₂₇ H ₄₆ NO ₈	1950	100	

Fig. 3 Bile acids annotated in SRM 1950 human reference plasma using MetaboScape




Fig. 4 Discovery and assignment of Bile acids using interactive visualizations in MetaboScape

Figure 3 shows that in total 25 bile acids were annotated in the human plasma extract using the Target List for 71 bile acids (containing, name, molecular formula, retention time and CCS value). Annotation quality (AQ) was assessed based on AQ score, a visual tool in MetaboScape incorporating multiple molecular identifiers (e.g., CCS and *m/z*).

The use of the Kendrick Mass Defect plot within MetaboScape readily enabled mining of the data for possible bile acids that were not annotated by the Target List. Figure 4 shows a zoom into the mass and mobility region for Glycocholic acid and Glycocholic acid. A feature was spotted with a similar mass and mobility (in the plot below the two annotated bile acids). The retention time for this potential bile acid is lower compared to the other two Glycine conjugated bile acids. A likely annotation for this bile acid is Isoglycocholic acid. Isoglycocholic acid contains a 3 beta -OH orientation compared to Glycocholic acid. This difference in orientation renders Isoglycocholic acid more polar and can lead to an earlier reversed phase elution.

Conclusion

- LC-TIMS-MS based method for 71 bile acid leverages CCS as a comprehensive criterion for annotation confidence.
- The LC-TIMS-MS method provides quantitative results covering up to 3.8 orders linear dynamic range.
- Determined bile acid concentrations in SRM 1950 human reference plasma are within the range of published reference values.
- Untargeted Profiling using the MetaboScape solution enabled to readily pinpoint and annotated a bile acid not contained in the custom Target List of 71 standards.

LC-TIMS-MS