

The SciLS™ autopilot – a new workflow facilitating MALDI Imaging measurement set-ups




¹Janina Oetjen, ¹Arne Fuetterer, ¹Christian Tabeling, ¹Arne Brix, ²Michael Easterling.

¹. Bruker Daltonics GmbH & Co. KG, Bremen, Germany. ². Bruker Daltonics Inc., Billerica, MA, USA


Introduction

The SciLS autopilot is a new software workflow that simplifies MALDI Imaging measurement set-up on the timsTOF fleX through automation for reliable and robust results.

Methods



Results




Slide adapter with two prepared IntelliSlides indicating the areas for automatic target profile determination and laser focus tuning (red circles).


m/z calibration

Three lipid standards were sprayed on IntelliSlides and the mass error determined in four independent measurements. The mass error for the protonated and sodiated forms was below 2 ppm.

We compared the outcome of an imaging run using the automated workflow with the results of a run starting the acquisition manually by an experience user.



Representative ion images of a rat testis dataset measured on the timsTOF fleX in the lipid range with 20 μm pixel size using either the automatic or the manual workflow to set-up the acquisition. The ion intensities and distributions are similar in both cases.



Segmentation map of the lipid imaging dataset after automatic and manual acquisition set-up demonstrating that the tissue type is the main discriminating factor rather than the acquisition set-up method.



Conclusions

The SciLS autopilot:

- Allows sample tracking based on permanently inscribed barcode.
- Enables easy operation that limits user input to a few key parameters.
- Stands for increased reliability and robustness due to less user introduced errors.
- Includes an integrated IntelliSlide workflow for automatic registration of reference fiducials.
- Contains automated algorithms for controlled instrument performance.
- Allows automatic tissue detection.