

4D-Lipidomics profiling of cardiolipins on heart tissue from Barth syndrome patients using timsTOF mass spectrometry

Yorrick R.J. Jaspers¹, Sven W. Meyer², Jill Hermans^{1,3}, Bauke V. Schomakers^{1,3}, Jan Bert van Klinken^{1,3}, Eric Wever^{1,3}, Stephan Kemp¹, Rieckelt H. Houtkooper¹, Frédéric M. Vaz^{1,3}, Michel van Weeghel^{1,3}, Victor Fursey², Michel Boisvert²


¹ Laboratory Genetic Metabolic Diseases, Amsterdam UMC, University of Amsterdam, Amsterdam Gastroenterology Endocrinology Metabolism, Amsterdam, Netherlands

² Bruker Daltonics GmbH & Co. KG, Bremen, Germany

³ Core Facility Metabolomics, Amsterdam UMC, Amsterdam, Netherlands

Introduction

- Barth Syndrome is caused by defects in the gene *TAFAZZIN* which is involved in the remodeling of cardiolipins (CL).
- CL are involved in the mitochondrial respiratory chain, mitochondrial dynamics, mitophagy and apoptosis.
- Defects in *TAFAZZIN* results in decreased levels of CL and increased levels of monolysocardiolipins (MLCL)



Method

- We performed 4D-Lipidomics (retention time, *m/z*, ion mobility and fragmentation) using a reversed phase HPLC-MS platform using a timsTOF Pro (Bruker) on heart tissue from Barth syndrome patients and controls.
- A vacuum-insulated probe heated ESI source was used and ions were separated and fragmented using Parallel Accumulation Serial Fragmentation (PASEF).
- Lipid annotation was done using an in-house bioinformatics pipeline in combination with MetaboScape[®] 2023 (Bruker).






Results

- Using 4D-Lipidomics we were able to generate an in-depth characterization of cardiolipins in heart tissue from Barth syndrome patients and controls.
- Our results confirm earlier findings that Barth syndrome results in elevated MLCLs and decreased CLs.
- Barth syndrome results in abnormal CL fatty acid composition which is characterized by less incorporation of unsaturated fatty acids (FA).
- PASEF enabled high speed MS/MS acquisition which in combination with MetaboScape[®] allowed for the characterization of the fatty acid composition of CL and MLCL.

4D-Lipidomics

MetaboScape 2023

RT [min]	<i>m/z</i> meas.	M meas.	$\Delta m/z$ [ppm]	CCS (Å ²)	mSigma	Ions	MS/MS	Name	Molecular Formula	Annotations	AQ
70	13.73	1258.88166	1240.84769	-0.042	349.7	13.8	□	CL 28:0/28:0	C ₆₅ H ₁₂₆ O ₁₇ P ₂	IS	■
71	13.95	1386.94370	1368.90987	0.417	364.2	32.0	□	CL 30:2_36:4	C ₇₅ H ₁₃₄ O ₁₇ P ₂	IS	■
72	14.12	1388.96006	1370.92624	0.119	365.4	21.9	□	CL 32:2_34:3	C ₇₅ H ₁₃₆ O ₁₇ P ₂	IS	■
73	14.38	1416.99697	1398.96315	4.933	369.5	28.4	□	CL 32:2_36:3	C ₇₇ H ₁₄₀ O ₁₇ P ₂	IS	■
74	14.03	1412.95968	1394.92585	0.075	367.7	22.4	□	CL 32:3_36:4	C ₇₇ H ₁₃₆ O ₁₇ P ₂	IS	■
75	14.82	1423.03840	1405.00458	2.016	375.1	32.4	□	CL 34:1/34:1	C ₇₇ H ₁₄₆ O ₁₇ P ₂	IS	■
76	14.88	1449.05388	1431.02005	0.772	378.3	23.7	□	CL 34:1_36:2	C ₇₉ H ₁₄₈ O ₁₇ P ₂	IS	■
77	14.73	1447.03978	1429.00703	1.441	376.9	13.0	□	CL 34:1_36:3	C ₇₉ H ₁₄₆ O ₁₇ P ₂	IS	■
78	14.58	1445.02456	1426.99273	1.862	374.8	35.9	□	CL 34:2_36:3	C ₇₉ H ₁₄₄ O ₁₇ P ₂	IS	■
79	14.43	1443.01095	1424.97599	3.899	373.0	41.8	□	CL 34:3_36:3	C ₇₉ H ₁₄₂ O ₁₇ P ₂	IS	■
80	14.26	1440.99529	1422.96098	4.718	372.4	7.5	□	CL 34:3_36:4	C ₇₉ H ₁₄₀ O ₁₇ P ₂	IS	■
81	14.12	1438.97637	1420.94254	1.424	371.2	38.7	□	CL 34:3_36:5	C ₇₉ H ₁₃₈ O ₁₇ P ₂	IS	■
82	14.94	1475.06836	1457.03454	-0.417	381.3	72.9	□	CL 36:2/36:2	C ₈₁ H ₁₅₀ O ₁₇ P ₂	IS	■
83	14.78	1473.05516	1455.02175	1.529	380.1	25.3	□	CL 36:2_36:3	C ₈₁ H ₁₄₈ O ₁₇ P ₂	IS	■
84	14.64	1471.04133	1453.00677	2.704	378.9	28.5	□	CL 36:3/36:3	C ₈₁ H ₁₄₆ O ₁₇ P ₂	IS	■



Highlights

- 4D-lipidomics using PASEF is a promising technique for in-depth characterization of cardiolipins in complex samples.
- PASEF enables high speed acquisition of MSMS spectra.
- Barth syndrome results in significant changes in cardiolipin levels. Further investigations are needed to analyze whether these systemic lipidomic changes can be tied to the pathology of Barth syndrome.

4D-Lipidomics