



 Rapid and simple determination of low level fipronil and amitraz in egg samples using the Bruker EVOQ<sup>™</sup> LC-TQ Elite MS/MS system

## Abstract

We have developed a fast, sensitive, and reliable method to determine the presence of fipronil and amitraz in eggs using a simple extraction process and the Bruker EVOQ LC-TQ Elite mass spectrometry system. The proposed assay uses a liquid-liquid extraction (LLE) of a homogenized egg sample with acetonitrile (ACN). After filtering, the organic extract is injected into the system autosampler for unattended analysis. A method reporting limit (MRL) of 5 µg/kg (ppb) was established, meeting the quality criteria under ISO standard 17025 and the European validation procedures described in the SANTE Guidelines [1].

The exceptional sensitivity of this method achieves sub-ppb limits of detection (LOD) of both fipronil and amitraz.

Keywords: Fipronil, Amitraz, Egg sample, Food Safety Alert, ISO 17025 SANTE Guidelines Introduction

In August 2017, the Netherlands Food and Consumer Product Safety Authority (known as the NVWA in Dutch) and various media outlets worldwide warned of unusually high levels of fipronil and amitraz in eggs distributed in at least 19 countries, primarily in Europe [2,3].

Fipronil is an insecticide that belongs to the phenylpyrazole family (CAS 120068-37-3) and is commonly used in pest management across Europe. It is banned in products destined for the food chain. Amitraz is an amidine used as an acaricide-insecticide (CAS 33089-61-1) [4]. The maximum residue limit established by the EU [5] for both compounds and their metabolites in eggs is 0.005 mg/kg forfipronil [6] and 0.01 mg/kg for amitraz [7].

The extraordinary sensitivity of the Bruker EVOQ LC-TQ Elite mass spectrometry system enables analysis with a minimal sample quantity and simple sample preparation, avoiding complex cleaning stages, to achieve reliable detection and quantitation below the established maximum residue limits.

The entire methodology was developed in accordance with the Quality Requirements for Testing and Calibration Laboratories (ISO standard 17025).

## **Experimental**

#### **Sample Preparation**

2g of homogenized egg (organic pesticide-free chicken eggs) are weighed in a 20mL conical centrifuge tube and 2mL of acetonitrile (ACN) are added to the sample. The mixture is shaken in an orbital shaker for 30 seconds. After this time, the samples are centrifuged at 4000 rpm and 1 mL of the organic supernatant is collected. The aliquot is then filtered (0.22  $\mu$ m nylon filter) using a syringe and the extract is placed in the system autosampler vial for unattended analysis.

Working (stock) solutions of the target analytes were prepared in ACN using individual standards supplied by Accu-Standard, Inc. (New Haven, CT, USA). The compounds analyzed are shown in Figure 1. Spiked egg samples were prepared with the addition of all target analytes (at the indicated concentrations)

prior to the addition of the ACN. Triphenylphosphate (TPP) was added as an internal standard for quantitation calculations.

#### Methodology

Instrument conditions are summarized in Table 1 and Figure 2. The MRM conditions for the analyzed compounds and their metabolites are shown in Table 2.

#### **Analytical method**

The acquisition method is rapid and simple to create due to the intuitive method editor built into the Bruker MS Workstation data acquisition software.

Bruker's TASQ<sup>™</sup> was used for the processing and statistical calculations of the quality parameters. This software allows for rapid and accurate review of results. Figure 3 provides an example of the detailed results display. Figure 4 shows the MRM chromatogram of an egg sample spiked with 5 µg/kg (ppb) of the compounds analyzed. Each com-




Figure 1: Compounds analyzed. Left to right: fipronil, fipronil sulfone, amitraz, and N-(2,4-Dimethylphenyl)-N'methylformamidine [DMPF, an amitraz metabolite]

#### Table 1: Mass Spectrometry Method Conditions

Mass Spectrometer	Bruker EVOQ™ LC-TQ Elite MS/MS system					
HESI	4000 (pos. mode), 3800 (neg. mode)					
Cone Temperature	250 °C					
Cone Gas Flow	20 psi					
Nebulizer Gas Flow	60 psi					
Heated Probe Temperature	400 °C					
Probe Gas Flow	50 psi					
CID Gas	Ar, 2.0 mtorr					
Detector Mode	EDR					
Polarity	Automatic polarity switching ( +/-)					
Liquid Chromatography	Bruker Elute™ UHPLC system					
LC Column	Bruker Intensity Solo C18 100 x 2.1 mm (P/N:BRKHSC18022100)					
Mobile Phase A	Water + 0.05% formic acid + 2 mM ammonium formate					
Mobile Phase B	Acetonitrile + 0.05% formic acid + 2 mM ammonium formate					
Flow Rate	400 µL/min					
Injection Volume	2 µL					
Column Oven Temperature	40 °C					
Total Run Time	14 min					
Software	HyStar 4.1/Bruker MSWS 8.2.1/TASQ 1.4 processing software					

### LC Gradient

Time table

Line:	Time(min):	Flow (mL/min):	%A:	%B:
01	0.00	0.400	98	2.0
02	0.40	0.400	98	2.0
03	0.50	0.400	65	35.0
04	7	0.400	2	98.0
05	11	0.400	2	98.0
06	11.1	0.400	98	2.0
07	13.5	0.400	98	2.0

💽 📄 👔 🚺 🗉 Advanced




Figure 2: LC gradient overview, as shown in the Elute series plug-in for the HyStar control software

pound presents two transitions (MRM) in order to meet the validation criteria.

## **Results and discussion**

For validation of the method, a control and quality assurance procedure was

applied using samples from organic, pesticide-free eggs spiked with the compounds to be analyzed and covering a wide range of concentrations.

#### Table 2: MRM transitions for the compounds analyzed

Compounds	RT (min.)	Precursor ion	Polarity	Quantita- tion ion	CE (V)	Confirma- tion ion	CE (V)
Amitraz metabolite DMPF	1.73	163	Positive	107	-26	122	-17
Fipronil	4.68	435	Negative	330	+12	250	+25
Fipronil Sulfone	5.10	451	Negative	282	+27	415	+14
Amitraz	6.48	294.3	Positive	163	-12	122	-28
TPP (IS)	4.90	327	Positive	152	-29	215	-20

The analytical quality criteria included evaluations of linearity, precision and accuracy.


Linearity

The linearity of the method was demonstrated using a calibration curve with egg samples spiked across the concentration range under study. Four solutions were prepared with the following concentrations: 5 ppb, 10 ppb, 25 ppb and 100 ppb. The applied quality criteria for linearity are:  $R^2 \ge 0.990$  with a variation (RSD %) of the curve's response factor  $\le 30\%$  ( $\sigma$ RF) for the entire concentration range under study.

Figure 5 shows the calibration curves for all compounds analyzed. A summary of the calibration results is shown in Table 3, which demonstrates that all of the compounds analyzed present an  $R^2 > 0.99$  and RSD < 15%, securely meeting the pre-established quality criteria.

#### Precision

To determine precision, expressed as repeatability, the RSD was calculated for five replicates (n=5) of the egg samples spiked at corresponding low concentration levels of  $5 \mu g/kg$  (ppb) and  $10 \mu g/kg$  (ppb).






Figure 4: MRM chromatogram of an egg sample spiked with 5 µg/kg (ppb) of the analyzed compounds, as shown in the TASQ software




Figure 5: Calibration curves for all compounds analyzed showing results for R<sup>2</sup> and RSD, as shown in the TASQ software

#### Table 3: Summary of linearity results for compounds analyzed

Table 4: Summary of repeatability results for egg samples spiked with the targeted analytes

Summar	y of linearity res	sults	Summary of repeatability results				
				RSD %			
Compounds	R <sup>2</sup>	RSD %	Compounds	Sample spiked with 5 ppb	Sample spiked with 10 ppb		
Amitraz metabolite	0.99986	4.25	Amitraz metabolite	2.06%	2.59%		
Fipronil	0.99881	5.27	Fipronil	4.84%	5.23%		
Fipronil sulfone	0.99982	3.90	Fipronil sulfone	4.28%	3.61%		
Amitraz	0.99639	14.99	Amitraz	5.64%	6.60%		

The TASQ software automatically performs the statistical calculation, as can be seen in Figure 6 for fipronil.

A summary of the repeatability results is shown in Table 4, demonstrating that all of the compounds analyzed present repeatability of 2-7% for low concentration levels.

#### Accuracy

Recoveries (in %) were calculated on the basis of five independent extractions (n=5) from egg samples spiked with all analytes at  $5 \mu g/kg$  and  $10 \mu g/kg$  concentration levels. All of the results are provided in Figure 7: the "Quantity" column indicates the measured concentration value and the "Accuracy" column shows the recovery percentage automatically calculated by the TASQ software.

Recoveries at 5 and  $10\,\mu g/kg$  concentration levels are between 80-102% for all analytes.

#### **Sensitivity and detection limits**

The high sensitivity of the method can be seen in Figure 8, which shows the MRM chromatograms for all compounds analyzed for the first calibration level of 5  $\mu$ g/kg in the extracted egg samples.

It has been demonstrated that both quantitation and confirmation ions have an excellent signal-to-noise (S/N) ratio at 5 ppb for all compounds analyzed, leading to the conclusion that sub-ppb limits of detection (LOD) can be easily achieved for all compounds analyzed in the egg matrix.

•••• RSD	Statistics Graph	RSD Statistics Ta	able 🛛				
Row	Concentration	rel.Std.Dev	Mean	Median	Std.Dev.	Minimum	Maximum
1	10.0	5.23	1225554	1224894	64153	1165218	1324668
2	5.0	4.84	635391	625135	30729	606549	687155

Figure 6: Evaluation of precision of fipronil replicate quantitation at 5 and 10 ppb, as shown in the TASQ software

Row	Data Set	Analyte	Residual [%]	∆RT [min]	Area	Quantity [ug	Height	RT [min] exp.	S/N	Quantity	Accuracy [%]
1	egg_5ppb_2_5	Amitraz	-18.6	0.02	434273	4.1	122444	6.47	1813	4.1 ug/kg	81.4
2	egg_5ppb_2_5	Amitraz deg	-0.7	0.01	16235	5.0	5351	1.73	152	5.0 ug/kg	99.3
3	egg_5ppb_2_5	Fipronil	-10.5	0.02	635388	4.5	186299	4.67	6287	4.5 ug/kg	89.5
4	egg_5ppb_2_5	Fipronil sulfona	-4.3	0.02	605302	4.8	180162	5.08	7395	4.8 ug/kg	95.7
Row	Data Set	Analyte	Residual [%]	ΔRT [min]	Area	Quantity [ug	Height	RT [min] exp.	S/N	Quantity	Accuracy [%
1	egg_10ppb_2_5	Amitraz	-9.0	0.01	607831	9.1	173780	6.47	2874	9.1 ug/kg	91.0
2	egg_10ppb_2_5	Amitraz deg	-2.5	0.01	28379	9.8	9254	1.73	281	9.8 ug/kg	97.5
3	egg_10ppb_2_5	Fipronil	1.8	0.01	1240197	10.2	356130	4.67	12274	10.2 ug/kg	101.8
4	egg_10ppb_2_5	Fipronil sulfona	-0.6	0.02	1118584	9.9	329380	5.08	15467	9.9 ug/kg	99.4

Figure 7: Recovery for egg samples spiked with all compounds at 5 µg/kg and 10 µg/kg concentration levels, as shown in the TASO software




Figure 8: MRM chromatogram of an egg sample spiked with all analytes at 5 µg/kg. Quantitation ions and confirmation ions are shown, as well as signal-to-noise (S/N), for each compound, as shown in the TASQ software.

# Conclusion

A rapid and ultra-sensitive method with a very simple sample preparation process has been developed for the detection of fipronil and amitraz in egg samples using the Bruker EVOQ<sup>™</sup> LC-TQ Elite mass spectrometry system. The method provides excellent analytical quality parameters for linearity, accuracy and precision in accordance with the ISO 17025 standard. This method has been validated and is ready for implementation in food quality control process laboratories.





# Learn More

You are looking for further Information? Check out the Link or scan the  $\ensuremath{\mathsf{QR}}$  Code.

www.bruker.com/evoq-lc

#### References

[1] European Commission, SANTE/11945/2015
[2] http://www.bbc.com/news/world-europe-40896899
[3] https://english.nvwa.nl/topics/fipronil-in-eggs
[4] http://www.eurl-pesticides-datapool.eu/
[5] http://ec.europa.eu/food/plant/pesticides/eu-pesticides-database
[6] Commission Regulation (EU) No 1127/2014
[7] Commission Regulation (EU) No 623/2017

For research use only. Not for use in diagnostic procedures.

#### Bruker Daltonics GmbH & Co. KG

## Bruker Scientific LLC

Bremen · Germany Phone +49 (0)421-2205-0 Billerica, MA · USA Phone +1 (978) 663-3660

ms.sales.bdal@bruker.com - www.bruker.com