



## Automatic Concentration Determination Of Compound Libraries by NMR

### ● CMC-q<sup>™</sup> Complete Molecular Confidence - Quantification

Complete Molecular Confidence for quantification (CMC-q) is a complete workflow solution that facilitates automatic NMR-based quality assurance in batches.

- Controls experiment setup
  - Manages data acquisition and processing
  - Automatically analyzes the data
  - Easy-to-use viewer for browsing and editing results
  - Export results to Excel or text files
- CMC-q provides quick access to automated NMR quality assurance and quantification of larger batches of samples. Delivering accurate, precise information on sample concentration and water content in typical screening samples, CMC-q also marks questionable structures and provides a suggestion for spectral assignment. This is ideally complemented with LC-MS information such as derived from Bruker's SmartFormula program.
  - Operating on a file-in, file-out basis, the user supplies an input file describing the samples to be measured, and receives an output file of results.
  - The spectra interpretation function can also analyze individual datasets and prepare spectra for publication or for further analysis.

# CMC-q Workflow



## Setup

- Sample preparation from well plates
- Samples accompanied by SD Files
- SD Files used for experiment set up

## Acquisition


- Experimental parameters for standard NMR experiments in deuterated solvents, and typical screening samples in non-deuterated DMSO

## Analysis

- Automatic NMR analysis of concentration, water content (for typical liquid store samples in non-deuterated DMSO), and plausibility with suggestion for assignment

## Report



- Automatic report generation with easy viewing and editing, exportable to Excel or text file formats



## Quality of results

- CMC-q analysis assigns concentrations to over 95% of all datasets
- CMC-q also finds over 85% of the structures plausible, which the human expert found consistent.
- On these, more than 90% of the automatically assigned concentrations deviate only 20% or less from the manual expert result.
- This data was collected on ca. 2500 expert-interpreted spectra
- The compounds are commercial Lipinsky rule compounds from TimTec's APEX screening library and compounds from Actelion's own compound library
- Typically, 6-8ul samples in non-deuterated DMSO where used and data was acquired on 400 MHz instruments equipped with 1 mm MicroProbes and SampleJets. The turnover was a few minutes per sample.

## Single and multiple result view



## Single and multiple result view

- Overview window of scaled dots representing concentration levels
- Dot color indicates structural plausibility status
- Overview of water content
- Concentration values that fall outside defined thresholds (low and high) can be highlighted
- Modify integral ranges, assignments and other results easily in the interactive CMC-q viewer
- User interaction can be logged

A	B	C	D	E	F	G	H	I	J	K	L
Index	Status	Structure	Sum formula	Mass	Structure used	Spectrum	Water [v/v]	Conc [mM]	Proton	Shift [ppm]	
2	1	<chem>C19Q1M02S</chem>	255.1293	yes	CMCQ_LibCompound1_1_1	34	74.3	1	1.51		
3	2	<chem>C10H14N02SF</chem>	231.0794	yes	CMCQ_LibCompound2_1_1	34	150.1	6	0.73		
4	3	<chem>C10H13N03SC0</chem>	286.9993	yes	CMCQ_LibCompound3_1_1	34	242.3	2	1.55		
5	4	<chem>C8H6N02SC0</chem>	263.9527	yes	CMCQ_LibCompound4_1_1	34	92.3	1	7.95		
6	5	<chem>C12H19N02S</chem>	241.1366	yes	CMCQ_LibCompound5_1_1	34	150.1	6	0.80		

## ● Bruker BioSpin

info@bruker-biospin.com  
www.bruker.com/cm-c-q

We would like to thank Actelion Pharmaceuticals Ltd.  
for the very fruitful collaboration.

