

New Capabilities with X-Detection

Cooled ¹³**C and** ¹⁵**N Preamplifiers and TXO CryoProbes**

High-field high-resolution NMR in structural biology has been very successful using inverse detection schemes. The most popular CryoProbe models for macromolecular work are therefore triple resonance inverse (TCI) or quadruple resonance inverse (QCI) probes. Active (cryogenically cooled) preamplifier electronics for ¹³C have been standard for a number of years, for almost all CryoProbes.

Over the past decade, this increased sensitivity has supported the development of direct ¹³C detection.¹ Direct ¹³C detection is now an established, powerful, complementary technique, especially for the study of intrinsically disordered proteins (IDP) and proteins with paramagnetic centers.² For very large globular proteins, ¹H and ¹³C approach their natural limits, while the sharp ¹⁵N lines of a few Hz line width still offer long magnetic coherence lifetimes and high resolution.



Recently, a collaboration between groups in Harvard and Tokyo³ has demonstrated great potential for using ¹⁵N direct detection in combination with TROSY. Additionally, the need for perdeuteration is reduced when using the ¹⁵N detected approach, thus facilitating protein expression and detection of amide protons involved in hydrogen bonds.





Fig. 2 1.8mM ¹³C/¹⁵N labeled Ubiquitin, ¹⁵N detected refocused INEPT at 950 MHz: Top: 298 K, bottom 278 K. Reduced HN-solvent exchange shows larger benefit for many resonances in ¹⁵N detection as compared to signal broadening due to slower tumbling rates. Typical linewidth 3-8 Hz. Adiabatic ¹³C (Chirp, 48 kHz sweep) CPD decoupling. NS = 256, experimental time 7 min 47 sec



Fig. 3 1.8mM ¹³C/¹⁵N labeled Ubiquitin. Left: ¹⁵N detected 2D C-N and CO-N projections of 3D (H)CCCON. NS = 16, TD = $2k \times 64 \times 96$; total experimental time 37 hours.

950 MHz 5mm TCI with cooled ¹⁵N preamplifier

Fig. 1 1.3mM ¹³C/¹⁵N labeled Liver Basic Fatty Acid Binding Protein (LB-FABP) dual receiver experiment. Two scans per increment, 12 hours total acquisition time. 1k x 40 x 32 complex points after IPAP processing. Structure ensemble from Vasile et al. J Biomol NMR; 25, 157-160, 2003.

Direct Observe CryoProbes offer several advantages over indirect probes. The increased X-nucleus sensitivity makes the development and use of ¹³C-¹H dual-receive experiments much more attractive, as the relative sensitivity between the two nuclei is much more balanced and hence the compromises to be taken are significantly reduced when concatenating two experiments. In addition, the effect of ionic strength on ¹H sensitivity is reduced, as the filling factor for ¹H is lower compared to inverse probes.

References

- 1. Bermel W, Bertini I, Felli IC, Kuemmerle R, Pierattelli R, 2003, J Am Chem Soc 125:16423-9; Bermel W, Bertini I, Felli IC, Piccioli M, Pierattelli R, 2006, Prog Nucl Magn Res Spec 48:25-45



Fig. 4 1.8mM ¹³C/¹⁵N labeled Ubiquitin. Left: ¹⁵N detected high resolution 2D refocused and **fully decoupled HX-INEPT** (298 K). NS = 4, TD = 4k x 128; total experimental time 11 minutes. Right: zoom in to a single cross-peak when acquiring a fully coupled HX-INEPT illustrating the TROSY effect. Data acquired at 278 K.

Summary

- ¹³C and ¹⁵N detection is attractive for IDR & IDP.
- ¹³C detection is a well established alternative to traditional indirect detection.
- ¹⁵N detection shows potential especially for high molecular weight globular macromolecules.
- Active preamplifier electronics for ¹³C and ¹⁵N



2. Bertini I, Lee YM, Luchinat C, Piccioli M, Poggi L, 2001, ChemBioChem2:550-558



field yields high resolution and sensitivity for protein NMR, J. Biomol. NMR, 63(4):323-31

are available for new TCI's and TXO's.