

# A novel supervised learning algorithm for real-time collision energy selection to optimize peptide fragmentation in ion mobility-mass spectrometry

# Introduction

- $\succ$  The ability to identify peptides, proteins, and their associated post-translational modifications using mass spectrometry is directly linked to the level of fragmentation of peptide precursor ions.
- $\succ$  Many precursor ion properties such as ion-mobility coefficient, mass-to-charge ratio and charge state affect the amount of collision energy required for optimal fragmentation.



### **Objective**

 $\succ$  To build an artificial neural network that selects collision energy that optimizes peptide fragmentation to improve peptide and protein identification sensitivity.

# Methodology

HeLa whole cell lysate digest dataset from Bruker timsTOF Pro: 1,354,136 peptides fragmented with different collision energy values (5-100eV)

# values



identification score.



Yun-En Chung<sup>1</sup>, Matthew Willetts<sup>2</sup>, Nagarjuna Nagaraj<sup>2</sup>, Jens Decker<sup>2</sup>, Jonathan Krieger<sup>2</sup>, Tharan Srikumar<sup>2</sup>, Mathieu Lavallée-Adam<sup>1</sup>

<sup>1</sup> Department of Biochemistry, Microbiology and Immunology and Ottawa Institute of Systems Biology, University of Ottawa, Ottawa, Ontario, Canada <sup>2</sup> Bruker Daltonics, Billerica, Massachusetts, United States

### > Train artificial neural network to predict peptide identification confidence score using precursor ion properties as input.



m/z and collision energy-dependence of peptide

### Predicted peptide identification scores correlate with observed

### > Use the trained artificial neural network to predict a peptide's identification score at different collision energy values.

### $\blacktriangleright$ Select optimal collision energy that maximizes identification score.

Artificial neural network-predicted v.s. observed identification scores at different collision energy values for peptide ion INEELESQYQQSMDSK (+2)

### Results

### Artificial neural network captures the relationship between

Bin-averaged PEAKS score

**Empirical** optimal CE: **51.0 eV** 

### **Collision energy values from the artificial neural network reflect** empirically-identified optimal collision energy values better than default instrument values.



Artificial neural network-selected optimal collision energy shows distinct patterns according to charge state and varies from default collision energy used.



Applying optimized collision energy in real-time improves fragmentation by reducing intensity of unfragmented precursor ion by 2.18-fold on average, without sacrificing peptide identification sensitivity.



Artificial neural network-prediction

Identified 9.7% more peptides with post translational modifications in phosphopeptide enriched samples when using collision energy values selected by the artificial neural network.



## **Conclusion and Future Directions**

- $\succ$  The peptide identification score of a given precursor ion can be accurately predicted by an artificial neural network using only its ion properties and a collision energy value, without any peptide sequence information.
- $\succ$  The artificial neural network enables selection of the optimal collision energy for a given peptide, which improves fragmentation and increases the number of phosphopeptides identified.
- $\succ$  We will assess how our approach improves characterization of samples that are challenging to characterize, such as post-translational modifications of different types (eg. acetylation) or cross-linked peptides.

## Acknowledgement

I would like to thank all the members of the Lavallée-Adam lab and our collaborators from Bruker Daltonics for their support and guidance on this project.



Contact: vchun060@uottawa.ca





