

## High-throughput UHPLC-TIMS-based lipidomics: Sub-5 min screening of plasma, cells, stool and tissues

## Introduction

- > In the last years ion mobility mass spectrometry has emerged as an interesting alternative to perform lipid profiling.
- Reversed-phase UHPLC-MS represents the golden standard in untargeted lipidomics.
- Conventional methods are characterized by long analyses times that are not suitable for large cohort screening.
- $\succ$  In this contribution we evaluated the implementation of trapped ion mobility mass spectrometry (TIMS) as tool to reduce analysis time while keeping accuracy and coverage.



## Results

- > The developed platform was first challenged against human plasma extract.
- > Following careful manual curation, the method showed average MS/MS scores,  $\triangle$ ppm and  $\triangle$ CCS (Å<sup>2</sup>) errors of: 876.5, -0.57 and 1.47, respectively.
- > We evaluated the coefficient of variation (CV%) of retention times and CCS values for more than 800 lipids over 6 consecutive days achieving average values of 0.54 % and 0.34 %, respectively.
- The method was extended to different biological specimens, such as AC-16 cell lines, mouse stool and brain homogenates, with 950, 518 and 645 lipids monitored respectively, and covering more than 30 subclasses.



<u>Erica Forsberg<sup>1</sup></u>, Eduardo Sommella<sup>2</sup>, Fabrizio Merciai<sup>2</sup>, Pietro Campiglia<sup>2</sup> <sup>1</sup>Bruker Scientific, San Jose, United States; <sup>2</sup>University of Salerno, Department of Pharmacy, Fisciano, Italy

