
# **Combining MALDI Imaging and Liquid Surface Extraction for Spatial Metabolomics**




IMSS II/Ourcon IV 2018

Jeremy Wolff<sup>1</sup>, Alain Creissen<sup>2</sup>, Matt Orcutt<sup>2</sup>, Jan Kobarg<sup>3</sup>, and Shannon Cornett<sup>1</sup>

<sup>1</sup>Bruker Daltonics Inc., 40 Manning Road, Manning Park, Billerica, MA 01821, USA <sup>2</sup>HTX Technologies, LLC, PO Box 16007 Chapel Hill, NC 27516 <sup>3</sup>Bruker Daltonik GmbH, Fahrenheitstraße 4, 28359 Bremen, Germany

### Introduction:


- While MALDI Imaging does an excellent job of providing localization of metabolites, lipids, and proteins, it can be challenging to ID and quantify these molecules.
- There is great interest in combining MALDI Imaging with Metabolomics to produce Spatial Metabolomics, where metabolites of interest can be both localized, identified, and accurately quantified.
- Current spatial metabolomics methods rely on excision of small regions of tissue, followed by homogenization/extraction and analysis by HPLC-MS.
- While this method can be used to accurately identify and quantify metabolites, it is not high throughput, cannot be easily automated, and has poor spatial resolution.
- Here, we combine liquid surface extraction of tissue slices using an HTX SepQuant with MALDI Imaging.



**Figure 1.** HTX SepQuant DropletProbe Workflow. **(A)** Drawing of surface sampling device using a PAL autosampler with a 10  $\mu$ L syringe. **(B)** Syringe is positioned 100 – 200  $\mu$ m above the tissue surface. **(C)** 0.5  $\mu$ L of extraction liquid is dispensed, making contact with the surface. **(D)** Extraction liquid is held on surface for ~2 seconds. **(E)** Extraction liquid is drawn back into syringe. Steps C – D are repeated 2 more times to optimize extraction. **(F)** Extraction liquid is injected into sample loop (LOAD) on a switching valve. **(G)** Switching valve switches (INJECT) and sample is separated on HPLC column.

## Methods:


- Fresh frozen rat brain was sagittally sectioned at 10  $\mu$ m and mounted on standard glass slides for SepQuant analysis, or ITO-coated glass slides for MALDI Imaging analysis.
- For SepQuant analysis, roomtemperature tissue sections were extracted with three,  $0.5-1~\mu L$  extractions of 50:50  $H_2O$ :ACN, and injected onto a 2.1 x 100 mm C18 column.
- The SepQuant workflow is shown in **Figure 1**. A zoom of the SepQuant performing liquid surface extraction is shown in **Figure 2**.
- A Bruker Elute pump performed a 5 minute separation gradient (10-90% B) at 500  $\mu$ L/min, and the eluted species were analyzed on the Bruker timsTOF and solariX XR instruments.
- For MALDI Imaging, tissue sections were coated with 9-AA using an HTX TM Sprayer. The coated tissue was analyzed on a Bruker solariX XR with 90  $\mu$ m pixel resolution.
- For MALDI Imaging, data was analyzed using FlexImaging and SCiLS. For liquid surface extraction, data was analyzed using DataAnalysis.
- Metabolites were identified by accurate mass, isotopic fine structure, MSMS fragmentation, and HMDB.




extraction interface is shown in the inset.

#### Results:


 MALDI Imaging of the rat brain produced typical negative ion results, as shown in Figure 3.






**Figure 3.** Negative ion MALDI Imaging results of rat brain, showing unique localization of lipids. **(A)** Localization of the lipid at m/z 885.55 and **(B)** localization of the lipid at m/z 888.630.


- Liquid Surface Extraction was used to probe various areas of a serial slice of rat brain. Extraction volumes of 0.5 1  $\mu$ L were used, which probed an area of 1 3 mm.
- Many small molecules were measured by both instruments and both ionization techniques. We observed a mix of ions that were detected by both ionization techniques, and were unique to each technique.
- One of the small molecules measured by the SepQuant with LCMS is the peak at m/z 309.1711, which eluted at ~3.35 minutes as shown in **Figure 4**.
- Using the high resolving power of the solariX XR to measure isotopic fine structure on an LC timescale, the formula for this peak was identified as  $C_{17}H_{26}O_5$ , as shown in **Figure 4**.



• Searching HMDB for either the accurate m/z or the formula generates the iridoid compound HMDB0040980 with the common name Valdiate. The structure of this compound is shown in **Figure 5**.




- Searching the MALDI imaging data for Valdiate at m/z 309.1696 produces the image shown in **Figure 6.** The compound is not strongly localized and is dispersed throughout the brain.
- The non-localization may be due to Valdiate associated with the membranes of the rat brain.




eluting at  $\sim$ 3.3 minutes from SepQuant using isotopic fine structure. Using SmartForula XR, the formula  $C_{17}H_{26}O_5$  was generated from the peak at m/z 309.17118.

- A similar extraction and LC analysis was performed on the **timsTOF Pro**, to generate ion mobility and MSMS data to aide in molecule identification.
- The peak at *m/z* 309.17 eluted at the same time with SepQuant LCMS on the timsTOF compared to SepQuant LCMS on the solariX XR, as show Figure 7.



• Using TIMS, the peak at m/z 309.17 produced a single mobility as shown in the Extracted Ion Mobilogram in **Figure 7**. This peak has a reduced mobility of 0.8 and a measured CCS of 167.8 Å<sup>2</sup>.



- Using the SepQuant and the timsTOF Pro, LCMSMS data for the peak at m/z 309.17 was measured. The MSMS spectrum of the peak at m/z 309.17 is shown in Figure 8.
- This MSMS spectrum did not match the theoretical MSMS spectrum for Valdiate from HMDB; further work is needed to confirm this MSMS data with Valdiate.

# **Conclusions:**

MALDI Imaging coupled with liquid surface extraction using the SepQuant provides a powerful tool to couple molecule identification and localization using isotopic fine structure, CCS measurements, and MSMS.

MALDI Imaging