

# Trapped Ion Mobility PASEF Based Lipidomics Highlights Potential Lipid Biomarkers of Covid-19 severity



Eduardo Sommella<sup>1</sup>, Fabrizio Merciai<sup>1</sup>, Emanuela Salviati<sup>1</sup>, Albino Carrizzo<sup>2</sup>, Paola di Pietro<sup>2</sup>, Michele Ciccarelli<sup>2</sup>, Aiko Barsch<sup>3</sup>, Sven W. Meyer<sup>3</sup>, Carmine Vecchione<sup>2</sup>, Pietro Campiglia<sup>1</sup>

<sup>1</sup> University of Salerno, Department of Pharmacy, Fisciano, Italy; <sup>2</sup> University of Salerno, Department of Medicine, Baronissi, Italy; <sup>3</sup> Bruker Daltonics GmbH & Co. KG, Bremen, Germany


## Introduction

- SARS-CoV-2** was first reported in late 2019 and to date more than 160 million people have contracted the virus and almost 4 millions confirmed deaths have been documented.
- Covid-19 disease** presents a **varied symptomatology**: In most cases it appears with **mild** symptoms but can evolve into a **severe** condition with an acute respiratory syndrome requiring **hospitalization** and **intensive cares**.
- The identification of possible **prognostic Biomarkers**, which may predict the course of the disease or the response to therapeutic treatments, becomes a **fundamental goal**.
- Lipids play a key role in the regulation of numerous pathways, in this study a **4D-Lipidomics** approach was performed for the analysis of **Covid-19 plasma samples** of patients with different degree of severity.

## Methods

- Lipids were extracted from plasma by the **Matyash<sup>11</sup> method** (MeOH/MTBE/H<sub>2</sub>O).
- RP-UHPLC-TIMS-MS/MS** was employed for untargeted lipidomics data acquisition.
- PASEF** mode was used for DDA, in both ESI(+) and ESI(-) ionization.
- MetaboScape 2021** and **MetaboAnalyst 5.0** were used for Identification and Statistical analysis.


[1] Matyash V., et al. 2008, J. Lipid. Res., <https://doi.org/10.1194/jlr.D700041-JLR200>



### Enrollment

- Patients were classified according to the clinical phenotype following World Health Organization (WHO) severity score as **Mild** (n= 44) patients or **Severe** (n= 54) patients.
- Healthy and SARS-CoV-2 negative subjects (n = 21) were recruited from healthcare workers.

## Study Design




### Blood sample processing

- From each patient, 5-8 mL of whole blood was drawn into EDTA vacutainers and centrifuged at 1,000 x g for 20 minutes at 25°C to separate blood cells and plasma. After collection, plasma samples were aliquoted and stored at - 80°C until analysis.

### Plasma lipidome extraction

- Lipids were extracted from 20 µl of plasma, spiked with a mix of deuterated standards (Splash Lipidomix®, Avanti Polar).
- A pooled QC was prepared to assess repeatability.

## Analytical Workflow



### Lipid annotation and statistical analysis


### Analytical platform

- Untargeted lipid profiling was performed by RP-UHPLC hyphenated to a timsTOF Pro (Bruker) Q-TOF mass spectrometer.
- Lipids were separated on a CSH 100 mm x 2.1 mm, 1.7 µm column (Waters), using as mobile phases: A) 10mM HCOONH<sub>4</sub> + 0.1% HCOOH H<sub>2</sub>O/ACN 60/40 (v/v) B) 10mM HCOONH<sub>4</sub> + 0.1% HCOOH IPA/ACN 90/10 (v/v). T=55°C.

- Identification was performed in **MetaboScape** (Bruker) by using LipidBlast spectral library and Rule-Based Annotation. MS/MS score, ΔCCS% values (vs predicted) and retention time linearity were considered for lipid annotation.
- Pre-processing**: Internal standard normalization, log transformation, auto-scaling.
- MetaboAnalyst**: Univariate (one-way ANOVA) and multivariate (PLS-DA) analysis were employed to identify statistically significant lipids and to visualize class separation. A Random Forest (RF) model was trained and tested to predict severity.


## Statistical Analysis

Figure 1




**Figure 1) left:** 2D PLS-DA model score plot showing the discrimination of different classes: Covid (-) red, mild: green, severe: blue;  $R^2$  and  $Q^2$  results estimated via cross validation and permutation test results based on 1000 iterations; **right:** The 15 highest scoring variable importance in projection (VIP) metabolites are shown. The number of VIPs was established by setting the VIP-score  $\geq 1.8$  as a cutoff value. The colored boxes on the right indicate the relative amount of the corresponding lipid compound in each group.

Figure 2



**Figure 2)** Comparison of normalized intensity (based on labelled internal standards mixture, log transformed and autoscaled dataset) of the selected lipid panels in mild and severe survivors compared to non-survivor patients, \*\*\*\*  $p < 0,0001$

Figure 3



**Figure 3)** ROC curves for severity (left) and outcome (right) obtained with the predictive model (RF) on the reduced lipid panel composed by LPC O-18:1, PC O-34-3, LPC 20:1, LPC O-16-1, LPC 18:0, LPC O-16:0.

## Results

- 348 lipids** were confidently annotated, belonging to different lipid classes. Avg. MS/MS score: 914.60, Δppm: 0.60, ΔCCS%: 1.30%, avg. CV% (rt and area)  $\leq 12\%$ .
- PLS-DA highlighted a distinct **class separation**, among 15 highest scoring VIP lipids **9 were LPCs**. 191 lipids were found significantly modulated ( $p < 0.01$ , FDR: 0.01%).
- A **progressive increase of TGs, Cer, PEs levels** was found from mild to severe patients, whereas a **marked decrease of SMs, HexCer and mostly LPCs, LPC-O and PC-O** was observed in severe patients.
- A **RF model** was built to evaluate the **predictive potential** of the lipidome signature.
- The model provided AUC values of 0.751 (95% CI: 0.599-0.887) for severity and 0.815 (95% CI: 0.662-0.944) for outcome. **Using a reduced panel of six lipids selected as the most significant**, the model showed the same **predictive ability**.

## Discussion

- The target of our study was to identify, at the time of hospitalization, a **potential lipid signature** able to **predict the progression of Covid-19**.
- The **RP-UHPLC-TIMS-PASEF** approach resulted in **high lipidome coverage** and showed **high repeatability** across the batches.
- LPC class** was found as the **most discriminant** across the different conditions. In addition, **non-survivor patients** (either mild and severe) **showed lower values** with respect to survivors.
- A **restricted panel of six lipids** composed by LPC O-18:1, LPC 20:1, LPC O-16-1, LPC 18:0, LPC O-16:0, PC O-34-3 showed **high predictive ability in discriminate severity and mortality**.
- These results could suggest the potential employment of the defined lipidome signature as **prognostic tool** in **targeted approaches** and to **evaluate the recovery** of patients after discharge from the hospital.