

Development of proteomic-biomarker landscape of AY9944-treated mouse brain tissues mimicking Smith-Lemli-Opitz syndrome (SLOS) by using DIA-PASEF and DIA-FAIMS mass spectrometry

Mahmud Hossain¹, Maria Stella Ritorto¹, Lin An¹, Tracie Treleaven², Guoqing Sheng², Elizabeth Gordon³, Diego M Assis³, Matthew Willetts³, and Bailin Zhang¹
¹Disease Profiling & Functional Genomics, PMCB, Sanofi, ²Neuroinflammation & MS, Sanofi, Framingham, USA, ³Bruker Scientific LLC, Billerica, USA

Goal: Comparison of two very recent DIA mass spectrometry techniques - DIA-PASEF & DIA-FAIMS - in analyzing proteomic biomarkers related to the compound AY9944-treated mouse tissues mimicking Smith-Lemli-Opitz syndrome (SLOS).

SLOS: significant unmet need

- **A Rare disease** - An autosomal recessive inherited disorders with an incidence of **1 in 10,000–70,000** and a carrier frequency as high as **1 in 30**, and is caused by 7-dehydrocholesterol reductase (DHCR7) deficiency in cholesterol biosynthesis pathway
- **Hypomyelination is one of hallmarks-CNS defects** in white matter, often involving corpus callosum absence or hypoplasia and reported to involve absence of myelin and demyelination
- **Pathogenesis** - Multiple congenital malformation, mental retardation with behavioral phenotype
- **Biological target** - Lost of function of DHCR7 enzyme activity
- **Biochemical biomarkers**- Accumulation of 7DHC (2.7–470 ug/ml; 10–2000-fold normal) and a deficiency of endogenous cholesterol





Fig. 2. (left-side) QC (sample pool) performances are satisfactory. **(right-side)** Ctrl- and treated-samples of 21- and 35-days are well segregated in 2D PCA plot.




Fig. 5. (left-side) PLS (partial least square regression) for biomarker selection of Bruker timsTOF Pro mass spectrometer acquired samples. Component 1: 21-days vs. 35-days; Component 2: Treated vs. Ctrl. **(right-side)** Top 10 proteins in component 2





Fig. 3. (left-side) PLS (partial least square regression) for biomarker selection of Thermo Exploris480 mass spectrometer acquired samples. Component 1: 21-days vs. 35-days; Component 2: Treated vs. Ctrl. **(right-side)** Top 10 proteins in component 2

Differential proteins (p<0.05)



Mom1 vs Mom2 (21-days) Mom3 vs Mom4 (35-days)

| Name | p-value | Overlap |
|--|----------|---------------|
| Spliceosomal Cycle | 9.04E-11 | 28.6 % 14/49 |
| Huntington's Disease Signaling | 2.40E-08 | 10.0 % 28/281 |
| Superpathway of Cholesterol Biosynthesis | 1.40E-06 | 27.6 % 8/29 |
| Synaptogenesis Signaling Pathway | 2.68E-06 | 8.3 % 26/314 |
| Cholesterol Biosynthesis I | 2.37E-05 | 38.5 % 5/13 |

Fig. 6. (top) t test on differential proteins. **(bottom)** Top canonical pathways derived from IPA (Integrated Pathway Analysis) of significant proteins from timsTOF Pro acquisition of 21-days samples




Fig. 1. Experimental workflow

Differential proteins (p<0.05)



Mom1 vs Mom2 (21-days) Mom3 vs Mom4 (35-days)

| Name | p-value | Overlap |
|--|----------|--------------|
| Synaptogenesis Signaling Pathway | 5.89E-07 | 5.1 % 16/314 |
| Mechanisms of Viral Exit from Host Cells | 9.21E-05 | 12.2 % 5/41 |
| Spliceosomal Cycle | 2.18E-04 | 10.2 % 5/49 |
| Superpathway of Cholesterol Biosynthesis | 2.95E-04 | 13.8 % 4/29 |
| Cholesterol Biosynthesis I | 3.67E-04 | 23.1 % 3/13 |

Fig. 4. (top) t test on differential proteins. **(bottom)** Top canonical pathways derived from IPA (Integrated Pathway Analysis) of significant proteins from Exploris480 acquisition of 21-days samples

CONCLUSIONS

- ❑ A total of 6,811 and 5,516 unique proteins have been identified from DIA-PASEF (Bruker nanoElute LC/Bruker timsTOF Pro MS) and DIA-FAIMS methods (Waters nanoAcquity UPLC/Thermo Exploris480 MS), respectively.
- ❑ Comparing the AY9944-treated samples versus control, a total of 204 and 460 differentially expressed proteins with a p<0.05 (an overlap of 68 proteins) were identified from 21- and 35-days samples, respectively, in DIA-FAIMS study. On the other hand, 511 and 512 proteins with an overlap of 215 proteins were identified from 21- and 35-days samples, respectively, with a significant p value (p<0.05) from DIA-PASEF study.
- ❑ Top canonical pathways related to differentially expressed proteins in both methodologies includes Superpathway of Cholesterol Biosynthesis, Synaptogenesis Signaling Pathway, Spliceosomal Cycle, and Cholesterol Biosynthesis I which are related to neurodegeneration.
- ❑ Overall, DIA-PASEF data provides better performance compared to DIA-FAIMS in this study.

ACKNOWLEDGEMENT

Matt Willetts, Bruker Scientific LLC, Billerica, MA 01821.