

Rat Liver and Kidney Contrast-Enhanced μ CT Imaging with the Albira CT

Stephanie Prince¹, Sarah Chapman², Todd A. Sasser³, W. Matthew Leevy^{1,2}

¹ Department of Biology, 236 Nieuwland Science Hall, University of Notre Dame, Notre Dame, IN 46556; ² Notre Dame Integrated Imaging Facility, University of Notre Dame, Notre Dame, IN 46556, ³ Bruker Preclinical Imaging, 44 Manning Rd, Billerica, MA, 01821.

Preclinical CT imaging is probably most commonly used in studies of bone anatomy and disease and as an anatomical reference for functional (i.e. PET and SPECT) modalities. Adipose tissue and lung tissue have innate contrast for CT imaging that can also be leveraged for preclinical studies of obesity and pulmonary disease (Sasser et al., 2012; Davison et al., 2013). Most soft tissue anatomy has limited innate CT contrast. However, some studies of soft tissue anatomy can be done using CT contrast agents (Lusic and Grinstaff, 2013).

There are a variety of commercial CT contrast agents. Exitron™ nano agents and Visipaque™ are used in small animal studies for liver and kidney contrast, respectively. Exitron™ nano 6000 and 12000 are alkaline earth metal nanoparticle agents that provide both liver and spleen CT contrast (Boll et al., 2011). Visipaque™ is an iodixanol based CT contrast agent that is approved for clinical arteriography/angiography, though in rodents Visipaque™ readily provides contrast in the kidneys and bladder (Heglund et al., 1995).

Contrast imaging, including liver contrast imaging using Exitron™ nano agents and kidney contrast imaging using Visipaque™, was demonstrated in mice and employing the Albira CT by Wathen et al. (2013). While mouse imaging techniques

are established with a variety of agents (Rothe et al., 2015), there is less information available on contrast imaging for rats, which are up to 10-fold larger by weight. Here, we have identified basic methods to image kidney and liver systems in rats using two (Exitron™ and Visipaque™) commercially available CT contrast agents. These methods may be applied in future studies of liver and kidney disease using rat models.

Methods

Contrast Agents and Imaging

Rats were anesthetized using a mixture of isoflurane gas (3 % in O₂ (v/v)) at a flow rate of 2 L/min for animal preparation and imaging. For liver contrast, three rats were injected via tail vein with a bolus injection of 250 μ L Exitron™ nano 12000 (Miltényi Biotec, Auburn, CA, USA) and imaged immediately after injection and 24 and 48 hours post injection. Albira CT acquisitions were made using the Albira "Good setting" (400 projections, 125 μ m voxels) at high dose (400 μ A) and high voltage (45 kVp).


For kidney imaging three rats were injected via tail vein with a single bolus injection of 750 μL of Visipaque™ (GE Healthcare, USA) 270 mg/ml and imaged immediately after injection, and at intervals for three hours thereafter. Albira CT acquisitions were made using the Albira “Standard setting” (200 projections, 125 μm voxels) setting at high dose (400 μA) and high voltage (45 kVp).

Image Presentation

All images were reconstructed with the Albira Reconstructor Software Advanced option at high resolution with a FBP algorithm. Images were processed in the PMOD (PMOD Technologies Ltd., Zurich, Switzerland) software and VolView (Kitware Inc.) software. 3D image presentation was performed according to the methods detailed by Wathen et al. (2012).


Results/Discussion

Preliminary studies (not shown) were made to identify a general working volume range of Exitron™ nano 12000 and Visipaque™ for contrast imaging in rats. Not surprisingly, a greater working volume is required for rat imaging than for mouse imaging. At 250 μL Exitron™ nano 12000 provided clear contrast for the liver and spleen for at least 48 hours after the initial bolus injection (Figure 1). This protocol provides a foundation for applied liver contrast imaging in rats. Previously Exitron™ agents have been used in μCT imaging mouse liver metastasis models (Boll et al., 2013), and the protocol described here could feasibly be translatable to imaging of rat liver tumor models.



Rat liver and spleen contrast enhanced μCT imaging using Exitron™ nano 12000. A) Liver and spleen are clearly visible in rat injected with Exitron™ nano 12000 and imaged by μCT , here showing sagittal, transverse, and coronal views. B) 3D view of rat liver and spleen contrast enhanced μCT imaging.

A working volume of 750 μL of Visipaque™ was found to provide good contrast for kidney and bladder contrast in rats (Figure 2). The contrast agent cleared from the kidneys completely within 3 hours post injection, though the bladder contrast was visible for the 3 hour duration of imaging studies. The results provide a basis for applied kidney imaging using the Albira CT system. One potential application for this protocol would be for enhanced anatomical CT reference for multi-modal PET and/or SPECT imaging. This could be particularly relevant in disease models where fine anatomical localization related to the kidney or bladder is required. Additionally, a method for evaluating kidney clearance rates in mice using Visipaque™ and X-ray imaging was reported (Orton et al., 2010). Similar analysis for kidney clearance in rats would likely be possible through adoption of the methods reported here, and using the Albira CT 2D dynamic scan setting.



Rat kidney contrast enhanced μCT imaging using Visipaque™. A) Kidneys are clearly visible in rats injected with Visipaque™ and imaged by μCT , here showing sagittal, transverse, and coronal views. In these view, the renal cortex and renal pelvis can be easily contrasted. B) 3D view of rat kidney contrast enhanced microCT imaging. In this view of the lower anatomy the ureter (orange) and bladder are clearly visible.

Conclusion

Rodent liver and kidney tissue has an innately low contrast with μCT imaging. Exitron™ and Visipaque™ contrast agents have been used in μCT studies of mouse liver and kidney disease. Workable volumes and useable timing for both agents to perform μCT contrast imaging in rats and employing the Albira CT could be identified. These methods should be applicable for studies in relevant rat liver, kidney, and bladder disease models.

References

- [1] Sasser TA, Chapman SE, Li S, Hudson C, Orton SP, Diener JM, Gammon ST, Correcher C, Leevy WM. (2012) Segmentation and measurement of fat volumes in murine obesity models using X-ray computed tomography. *J Vis Exp.* 4;(62):e3680.
- [2] Davison CA, Chapman SE, Sasser TA, Wathen C, Diener J, Schafer ZT, Leevy WM. (2013) Multi-modal optical, X-ray CT, and SPECT imaging of a mouse model of breast cancer lung metastasis. *Curr Mol Med.* 13(3):368-76.
- [3] Lusic H, Grinstaff M W (2013) X-ray-Computed Tomography Contrast Agents. *Chem. Rev.* 113(3): 1641-1666.
- [4] Boll H, Nittka S, Doyon F, Neumaier M, Marx A, Kramer M, Groden C, Brockmann MA (2011) Micro-CT Based Experimental Liver Imaging Using a Nanoparticulate Contrast Agent: A Longitudinal Study in Mice. *PLoS ONE* 6(9).
- [5] Heglund, IF, Michelet, AA, Blazak WF. Furuhashi, K, Holtz, E. (1995) Preclinical pharmacokinetics and general toxicology of iodixanol. *Acta Radiol Suppl.* 399:69-82.
- [6] Wathen CA, Foje N, van Avermaete T, Miramontes B, Chapman SE, Sasser TA, Kannan R, Gerstler S, Leevy WM. (2013) In vivo X-ray computed tomographic imaging of soft tissue with native, intravenous, or oral contrast. *Sensors (Basel).* 27;13(6):6957-80.
- [7] Rothe JH, Rudolph I, Rohwer N, Kupitz D, Gregor-Mamoudou B, Derlin T, Furth C, Amthauer H, Brenner W, Buchert R, Cramer T, (2015) Time course of contrast enhancement by micro-CT with dedicated contrast agents in normal mice and mice with hepatocellular carcinoma: comparison of one iodinated and two nanoparticle-based agents. *Acad Radiol.* 22(2):169-78.
- [8] Wathen C, Sasser TA, Chapman SE, Leevy MW (2012) Application protocol for small animal lung volume segmentation employing the Albira CT and PMOD software. Bruker AP0121. Bruker Molecular Imaging.
- [9] Boll H, Figueiredo G, Fiebig T, Nittka S, Doyon F, Kerl HU, Nölte I, Förster A, Kramer M, Brockmann MA. (2013) Comparison of Fenestra C, ExiTron nano 6000, and ExiTron nano 12000 for micro-CT imaging of liver and spleen in mice. *Acad Radiol.* 20(9):1137-43.
- [10] Orton S, Leevy WM, Feke G, McLaughlin W (2010) Dynamic X-ray imaging of iodinated contrast-agent clearance from a living mouse. AP0101. Bruker Molecular Imaging.