# **Reproducibility of MALDI Imaging Based Tissue Classifications – Results of a Multi-Center Study**

### **ASMS 2019 WP-373**

Sören-Oliver Deininger<sup>1</sup>, Rita Casadonte<sup>2</sup>, Petra Wandernoth<sup>2</sup>, Kristina Schwamborn<sup>3</sup>, Christine **Bollwein<sup>3</sup>, Christian Marsching<sup>4</sup>, Katharina** Kriegsmann<sup>5</sup>, Carsten Hopf<sup>4</sup>, Wilko Weichert<sup>3</sup>, Jörg Kriegsmann<sup>2</sup>, Peter Schirmacher<sup>6</sup>, Mark Kriegsmann<sup>6</sup>, Alice Ly<sup>1</sup>

<sup>1</sup>Bruker Daltonik GmbH, Bremen, Germany <sup>2</sup>Proteopath GmbH, Trier, Germany <sup>3</sup>Institute of Pathology, TU Munich, Germany <sup>4</sup>CeMOS, Mannheim University of Applied Sciences, Mannheim, German <sup>5</sup>Department of Hematology, Oncology and Rheumatology, University Hospital Heidelberg, Germany <sup>6</sup>Institute of Pathology, University Hospital Heidelberg, Germany

### Introduction

Classification of tissues based on label-free mass spectrometric phenotypes measured directly from sections is a promising tool for clinical research. However, reproducibly measuring mass spectra can be challenging and comprehensive studies assessing the variation across different sites are largely lacking. In this work we have compared the reproducibility of Matrix-Assisted-Laser-Desorption/Ionization MALDI mass spectrometric imaging (MALDI-MSI) based tissue classifications measured at three different sites.

### Methods

Mouse intestine and a tissue microarray (TMA) containing samples from 95 subjects (6 tumor types sampled at 3 sites) was used (layout shown in figure 2). Samples were sectioned onto conductive glass slides (Bruker Daltonik GmbH), and underwent deparaffinization and antigen retrieval. Serial sections were prepared and measured as previously published [1] at different sites.

Samples were sprayed with trypsin and incubated. After digestion, matrix was applied and sections were measured at 50 µm step size with a rapifleX MALDI Tissuetyper TOF mass spectrometer (Bruker). Data were analyzed using flexImaging (Bruker), SCiLS Lab Pro (Bruker) and R (R-project).

The classification was calculated by linear discriminant analysis (R-package MASS). The performance was estimated as the accuracy of the classifications with different crossvalidation-scenarios: Leave-one-samplingsite-out, Leave-one-TMA-out or a two-step Leave-one-TMA-out-leave-one-patient-out.

### Results

Spatial segmentation on mouse gut samples measured at five sites and two time points showed that the spectra were clustered



using annotations and feature selection is according to tissue type, not according to site or time of measurement (Figure 1). For shown in figure 3. In this case the classification result was 84.0%. the classification of the TMA data, the TMA was prepared and measured at three sites. A total of 407 monoisotopic mass spectral features were found in the average spectrum. A forward feature selection was Mantle Cell Lymphoma used to identify the 25 most relevant Squamous Cell Carcinoma of the Lung features for the classification. Breast Cancer

The histological analysis of the H&E stained TMAs showed that some cores were heterogenous with only small actual tumor areas (Figure 3). Some cores contained no tumor. The data were analyzed with and without histologic annotation of the tumor areas and feature selection taken into account. The results are shown in table 1.

A classification result of a leave-one-TMAout-leave-one-subject out classification



**Figure 1**. Segmentation analysis of mouse intestine samples measured at five sites over three time points. (Scale: 3mm)

Ac
92
73
79
78
84

 
 Table 1. Accuracy of
classification with and without cross validation and feature selection.

9.1% 8.1% 4.5%



#### **Figure 2**. Layout of the Multi-Tumor-TMA



Figure 3. Histological annotation of tumor regions for one TMA section and classification map of a Leave-One-TMA-Out-Leave-One-Subject-Out Cross validation

### Summary

We have shown that MALDI imaging of FFPE tissues can be performed reproducibly across different sites, operators and instruments. Although TMAs contain pre-selected tissue, a detailed histological annotation is necessary to obtain optimal results. Restricting the number of mass signals used for classification improves the performance of tissue based classifiers on new data. MALDI imaging based classifiers are able to generalize across sites.

References

(1) Ly et al.; Proteomics Clin. Appl. 2019, 13, 1800029

### Conclusions

- classifiers
- data

MALDI Imaging on FFPE-Tissue can be performed reproducibly

Tissue classifiers are able to generalize across sites

Feature selection improves the performance of

Histological annotations are necessary even on TMA

## MALDI Tissuetyper