

TXRF analysis of cell culture media and medical microsamples – saving time, money and lives



Bruker Nano Analytics, Berlin
Webinar, January 18th, 2018



Welcome



Speakers

Dr. Hagen Stosnach
Applications Scientist TXRF
Berlin, Germany



Dr. Armin Gross
Global Product Manager TXRF
Berlin, Germany



Part I: Cell cultures

- Cell cultures – background
- Method development for element analysis of cell culture media
- TXRF analysis of high performance cell culture media
- Summary and conclusion

Part II: Biological and medical microsamples

- Protein analysis
- Bioassay analysis
- Summary and outlook

Part I: Cell cultures

Background

Cell cultures

Background



Application of cell cultures

- Biopharmaceuticals
pharmaceutically active proteins and nucleic acids
- Clinical research: tissue cells for therapeutic purposes
- Study of elementary cellular processes without classical animal experiments



Cell cultures

Background



Cell culture media

- Natural environment of cells has to be imitated
- Standardized media, e.g. DMEM, 1 g/l or 4.5 g/l glucose
- Media contain amino acids for protein biosynthesis, vitamins and inorganic salts for metabolic processes, carbohydrates for energy production
- Nutrient (P, K, Ca, Zn, Fe) and trace elements (Se, Mn) nutrient sources, co-factors in enzymatic reactions
- Contaminations reduce yield of the biopharmaceuticals
may cause necrosis of the cells



Cell culture media

Objective



Objective 1

- Development of a TXRF method for cell culture media
 1. Optimize sample layer on disc, standard DMEM media
 2. Determination of LOD / LOQ
 3. Reproducibility

Objective 2

- Measurement of nutrient elements in cell culture media

Objective 3

- Apply method to modern high performance media
 1. QC: identify contaminants
 2. Spike experiments for method detection limits

Cell culture media

Objective



Samples

- Commercially available DMEM, 1 g/l or 4.5 g/l glucose
- Mammalian media from two German suppliers
Sugars conc. 8 g/l
 - Human Embryonic Kidney (HEK)
 - Chinese Hamster Ovary (CHO)
 - Baby Hamster Kidney (BHK)
 - Invertebrate (IKZ)
 - Feed media (20 g/l glucose)

Cell culture media TXRF spectrometer



S4 T-STAR - Unique benefits

- Three excitation modes to detect most elements of the PSE
- 60/100 mm² detectors improved sensitivity for lowest limits of detection
- New analytical capabilities angle scan for depth profiling and layer analysis
- Sample geometry flexibility measurement of discs, microscopy slides, wafers etc.
- Motorized beam path automatic beam adjustment and QC procedures
- Large sample capacity up to 90 sample discs, multi-user operation
- Most modern software instrument/measurement status display, statistical functions





Method development for element analysis of cell culture media

Cell culture media

Method development



1. Optimize sample layer on disc

Comparison of different modifier solutions
(smoothing agents)

- DMEM media with 50 – 200 ppm P, S, K, Ca
- Addition of **10% modifier**
(1% Triton X-100, Ethanol, PVA (0,3 g/l))
- Dilutions down to 1 : 100
- Internal standard Sc

Measurement parameter S4 T-STAR

- Mo excitation, 50 kV, 1000 μ A
- W-L excitation, 50 kV, 1000 μ A
- Measurement time 1000 s


Cell culture media

Method development



Results

- Ethanol leads to satisfying recovery rates and standard deviations < 5% (< 10% at 1:100)
- Overlap of IS Sc with Ca, following experiments with V
- Mo excitation did not provide good data for P and S




Cell culture media

Method development



Results

- W-L excitation delivers reproducible results for P and S, if Ethanol was used




Cell culture media

Method development



2. Reproducibility

- Two preparations
 - 700 μl sample + 10 μl IS + 290 μl EtOH
 - 100 μl sample + 100 ml EtOH + 790 μl H₂O + 10 μl IS
- Better than 10%, typically < 3%
- Dilution recommended due to high amount of macro nutrient
- EtOH contained trace impurities (Fe) -> ultrapure or HPLC grade required



Cell culture media

Method development



3. LOQ


- EtOH treated media

	LOQ (mg/l)
Magnesium*	122,7
Phosphor	4,165
Sulfur	1,205
Potassium	0,200
Calcium	0,122
Iron	0,016

*) W-L excitation

- Sample preparation

- 700 µl media
- 290 µl EtOH
- 10 µl V (IS)





TXRF analysis of high performance cell culture media

Modern high performance media

Quality control



Analysis of mammalian cell culture media

- Method developed for standard media
- 24 media in triplicate
= 72 discs

Measurement parameter S4 T-STAR

- Mo excitation, 50 kV, 1000 μ A
- W-L excitation, 50 kV, 1000 μ A
- W-Brems excitation, 50 kV, 1000 μ A
- Measurement time 1000 s

Ethanol	290 μ l
IS Vanadium	10 μ l
Sample	700 μ l

Volume on disc	10 μ l
Drying	30° C


Modern high performance media

Quality control



Results mammalian cell culture media

- Micro-nutrient and contamination test of different media batches
- Application of Mo and W-L excitation
- Typical concentration range 1 to 1000 ppb
- Distinct differences between cell lines



CHO	Chinese Hamster Ovary
HEK	Human Embryonic Kidney
IZK	Invertebrate cell cultures
Feed	Feed media


Modern high performance media

Quality control



Limit of detection

- LOD down to single digit ppb
- W-L significantly improves the detection of light elements (factor 2 - 4)
- Recommended measurement time for routine analysis = 300 s



Modern high performance media

Spike experiments



Preparation of spike media

- 1 ml of one high performance medium
- Two spike elements per sample concentrations see table
- Measurement in triplicate
15 samples = 45 discs

Spiked element concentrations ($\mu\text{g/l}$)					
Pb	Mn	Se	Ni	Cr	Cd
5	20	5	10	15	300
10	50	10	20	30	500
20	100	20	50	50	1000
50	200	50	100	100	2000
100	500	100	200	200	4000

Measurement parameter S4 T-STAR

- Mo excitation, 50 kV, 1000 μA
- W-L excitation, 50 kV, 1000 μA
- W-Brems excitation, 50 kV, 1000 μA
- Measurement time 1000 s


Modern high performance media

Spike experiments




Mo spectra

- Cr spike (200 ppb)



W-Brems spectra

- Cd spike (4 ppm)



Modern high performance media

Spike experiments





Measurements with Mo excitation

ICP	2,1 µg/l
TXRF spike	7,5 µg/l
TXRF direct	5,3 µg/l

- Very low concentrations of Ni, Se
- Values close to ICP

ICP	5,2 µg/l
TXRF spike	3,6 µg/l
TXRF direct	3,4 µg/l



Modern high performance media



Spike experiments



Measurements with Mo excitation

ICP	n.a.
TXRF spike	17 µg/l
TXRF direct	7,6 µg/l

ICP	2,6 µg/l
TXRF spike	4,4 µg/l
TXRF direct	15,8 µg/l




Modern high performance media

Spike experiments




Measurements with W-L and W-Brems excitation

ICP	1,9 µg/l
TXRF spike	24,7 µg/l
TXRF direct	14,7 µg/l



- ICP values for Cr questionable
- Cd quantification not successful

ICP	< 0,2 µg/l
TXRF spike	< 45 µg/l
TXRF direct	< 77,3 µg/l



Summary and conclusion



- A rapid method for TXRF measurements of cell culture media was developed
- EtOH has to be used as smooting agent for optimal layer formation on quartz sample discs
- For contamination control detection limits in the one digit ppb range can be achieved
- The measurement of nutrient elements after a 1:10 dilution provides a reproducibility < 10%



Part II: Biological and medical microsamples

Protein analysis


Introduction



Metal ions and enzymes

- Metal ions are important for the biological function of enzymes
- Various modes of metal-protein interaction: metal-, ligand-, enzyme-bridge complexes
- Metals serve as electron donors or acceptors, Lewis acids or structural regulators

(Riordan JF.: "The role of metals in enzyme activity.", Ann Clin Lab Sci. 1977 Mar-Apr; 7(2):119-29



C. Arnaud: CEN, January 7, 2008
Volume 86, Number 1, p. 8

Protein analysis

Introduction



Examples of metal ions in enzymes

Role	Metals	Protein
Oxygen transport and storage	Fe, Cu	Haemoglobin
Electron transport	Fe, Cu	Cytochromes
Nitrogen fixation	Fe, Mo, V	Nitrogenase
Oxygen atom transfer	Mo, W	Oxidases, Reductases
Alkyl group transfer	Co	Vitamin B ₁₂
Hydrolysis	Zn, Cu, Mn	Hydroxylases, Peptidases
Storage and transport	Fe, Cu, Zn	Ferritin, Metallothioneins

Protein analysis

Introduction



Commonly used analytical techniques


- Atomic Absorption Spectrometry
 - Ca, Co, Cu, Fe, Mg, Mo, Ni, Se, Zn
- Flame/Inductively Coupled Emission Spectrometry
 - for most metals occurring in proteins
- Neutron activation analysis
 - limited availability
 - most reliable technique for a number of elements
 - very sensitive for some elements
 - some critical interferes, e.g. Na
- Inductively Coupled Mass Spectrometry
 - detection of ppb levels of >40 elements in one minute

Protein analysis Introduction



Analytical issues

- High matrix (polypeptides, buffers, salts)
- Sample viscosity, turbidity
- Removal of buffer, salts by dialysis or gel filtration
total destruction of organic matter
- Significant sample amount needed
microanalysis impossible
- Methods are expensive and laborious



Mounicou et al. (2004), *Analyst*, (2), 116-123

Protein analysis

TXRF spectroscopy



Samples

Certified reference materials

- "BCR 273 – Single cell protein"
- "BCR 274 – Single cell protein"

500 mg powderous sample were suspended in 25 ml pure water



Protein analysis

TXRF spectroscopy




Sample preparation

First approach:
mixing of microsamples in **vials**

- 100 μ l sample + internal standard
+ 10 μ l Sc (10 mg/l)
+ 10 μ l Ga (10 mg/l)



- Preparation of
 - 10 μ l
 - 5 μ l
 - 3 μ l
 - 1 μ l
 - 0,5 μ l



Protein analysis


TXRF spectroscopy



Sample preparation

Second approach:
mixing of microsamples on **discs**

- 5 μl sample + 5 μl Sc/Ga solution (20 mg/l)
- 3 μl sample + 3 μl Sc/Ga solution (20 mg/l)
- 1 μl sample + 1 μl Sc/Ga solution (20 mg/l)



Protein analysis

TXRF spectroscopy



S4 T•STAR

- Mo tube, 50 kV/1000 μ A
- W-tube, 50 kV/1000 μ A
- 60 mm² XFlash SDD
- 90 position sample changer
- Mo-K excitation, 1000 s
- W-L excitation, 1000 s
- W-Brems, 1000 s




Protein analysis

TXRF spectroscopy




Results mixing in vials



Protein analysis TXRF spectroscopy



Results mixing in vials




Protein analysis

TXRF spectroscopy



Results mixing in vials

Other major and trace elements (K, Ca, Mn, Fe, Cu, Zn and Se) are not affected




Protein analysis

TXRF spectroscopy



Results mixing in vials

- The detection limits are directly related to the deposited amount




Protein analysis

TXRF spectroscopy




Results mixing on discs



Protein analysis TXRF spectroscopy



Results mixing on discs




Protein analysis

TXRF spectroscopy



Results mixing on discs

Other major and trace elements (K, Ca, Mn, Fe, Cu, Zn and Se) are not affected




Protein analysis

TXRF spectroscopy



Results mixing on discs

- The detection limits are directly related to the deposited amount



Bioassay analysis

Introduction



Saving lives...?

Not directly theirs...



Bioassay analysis

Introduction



Saving lifes...?

But theirs...



Bioassay analysis

TXRF spectroscopy




Sample preparation

Sampling and preparing non-lethal amounts of blood (urine, saliva ...)

20 μ l sample
+ 10 μ l Ga (4 mg/l)



Preparation of
• 10 μ l (duplicate)




Bioassay analysis

TXRF spectroscopy



Results

- Seronorm human serum CRM, S2 PICOFOX





Bioassay analysis

TXRF spectroscopy



Results

- Serum samples from lab mice (publication in preparation)



Summary and Outlook



The analysis of minute protein or serum samples by means of TXRF is possible

- Mixing of internal standard and sample in a vial is recommended
- Determination of detection limits for trace elements requires careful consideration of the sample amount

Outlook


- Intensive analysis of bioassays is ongoing
- Measurements with W-Brems excitation (S4 T•STAR) to detect elements like Mo, I
- Method development for light elements with optimized calibration factors


Q & A



Any Questions?

Please **type in** the questions you may have for our speakers in the **Questions Box** and click **Submit**





Innovation with Integrity

Thank you for your attention!

www.bruker.com / www.s4tstar.com