


Quickly analyse nutrients, additives, and contaminants in food, feed, supplements, and beverages with PXRF and TXRF



Kimberley Russell, Esa Nummi, Hagen Stosnach, Armin Gross
Bruker Nano GmbH, Berlin, Germany
June, 2020



Welcome



Speakers

Kimberley Russell
Food Safety, Agriculture and Environmental
Market Segment Manager, HMP/XMA



Esa Nummi
Director Product Management & Business
Development HMP/XMA



Dr. Hagen Stosnach
Application Scientist TXRF
Berlin, Germany




Dr. Armin Gross
Global Product Manager TXRF
Berlin, Germany



Itinerary



- Introduction
- PXRF and TXRF – How does it work?
- Measurements with portable XRF
 - Measurement of mineral nutrients in animal feed
 - Elements monitored in processing chocolate and edible oil
 - Physical contaminants in human and pet food
- Ultra-trace analysis with TXRF
 - micronutrients and toxic elements in solid food samples
 - dietary supplement pills
 - milk samples
- Q & A







PXRF and TXRF

How does it work?


PXRF Technology




PXRF Technology



PXRF Technology



PXRF Technology




24	52.00	26	55.85	28	58.69	42	95.94
Cr	7.15	Fe	7.87	Ni	8.91	Mo	10.22
Chromium		Iron		Nickel		Molybdenum	
K α 5.415		K α 6.405		K α 7.480		K α 17.480	
L α 0.572		L α 0.705		L α 0.849		L α 2.292	

316SS				
124 Match 9.8 11-24 17:26				
Time 46.0 HiZ SS Fe				
El	Min	%	Max	+/- [*2]
Fe	60.00	70.40	75.00	1.07
Cr	16.00	16.25	18.00	0.45
Ni	10.00	11.40	14.00	0.50
Mo	2.00	1.63	3.00	0.10
Sn		0.21		0.09
Nb		0.11		0.06
Ti	0.00	< LOD	0.00	0.14
Mn	0.00	< LOD	2.00	0.24



Spectrum

Total reflection X-ray fluorescence spectroscopy



- Samples must be prepared on a reflective media
- Dried to a thin layer
- Matrix effects are negligible
- Quantification by internal standardization

Total reflection X-ray fluorescence spectroscopy



- Samples must be prepared on a reflective media
- Dried to a thin layer
- Matrix effects are negligible
- Quantification by internal standardization

$$C_i = \frac{C_{IS} \cdot N_i \cdot S_{IS}}{N_{IS} \cdot S_i}$$

Elements of interest in food analysis



Quickly analyze nutrients, additives, and contaminants in food, feed, supplements, and beverages with portable XRF and T-XRF

- Macro mineral nutrients: Ca, P, Mg, Na, K, Cl, S
- Trace mineral nutrients: Fe, Mn, Cu, I, Zn, Co, F, Se
- Toxic elements: As, Cd, Hg, Pb
- Processing elements: P, Fe, and others
- Physical contaminants: bits of metal, ceramic, glass and plastic



S4 T-STAR TXRF for ultra low detection limits of elements from sodium to uranium. High volume sample analysis is possible with stackable trays.



TRACER 5i battery operated, handheld XRF for analysis of elements from sodium to uranium wherever and whenever needed.



Measurement of mineral nutrients in
animal feed with a portable XRF

Application




- Macromineral nutrients, such as calcium (Ca) and phosphorus (P), are just as important in animal food as they are in human food.
- The optimum ratio of Ca to P (Ca:P) in feed is particularly critical for laying hens because it helps balance support for skeletal growth, eggshell formation, and immune responses.
- Portable XRF enables the measurement of these and other elements where feed materials are received, at feeding sites, or in the lab.

Methods and equipment

Simple steps to measure elements in feed with a portable XRF

1. Create or purchase an elemental calibration for feed to be installed on the PXRF using reference samples with a similar base matrix and concentrations covering the range of interest.



2. Put crushed feed in sample cup and cover with film.



3. Place sample cup on XRF analysis window and close the cover.

4. Select the calibration and press start.



Results

View results on screen as elements of interest, full composition, P/F, or spectra. Results are stored in the instrument, but can be saved onto a USB stick or transferred to another data collection system.



Food QualityDual
424 03-25 08:59
Time 30.0

El	%	+/- [*2]
Ca	0.48	0.00
P	0.30	0.01

Use in Average

Food QualityDual
429 03-25 09:07
Time 24.0

El	%	+/- [*2]
Mg	0.67	0.42
Si	0.35	0.02
P	0.30	0.01
S	0.20	0.01
K	0.93	0.01
Ca	0.48	0.01
Cr	0.01	0.01
Mn	0.01	0.01
Fe	0.05	0.00



Application

- Salt is a critical component of a nutritionally balanced diet for animals.
- Insufficient amounts of sodium can lead to serious health and behavioral issues as well as to a decrease in feed utilization.
- Recommended amounts of salt vary depending on the animal, its location, and its activity.
- Portable XRF enables the measurement of sodium and other elements where feed materials are received, at feeding sites, or in the lab.



Methods and equipment

- Sodium is a very light (low energy) element which presents a challenge for PXRF analysis.
- The TRACER 5 handheld XRF can perform Na measurements with the use of a helium path and the removal of windows and films between the analyzer and the sample.
- These steps increase the transmission of light element X-rays which in turn increases the sensitivity of the measurement enabling the analysis of elements as light as Na.



Methods and equipment

- The TRACER 5 platform provides the ability to control the excitation conditions.
- Users can adjust the current and voltage directly on the interactive touch screen display.
- Users can select filters from the integrated filter wheel or insert user-designed filters.
- Users can select one of three measurement paths – air, vacuum, or helium.
- TRACER 5 control and live spectra analysis via PC is also possible.



Methods and equipment

- Good working detection limits of light elements can be achieved with the TRACER 5g (1 μm graphene detector window) or the TRACER 5i (8 μm beryllium detector window).

TRACER 5g		
Element	LOD (PPM)	Sensitivity (Counts/PPM)
Na	312	0.62
Mg	122	2.13
Al	134	4.41
Ca	24	14.04
Fe	50	25.87


TRACER 5i		
Element	LOD (PPM)	Sensitivity (Counts/PPM)
Na	828	0.20
Mg	185	1.17
Al	177	3.08
Ca	22	16.10
Fe	48	30.03

- Selection of a helium flow path and the removal of the instrument analysis window are required.
- Samples should be put in a cup without any film covering and placed under the tube (nose down) to prevent debris from falling inside.



Methods and equipment

- Users can create or purchase a custom calibration to be installed on the TRACER for Na analysis, using reference samples with a similar base matrix and with concentrations covering the range of interest.
- The measurement conditions must be the same for reference and unknown samples, including the use of helium and the preparation and presentation of the samples.
- Qualitative identification or screening tests for Na in feed are also possible without the need of an installed calibration curve.



Correlation for measuring Na using balloon He atmosphere with a handheld XRF



Elements monitored in processing
chocolate and edible oil

Application




- Iron (Fe) is an important mineral nutrient measured and reported in cocoa.
- It's also important to monitor Fe during production for taste consistency. In oxidative reactions, some iron compounds can result in off-flavors, colors, or odors.
- Iron is also monitored during production to prevent accidental Fe contamination from process equipment.
- Portable XRF enables the measurement of Fe and other elements where starting materials are received, on the production line, or in the lab.

Methods and equipment

Simple steps to measure Fe in cocoa with a portable XRF

1. Create or purchase a calibration for Fe in cocoa to be installed on the PXRF using reference samples with a similar base matrix and concentrations covering the range of interest.
2. Put cocoa powder in sample cup and cover with film.
3. Place sample cup on XRF analysis window and close the cover.
4. Select the calibration and press start.



Results

View results on screen as elements of interest, full composition, P/F, or spectra. Results are stored in the instrument, but can be saved onto a USB stick or transferred to another data collection system.



83 BChoc PASS (9.9)				
Time 20.0				
El	Min	PPM	Max	+/- [...]
Fe	0	60	100	4
Zn		37		2
Mn		43		24
Cu		23		3
Pb		< LOD		1

<	<input type="checkbox"/> Use in Average	>
Averaging		Calculate Average
Spectrum	Edit Info	Back

82 Cocoa PASS (9.9)				
Time 20.0				
El	Min	PPM	Max	+/- [...]
Fe	0	81	100	6
Zn		77		2
Mn		57		24
Cu		35		3
Pb		2		1

<	<input type="checkbox"/> Use in Average	>
Averaging		Calculate Average
Spectrum	Edit Info	Back

Application





- Phosphorus is not an edible oil nutrient; however, its presence is important to monitor during processing.
- Measuring phosphorus helps monitor the refining of edible oils, especially important during the degumming process of seed oils, the settling of fruit oils, and the neutralization of all edible oils.
- Phosphorus measurements are indicative of phosphatide content which ultimately helps determine edible oil quality.
- Portable XRF enables the measurement of phosphorus and other elements where starting materials are received, on the production line, or in the lab.

Methods and equipment

Simple steps to measure Phosphorus in edible oil with a portable XRF

1. Create or purchase a calibration for Phosphorus in oil to be installed on the PXRf using reference samples with a similar base matrix and concentrations covering the range of interest.
2. Put oil in sample cup and cover with film.
3. Place sample cup on XRF analysis window and close the cover.
4. Select the calibration and press start.




Results

View results on screen as elements of interest, full composition, P/F, or spectra. Results are stored in the instrument, but can be saved onto a USB stick or transferred to another data collection system.



Food QualityDual			
465 03-26 12:50			
Time 13.0			
EI	PPM	+/- [*2]	
P	231	12	
<input type="checkbox"/> Use in Average			
Averaging		Calculate Average	
Spectrum	Edit Info	Back	

Food QualityDual			
462 03-26 12:46			
Time 18.0			
EI	PPM	+/- [*2]	
P	232	16	
K	87	1	
Ca	53	1	
Fe	44	4	
S	19	1	
Cr	15	7	
Mn	10	5	
Cu	2	1	
			
Spectrum			

Physical contaminants in human and pet food

Application

Illustration of baking production line equipment




Physical contaminants can enter food products from wear and tear of equipment



Metal



Stone and ceramic



Glass



Plastic and rubber

Application

Typical prevention



- 1. **Inspect** food material with X-ray or metal detector
- 2. **Detect** foreign body
- 3. **Remove** product containing foreign body

HHXRF for enhanced prevention



- 1. **Identify** found contaminant to determine its source
- 2. **Correct** issue based on source of the contaminant




Auger



Conveyor



Grinder/cutter



Mixer



Roller mill



Sorter/packager

Methods and equipment



- Small physical contaminants found in food products can be identified on the production floor or in the lab with handheld XRF
- A found metal object can be ID'd with handheld XRF by quickly determining its composition, the type of metal alloy, and the alloy grade name.

Methods, equipment results

Simple steps to ID small metal contaminants with Handheld XRF

1. Prepare found object



2. Select "Alloys" Cal



3. Position found object



4. Test & View Results





TXRF analysis of micronutrients and toxic elements in solid food samples

Nutrition-relevant elements, also called "minerals"

- Minerals are inorganic substances required by the body in small amounts for a variety of different functions.
 - involved in the formation of bones and teeth
 - essential constituents of body fluids and tissues
 - components of enzyme systems
 - involved in normal nerve function
- The body requires different amounts of each mineral; requirements depend on age, sex, physiological state (e.g. pregnancy) and state of health
- Nutrition-relevant elements are calcium, phosphorus, magnesium, sodium, potassium, chlorine, iron, zinc, copper, selenium and iodine
(Source: British Nutrition Foundation)



Toxic elements

- Because of the high consumption of food the control of toxic metals is also of high importance
- Limit values for As, Cd, Pb and Hg are applied

Heavy metals	Stated limit	Calculated daily
	(PTWI, weekly)	limit (adult, 70 kg)
Arsenic	15 µg inorganic arsenic/kg bw	150 µg
Cadmium	7 µg cadmium/kg bw	70 µg
Lead	25 µg lead/kg bw	250 µg
Mercury	1,6 µg methylmercury/kg bw	16 µg

Source: JEFCA: The Joint FAO/WHO Expert Committee on Food Additives;
PTWI: provisional tolerable weekly intake

- The control of other toxic metals like Cr, Ni, Co and U would be beneficial

Samples



- Round robin Cereal sample PTNATIAEA16 (<http://www.pt-nsil.com/>)
 - Nominal values from round robin test, applying different analytical techniques
 - Partly high number of outliers for the nominal values
- Certified reference standard NIST 1568a (rice flour)
 - Reference values based on measurements with high number of different analytical methods and laboratories

Methods and Equipment

Sample preparation



Sample preparation



Methods and Equipment Measurements



S2 PICOFOX

- Mo tube, 50 kV/1000 μ A
- 60 mm² XFlash SDD
- 25 position sample changer

Measurements

- Mo-K excitation, 1000 s



Methods and Equipment Measurements



S2 PICOFOX - Measurement program

1 H Hydrogen																	2 He Helium				
3 Li Lithium	4 Be Beryllium															5 B Boron	6 C Carbon	7 N Nitrogen	8 O Oxygen	9 F Fluorine	10 Ne Neon
11 Na Sodium	12 Mg Magnesium															13 Al Aluminium	14 Si Silicon	15 P Phosphorus	16 S Sulphur	17 Cl Chlorine	18 Ar Argon
19 K Potassium	20 Ca Calcium	21 Sc Scandium	22 Ti Titanium	23 V Vanadium	24 Cr Chromium	25 Mn Manganese	26 Fe Iron	27 Co Cobalt	28 Ni Nickel	29 Cu Copper	30 Zn Zinc	31 Ga Gallium	32 Ge Germanium	33 As Arsenic	34 Se Selenium	35 Br Bromine	36 Kr Krypton				
37 Rb Rubidium	38 Sr Strontium	39 Y Yttrium	40 Zr Zirconium	41 Nb Niobium	42 Mo Molybdenum	43 Tc Technetium	44 Ru Ruthenium	45 Rh Rhodium	46 Pd Palladium	47 Ag Silver	48 Cd Cadmium	49 In Indium	50 Sb Antimony	51 Sn Tin	52 Te Tellurium	53 I Iodine	54 Xe Xenon				
55 Cs Cesium	56 Ba Barium	57 La Lanthanum	72 Hf Hafnium	73 Ta Tantalum	74 W Tungsten	75 Re Rhenium	76 Os Osmium	77 Ir Iridium	78 Pt Platinum	79 Au Gold	80 Hg Mercury	81 Tl Thallium	82 Pb Lead	83 Bi Bismuth	84 Po Polonium	85 At Astatine	86 Rn Radon				
87 Fr Francium	88 Ra Radium	89 Ac Actinium																			
L Lanthanides	58 Ce Cerium	59 Pr Praseodymium	60 Nd Neodymium	61 Pm Promethium	62 Sm Samarium	63 Eu Europium	64 Gd Gadolinium	65 Tb Terbium	66 Dy Dysprosium	67 Ho Holmium	68 Er Erbium	69 Th Thulium	70 Yb Ytterbium	71 Lu Lutetium							
Ac Actinides	90 Th Thorium	91 Pa Protactinium	92 U Uranium	93 Np Neptunium	94 Pu Plutonium	95 Am Americium	96 Cm Curium	97 Bk Berkelium	98 Cf Californium	99 Es Einsteinium	100 Fm Fermium	101 Md Mendelevium	102 No Nobelium	103 Lr Lawrencium							



Analysed using K-lines



Internal standard



Analysed using L-lines

Methods and Equipment Measurements



S4 T-STAR

- Mo tube, 50 kV/1000 μ A
- W-tube, 50 kV/1000 μ A
- Monochromator system for Mo-K, W-L and W-Brems monochromatization
- 60 mm² XFlash SDD
- 90 position sample changer

Measurements

- Mo-K excitation, 1000 s
- W-L excitation, 1000 s
- W-Brems excitation, 1000 s



Methods and Equipment Measurements



S4 T-STAR – Measurement program

1 H Hydrogen																	2 He Helium
3 Li Lithium	4 Be Beryllium											5 B Boron	6 C Carbon	7 N Nitrogen	8 O Oxygen	9 F Fluorine	10 Ne Neon
11 Na Sodium	12 Mg Magnesium											13 Al Aluminium	14 Si Silicon	15 P Phosphorus	16 S Sulphur	17 Cl Chlorine	18 Ar Argon
19 K Potassium	20 Ca Calcium	21 Sc Scandium	22 Ti Titanium	23 V Vanadium	24 Cr Chromium	25 Mn Manganese	26 Fe Iron	27 Co Cobalt	28 Ni Nickel	29 Cu Copper	30 Zn Zinc	31 Ga Gallium	32 Ge Germanium	33 As Arsenic	34 Se Selenium	35 Br Bromine	36 Kr Krypton
37 Rb Rubidium	38 Sr Strontium	39 Y Yttrium	40 Zr Zirconium	41 Nb Niobium	42 Mo Molybdenum	43 Tc Technetium	44 Ru Ruthenium	45 Rh Rhodium	46 Pd Palladium	47 Ag Silver	48 Cd Cadmium	49 In Indium	50 Sn Tin	51 Sb Antimony	52 Te Tellurium	53 I Iodine	54 Xe Xenon
55 Cs Cesium	56 Ba Barium	57 La Lanthanum	72 Hf Hafnium	73 Ta Tantalum	74 W Tungsten	75 Re Rhenium	76 Os Osmium	77 Ir Iridium	78 Pt Platinum	79 Au Gold	80 Hg Mercury	81 Tl Thallium	82 Pb Lead	83 Bi Bismuth	84 Po Polonium	85 At Astatine	86 Rn Radon
87 Fr Francium	88 Ra Radium	89 Ac Actinium															
L Lanthanides			58 Ce Cerium	59 Pr Praseodymium	60 Nd Neodymium	61 Pm Promethium	62 Sm Samarium	63 Eu Europium	64 Gd Gadolinium	65 Tb Terbium	66 Dy Dysprosium	67 Ho Holmium	68 Er Erbium	69 Th Thulium	70 Yb Ytterbium	71 Lu Lutetium	
Ac Actinides			90 Th Thorium	91 Pa Protactinium	92 U Uranium	93 Np Neptunium	94 Pu Plutonium	95 Am Americium	96 Cm Curium	97 Bk Berkelium	98 Cf Californium	99 Es Einsteinium	100 Fm Fermium	101 Md Mendelevium	102 No Nobelium	103 Lr Lawrencium	

Impossible to analyse

Analysed using K-lines (Mo-K excitation)

Analysed using K-lines (W-Brems excitation)

Internal standards

Analysed using L-lines (Mo-K excitation)


Analysed using K-lines (W-L excitation)

Results

Micronutrients



Micronutrients




Results

Micronutrients



Micronutrients




- Fe, Cu and Zn could be analyzed with good accuracy
- Se could be quantified accurately in rice flour, but was below the detection limit in cereals (e.g. 110 $\mu\text{g}/\text{kg}$ for the S4 T-STAR)
- Iodine was present in concentrations distinctly below the detection limit of 20 mg/kg

Results

Toxic metals



Toxic heavy metals




Results

Toxic metals



Toxic heavy metals




- As could be analyzed in both samples
- Cadmium was distinctly below the detection limit of 2,8 mg/kg
- Hg was slightly below the detection limits (70 – 200 $\mu\text{g}/\text{kg}$)
- Pb could be analyzed in the cereal sample, but was below the detection limit in rice flour samples ($\sim 130 \mu\text{g}/\text{kg}$)

Results

Detection limits



3-sigma detection limits



Summary



- All food sample types, which can be dried and ground to a fine powder, can be analyzed with TXRF applying the described easy suspension preparation
- The quantification is done by means of internal standardization; no external calibrations are necessary
- Elements are detected simultaneously and independent from the concentration ranges
- It was demonstrated that macro nutrients Mg, P, Cl, K and Ca as well as micronutrients Fe, Cu, Zn and Se could be analyzed
- The control for the toxic heavy metals As, Pb and with restrictions Hg is possible
- The S4 T-STAR offers best detection limits and a high analytical flexibility
- The S2 PICOFOX is optimally suited for rapid onsite screening



TXRF analysis of dietary supplement pills

Objective



- The purpose of food supplements is to supplement the normal diet
- These supplements are concentrated sources of nutrients or other substances with a nutritional or physiological effect
- Relevant regulations are the European Directive 2002/46/EC and the US-FDA Dietary Supplement Health and Education Act of 1994



Methods and Equipment

Samples and sample preparation



Samples

2 multivitamin/mineral pills

Analytical task

Analysis of Mg, P, Ca, Cr, Zn, Se, Mo and I

Sample preparation

- Grinding in agate mortar
- Weighing of 100 ± 5 mg sample into vials
- Suspension in 5 ml aqueous Triton X 100 solution (1 vol.-%)
- Addition of 10 μ l Ni (1 g/l) and 10 μ l Pd (1 g/l)
- Homogenization
- Pipetting of 10 μ l sample onto quartz glass carriers
- Drying in vacuum



Methods and Equipment Measurements



S4 T-STAR


- Mo tube, 50 kV/1000 μ A
- W-tube, 50 kV/1000 μ A
- 100 mm² XFlash SDD

Measurements

- Mo-K excitation, 1000 s
- W-Brems excitation, 1000 s



1 H Hydrogen																	2 He Helium						
3 Li Lithium	4 Be Beryllium																	5 B Boron	6 C Carbon	7 N Nitrogen	8 O Oxygen	9 F Fluorine	10 Ne Neon
11 Na Sodium	12 Mg Magnesium																	13 Al Aluminium	14 Si Silicon	15 P Phosphorus	16 S Sulphur	17 Cl Chlorine	18 Ar Argon
19 K Potassium	20 Ca Calcium	21 Sc Scandium	22 Ti Titanium	23 V Vanadium	24 Cr Chromium	25 Mn Manganese	26 Fe Iron	27 Co Cobalt	28 Ni Nickel	29 Cu Copper	30 Zn Zinc	31 Ga Gallium	32 Ge Germanium	33 As Arsenic	34 Se Selenium	35 Br Bromine	36 Kr Krypton						
37 Rb Rubidium	38 Sr Strontium	39 Y Yttrium	40 Zr Zirconium	41 Nb Niobium	42 Mo Molybdenum	43 Tc Technetium	44 Ru Ruthenium	45 Rh Rhodium	46 Pd Palladium	47 Ag Silver	48 Cd Cadmium	49 In Indium	50 Sn Tin	51 Sb Antimony	52 Te Tellurium	53 I Iodine	54 Xe Xenon						
55 Cs Cesium	56 Ba Barium	57 La Lanthanum	72 Hf Hafnium	73 Ta Tantalum	74 W Tungsten	75 Re Rhenium	76 Os Osmium	77 Ir Iridium	78 Pt Platinum	79 Au Gold	80 Hg Mercury	81 Tl Thallium	82 Pb Lead	83 Bi Bismuth	84 Po Polonium	85 At Astatine	86 Rn Radon						
87 Fr Francium	88 Ra Radium	89 Ac Actinium																					




Results

Element concentrations



Element concentrations




Results

Detection limits





Detection limits



Summary



- All elements of interest can easily be analyzed with the **S4 T·STAR**
- The detection limits are sufficient for the accurate analysis of major and trace nutrients
- The analysis of Cr, Se and Mo according to several AOAC and EN norms is possible (AOAC 2011.19, AOAC 986.15, AOAC 996.17, AOAC 2006.03, EN14627, EN 14082, EN 14083)
- In addition, iodine can be monitored in dietary supplements without laborious sample treatment





TXRF analysis of milk samples

Objective

Samples and sample preparation



Samples

- 8 milk samples

Analytical task

- Analysis of Na, Mg, I and transition metals

Sample preparation

- Weighing of 1 ml milk into reaction vials
- Dilution with 1,5 ml ultrapure water
- Addition of 10 μl V (1 g/l), 10 μl Ga (100 mg/l) and 5 μl Pd (1 g/l)
- Homogenization
- Pipetting of 10 μl sample onto quartz glass carriers
- Drying in vacuum



Measurements



S4 T-STAR

- Mo tube, 50 kV/1000 μ A
- W-tube, 50 kV/1000 μ A
- 60 mm² XFlash SDD

Measurement conditions

- Mo-K excitation, 1000 s
- W-L excitation, 1000 s
- W-Brems excitation, 1000 s





Measurements



Spectra

- Typical spectra of milk measured with W-L (left) and Mo (right) excitation




Results



Element concentrations

- Different colored bars represent 8 milk samples
- Nutrient element concentration range from 50 ppm to 0,1%
- Trace elements can be detected below 10 ppb




Results & Discussion



Detection limits

- Untagged elements were measured with Mo excitation
- W-L provides 3 to 10 times better LLD for elements lighter than Cr
- Identification of toxic elements As, Pb, Cd, Hg at low levels is possible




Results & Discussion




Challenging elements

- Na: detectable with W-L excitation
 - Detection limit about 50 mg/kg
- I: improved detection limit of 0,7 mg/kg without sample dilution

W-L spectrum with 200 ppm sodium




W-Brems spectrum with 42 ppm iodine



Summary



- Applying the S4 T-STAR the entire element range from Na to U can be detected in milk samples
- Detection limits are in the ppb range for almost all elements





Q & A



Any Questions?

Please **type in** the questions you may have for our speakers in the **Questions Box** and click **Submit**





Innovation with Integrity

Thank you for your attention!

www.bruker.com

www.bruker.com/hhxrf

www.s4tstar.com