



Micro-XRF Back to the Roots – Part III



The speakers



- Falk Reinhardt
- Senior Application Scientist,
Bruker Nano Analytics, Berlin, Germany



- Dr. Roald Tagle
- Head of XMP Application,
Bruker Nano Analytics, Berlin, Germany



Overview

- **Part I** of this series aimed to introduce micro-XRF as a technique and give an overview of the individual components of a micro-XRF instrument and why they are the way they are, today.

- **Part II** focused on qualitative micro-XRF analysis.

- **Part III** will discuss quantitative XRF analysis.
 - What is quantification?
 - What are prerequisites for quantitative analysis?
 - Which samples can be quantified and how?

What is XRF quantification? The very basics




Element	AN	Series	Net intensity [cps]
O	8	K series	0.00
Na	11	K series	340.84
Mg	12	K series	172.30
Al	13	K series	153.60
Si	14	K series	12481.62
S	16	K series	49.57
K	19	K series	184.71
Ca	20	K series	7147.86
Ti	22	K series	56.83
Mn	25	K series	10.70
Fe	26	K series	327.96
As	33	K series	6.21
Rb	37	K series	17.30
Sr	38	K series	32.57
Zr	40	K series	23.02

Element	AN	Series	Compound	norm. stoich. C [wt.%]
O	8	K series		0.000
Na	11	K series	Na2O	14.282
Mg	12	K series	MgO	3.870
Al	13	K series	Al2O3	1.521
Si	14	K series	SiO2	71.338
S	16	K series	SO3	0.175
K	19	K series	K2O	0.351
Ca	20	K series	CaO	8.331
Ti	22	K series	TiO2	0.045
Mn	25	K series	MnO	0.003
Fe	26	K series	Fe2O3	0.076
As	33	K series	As2O3	0.001
Rb	37	K series	Rb2O	0.002
Sr	38	K series	SrO	0.003
Zr	40	K series	ZrO2	0.003

What is XRF quantification?

Obtaining net peak intensities

- The spectrum is composed of fluorescence lines, some background, and different artifacts.
- The aim is to determine the net peak intensities of the fluorescence lines.
- There are two fundamentally different approaches:
 1. “deconvolve” the spectrum, i.e. fit peaks into the spectrum.
 - Line overlap needs some thinking-about.
 2. Assume a concentration and forward-calculate what the spectrum of this sample would look like.
 - The net peak intensities are the result of the quantification.



What are prerequisites for quantitative analysis?

What is intensity?

- DIN ISO 22309 “Microbeam analysis - Quantitative analysis using energy-dispersive spectrometry (EDS) for elements with an atomic number of 11 (Na) or above” states:

3.21**peak intensity**

total number of X-rays (counts) under the profile of a characteristic X-ray peak after background subtraction

NOTE This is sometimes referred to as the peak integral.

- Usually in physics, intensity is a time-normalized measure, i.e. a candle’s flame does not become more intense just because you look at it longer.
- We will be using **counts per second** as a measure for intensity, not counts.


What are prerequisites for quantitative analysis?

- In order to obtain "good results" the expectations must be clear :
 - What is the analytical question? Can (and must) it be answered quantitatively?
 - What is the instrument's capability to solve the task?
 - What measurement conditions are ideal?
 - What sort of sample can be quantified with the methods employed?

Homogeneity, surface structure, dark matrix, crystallinity, concentration ranges, element combinations, orientation

What are prerequisites for quantitative analysis? Accuracy, Trueness and Precision

- Since 1994 ISO 5725-1 defines accuracy as a convolution of trueness and precision.




- This definition only works under the assumption that the reference value is the true value.

What are prerequisites for quantitative analysis?

Uncertainty of the net peak intensity



- XRF is a statistical measurement method, i.e. whether a photon is detected at a specific time is pure chance.
- With long enough time, several photons are detected and a peak becomes visible in the spectrum.
- This peak contains N counts and the statistical nature of this measurement dictates an uncertainty of \sqrt{N} .
- The relative error then is $\frac{\sqrt{N}}{N} = \frac{1}{\sqrt{N}}$
- The more counts, the smaller the relative uncertainty of the measurement → since the count rate is fixed, it needs time!
- This uncertainty is NOT reflected in the width of the fluorescence peak, it's how well we can determine it's area!



What are prerequisites for quantitative analysis?


Net peak intensity and quantification

- Since all quantification is based on obtaining the net intensities of the fluorescence peaks, it is safe to assume that longer measurement time improves the precision of the quantification ...



- ... but there is a limit.

Which samples can be quantified? Homogeneity




Instrument and analytical method performance

20 times same position
30 s each

Spectrum	Al	Mn	Fe	Cu	Zn	Pb	Bi
Al 2011 fixedPos 1	92.71	0.118	0.49	5.69	0.030	0.46	0.51
Al 2011 fixedPos 2	92.72	0.116	0.49	5.68	0.031	0.46	0.51
Al 2011 fixedPos 3	92.70	0.118	0.49	5.69	0.032	0.46	0.51
Al 2011 fixedPos 4	92.69	0.118	0.49	5.70	0.031	0.46	0.51
Al 2011 fixedPos 5	92.70	0.120	0.48	5.69	0.031	0.46	0.51
Al 2011 fixedPos 6	92.72	0.118	0.49	5.68	0.031	0.46	0.50
Al 2011 fixedPos 7	92.71	0.120	0.49	5.69	0.031	0.46	0.51
Al 2011 fixedPos 8	92.70	0.117	0.49	5.69	0.032	0.47	0.50
Al 2011 fixedPos 9	92.70	0.118	0.49	5.69	0.032	0.47	0.51
Al 2011 fixedPos 10	92.68	0.121	0.49	5.70	0.031	0.47	0.51
Al 2011 fixedPos 11	92.71	0.121	0.48	5.68	0.031	0.47	0.51
Al 2011 fixedPos 12	92.69	0.120	0.49	5.70	0.031	0.46	0.51
Al 2011 fixedPos 13	92.69	0.121	0.48	5.70	0.030	0.47	0.51
Al 2011 fixedPos 14	92.71	0.122	0.48	5.69	0.031	0.46	0.50
Al 2011 fixedPos 15	92.71	0.123	0.48	5.68	0.031	0.47	0.51
Al 2011 fixedPos 16	92.70	0.123	0.48	5.70	0.031	0.47	0.51
Al 2011 fixedPos 17	92.69	0.122	0.48	5.70	0.031	0.46	0.51
Al 2011 fixedPos 18	92.69	0.122	0.48	5.70	0.031	0.47	0.51
Al 2011 fixedPos 19	92.69	0.125	0.49	5.69	0.031	0.47	0.50
Al 2011 fixedPos 20	92.69	0.125	0.48	5.70	0.030	0.47	0.51
Mean value / wt. %:	92.70	0.120	0.48	5.69	0.031	0.47	0.51
Std. dev. / wt. %:	0.010	0.003	0.002	0.007	0.001	0.002	0.002
Std. dev. rel. / %:	0.01	2.17	0.44	0.13	1.78	0.53	0.43

Same as before plus sample homogeneity

20 different positions
30 s each




Spectrum	Al	Mn	Fe	Cu	Zn	Pb	Bi
Al 2011 sample pos 1	93.20	0.017	0.46	5.58	0.028	0.34	0.37
Al 2011 sample pos 2	92.52	0.065	0.47	5.98	0.034	0.44	0.49
Al 2011 sample pos 3	93.06	0.043	0.49	5.59	0.032	0.38	0.40
Al 2011 sample pos 4	92.83	0.066	0.46	5.61	0.032	0.49	0.52
Al 2011 sample pos 5	92.90	0.040	0.44	5.56	0.049	0.49	0.52
Al 2011 sample pos 6	92.44	0.030	0.51	6.11	0.035	0.42	0.46
Al 2011 sample pos 7	92.93	0.036	0.48	5.69	0.033	0.42	0.42
Al 2011 sample pos 8	92.85	0.019	0.49	5.66	0.046	0.47	0.48
Al 2011 sample pos 9	92.80	0.042	0.45	5.67	0.030	0.49	0.52
Al 2011 sample pos 10	93.20	0.144	0.44	5.58	0.032	0.27	0.32
Al 2011 sample pos 11	92.80	0.045	0.49	5.67	0.031	0.47	0.50
Al 2011 sample pos 12	92.96	0.092	0.43	5.52	0.026	0.48	0.49
Al 2011 sample pos 13	92.96	0.147	0.44	5.54	0.044	0.42	0.45
Al 2011 sample pos 14	93.04	0.034	0.47	5.63	0.027	0.38	0.42
Al 2011 sample pos 15	92.42	0.025	0.50	5.89	0.031	0.56	0.58
Al 2011 sample pos 16	92.84	0.021	0.49	5.60	0.027	0.49	0.53
Al 2011 sample pos 17	93.05	0.030	0.45	5.54	0.029	0.44	0.46
Al 2011 sample pos 18	92.92	0.035	0.44	5.51	0.028	0.53	0.55
Al 2011 sample pos 19	92.97	0.036	0.47	5.58	0.029	0.44	0.47
Al 2011 sample pos 20	92.91	0.017	0.48	5.60	0.027	0.47	0.49
Mean value / wt. %:	92.88	0.049	0.47	5.65	0.033	0.45	0.47
Std. dev. / wt. %:	0.214	0.038	0.023	0.158	0.007	0.066	0.061
Std. dev. rel. / %:	0.23	76.59	4.87	2.79	19.96	14.72	13.02

Which samples can be quantified?

What “model” is implemented in the quantification algorithm?



In XRF, samples are classified into 3 groups:

- **Infinitely thick samples**
 - A thicker sample would not change the XRF signal.
 - Bulk-XRF
- **Thin samples**
 - Attenuation and self-absorption effects can be neglected.
 - TXRF, some synchrotron science.
- **Intermediate-thick samples**
 - Not easy!
 - Almost all **layer thickness analysis**.
 - Makes a sample a “non-ideal” sample.



Which samples can be quantified? What is an ideal sample? What is not?

- Actually, anything that is not flat and homogeneous is not ideal.
- Homogenous, means **one composition** within the analytical volume!
- For bulk-XRF, in addition, the sample needs to be infinitely thick (for the analyzed elements).
- Defined sample geometry.




Which samples can be quantified?

What “model” is implemented in the quantification algorithm?

- What is “infinitely thick”?
- The M4 TORNADO software comes with a Fundamental Parameter (FP) bulk quantification algorithm.
- For quantification of “non-ideal” samples XMethod software package can be used.

Information depths of selected element fluorescence lines in different matrices



xrfcheck.bruker.com

Which samples can be quantified? ...and with what algorithm?

What can be quantified?

Ideal samples

FP

- Flat, homogeneous bulk
- Flat layers

Non-ideal samples


FP with larger uncertainties, i.e.
lower accuracy
Empirical, if the standard
matches the samples

Samples difficult to describe

empirical

- Wherever the standard
matches the sample

Bulk-FP quantification of “ideal” samples Steels








Grade-IARM	Al	Si	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Nb	Mo	W
AISI 316-5D	0.0050	0.5	0.0	0.0	16.6	1.8	68.2	0.1	10.4	0.17	0.004	2.1	0.0
AISI 321-6D	0.1100	0.3	0.6	0.1	17.5	1.5	69.4	0.2	9.4	0.30	0.039	0.4	0.1
AISI 330-7B	0.0230	1.4	0.0	0.0	19.3	1.5	41.3	0.1	35.8	0.21	0.023	0.2	0.0
AISI 347-8D	0.0040	0.4	0.0	0.1	17.3	1.8	69.3	0.1	9.2	0.47	0.720	0.4	0.1
AISI 410-9C	0.0140	0.4	0.0	0.1	12.0	0.4	86.2	0.0	0.3	0.06	0.005	0.2	0.1
AISI 416-10C	0.0030	0.4	0.0	0.0	12.3	0.4	86.0	0.0	0.2	0.16	0.003	0.1	0.0
AISI 420-154B	0.0020	0.5	0.0	0.1	12.2	0.4	86.1	0.0	0.2	0.09	0.003	0.1	0.0
AISI 422-205B	0.0090	0.4	0.0	0.3	11.7	0.7	83.7	0.0	0.7	0.15	0.018	1.0	1.1
AISI 430-11C	0.0100	0.5		0.0	17.7	0.5	80.8	0.0	0.2	0.07	0.005	0.1	
AISI 431-12B	0.0030	0.6	0.0	0.0	16.0	0.6	80.1	0.0	2.2	0.14	0.011	0.1	0.0
AISI 440C-13C	0.0030	0.7	0.0	0.1	16.8	0.4	80.2	0.0	0.1	0.03	0.004	0.5	
AISI 446-14B	0.0040	0.5	0.0	0.1	23.6	0.4	74.8	0.0	0.3	0.07	0.006	0.1	0.0
Nitronic 40-19B	0.0100	0.5	0.0	0.1	20.0	9.3	62.3	0.1	6.8	0.17	0.057	0.3	0.0
Nitronic 50-17B	0.0030	0.4	0.0	0.2	21.3	5.1	56.3	0.1	13.4	0.17	0.220	2.3	0.1
17-7PH-152B	1.1600	0.4	0.1	0.1	16.9	0.8	72.3		7.2	0.31	0.033	0.5	0.1

- Steels are usually homogeneous, they very quickly become infinitely thick, and they are easily made flat.
- Out of 51 ARMI reference alloys, we measured 15 stainless steels of certified composition.
- Average of 10 positions per sample with 10 s measurement time per point.


Bulk-FP quantification of “ideal” samples Steels

- Due to systematic physical effects (tertiary excitation in stainless steel samples) the quantification is not “spot-on”.
- The linearity of the regression shows a very stable recovery rate over the whole concentration range. → high precision
- The lack of trueness can be corrected for by calibration with one sample



Bulk-FP quantification of "ideal" samples Steels

- Type-calibration: adapting the sensitivities of the FP algorithm for the elements to correct the recovery rate to 1



Configuration - steel_TC

Elements

Use spectrum elements
 Use list elements
 Search additional elements

Double click an element to open element editor Clear all

Special properties of selected elements

Compound	Fix %	Dec. Dif.	Fact.
			1,00
Al			1,15
Si			1,50
V			1,16
Cr			0,97
Mn			0,98
Ni			1,07
Cu			1,15
Nb			1,15
Mo			1,07
Rh		<input checked="" type="checkbox"/>	1,00
W			1,10

Legend

- Fe Fixed list
- Compound
- Stoichiom. elements
- Fix concentration
- Deconvolution only
- Excluded element
- Difference element





Global options

Background cycles: Default Manual (60)


Minimum concentration: 0,00 %

Description

Load... Save... OK Cancel




Layer-FP quantification of “ideal” samples Metallic layers



- To quantify layered samples, the quantification algorithm has to know the concentration of all the elements that are in each layer → appropriate method editor needed

Calibration with two reference samples




- Even though FP methods work without calibration, sometimes (f.e. industrial applications) a calibration is required for highest accuracy (or proof thereof)


Layer-FP quantification of “ideal” samples Metallic layers

- The calibration corrects the trueness.
- The precision – if the sample allows – can be changed by different measurement times

Short measurements



long measurements



- For this sample system 20 s measurement time per point yielded sufficient stability.

10 points	Ni layer thickness	Pd layer thickness	Au layer thickness
1 s real time	/ μm	/ nm	/ nm
mean	2.09	50	28
sigma	0.05	14	3
sigma rel.	2.2	27.3	9.0

10 points	Ni layer thickness	Pd layer thickness	Au layer thickness
10 s real time	/ μm	/ nm	/ nm
mean	2.09	55	27
sigma	0.02	1.3	0.8
sigma rel.	0.7	2.3	2.8

10 points	Ni layer thickness	Pd layer thickness	Au layer thickness
60 s real time	/ μm	/ nm	/ nm
mean	2.10	54	27
sigma	0.01	0.7	0.3
sigma rel.	0.5	1.3	1.1

10 points	Ni layer thickness	Pd layer thickness	Au layer thickness
200 s real time	/ μm	/ nm	/ nm
mean	2.09	54	28
sigma	0.01	0.8	0.2
sigma rel.	0.4	1.5	0.8


Layer-FP quantification of “ideal” samples

Metallic layers

- Verification on an independent layer thickness standard:

Au: 53 nm
 Pd: 105 nm
 Ni: 4.02 μm




10 points	Ni layer thickness	Pd layer thickness	Au layer thickness
30 s real time	/ μm	/ nm	/ nm
mean	4.01	101	52
sigma	0.03	2	1
sigma rel.	0.7	2.1	1.2



	Substrate	Layer 1		Layer 2		Layer 3	
Spectrum	Cu [%]	Thickn. [μm]	Ni [%]	Thickn. [μm]	Pd [%]	Thickn. [μm]	Au [%]
T_H_10	100.00	4.251	100.00	0.233	100.00	0.061	100.00
T_H_9	100.00	3.644	100.00	0.239	100.00	0.061	100.00
T_H_8	100.00	4.148	100.00	0.235	100.00	0.060	100.00
T_H_7	100.00	4.173	100.00	0.245	100.00	0.060	100.00
T_H_6	100.00	3.742	100.00	0.238	100.00	0.062	100.00
T_H_5	100.00	3.715	100.00	0.235	100.00	0.062	100.00
T_H_4	100.00	3.991	100.00	0.224	100.00	0.060	100.00

Bulk-FP of “ideal” samples

Nano-milled geological samples – major and minor elements





- Milling geo SRMs to < 200 nm gives the resulting powder a glass-like quality.
- There is no measurable inhomogeneity anymore. When pressed, they stick without binder. The grains are smaller than the information depth of the elements!

Bulk-FP of “ideal” samples

Nano-milled geological samples – major and minor elements

- Again, the recovery rate is very good over a wide concentration range. (especially when considering ~ 50 % “dark matrix”)
- Remember, the quantification assumes a sample composition, based on the selected elements and iteratively calculates all physical effects, like absorption and secondary excitation, and adapts the concentrations until the theoretical spectrum matches the measured one.



Empirical bulk-quantification Nano-milled geological samples – trace elements

- When looking for trace elements, the signal-to-noise ratio of the measurement becomes crucial.
- To optimize the SNR for element like Rb, Sr, Y, Zr, Nb, a strong primary filter can be applied.
- BUT: using a strong filter reduces the excitation intensity for the light elements (which are major elements), therewith strongly afflicting the precision of the net peak intensity and, hence, the quantification.
- It is very ambitious to try trace element quantification like this.
- Solution: de-couple the quantification of the traces from the major elements.



http://georem.mpch-mainz.gwdg.de/sample_query.asp

Elements Analyzed	Concentration Range / ppm
Rb	6 – 390
Sr	3 – 246
Y	17 – 75 (184*)
Zr	48 – 780
Nb	9 – 110
Pb	5 – 45
Th	5 – 87
U	2 – 18

Reference samples utilized:

GH, JR-2, AC-E, RGM-1, JA-2, JB-2, BHVO-2, SARM-1, NIST 620

Empirical bulk-quantification Nano-milled geological samples – trace elements

- The empirical method focuses solely on the elements of interest.
- This is a valid approach, because in oxidic matrix there are no relevant inter-element effects for these trace elements.
- So, ignoring 99.9 % of the sample, we basically determine the instrument's sensitivity (cps per ppm) for these elements.

The screenshot shows the 'METHOD EDITOR' window for 'Oxide trace test 1'. It is divided into several sections:

- Method data:** Description: Oxide trace test 1; Type: polynomial calibration; Spectrum deconvolution: Bayes deconvolution.
- Measurement parameters:** HV / kV: 50; Collimator / mm: 15 μm LENS; Atmosphere: Vacuum; Current / μA: 600; Measure time / s: 20.
- Layer parameters:** Name: base; Compound: (empty); Start thickness: 99,999.9; Unit: μm (1); Fixed thickness: checked; Calc. mode: Emission; Normalize sample/standard: checked; Target value (%): 100.00; Density: Default (1.00) selected; Use tolerances?: unchecked.
- Structure:** A table showing the layer 'base' with chemical elements Rb, Sr, Y, Zr, Nb.
- Element overview:** A table listing elements and their start concentrations.

Element	Z	Main line	Start conc.
Rb	37	KA	0 ppm
Sr	38	KA	0 ppm
Y	39	KA	0 ppm
Zr	40	KA	0 ppm
Nb	41	KA	0 ppm

Empirical bulk-quantification Nano-milled geological samples – trace elements

- Used calibration curve: linear correlation with offset

Calibration samples

Layer	Parameter	W.	Std	calib.	Std	calib.	Y(ppm)	Std	calib.	Std	calib.	Nb(ppm)	Std	calib.
13	NANO CRM GH T Yuan	1	390	385	10	7	75	74	150	168	85	91		
20	NANO CRM SARM-1 T Yuan	1	325	317	10	12	143	95	300	297	53	52		
16	NANO CRM JR-2 T Yuan	1	297	313	8	11	51	52	97	102	19	18		
14	NANO CRM JR-1 T Yuan	1	257	263	30	41	45	45	101	107	16	18		
11	NANO CRM AC_E T Yuan	1	152	145	3	1	184	101	780	778	110	105		
15	NANO CRM RGM-1 T Yuan	1	150	147	105	129	23	24	228	224	9	8		
19	NANO CRM JA-2 T Yuan	1	70	59	246	256	17	15	109	104	9	13		
18	NANO CRM BHVO-2 T Yuan	1	9	14	294	294	26	12	171	111	18	14		
17	NANO CRM JB-2 T Yuan	1	6	14	178	134	24	12	48	34	6	6		

Deviation (σ)

Rb	9	19	1	10	4
----	---	----	---	----	---


Normation type
 Sample to 100.00 %
 Norm target value: 100

Calibration curve offset
 Use offset:

Calibration ranges
 Order of polynomial: 1

1000000 ppm Active

Rb average deviation = 9 ppm
 Deviation of Rb in all samples ± 15 ppm





Empirical bulk-quantification Nano-milled geological samples – trace elements

- Detection limits in micro-XRF depend on instrument performance, samples used, and measurement conditions. In combination with a simple mathematical model allows fast and accurate quantification of geological materials or their derivatives such as ceramics.
- The major uncertainty is the “known” values and (in)homogeneity for the chosen reference samples. Issues with sample inhomogeneity can be overcome in “ideal” samples such as the nano-powders described previously (or see presentation on nano-powders).

Element	Avg dev
Rb	9
Sr	19
Y	14
Zr	10
Nb	4


Values in ppm mg/kg sample	Rb	Sr	Y	Zr	Nb	Pb	Th	U
Estimated limit of detection*	1	1	1	3	2	1	1	1
Estimated limit of quantification**	3	3	3	5	5	3	3	4
* LOD= 3*concentration/SNR								
* LOQ= 9*concentration/SNR								

These limits represent method limits for oxide sample analysis with the M4 Tornado under the previously described condition.

Empirical bulk-quantification

Really adverse geometry

- Measuring sulfur in oil can be challenging, especially when only limited amounts of oil are available.
- Small amounts of oil for a droplet, which really does not fulfill the criteria of “ideal samples”.
- It's round, most of the X-rays go through, whether the S is at the surface or homogeneously distributed is difficult to assess, ...
- But: There is limited information depth for S and the self-absorption effects are negligible.
- Assumption: concentration of S oil is linear with measured intensity





Oil droplet



Empirical bulk-quantification

Really adverse geometry

- 10 µl of each reference sample (CONOSTAN) were deposited on a clean SiO₂ disc
 - Concentration of S from 0 to 10000 ppm (1 %)
 - The sample description merely accounts for the elements that interfere with the fitted S intensity
- base S Rh Ar
- Still the approach yields very good quantification.







Calibration samples

No.	Layer Parameter Sample	W.	S(ppm)		Rh(%)		Ar(%)	
			Std	calib.	Std	calib.	Std	calib.
124	S_OIL_10000 ppm	1	10000	9993	0.00	0.00	0.00	0.00
123	S_OIL_5000 ppm	1	5000	5014	0.00	0.00	0.00	0.00
122	S_OIL_1000 ppm	1	1000	995	0.00	0.00	0.00	0.00
120	S_OIL_500 ppm	0	500	458	0.00	0.00	0.00	0.00
119	S_OIL_100 ppm	1	100	93	0.00	0.00	0.00	0.00
118	S_OIL_25 ppm	1	25	23	0.00	0.00	0.00	0.00
116	S_OIL_0 ppm	1	0	6	0.00	0.00	0.00	0.00

Empirical bulk-quantification

Soaked-up liquids


- Deposition of a defined amount of solution in a defined area, scanning, and integration the detected fluorescence intensity.
- For the analysis 5 μ l-droplets were deposited on a $\sim \varnothing$ 6 mm absorbent tissue.
- To evaluate impact of poor sample preparation, two samples were “folded”.



Empirical bulk-quantification


Soaked-up liquids

- There is a very good linear correlation between the extracted intensity and the concentration. From the slope (sensitivity) the concentration in the sample can be derived. The main uncertainty results from the droplet preparation (see the deviation in the 3 dots for the same concentration).
- The folded samples also results in small deviations (see arrow).



Summary

- Micro-XRF is versatile quantification tool for the analysis of multiple sample types from solid to liquids or powders.
- The combination of high element sensitivity from major to trace elements combined with the easy operation allows the use of multiple empirical and fundamental parameter-based models.
- Every sample with a defined composition within the analytical volume can be quantified as bulk.
- If the sample is not homogenous you can either make compromises in terms of trueness as the sample does not have a real or true “composition”, or you move to empirical approaches
- Layer samples can be quantified with high accuracy as long as the layer succession and composition is known.
- Liquid samples can be measured depending on the element of interest (low-Z, high-Z) or nature of the liquid (drying easily or not).



Innovation with Integrity