

IVD



# **DNA**•STRIP Technology

Convenient diagnostic assays

## **DNA•**STRIP technology for dependable diagnostics

- Reliable
- Highly sensitive and specific
- Cost-efficient
- User-friendly

Innovation with Integrity

PCR

## **DNA**•STRIP Technology

#### **Assay principle**

DNA•STRIPs are coated with specific probes which are complementary to the amplified nucleic acid (amplicon). After denaturation the single-stranded amplicon specifically binds to the probes (hybridization) and is visualized in a subsequent enzymatic colour reaction. As a result, a specific banding pattern develops on the DNA•STRIP. This procedure can be performed manually or can be automated.





## Benefits of using the DNA•STRIP technology

Reliable: Internal controls document valid results, and secure safe and impeccable test procedures. Thus, a high diagnostic reliability is guaranteed. **Highly sensitive and specific:** The chosen targets of the different test systems ensure high specificity combined with maximum sensitivity. Cost-efficient: Only minimum technical equipment is needed for processing. This allows for a cost-effective implementation in all laboratories. User-friendly: In contrast to conventional methods, the DNA•STRIP technology saves valuable time and can easily be integrated in your laboratory routine.

#### Internal controls ensure valid results

The evaluation of the DNA•STRIP can easily be The combination of specific amplification and hybridperformed by aligning it to a template. With the ization guarantees a high level of diagnostic reliability. corresponding IFU the developed banding pattern can Internal controls ensure valid results: be interpreted accurately and guickly. This allows a • All DNA•STRIPs: The Conjugate Control documents reliable detection of the bacteria or genotype present the efficiency of the colour reaction. in the sample.

- Microbiological DNA•STRIPs: The additionally integrated Universal or Amplification Control shows that the test was performed correctly.
- Human genetic DNA•STRIPs: Gene site-specific Sensitivity Controls confirm the sensitivity of the hybridization reaction. If hybridization took place under unspecific test conditions, this is documented by the Specificity Control.

#### Sample preparation

Starting point of the analysis is a specimen from which the nucleic acid is extracted. The nucleic acid is selectively amplified in a subsequent PCR reaction. Since single-stranded nucleic acid is required for the next step, amplification is followed by a denaturation step.

#### **Evaluation**

## Tests based on the DNA•STRIP technology

#### **Human Genetics**

ThromboType®	Factor V Leiden, Factor II G20210A
ThromboType <sup>®</sup> plus	Factor V Leiden, Factor II G20210A, MTHFR C677T, A1298C
GenoType CVD	Eight different thrombophilia-associated mutations
GenoType MTHFR	Most important MTHFR polymorphisms
GenoType ApoE	Alleles $\epsilon 2$ , $\epsilon 3$ , $\epsilon 4$ of the ApoE gene
GenoType PAI-1	Most important PAI-1 polymorphisms
GenoType HH	Hereditary hemochromatosis
GenoType AAT	Alpha-1-antitrypsin deficiency allele
GenoType LCT	Most important polymorphisms of the lactase gene
GenoType SugarTol	Polymorphism in lactase gene C-13910T and 3 polymorphisms
	in aldolase B gene
Microbiology	
Genolype MIBC	MTB complex differentiation from culture
Genolype CM <i>direct</i>	MTB complex, NTM differentiation from clinical specimens
GenoType Mycobacterium CM	MIB complex, NIM differentiation from culture
GenoType Mycobacterium AS	Further NTM differentiation from culture
GenoType NTM-DR	NTM differentiation, resistance to microlides and
	aminoglycosides from culture
GenoType MTBDR <i>plus</i>	MTB complex, resistance to rifampicin/isoniazid from clinical
	specimens and culture
GenoType MTBDR <i>sl</i> VER1.0	MTB complex, resistance to fluoroquinolones/aminoglycosides/
	cycl. peptides/ethambutol from clinical specimens and culture
GenoType MTBDR <i>sl</i> VER2.0	MTB complex, resistance to fluoroquinolones/aminoglycosides/
	cycl. peptides from clinical specimens and culture
GenoType LepraeDR	<i>M. leprae</i> , resistance to rifampicin/ofloxacin/dapsone
	from clinical specimens
GenoType MRSA	S. aureus, S. epidermidis, mecA, mecC, PVL from culture
GenoType HelicoDR	<i>H. pylori</i> , resistance to fluoroquinolones/clarithromycin
	from clinical specimens and culture
GenoType EHEC	Shiga toxins, virulence factors from culture
GenoType Enterococcus	Species differentiation, resistance to vancomycin from culture
micro-IDent <sup>®</sup>	5 periodontopathogenic bacterial species
micro-IDent®plus11	11 periodontopathogenic bacterial species

Please contact your local representative for availability in your country. Not for sale in the USA. ThromboType® and microIDent® are registered trademarks of the Bruker corporate group.

#### Hain Lifescience GmbH

Hardwiesenstrasse 1 72147 Nehren - Germany Phone +49 (0) 7473-9451-0



info.mdx.de@bruker.com - www.hain-lifescience.de