"Deep Matter" annotation using 4D-MetabolomicsTM in *Glycyrrhiza uralensis* (東北甘草, 西北甘草, 炙甘草) used for Traditional Chinese Medicine



Bruker Japan K.K., 3-9 B6F, Moriya-Kanagawa-ku, Yokohama-city cho, KANAGAWA, 221-0022, Japan

² Bruker Scientific LLC 101 Daggett Drive San Jose, CA ,95134 USA

Novel Aspect

LC-TIMS-MS separation and annotation workflow for isomeric and isobaric

Methods

Root samples of *Glycyrrhiza uralensis* three different regions were from purchased from Uchida Wakanyaku Ltd.: Tohoku (東北甘草), Ningxia Hui (西北甘 草), and Gansu - roasted (炙甘草). Samples were milled and extracted with methanol using a BeatBox 80% homogenizer (PreOmics) and centrifuged. Supernatants were analyzed in triplicate using a Bruker timsTOF HT with a VIP-HESI source and an Agilent 1290 liquid chromatograph on a reverse phase gradient. Data was acquired with default 4D-Metabolomics parameters which employs DDA-PASEF, in both positive and negative polarities. Data analysis was performed in MetaboScape® 2024b.

Summary of annotations										
	MS #	MS/MS #	MOMA #	MS ≤ 5 ppm	MS, MS/MS ≤ 5 ppm, ≥ 500 score	MS, CCS ≤ 5 ppm, ≤ 4%	MS, MS/MS, CCS \leq 5 ppm, \geq 500 score, \leq 4%			
Positive	9481	6943	900	1643	1483	1108	1054			
Negative	10355	7573	1230	971	852	331	317			

MOMA: easy to find isomeric and isobaric metabolites tolerance, RT $\Delta 1$ s; *m/z*, $\Delta 3$; CCS, $\Delta 10$ Å





secondary metabolites from natural products

Introduction

Plant derived natural products have been traditional medicines for used as They contain secondary centuries. metabolites that often possess bioactivity that may be beneficial as therapeutics. Characterization of natural products necessitates in-depth mass spectrometry (MS) analysis. In addition to accurate mass, isotope pattern and fragmentation pattern, ion mobility adds an orthogonal measurement, collisional cross section (CCS) that is unique to the size and shape of a molecule. Trapped ion mobility spectrometry (TIMS) provides highly accurate CCS values for more confident annotations. Here, we demonstrate characterization Glycyrrhiza of metabolites using TIMS-MS.

LC conditions

Mobile phase	: Solvent A, H ₂ O inclu : Solvent B, Acetonitr	0 including 0.1 % formic acid etonitrile including 0.1 % formic acid						
Column	: HSS 2 mm x 100 mm 1.7 μm (Waters)							
Column oven	: 40 °C							
Gradient	: Right		O a h a A					
Inidation	· 5	min	SOIV. A	SOIV. B				
Injection	. υμί	0.0	95	5				
Flow rate	: 0.4 mL/min	10.0	0	100				
		12.5	0	100				
		12.6	95	5				
		15.0	95	5				





Results

Data processing resulted in 9481 and 10355 features and 6943 7573 and MS/MS spectra, in positive and negative mode, respectively. Finally, 2439 features were annotated by merging data from polarities. Principal component both analysis (PCA) showed clear separation between 東北甘草 and 西北甘草, and 炙甘 草, suggesting differences in metabolites depended on roasting. Features were with open-source annotated natural product database, KNApSAcK, containing InChI or SMILES to generate CCS values using CCS-Predict Pro, and in silico fragmentation patterns using MetFrag. Annotations were inspected manually to confirm accuracy. The separation of isomeric and isobaric metabolites using TIMS the number of increases metabolites detected. MOMA (Mobility-Offset Mass-Aligned) characterized pairs of isomer and isobaric metabolites. These metabolites were detected at the same retention times, but were completely separated in TIMS regardless of their size.

Fig. 3 Summary of annotations using MetaboScape and examples of characterized isomer/isobaric metabolites using MOMA function.

Summary

- TIMS separates small molecules regardless of their size and characterizes isomeric and isobaric metabolites, which are not separated under the LC conditions.
- MetaboScape using Target List including InChl SMILES Or specialized metabolites annotated including 18α -glycyrrhizin, а metabolites representative ÍN Glycyrrhiza uralensis.





SmartFormula Element-based

Χ

Preparing InChI or SMELIS of G. uralensis in Excel using KNApSAcK database. Importing the list to MetaboScape Annotating specialized metabolites with *in silico* fragmentation and predicted CCS value.



Fig. 1. Introduction of MS instrument, software, and how to annotate metabolites using MetaboScape functions.



Fig. 2. Result of PCA in positive (left) and negative ion mode (right). For Research Use Only. Not for use in diagnostic procedures.

Conclusion

4D-Metabolomics[™] using timsTOF is a strong tool for unearthing Deep Matters. TimsTOF is useful for metabolomics in plants

including Traditional Chinese Medicine research.

4D-MetabolomicsTM

Innovation with Integrity

© 2024 Bruker