

A spatial multiomics workflow on a new benchtop MALDI-TOF instrument deciphering the lipid and protein landscape of tissues

Janina Oetjen¹, Corinna Henkel¹, Dagmar Niemeyer¹, Ulrike Schweiger-Hufnagel¹, Tanja Bien¹, Sebastian Böhm¹, Christoph Nordmann¹, Andree Schuster¹, Katherine Stumpo², Michael Easterling², Andras Kiss⁴, Andrea Tannapfel³, Jens Christmann³

¹Bruker Daltonics GmbH & Co. KG, Fahrenheitstraße 4, Bremen, Germany; ²Bruker Scientific, 40 Manning Road, Billerica, MS, USA; ³Institute for Pathology, Ruhr-University Bochum, Bürkle-de-la-Camp-Platz 1, 44789 Bochum, Germany; ⁴Bruker Austria GmbH, Lemböckgasse 47b, 1230 Wien, Austria;

Introduction

Spatial biology is an emerging field combining spatial cellular information with molecular information of biological tissues, and as such is perfectly complementary to MALDI mass spectrometry imaging (MALDI Imaging). As a non-destructive method, MALDI Imaging can tie various modalities into one dataset by characterizing tissue phenotypes according to lipid or metabolic profiles and integrating spatial targeted protein expression.

Here we present a spatial multiomics workflow combining MALDI Imaging of lipids and MALDI HiPLEX-IHC on a new benchtop axial MALDI-TOF mass spectrometer using a clinically relevant research sample.

Methods

Formalin-fixed and paraffin-embedded (FFPE) colorectal cancer sections were processed for MALDI Imaging of lipids. Briefly, paraffin was removed by xylene washes and DHB matrix was sprayed using a M3+ sprayer (HTX-Technologies). MALDI Imaging data was acquired on a neofleX Imaging Profiler in positive reflector mode with 20 μm raster width at m/z 500-1100. Initial mass calibration was performed on a spot of red phosphorous. During measurement, an online calibration was conducted on endogenous lipid peaks. A consecutive tissue section was measured on the timsTOF fleX to obtain a reference list of annotated lipids for this sample. A MALDI HiPLEX-IHC experiment, using technology from AmberGen, Inc. was conducted on the same section as used before for lipid imaging on the neofleX. The 14 mass tagged antibodies used are listed in (3). After HCCA matrix application, data was acquired on the neofleX using a default acquisition method with 20 μm raster width and SCiLS autopilot for automated setup ensuring controlled instrument conditions for stable measurements. An OME-TIFF file was generated automatically for visualizations in SCiLS Scope. The section was stained with H&E post MALDI, and microscopy images scanned on Hamamatsu nanozoomer. Multimodal data integration, statistical analysis and visualizations were done in SCiLSTM Lab 2024b.




Fig. 1 Spatial multiomics workflow

Results




Fig. 2 MALDI Imaging of lipids of a colorectal cancer sample and confirmation of annotations.




Fig. 3 MALDI HiPLEX-IHC of the same sample

A) A 14-plex HiPLEX experiment with the antibodies indicated on the right of the figure was conducted on the same colorectal cancer tissue section used for MALDI Imaging of lipids. An OME-TIFF file was automatically generated and visualized in the new software, SCiLS Scope. The ion images of the mass tags from the selected antibodies were visualized by false-color coding in a multichannel image. High intensity of tumor marker PanCK (pan cytokeratin) was observed in a specific region (red) where also immune cells infiltrated (CD68, macrophage marker, yellow). B) SCiLS Scope enables measurement of distances as shown here for two blood vessels.




Fig. 4 MALDI Imaging of lipids of a colorectal cancer sample and confirmation of annotations.

Conclusion

- Our new spatial multiomics workflow on the neofleX Imaging Profiler provides insight into the lipid and targeted protein expression landscapes of tissues
- SCiLS Lab 2024b enables multimodal data integration and high-resolution optical image support
- Automatically generated OME-TIFF files can readily be visualized with SCiLS Scope

Imaging: Spatially-Resolved Omics