




Disentangling distorted distributions - improving z' ion assignment confidence and structural characterization of proteins by top-down MALDI-in-source decay MS

Introduction




Protein structural analysis by mass spectrometry relies on the use of different fragmentation techniques to provide comprehensive determination of protein sequence, including any post-translational modifications. MALDI, combined with in-source decay (ISD) fragmentation (Figure 1 and 2) and FT-ICR-MS, enables top-down C- and N-terminal sequencing of proteins. Here, we evaluate why isotopic distributions of z' ions, in top-down protein MALDI-FT-ICR mass spectra, deviate from theoretical profiles: a result of interference from the overlap with isotopic distributions of both z and y-NH₃ ions (which can co-occur and are isomeric with each other). Recognizing this interference can result in higher confidence structural characterization of proteins, which is important for manufacturing of pharmaceutical products such as therapeutic monoclonal antibodies (mAbs).

Method

MALDI-(ISD) MS measurements were performed on a 15T solariX XR FT-ICR MS system (Bruker Daltonics). 1,5-DAN and PPDA were used as MALDI matrices. NIST mAb was used as model compound. MALDI-MS mass spectra of NIST mAb were processed in both magnitude (mFT) and absorption (aFT) modes, using AutoVectis (Spectroswiss). A figure-of-merit (FoM) was calculated to provide a measure of how closely an isotopic peak distribution matched the theoretical distribution. Two synthetic peptides that contain either normal or deuterated alanine residues were used to investigate the presence and unravel the identity of isomeric z and y-NH₃ fragment ions.




The isotopic distributions of z' ions are distorted due to a partial overlap with isomeric z and y-NH₃ ions



Results

- Isotopic distributions of z' and z/y-NH₃ ions detected in mFT (A and B) and aFT (C and D) MALDI-MS FT-ICR MS spectra of the heavy chain of NIST mAb. The colored marks indicate theoretical isotopic distributions. The FoM values obtained for aFT isotopic distributions of z'16 and z'46 ions were 0.95 and 0.81, respectively.
- Evaluation of the apparent distortions of the isotopic distribution of z' ions induced by isomeric z and y-NH₃ ions on theoretical MALDI-MS FT-ICR mass spectra of the heavy chain of NIST mAb. The goodness of the isotopic distribution of each z' ions, expressed as FoM, was calculated considering different z/y-NH₃-to-z' ratios. FoM values lower than 1 indicate a deviation of the isotopic distributions induced by the presence of z/y-NH₃ ions.
- Evaluation of isomeric z and y-NH₃ fragment ions in MALDI-MS FT-ICR MS spectra of human insulin analyzed using A) 1,5-DAN and B) PPDA as MALDI matrices.



(4) Evaluation of isomeric z and y-NH₃ ions in MALDI-MS FT-ICR mass spectra obtained from the analysis of a pair of peptides, A) one with normal alanine (m/z 2082.0814) analyzed with 1,5-DAN and one with heavy alanine (m/z 2094.1591) analyzed with both B) 1,5-DAN and C) PPDA. The abstraction of one deuterium instead of one hydrogen from alanine-6 led to the detection of an additional peak at m/z 1684.9307 (or m/z 1684.9299). The isotopic distribution on z15*-D resulted distorted as a consequence of the presence of y15-NH₃ ions.

Conclusion

In this study, we showed that the detected isotopic distribution of z' fragment ions, generated by MALDI-MS FT-ICR MS of proteins, can be affected by the presence of interfering z and y-NH₃ fragment ion isotopic distributions. The z and y-NH₃ fragment ions are isomeric but are generated by different fragmentation paths. The evaluation of the presence of these ions increases the confidence of the identifications and explains decreased FoM values of certain z' ions.