Dimension Industrial AFMs

● Automated Nanometrology and Characterization

Cost-Effective Precision Metrology for R&D/QA/QC/FA

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● Exclusive Bruker AFM Technology for Industry

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Revolutionary PeakForce Tapping Enables Exceptional Surface Characterization

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FastScan Technology Delivers Highest Throughput

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

Proprietary FastScan AFM Probes Provide Lowest Cost per Measurement

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An alkali-aluminosilicate sheet glass main property is characterized by its strength and scratch resistance, making it challenging to measure its nanometer topography consistently and reliably. Using FastScan AFM and a single Bruker HPI probe, a first image (A) rendered roughness $R_q=0.99$ nm, $R_a=0.79$. Image no. 615 (B) rendered the same roughness; $R_q=0.99$ nm, $R_a=0.79$. Insets are digitally zoomed images of scan no. 1 and no. 615 revealing insignificant topographical changes.

● Easy-to-Use Software Makes Every User an AFM Expert

ScanAsyst Eliminates the Complexity of AFM Operation

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Site Assignments

- ▲ Sites
 - ▲ Trench
 -  30um Overview
 -  7um Depth
 -  1um Detail

Easy measurement recipe creation, allows engineers to define location by name, assign any type, and number measurements at each location.

Automation Software Ensures Exceptional Ease of Use

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Export


Edit

NanoScope Provides Fast, Intuitive Sample Navigation

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Data Analysis Flexibility for Varied Measurement Needs

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Consistently better data and faster metrology results, regardless of the operator's AFM experience level.

● Precision, Repeatability, and Flexibility for Production

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Highest Accuracy Roughness, Depth, and Height Measurements

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
Most Comprehensive Defect Review Characterization

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Fast Nanoelectrical Measurements

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Trenches are difficult to image with typical imaging modes due to excessive damping of the probe oscillation. ScanAsyst easily reaches the bottom of trenches. In this example, trenches with a depth of ~65 nm and a width of ~50 nm are routinely measured with standard Bruker probes.

Precise Nanoscale Mechanical Mapping

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● Unlimited Industrial Metrology Applications

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
High-Brightness LED and Solar Materials

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

Semi-automated 24/7 operations using a Dimension AFM to characterize 3D geometry for process control.

Modulus map of a multi-component polymer blend imaged using PeakForce QNM. There are three different components clearly present, the light blue component (A), the darker blue component (B), and the red/black component (C). (7 μm scan).

Bruker's AFMs and high-performance industrial probes provide consistent and reliable measurements for the most demanding critical-to-quality nanoscale measurements.

Dimension Industrial AFMs

Most Reliable Platform and Unrivalled Support

Bruker has the largest service and applications support network of any AFM supplier, with regional customer support and programs for extended assistance. Major industrial companies and leading data storage and semiconductor manufacturers throughout the world have chosen Bruker systems for our reputation as the most reliable precision metrology platforms and quick, knowledgeable support.

		Dimension HPI		Dimension Pro (300 mm)	
		Icon	FastScan	Icon	FastScan
STAGES	XY accessible area, without manual interaction (mm)	150x105	150x150	210x200 300x300	210x230 300x300
	Sample size (mm)	210x210		200x200 or 300x300	
	XY stage repeatability, bi-directional, (µm)	3		6	
MODES	SCM, CAFM, TUNA, SSRM, sMIM, nano-TA, SThM	●	○	●	○
	PF-TUNA, PF-KPFM, PF-sMIM	●	○	●	○
	AutoMET software	●	●	●	●
	PF-QNM	●	●	●	●
	Contact, Tapping, PeakForce Tapping, EFM, MFM	●	●	●	●
SPECIFICATIONS*	X-Y Scan Range (µm)	≥ 90x90	≥ 35x35	≥ 90x90	≥ 35x35
	Z Range (µm)	≥ 12	≥ 3	≥ 12	≥ 3
	Vertical Noise Floor - RMS/Adev (pm)	≤ 30/25 <small>Height in appropriate environment, typical imaging BW (up to 625Hz)</small>	≤ 40/32 <small>Sensor in appropriate environment, typical imaging BW (up to 625Hz)</small>	≤ 35/28 <small>Height in appropriate environment, typical imaging BW (up to 625Hz)</small>	≤ 45/36 <small>Sensor in appropriate environment, typical imaging BW (up to 625Hz)</small>
	X-Y Position Noise, (closed-loop) (pm)	≤ 125/100	≤ 145/116	≤ 125/100	≤ 145/116
	Z Sensor Noise, - RMS/Adev (pm)	≤ 20/15 <small>Typical imaging BW (up to 625Hz)</small>		≤ 20/15 <small>Typical imaging BW (up to 625Hz)</small>	
	Laser and PD alignment	Manual	Automated	Manual	Automated

Legend: Included/Optional ● Not Available ○
* Specifications for Dimension PRO are for 300 mm stage option.

Bruker Nano Surfaces Division is continually improving its products and reserves the right to change specifications without notice. Dimension, FastScan, Icon, NanoScope, PeakForce KPFM, PeakForce QNM, PeakForce Tapping, PeakForce TUNA, and ScanAsyst are trademarks of Bruker Corporation. All other trademarks are the property of their respective companies. © 2019 Bruker Corporation. All rights reserved. B088; Rev. A0

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