



# Elemental Analysis for the Iron, Steel, and Metal Industry by OES, XRF and Combustion/Fusion

Bruker Webinar  
July 09, 2020



# Elemental Analysis for Iron, Steel and Metal Industry



## Topics

- Introduction to steel and steel making process  
*Iron & steel treated as representatives for all other metals*
  - What is steel?
  - Workflow in an integrated steel plant
  - Individual strengths of the 3 complementary elemental analysis techniques: OES, XRF, Combustion/Fusion
- Full elemental analysis of metals by Optical Emission Spectrometry (OES)
  - Sample taking & preparation
  - OES principles
  - Examples with Q8 MAGELLAN
  - Automation with OES – Q8 MAGELLAN *online*
- Material and slag analysis by XRF
  - Principles & differentiation (ED vs WD)
  - Analysis of Raw-, Intermediate and final materials, online XRF
  - Optimized process control by slag analysis
  - Final material inspection

# Elemental Analysis for Iron, Steel and Metal Industry



## Topics

- Raw material- and slag-analysis by XRF
  - Principles & benefits
  - Raw materials analysis
  - Optimized process control by slag analysis
- Fast and precise analysis of C, S, O, N, H by combustion / fusion
  - Principles of carbon and sulfur analysis by combustion
  - Principles of oxygen, nitrogen and hydrogen analysis by fusion
  - Applications in the steel making process
  - Examples
  - Analysis of finished materials (e.g. welding rods) for diffusible hydrogen

# Meet your speakers



## Optical Emission Spectrometry OES



**Peter Papelewski**  
Product Manager OES, CS/ONH  
Karlsruhe, Germany

## X-ray Fluorescence Spectrometry XRF



**Adrian Fiege**  
Product Manager XRF  
Karlsruhe, Germany


## Combustion/Fusion Analysis CS/ONH



**Kristin Odegaard**  
Technical Sales Specialist  
CS/ONH, OES  
Madison, WI, USA

# Introduction to steel and steel making

## What is steel?



- **Steel**  
General name given to a large family of **alloys of iron with carbon**, containing less than 2.1% carbon and other elements
- Steel is such a versatile material because we can adjust its **composition** and internal **structure** to tailor its properties for certain applications
- More than 3.500 different grades of steel, of which **60% are younger than 4 years!**
- If the Eiffel Tower were to be rebuilt today: only one-third of the steel originally used would be needed

# Introduction to steel and steel making

## Alloys of iron with carbon




The final properties are influenced by

- Alloying: Adding other elements such as: **Si, Mn, C, P, S, Al, Ti, Mo, Cu, Cr, Ni, etc.** during furnace tapping (BOF, AOD, EAF) and/or in secondary metallurgy (Ladle furnace, RH Degasser, VD/VOD-Plant, CAS-OB , etc.)
- By rolling process (incl. heat treatments)

■ Ultra High Strength Steel  
■ Extra High Strength Steel  
■ Very High Strength Steel  
■ High Strength Steel  
■ Mild Steel / Forming Grades  
■ Aluminium  
■ Magnesium



# Introduction to steel and steel making Workflow in an integrated steel plant



# Workflow in an integrated steel plant

## Raw materials



### 1. Raw materials

- Deposited in massive piles, ready to be processed
- Iron ore is analyzed immediately as it arrives to ensure it is the grade and quality ordered (and paid for)
- All input materials are analyzed for quality and purity
- Coke and associated waste are analyzed for quality and discharge purposes
- Sinter is analyzed prior going into the Blast Furnace
- Scrap is analyzed before going into Electric Arc Furnace




Iron Ore ( $\text{Fe}_x\text{O}_y$ )



Tata Scunthorpe, UK



# Introduction to steel and steel making Workflow in an integrated steel plant



# Workflow in an integrated steel plant

## Blast Furnace



### 2. Blast Furnace

- Aim: Reduce iron oxides to metallic iron
- The liquid iron is called “hot metal” (“pig iron” when solidified)
- Hot metal temperature: 1340 – 1380 deg C
- Hot metal transported by torpedo ladles to the BOF converter
- Important elements in hot metal are:

Element	Concentration Level
Carbon, <b>C</b>	4.2 - 4.5 %
Silicon, <b>Si</b>	0.2 - 1.2 %
Manganese, <b>Mn</b>	0.2 - 0.5%
Sulfur, <b>S</b>	0.01 - 0.0.7%
Phosphorus, <b>P</b>	0.05 – 0.13 %





Figure 1 - Blast Furnace Operation

<https://www.metalfurnace.ooo/2018/09/blast-furnace.html>



# Introduction to steel and steel making Workflow in an integrated steel plant




# Workflow in an integrated steel plant

## Primary & Secondary Metallurgy



### 3. Primary & Secondary Metallurgy

- Carbon reduction and impurity removal:
  - BOF converter: loaded with scrap & hot metal, addition of fluxes to form slags, a lance blows oxygen to lower concentrations of: **C, Si, P, Mn**
  - Electric Arc Furnace (EAF): loaded with scrap & direct reduced iron, addition of fluxes
- Fluxes (burnt-lime, dolomite, etc.)
  - Remove impurities by formation of slags, floating on top of the melt, slag reduces energy consumption and protects furnace lining refractories
  - Slag and steel are separated during tapping
- **Secondary metallurgy (ladle furnace refining)**
  - Goal: Deliver liquid steel to the caster...
  - **On time**, with **right temperature**, meeting total **chemical specifications**



# Workflow in an integrated steel plant

## Secondary Metallurgy



### 3. Secondary Metallurgy – Ladle Furnace


- Deoxidation of steel and slags (add Si, Al)
- Addition of alloying elements (Cr, Ni, Mo, V, ...)
- Removal of impurities: refining (removal of S)
- Homogenization (temperature and chem. composition)
- Control of (non-metallic) inclusions
- Determines timing & productivity
- *Slag-Metal-Refractories Interactions*
- *Mass balance calculations*



# Introduction to steel and steel making Analytical Techniques



- Integrated Plants (BOF), larger producers:  
typical minimum requirements: 1 XRD, 2 XRF, 3 OES, 2 CFA (CFA: combustion/fusion analysis)
- Mini-mills Plants (EAF), medium & small producers (incl. stainless steel\*)  
typical minimum analytical requirements: 1 XRF, 2 OES, 2 CFA, 1 XRD\*




# Elemental Analysis for the Iron, Steel and Metal Industry

## 3 Complementary Techniques



The techniques are complementary and complete each other

Metal analysis for high alloys,  
slag analysis



Fast and reliable, complete  
chemical composition in seconds



C, S, O, N, H

Bruker has a complete portfolio of techniques required by the steel and metals industry to provide a complete solution in elemental analysis

# Elemental Analysis for the Iron, Steel and Metal Industry

## Characteristics of each technique



	OES	XRF	CFA (CS/ONH)
Element range	Wide element range	Wide element range	Limited (but some elements are not accessible by other techniques)
Concentration range	Perfect for metals: sub-ppm to high % (difficult @ high concentrations)	Very good: Sub-ppm to 100 %	Full: Sub-ppm to 100 %
<b>Speed</b>	Very fast (~30 s)	Medium to slow (depends on class and technique)	Fast (< 1 min)
<b>Sample Types</b>	Solid pieces, conductive material ⇒ metals	Any liquid or solid (no gases)	Inorganic Solids (no organics, no liquids)
Sample Forms	Solid pieces	any	any
Sample Prep. Requirements	Low - medium	Medium - high (depends on application, low compared to ICP)	Low - medium
Suitable for trace / impurities	yes / yes	yes / (yes) (depends on application and class)	Yes / yes
<b>Statistical Representation</b>	Medium (low ablation depth)  Calibration & method matrix dependent	Low (low penetration depth in metals)  Calibration matrix dependent	Complete (entire sample)  Calibration matrix independent
Operation / Acquisition costs	Low / low - medium	Low / medium - high	High / low - medium



# Audience Poll



Which of the following complementary analytical techniques for metals were you aware that Bruker offers? (Check all that apply.)

- Optical Emission Spectrometry (OES)
- X-ray Fluorescence (XRF)
- Combustion & Fusion Analysis for CS/ONH
- Handheld XRF ("guns")
- X-ray Diffractometry (XRD)
- Other (e.g. NMR, FTIR/NIR, ...)





# Elemental Analysis for the Iron, Steel & Metal Industry

## **Optical Emission Spectrometry (OES)**



# Elemental Analysis for the Iron, Steel and Metal Industry

## Producing samples for elemental analysis



### **Sample taking & sample preparation are both essential steps in the analytical process!**

- Must be homogeneous and representative
- Should be free from inclusions, cavities, fissures, burrs and ridges
- If taken from melt, cooling process should be controlled to obtain samples with defined structures; rapid cooling is usually advantageous



# Elemental Analysis for the Iron, Steel and Metal Industry

## Sample shapes in metal industry



<b>Non-Ferrous Metals</b>	<b>Metals – Fe, Ni, Co Bases</b>	<b>Cast Irons</b>	<b>Traces</b>	<b>Steel</b>
<p>The mushroom form sample taking is the most used for non-ferrous metals.</p>	<p>Very simple mold, more used for alloys of Fe, Ni and Co bases.</p>	<p>The mold is specific adapted for general casting irons, and pig iron in particular. Very high cooling speed.</p>	<p>The ring mold type is recommended for casting pure metals (trace analysis) or metals that have low tendency to segregation. The sample is homogeneous on its surface, but only on some few millimeters in depth.</p>	<p>The sample taking can be greatly simplified with a <i>SPEMIS</i> probe that can take the liquid metal directly from the casting.</p>

# Elemental Analysis for the Iron, Steel and Metal Industry

## Sample shapes for steel




Figure 4.4: Different sample shapes (printed with friendly permission of the company MINKON GmbH, Heinrich-Hertz-Str. 30–32, D-40699 Erkrath, Germany).

# Elemental Analysis for the Iron, Steel and Metal Industry

## Sample preparation



### Requirements


- Sample as homogenous as possible
- Clean & flat surface for reliable and reproducible measurements

### Grinding (disc or belt grinder)

- For Fe, Co, Ni, Ti, Cu alloys
- Grid 40-80 of
  - corundum paper (most common & economic)
  - silicon carbide (beware of Si, C content)
  - zirconium oxide (long lifetime, more expensive)

### Milling

- For Al, Mg, Pb, Sn, Zn, Cu & Fe alloys
- Milling machine, clamping unit, cutting plates must suit the application (e.g. hard skin layer on cast iron)






# Full Alloy Analysis by OES

## Physical Principles




- A metallic (conductive) sample is placed above an electrode; a high frequency (AC), high voltage is applied between electrode and sample  
 $\Rightarrow$  high energetic sparks are discharged onto the sample
- The sparking happens in a pure atmosphere of argon (Ar), Ar gets ionized and a stable Ar plasma is formed
- The spark hitting the sample surface ablates/atomizes parts of the sample
- These generated atoms are not in the electronic ground state but in an excited state
- When an electron "falls" back to the ground state, a photon of specific energy (wavelength) is emitted



# Full Alloy Analysis by OES

## Physical Principles



- Emitted light is guided into the optical system and dispersed by a diffraction grating.
- Detectors on specific positions on the Rowland circle are quantifying emitted light.

Useful wavelength in OES: 120 – 800 nm



- The intensity of a specific line correlates directly to the abundance (concentration) of the element in the sample.



# Full alloy analysis by OES Calibration



- The calibration process involves determining the exact mathematical relationship (function) between the intensity ( $I$ ) of the light emitted at the certain wavelength ( $\lambda$ ) and the concentration ( $c$ ) of the element producing the light.
- In fact, OES uses intensity ratios (elemental / matrix reference) to improve calibration and long term stability.
- Different emission lines with different sensitivities might be used for trace vs. bulk analysis.
- Factory calibration is built with certified reference material, but possible to include e.g. your own standards



# Full Alloy Analysis by OES Product Overview



## **Q2 ION**

*Ultra compact metals analyzer*

- Patented optics with AAC for instant operation readiness
- Incoming Material Inspection, PMI and Quality Assurance
- Ultra compact, portable, affordable

## **Q4 TASMAR**

*Advanced multi-optics benchtop OES*

- 3 models cover every QA/QC task
- Outstanding analytical performance with low cost of ownership
- Pneumatic sample clamp for productivity and ease-of-use

## **Q8 MAGELLAN**

*High-end PMT-based vacuum spectrometer*

- Best LODs and stability
- One-button operation (automatic clamp)
- Q8 MAGELLAN *online* ready for full automation

# Full Alloy Analysis by OES At a glance: Q8 MAGELLAN






## Designed for demanding applications

- Paschen-Runge mount, 750 mm focal lengths
- Photomultiplier (PMT) technology, up to 80 channels
- High vacuum optic (no purge)
- Single-Spark reading and evaluation (SSE)

## Figures for Fe-base

- Analysis time:  $\sim 25$  s
- Burn-spot size:  $\sim 4$  mm (spark-stand plate opening: 12 mm)
- Penetration depth:  $\sim 18$   $\mu\text{m}$
- Ablation rate:  $\sim 70$  ng/spark (average)
- Total sparks for 1 "burn": 4800-6000
- Total sample mass:  $\sim 0.4$  mg / burn



# Full Alloy Analysis by OES Fe 110 (low alloy steel)



ELEMENTAL ANALYSIS - User: AppDev - Instrument: Q8S\_B

File View Analysis Help

METHODS EDITOR ANALYSIS

Flush

Sample: St 52 / 1.0580 /

	C [%]	Si [%]	Mn [%]	P [%]	S [%]	Cr [%]	Mo [%]	Ni [%]	Al [%]	Co [%]
2	0.152	0.207	1.256	0.012	0.012	0.087	0.015	0.092	0.027	0.0081
3	0.153	0.205	1.256	0.012	0.012	0.086	0.014	0.090	0.027	0.0081
4	0.153	0.207	1.259	0.012	0.012	0.087	0.015	0.091	0.027	0.0082
5	0.152	0.206	1.254	0.012	0.012	0.089	0.016	0.094	0.028	0.0084
∅	0.22	0.55	1.6	0.045	0.045					
∅	<b>0.152</b>	<b>0.206</b>	<b>1.255</b>	<b>0.012</b>	<b>0.012</b>	<b>0.087</b>	<b>0.015</b>	<b>0.092</b>	<b>0.027</b>	<b>0.0082</b>
∅										
σ	0.0010	0.0008	0.0029	0.00010	0.0003	0.0015	0.0008	0.0018	0.0003	0.0001
∪	0.67	0.38	0.23	0.83	2.86	1.72	5.17	1.92	1.22	1.71

	Cu [%]	Nb [%]	Ti [%]	V [%]	W [%]	Pb [%]	Sn [%]	Mg [%]	As [%]	Zr [%]
2	0.202	0.0043	0.0013	0.0013	0.0027	0.0027	0.017	0.0003	0.0064	0.0006
3	0.201	0.0043	0.0013	0.0013	0.0027	0.0028	0.016	0.0003	0.0063	0.0006
4	0.202	0.0044	0.0013	0.0013	0.0031	0.0028	0.017	0.0003	0.0066	0.0006
5	0.201	0.0043	0.0013	0.0013	0.0032	0.0030	0.017	0.0004	0.0067	0.0007
∅	<b>0.201</b>	<b>0.0043</b>	<b>0.0013</b>	<b>0.0013</b>	<b>0.0029</b>	<b>0.0028</b>	<b>0.016</b>	<b>0.0003</b>	<b>0.0065</b>	<b>0.0006</b>
∅										
σ	0.0006	0.00004	0.00002	0.00005	0.0002	0.00009	0.0003	0.00002	0.0002	0.00007
∪	0.31	0.99	1.26	3.84	8.57	3.29	1.77	5.92	2.34	11.17

	Bi [%]	Ca [%]	Ta [%]	B [%]	La [%]	N [%]	O [%]	Fe [%]	Ceq [%]
2	<0.0005	0.0020	0.0097	0.0009	<0.0005	0.0089	<0.0020	97.87	0.540
3	0.0005	0.0022	0.0091	0.0010	<0.0005	0.0059	<0.0020	97.88	0.538
4	0.0006	0.0022	0.016	0.0011	<0.0005	0.010	<0.0020	97.86	0.541
5	0.0007	0.0025	0.010	0.0013	<0.0005	0.011	<0.0020	97.86	0.545
∅	<b>0.0005</b>	<b>0.0023</b>	<b>0.011</b>	<b>0.0010</b>	<b>&lt;0.0005</b>	<b>0.0085</b>	<b>&lt;0.0020</b>	<b>97.87</b>	<b>0.540</b>
∅									
σ	0.0001	0.0002	0.0027	0.0001	0.00003	0.0022	0.0016	0.012	0.0039
∪	21.67	8.16	23.84	12.39	10.63	26.17	3.04	0.01	0.72

Activate Method Complete Analysis Report Analysis Reject Analysis

View History View Analyses

PMI STDZ Stop Start

Fe Fe110 Routine Element Concentrations Mask: -500 Waiting for Start

Bruker AXS Analysis Report

Sample: St 52 Material: 1.0580 Method: Fe 110

Function No: 100239 Head No: User: [ ]

Analysis Time: 2020-07-01 10:22:39

	C	Si	Mn	P	S	Cr	Mo
∅	0.225	0.206	1.255	0.045	0.045	0.087	0.015
∅	0.225	0.206	1.255	0.045	0.045	0.087	0.015
w	1.47	0.38	1.26	0.83	2.86	1.72	5.17
∅	0.152	0.206	1.254	0.012	0.012	0.087	0.015
∅	0.152	0.207	1.259	0.012	0.012	0.087	0.015
∅	0.152	0.205	1.254	0.012	0.012	0.089	0.016
∅	0.152	0.207	1.259	0.012	0.012	0.087	0.015
∅	0.152	0.206	1.254	0.012	0.012	0.089	0.016

	Ni	Al	Co	Nb	Ti	V
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013
w	1.92	1.22	1.71	0.99	1.26	3.84
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013
∅	0.092	0.027	0.0082	0.0043	0.0013	0.0013

	W	Pb	Sn	Mg	As	Zr	Bi
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007
w	8.57	0.29	1.77	0.92	2.34	11.17	21.67
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007
∅	0.0029	0.0003	0.0016	0.0003	0.0065	0.0006	0.0007

	Ca	Ta	B	La	N	O	Fe
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87
w	8.16	23.84	12.39	10.63	26.17	3.04	97.87
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87
∅	0.0023	0.011	0.0010	<0.0005	0.0085	<0.0020	97.87

2020-07-01 10:44:27 1/2

# Full Alloy Analysis by OES Fe 130 (Cr/CrNi)



ELEMENTAL ANALYSIS - User: AppDev - Instrument: Q8S\_B

File View Analysis Help

METHODS EDITOR ANALYSIS

Flush Sample: X2 CrNi 19 11 - 304L / 1.4306 /

Name	Unit	Low (Alarm)	Average	High (Alarm)	Abs. Std. Dev.	Rel. Std. Dev.	1	2	3	4	5
C	%	0	0.017	0.03	0.0013	7.69	0.015	0.016	0.017	0.017	0.018
Si	%	0	0.343	1	0.0016	0.46	0.343	0.344	0.344	0.344	0.340
Mn	%	0	1.102	2	0.0070	0.64	1.096	1.094	1.102	1.112	1.104
P	%	0	0.023	0.045	0.0006	2.60	0.023	0.023	0.024	0.024	0.023
S	%	0	0.030	0.03	0.0005	1.52	0.030	0.029	0.030	0.030	0.030
Cr	%	18	18.18	20	0.036	0.20	18.20	18.14	18.15	18.20	18.23
Mo	%	0	0.191	0.2	0.0015	0.79	0.190	0.190	0.193	0.193	0.192
Ni	%	10	10.08	12	0.050	0.50	10.05	10.05	10.11	10.16	10.04
Al	%		0.0068		0.0004	5.24	0.0068	0.0062	0.0071	0.0069	0.0068
Co	%		0.149		0.0008	0.50	0.149	0.150	0.150	0.149	0.149
Cu	%		0.467		0.0037	0.79	0.468	0.462	0.468	0.472	0.464
Nb	%		0.014		0.0001	1.04	0.014	0.014	0.014	0.014	0.014
Ti	%		0.0050		0.00006	1.15	0.0050	0.0050	0.0051	0.0050	0.0051
V	%		0.081		0.0004	0.48	0.080	0.081	0.081	0.081	0.081
W	%		0.010		0.0003	2.67	0.010	0.0099	0.0099	0.0098	0.010
Sn	%		0.0067		0.0001	1.84	0.0066	0.0066	0.0068	0.0066	0.0068
As	%		0.0053		0.0002	3.00	0.0051	0.0054	0.0053	0.0052	0.0055
Se	%		0.0062		0.0010	16.19	0.0067	0.0066	0.0044	0.0069	0.0063
B	%		0.0018		0.0002	8.45	0.0017	0.0016	0.0018	0.0019	0.0020
N	%	0	0.038	0.11	0.0018	4.70	0.038	0.041	0.038	0.037	0.036
Fe	%	64.3	69.24	71.5	0.073	0.11	69.26	69.33	69.24	69.12	69.23

Activate Method Complete Analysis Repo View History View Analyses PMI STDZ Stop Start

Fe Fe130 Routine Element Concentrations Mask: -500 Waiting for Start

Braker AXS Analysis Report

Sample: X2 CrNi 19 11 - 304L Material: 1.4306 Method: Fe130  
 Furnace No: 8000 Heat No: 0000 Order No:  
 Analysis Time: 2020-07-01 09:33:24 User: Unit: [%]

	C	Si	Mn	P	S	Cr	Mo
11	0.030	1.000	2.000	0.045	0.030	20.00	0.200
⊙	0.017	0.343	1.102	0.023	0.030	18.18	0.191
11	0.000	0.000	0.000	0.000	0.000	18.00	0.000
σ	0.0013	0.0016	0.0070	0.0006	0.0005	0.036	0.0015
u	7.69	0.46	0.64	2.60	1.52	0.20	0.79
1	0.015	0.343	1.096	0.023	0.030	18.20	0.190
2	0.016	0.344	1.094	0.023	0.029	18.14	0.190
3	0.017	0.344	1.102	0.024	0.030	18.15	0.193
4	0.017	0.344	1.112	0.024	0.030	18.20	0.193
5	0.018	0.340	1.104	0.023	0.030	18.23	0.192
11	10.08						
⊙	10.08	0.0068	0.149	0.467	0.014	0.0050	0.081
11	10.00						
σ	0.050	0.0004	0.0000	0.0037	0.0001	0.00006	0.0004
u	5.24	5.24	0.50	0.79	1.04	1.15	0.48
1	10.05	0.0068	0.149	0.468	0.014	0.0050	0.080
2	10.05	0.0062	0.150	0.462	0.014	0.0050	0.081
3	10.11	0.0071	0.150	0.468	0.014	0.0051	0.081
4	10.16	0.0069	0.149	0.472	0.014	0.0050	0.081
5	10.04	0.0068	0.149	0.464	0.014	0.0051	0.081
11							
⊙	0.010	0.0067	0.0053	0.0062	0.0018	0.038	0.041
11						0.000	0.037
σ	0.0003	0.0001	0.0002	0.0010	0.0002	0.0018	0.0018
u	2.67	1.84	3.00	16.19	8.45	4.70	4.70
1	0.010	0.0066	0.0051	0.0067	0.0017	0.038	0.041
2	0.0099	0.0066	0.0054	0.0066	0.0016	0.041	0.041
3	0.0098	0.0068	0.0053	0.0066	0.0018	0.038	0.038
4	0.0098	0.0066	0.0052	0.0069	0.0019	0.037	0.037
5	0.010	0.0068	0.0055	0.0063	0.0020	0.036	0.036

2020-07-01 10:59:39 1/1





# Elemental Analysis for the Iron, Steel & Metal Industry

## **X-ray Fluorescence (XRF)**



# Audience Poll




What are your main goals with elemental analysis – What are you trying to achieve? (Check all that apply.)

- Check quality of incoming material
- Increase percentage of recycled material in production
- Optimize production costs
- Increase throughput
- Optimize final product quality
- Other (specify)



# EDXRF or WDXRF?

## Depends on the requirements!



### WDXRF

- Sophisticated setup
- Best resolution and sensitivity
- Detects intensity at a given wavelength


### EDXRF

- Simple, robust, compact
- Detects energy and intensity





# XRF for Raw Materials



# S2 KODIAK – Online-Analysis of Scrap Sorting and Certification



- Steel producers aim toward **circular economy**
  - Increase the use of recycled steel to reduce CO<sub>2</sub> emission
- **Tramp elements** (Cu, Cr, Ni, Mo...) need to be monitored
  - The thresholds depend on the production process
    - Blast furnace / converter
    - Electric arc furnace
- **Online elemental analysis** to
  - Enhance the sorting procedure
  - Determine the average composition of a truck load
- **Integration:** Sicon Germany/America



time	t/h	Fe %	Cu %	Mn %	Mo %	Ti %	Zn %
6 AM	72.20	97.14	0.406	0.301	0.004	0.004	0.045
7 AM	79.41	97.11	0.432	0.292	0.004	0.004	0.053
8 AM	62.45	97.06	0.499	0.295	0.004	0.004	0.049
9 AM	71.71	97.06	0.478	0.300	0.004	0.004	0.052
...	...	...	...	...	...	...	...
Sum	652.0	--	--	--	--	--	--
AVG	--	97.10	0.46	0.30	0.01	0.01	0.05

*Exemplary report*

# S2 KODIAK

## The perfect fit for online applications



Online EDXRF multi-element analyzer for real-time results for Mining, **Material Sorting / Recycling**, Coatings and more.

- Encapsulated, ruggedized design; *tested*:
  - Shock & Vibration
  - Temperature
  - Water spray and Condensing water
- Enhanced system cooling
- Automatic system startup and self alignment
- Autonomous operation 24/7
- Integrated UPS, camera, distance sensor
- TCP/IP data transfer to plant control SW
- Maintenance-free operation
- No radioactive sources



# S8 TIGER for Incoming Material Inspection Iron Ore/Sinter




- According to ISO 9516 (WDXRF)
- Fused beads
  - At commercial and service labs
  - At central labs
  - At customs
- Maximum 60 s measurement time
- Excellent long-term precision for 200 measurements

Sample	Average (%)	Abs. Std. Dev. (%)	Rel. Std. Dev. (%)
Fe (%)	47.065	0.013	0.028
SiO <sub>2</sub> (%)	0.673	0.006	0.869
P (%)	0.031	0.001	0.746
Al <sub>2</sub> O <sub>3</sub> (%)	4.081	0.014	0.334
Mn (%)	15.348	0.016	0.102
CaO (%)	0.032	0.001	1.867
MgO (%)	0.032	0.002	5.401
TiO <sub>2</sub> (%)	0.023	0.003	11.705
K <sub>2</sub> O (%)	0.160	0.001	0.484



# GEO-QUANT Iron Ore Solution



- **GEO-QUANT Iron Ore contains**
  - 16 certified reference materials (no synthetic standards)
  - Evaluation samples
  - User Manual
  - Preparation Manual
- **GEO-QUANT Iron Ore requires**
  - LIF 220 crystal
  - XS 100 crystal
    - Covers F to Cl
    - Increased sensitivity compared to the XS-CEM
    - Enhanced long term stability compared to the PET
    - Not temperature sensitive
    - **Shorter measurement time for elements from F to Cl**


# S6 JAGUAR - Benchtop WDXRF Coal and Coke



Coal is widely used in metal making (aluminum, steel) and power generation (electricity)

Analysis of coal, coke and carbon products is vital:

- Prevent contamination of metal products
- Inhibit (steel) corrosion (monitoring of Cl)
- Reduce environmental impact (reducing S content)



# S6 JAGUAR - Benchtop WDXRF Coal and Coke




Element	XRF conc.	Cert. Conc.	Abs Std Dev.	Rel Std Dev.
S [%]	3.40	3.30	0.02	0.72
Ni [ppm]	128.00	124.00	4.48	3.50
Si [ppm]	24.00	28.00	1.15	4.80
Fe [ppm]	276.00	266.00	2.65	0.96
Na [ppm]	623.00	645.00	43.61	7.00
Al [ppm]	153.00	150.00	5.66	3.70
Ca [ppm]	112.00	107.00	2.35	2.10
K [ppm]	17.00	17.00	1.84	10.80
Cl [ppm]	100.00	n.a.		
Ti [ppm]	5.00	4.00	0.31	6.10
Zn [ppm]	40.00	41.00	0.10	0.25
V [ppm]	302.00	300.00	2.33	0.77

Optimal accuracy for:

- Minor elements, such as S
- Traces, e.g. Si, Cl, Fe

# XRF for Production Monitoring





# S2 PUMA Series 2

## Slag analysis made fast and reliable



Slag composition needs to be monitored to achieve e.g. optimal metal purity and least attack on refractory lining.

There are different types for slags in steel production, e.g.:

- Electric Arc Furnace (EAF) Slag
- Ladle Metallurgy Furnace (LMF) Slag

The difference in steel production procedure results in different slag compositions and, thus, in different requirements for analytical solutions.

Reliable results are required within minutes!



# S2 PUMA Series 2

## Slag analysis made fast and reliable



### Optimal excitation of the sample is ensured by:

- High power 50-Watt X-ray tube
- Closely coupled optics
- 10-position primary beam filter
- The Next generation silicon drift detectors (SDD)


### S2 PUMA Series 2 Single

- Convenient and fail-safe sample handling
  - Load the sample, touch the button and get results within minutes!
- Ease-of-use thanks to TouchControl™
- Sample Rotation for heterogeneous samples
- SampleCare™ protects vital system components
- Dedicated Light Element configuration (e.g. for F)
  - Optional vacuum mode for best light element performance and lowest cost of operation (no He)




# S2 PUMA Series 2

## Slag analysis made fast and reliable



- Calibration curve for MgO



- Calibration peaks for Mg K $\alpha$ 1, Al K $\alpha$ 1 and Si K $\alpha$ 1

# S2 PUMA Series 2

## Slag analysis made fast and reliable




- Equipped with the HighSense™ LE detector, the S2 PUMA Series 2 achieves outstanding performance – even for light elements and with a single analytical range – in just **2 minutes** counting time!
- Precision is typically < 1 % RSD

Compound	LMF Slags [wt%]
F	0.03 – 7.9
MgO	0.2 – 21.2
Al <sub>2</sub> O <sub>3</sub>	1.8 – 38.6
SiO <sub>2</sub>	7.4 – 51.4
P <sub>2</sub> O <sub>5</sub>	0 – 1.6
S	0 – 1.2
CaO	0.6 – 60.4
TiO <sub>2</sub>	0.01 – 2.2
MnO	0.06 – 14.9
FeO	0.1 – 17.2


Elements	Voltage [kV]	Current [mA]	Measurement time [s]	Beam Filter	Mode
F, Mg, Al, Si, P, S, Ca, Ti, Mn, Fe	20	automatic*	120	none	Vacuum

\*Current is maximized automatically for best count statistic.

# XRF for Quality Control of Steel Products





# S6 JAGUAR HighSense™ Goniometer: High Resolution



The S6 JAGUAR with WDXRF HighSense Goniometer exceeds EDXRF-based systems in resolution and analytical precision

# S6 JAGUAR versus EDXRF F in Foundry Products



Strong overlap of F Ka and Fe La with **EDXRF** leads to medium accuracy and precision: **Min 3.59 % -> 3.78 % <- Max 4.07 %**

**S6 JAGUAR**: Optimal resolution, clear separation of both lines, high sensitivity with 400 W power: **Min 3.97 % -> 4.03 % <- Max 4.07 %**

# S6 JAGUAR

## Low Alloy Steels



Quality Control of low alloy and mild steels:

Analysis of Fe, Ti, V, Cr, Mn, Co, Ni, Cu, Mo, W


Impurities with negative impact on steel quality:

Al, Si, P, S, Cu, As, Sb, Pb

- Incoming inspection of raw material for manufacturing
- Specifying alloy types



Current WDXRF spectrometers are suffering for major elements due to saturation of the scintillation counter

Sum, escape peaks and crystal artefacts will even worsen the situation (Pulse Height Distribution)





# S6 JAGUAR Low Alloy Steels



- Separation of S Ka from Mo La with high resolution WDXRF
- Handling of line overlaps and FP matrix correction with variable alphas

**New solid state HighSense XE detector** replaces conventional scintillation detector:

- State-of-the-art technology
- Medium and heavy element detection
- Linear range of 2 Mcps
- 2 times better energy resolution

# S8 TIGER for highest throughput Cast iron



**Process Control (CAST IRON)** with the  
S8 TIGER, 4 kW:

**15 elements:**

C, Si, P, S, Ti, V, Cr, Mn, Fe, Ni,  
Cu, Zn, As, Sn, Pb

Time-to-result:

**1 min 47 sec.** per sample

Sample throughput:

**30 samples/h**



Fast, accurate and precise process and  
quality control of cast iron, including  
carbon analysis

**HIGHEST SPEED ANALYSIS**

# S8 TIGER for highest throughput Cast iron



Sample	Si (%)	Mn (%)	P (%)	S (%)	V (%)	Ti (%)	C (%)
cast	1.488	3.08	0.321	0.175	0.036	0.265	3.74
cast/1	1.488	3.08	0.320	0.175	0.035	0.264	3.78
cast/2	1.488	3.08	0.323	0.174	0.036	0.265	3.78
cast/3	1.488	3.08	0.323	0.175	0.035	0.266	3.78
cast/4	1.491	3.08	0.324	0.175	0.035	0.265	3.79
cast/5	1.483	3.08	0.323	0.175	0.035	0.265	3.80
cast/6	1.487	3.08	0.323	0.174	0.036	0.264	3.80
cast/7	1.490	3.08	0.323	0.175	0.036	0.266	3.80
cast/8	1.490	3.08	0.323	0.175	0.036	0.264	3.82
cast/9	1.490	3.08	0.323	0.175	0.036	0.264	3.83
cast/10	1.490	3.08	0.322	0.175	0.036	0.264	3.82
<b>Mean</b>	<b>1.488</b>	<b>3.08</b>	<b>0.323</b>	<b>0.175</b>	<b>0.036</b>	<b>0.265</b>	<b>3.79</b>
Std. Dev.	0.002	0.001	0.001	0.0002	0.001	0.001	0.026
RSD	0.16	0.04	0.36	0.24	1.45	0.31	0.69





# Elemental Analysis for the Iron, Steel & Metal Industry

## **Combustion/Fusion Analysis**




# Audience Poll



Prior to this webinar, were you aware that Bruker offers C/S and O/N/H analyzers?

- Yes
- Somewhat
- Not really
- No



# Elemental Analysis of Light Elements (C,S / O,N,H) Topics




- Introduction to the term Elemental Analysis
  - C/S & O/N/H for *inorganic* materials
    - Basic Principals, Definitions, and Benefits
- Overview of BAXS CGA Product Line
  - C/S (by HF-Induction, NDIR-Detection)
    - ➔ G4 ICARUS HF
  - O/N/H (by Inert Gas Fusion)
    - ➔ G6 LEONARDO, G8 GALILEO

# Elemental Analysis Definitions




**Elemental Analysis: Qualitative & Quantitative** determination of elemental (and isotopic) composition of a material (major, trace and ultra-trace). Not limited to any specific matrix, target analyte or method.



**Qualitative:** What elements are present?

Calibration Curve



**Quantitative:** How much of those elements are present?

# Elemental Analysis Definitions



- **Combustion Gas Analysis:** Uses oxygen as a carrier gas (**Carbon/Sulfur**) and a combustion occurs.
- **Fusion Gas Analysis:** Sample is fused in the absence of Oxygen (but using inert carrier gases like He, Ar, N<sub>2</sub>). No combustion occurs but decomposition/release of gaseous reaction products (**Oxygen, Hydrogen, Nitrogen**).



# The basic principle

## Light element analysis in inorganic solids



- Weighing of dry solid sample



Turnings



Billets



Bars




Powders

- The sample is placed in a crucible and then heated HF-Induction furnace (C/S) or Electrode furnace (O/N/H)



Combustion




Fusion

# Elemental Analysis of Light Elements (C,S / O,N,H) Gas Release and Detection



- The elements: C, S, N, O, H are released from the sample in a stable gaseous form (= molecules) and transported by the carrier gas stream to detectors for quantification. Eventually disturbing by-product are removed upfront
- Detection Principles: IR-absorption (NDIR), Thermal Conductivity (TC) or Mass Spectrometry (MS)



# ONH – why analyze them?



## Oxygen

- Forms **oxidic inclusion** with Al, Si, Ca... in steel
- Inclusions degrade mechanical properties (e.g. impact strength, brittleness)
- Oxygen entering steel while forging causes red shortness



## Nitrogen

- Forms **inclusions** which reduce toughness, aging stability and increase risk due to segregation (> 120 ppm)
- **Alloying** of N (e.g. 0.4%) for stabilizing of austenitic grades (CrNiMn steel) and improving the hardness of Mn steels



## Hydrogen

- Strong influence on mechanical properties of steel: forming pores and cracks
- Can lead to most diverse form of damage: flake formation, pickling, blistering, hydrogen-induced embrittlement



# The importance of C and S analysis



## Carbon

- Most alloying element in steel
- Influences hardness, wear resistance, workability, etc.
- Typical concentrations:
  - Carbon steel:  $\sim 0.05\text{-}2\%$  C
  - Cast iron:  $\sim 2\text{-}4\%$  C
  - Tungsten carbide:  $\sim 6.1\%$  C

## Sulfur

- Undesirable in steel
- Influences brittleness, conductivity, workability, formability, etc.
- Typical concentration in steel:  $\sim 100\text{-}200$  ppm S



# Typical Applications

## C/S, O/N/H



**C/S**

**Application Note**  
Carbon and Sulfur in Coal and Coke  
Analyzer G4 ICARUS HF

**I. Measuring Principle**

For rapid and precise carbon and sulfur measurements the G4 ICARUS HF, using the combustion method with high frequency induction furnace and infrared detection proves highly effective especially with metallic materials and also with inorganic materials like cement, clays and many others.


The solid sample, placed in a ceramic crucible together with accelerator material, is combusted in the high-frequency furnace in an oxygen stream. The carbon and sulfur components in the sample are oxidized to release CO<sub>2</sub> and SO<sub>2</sub> respectively. These measuring components are swept by the carrier gas O<sub>2</sub> to the solid-state NDIR detector system of maximum selectivity and stability.



Coke Oven Plant



Sinter Plant



Blast Furnace

**C/S**



Lab Report CS/ONH 21  
**G4 ICARUS Series 2**  
Fast and Reliable Carbon and Sulfur Determination in Limestone, Dolomite and Lime



Converter



Ladle Furnace



Ladle

**C/S**



Lab Report CS/ONH 22  
**G4 ICARUS Series 2**  
Fast and reliable Carbon and Sulfur Determination in Cast Iron

Continuous Casting



Scrap Yard



Pig Iron



Pellet Plant



Drawing Unit



Wire Rod



Electric Arc Furnace

**O/N**



Lab Report CS/ONH 24  
**G6 LEONARDO**  
Fast and reliable Oxygen and Nitrogen Determination in Steel, Iron, Nickel, and Cobalt Alloys



Reheating Furnace

Reheating Furnace



Cool


**O/N**

**Application Note**  
Oxygen and Nitrogen in Steel  
Analyzer G6 GALILEO

**I. Measuring Principle**

The determination of the oxygen and nitrogen content in steel with the G6 GALILEO is made using the carrier gas method with melt extraction. The sample is heated in a graphite crucible under a flowing inert gas stream. The finely programmable temperature of the electrode furnace is monitored and controlled by a contact-free optical sensor.

The oxygen in the sample, which in general is present in form of oxides, reacts at the applied high temperature and the excess of carbon of the graphite crucible, accordingly by a major reaction to carbon monoxide (CO). Carbonic oxide (CO) in the sample is decomposed evolving N<sub>2</sub> gas.



Rebar, Annealed Wire, Ribbed reinforced metal


**dH**

**Application Note**  
Diffusible Hydrogen in Weld Seams  
G4 PHOENIX DH

**I. Measuring Principle**

The G4 PHOENIX DH has been developed for the determination of the diffusible hydrogen in different sample matrices by means of the carrier gas hot extraction method. The analysis system comprises a rapid heating and cooling infrared cleanshell furnace and/or a wire-heated tube furnace, both equipped with a quartz tube.

For the determination of diffusible hydrogen in weld seams and welding materials according to ENISO 3690 and AWS A4.3 the infrared heated furnace with a tube diameter of 50 mm is used.



Coiled wire, galvanized wire

**H**

**Application Note**  
Hydrogen in Steel  
Analyzer G6 GALILEO


**I. Measuring Principle**

As a result of the harmful influence of hydrogen on the mechanical properties of metal it is due to embrittlement and to the formation of cracks or pores, great quality content in such materials not only for the process control and the product quality but also for the development of new materials.



Structural Profiles

# G4 ICARUS Series 2 Carbon and Sulfur



# Overview of Product Line ON/H



## G8 GALILEO

Most flexible system offered, able to quantify ON/H



## G6 LEONARDO

Most economic system offered, able to quantify O or N or H, ON or OH



# G6 LEONARDO

## Results O/N



Standard B.S. HON T

Oxygen 0.0044 Wt % (+/- 0.0004)

Nitrogen 0.0365 Wt % (+/- 0.0008)

<b>BS-HONT</b> (1g steel pin)		
certified values:		
	O: 0.0044 ( $\pm$ 0.0004)%	
	N: 0.0365 ( $\pm$ 0.0008)%	
Mass / g	Oxygen / %	Nitrogen / %
1.0188	0.0045	0.0364
1.0178	0.0045	0.0366
1.0193	0.0042	0.0359
1.0191	0.0044	0.0369
1.0196	0.0043	0.0366
Mean <sup>3)</sup>	<b>0.0044</b>	<b>0.0365</b>
STD <sup>3)</sup>	0.0001	0.0004



# G6 LEONARDO

## Results O/N



### Standard YSB C 41340B-2011

Oxygen 0.0008 Wt % (+/- 0.00008)

Nitrogen 0.0084 Wt % (+/- 0.0003)

#### YSB C 41340B-2011 (1g ball)

certified values:

O: 0.0008 ( $\pm$  0.00008)%

N: 0.0084 ( $\pm$  0.0003)%

Mass / g	Oxygen / %	Nitrogen / %
0.9543	0.00083	0.0088
0.9542	0.00095	0.0084
0.9535	0.00085	0.0080
0.9549	0.00087	0.0082
0.9539	0.00081	0.0086
Mean	<b>0.00086</b>	<b>0.0084</b>
STD	0.00005	0.0003

# G4 PHOENIX

## Diffusible Hydrogen for Welding



- Diffusible Hydrogen is responsible for material failure in steel and other metals:
  - Hydrogen diffuses through the material and is absorbed into so-called "traps" (voids, pores, grain boundaries, micro-cracks, substituted atoms)
  - Can cause embrittlement of the material
  - **Extremely important for automotive, welding and aerospace!**
- G4 PHOENIX and G8 GALILEO + IR07 can measure the amount of diffusible H in a sample using an IR furnace coupled to either TCD or MS



# Lab Reports Available



Lab Report CS/ONH 23

## G4 ICARUS Series 2

- Fast and reliable Carbon and Sulfur Determination in Ore Concentrates, Metal-Bearing Ores and related materials


Many metals occur naturally in form of sulfide minerals like Pyrite ( $\text{FeS}_2$ ), Chalcocite ( $\text{Cu}_2\text{S}$ ), Argentite ( $\text{Ag}_2\text{S}$ ), Galena ( $\text{PbS}$ ), Sphalerite ( $\text{ZnS}$ ) and many more. Since the amount of economically important metal relative to the ore is small, a first step in the metal extraction process is the separation of target minerals from unwanted rock, gangue and waste minerals. The result, called ore concentrate, is a mixture of the target metal sulfide ores, other sulfide minerals and some gangue.

The sulfur level of the concentrate provides an estimation of the metal-in-concentrate and is used as a quality control measure, determining the efficiency of the concentration process.

In the metal recovery process, determining the carbon and sulfur level in metal-bearing ores is a necessary step for the mine operation to control the metallurgical process kinetics and to assess environmental aspects of the process waste. This is possible due to the fact that most metal-bearing

Elemental Analysis

Innovation with Integrity



Lab Report CS/ONH 22

## G4 ICARUS Series 2

- Fast and reliable Carbon and Sulfur Determination in Cast Iron

### Introduction

The importance of the element carbon in iron and steel making is well known. Properties like hardness, hardenability, ductility and brittleness depend to a large extent on the content of carbon in the material. That makes carbon the main alloying element in cast iron with typical levels greater than 2%. In the iron matrix, carbon can be present in different forms, e.g. dissolved in the matrix, bound to other elements as carbides, or elementary as free graphite.

Sulfur in contrast is typically an unwanted, contaminant element in cast iron. It can e.g. prevent the graphite formation in cast iron and thus impact the mechanical properties and its workability. Quality- and Process-Control require the fast and accurate analysis of the carbon and sulfur content in cast iron and casted products. This Lab Report shows the simplicity, speed and reliability of carbon and sulfur determination by the G4 ICARUS Series 2 using a high-frequency (HF) induction furnace.

Elemental Analysis

Innovation with Integrity

# Lab Reports Available



**G6 LEONARDO**

- Inert Gas Fusion Analysis of Inorganics

Innovation with Integrity

Elemental Analysis



**Lab Report CS/ONH 24**  
**G6 LEONARDO**

- Fast and reliable Oxygen and Nitrogen Determination in Steel, Iron, Nickel, and Cobalt Alloys

**Introduction**  
Oxygen is first used in the steel making process, which is a controlled oxidation process, to remove excess carbon and impurities such as P, Si, and Mn from the hot metal. The level of dissolved oxygen in liquid steel must be lowered because oxygen reacts with carbon, forming CO, which results in excessive porosity during solidification. After the refining, a deoxidation process for removal of residual oxides (inclusions) follows. The nitrogen level present in the alloy has an influence on its mechanical and corrosion properties.

Increasing nitrogen levels can improve hardness and corrosion resistance of certain alloys but too high nitrogen levels can cause embrittlement, low ductility and inconsistent mechanical properties.

Thus, determining the amount of oxygen and nitrogen is critical to process and quality control when making steel but also nickel, and cobalt base alloys. The G6 LEONARDO provides a cost-effective, efficient and reliable means of determining oxygen and nitrogen in these materials.

Innovation with Integrity

Elemental Analysis



# Questions and Answers



Any questions?

Please type any questions you may have for our speakers in the [Q&A panel](#) and click Send.

Thank you!





Like what you learned in this webinar?

Subscribe to Bruker's ***FIRST Newsletter*** to get webinar announcements, technical articles, and X-ray crystallography news delivered right to your inbox.

**Subscribe at:**

<https://www.bruker.com/news/newsletter/first-newsletter-registration.html>



# Webinars – Live and On-Demand

Register for future webinars and view webinar recordings at:

<https://my.bruker.com/acton/media/2655/virtual-events-2020>



## Virtual Events for X-ray Diffraction and Elemental Analysis

There has never been a better time to learn something new.

You've come to the right place. In these trying times of social distancing and working from home, Bruker wants to assure everyone that we are here for you and are committed to keeping you up to date on the latest in our solutions and applications. This webpage will serve as your resource for all of our X-ray Diffraction and Elemental Analysis digital events, including webinars, online trainings, and virtual meetings. We look forward to seeing you online soon!

### Live Webinars


*Jul 15 - Introduction to Analysis of Light Elements (C,S - O,N,H) in Inorganic Materials*



The analysis of light elements can be particularly difficult, and with some analytical techniques, almost impossible. Join us for this 45-minute educational webinar to learn how combustion and fusion gas analysis can obtain accurate and precise results for inorganic materials. Whether you're in the steel, cement, mining, aerospace, or additive manufacturing industry, discover how Bruker AXS gas analyzers make your lab life a lot easier.

[Register now >](#)

*Jul 9 - Elemental Analysis for the Iron, Steel and Metal Industry by OES, XRF and Combustion/Fusion*



Working in the steel, copper, aluminum, or metal recycling industry? Join us for this 1-hour webinar as our specialists introduce 3 major techniques for elemental analysis:

- Optical Emission Spectrometry (OES)
- X-ray Fluorescence Spectrometry (XRF)
- Combustion/Fusion Analysis (CS/ONH)

The strengths and applications of each technique will also be covered.

[Register now >](#)

*Jul 8 - Ways to Improve the Productivity of Your Single Crystal Diffractometer*



In chemical crystallography or structural biology, the best X-ray diffractometers are those that are always in use, consistently producing publication-quality data. Whether you are running an older system or a state-of-the-art instrument, join us for this 45-minute webinar as we show you ways to upgrade to the latest technologies or add accessories that will enhance your diffractometer's versatility and attract more users to your facility.

[Register now >](#)

# Upcoming CS/ONH Webinar



Webinar



Join us for an educational webinar:

## Introduction to Analysis of Light Elements (C,S – O,N,H) in Inorganic Materials

The analysis of light elements can be particularly difficult, and with some analytical techniques, almost impossible. Join us for this educational webinar to learn about how combustion and fusion gas analysis can obtain accurate and precise results for inorganic materials.

Whether you're in the steel, cement, mining, additive manufacturing, or aerospace industry, our gas analyzer product line is sure to make your lab life a lot easier. With our dedicated Application Scientists, Service and Sales Engineers, Bruker is able to help your company meet their analytical needs in the most efficient and affordable way.

During this 40 minute complimentary webinar, specialists from Bruker will:

- Introduce Theory of Elemental Analysis
- Instrumentation for Light Element Analysis
- Bruker's Product Line

The webinar will conclude with an interactive Q&A session, where attendees may ask questions to the panelists.

Register now for the session that best fits your schedule.

[Register for Session 1 \(9 am CDT / 10 pm SGT\)](#)

[Register for Session 2 \(3 pm CDT / 10 pm CEST\)](#)

If you cannot attend, feel free to register anyway, and we will send you a link to the recording to watch at your convenience.

We look forward to seeing you there!

Wednesday, July 15, 2020

### Session 1:

[9:00 am CDT \(Chicago\)](#)  
[4:00 pm CEST \(Berlin\)](#)  
[10:00 pm SGT \(Singapore\)](#)

### Session 2:

[3:00 pm CDT \(Chicago\)](#)  
[10:00 pm CEST \(Berlin\)](#)

### Speakers:



**Kristin Odegaard**

Technical Sales Specialist, CGA & OES  
Madison, Wisconsin, USA



**Christian Zuehke**

Business Development Manager  
Combustion/Fusion Analysis (CS/ONH)  
Karlsruhe, Germany

## Introduction to Analysis of Light Elements (C,S – O,N,H) in Inorganic Materials


Wednesday, July 15, 2020

[Register for Session 1 \(9 am CDT / 10 pm SGT\)](#)

[Register for Session 2 \(3 pm CDT / 10 pm CEST\)](#)

If you cannot attend the live webinar, please register anyway, and we'll send you a link to the recording to watch at your convenience.





Innovation with Integrity