





GC-APCI-QTOF as an alternative to GC-MS/MS for the wide-scope screening of pesticides in orange

AACD 2018

S.E. Kokolakis¹, D.E. Damalas¹, C. Baessmann², N.S. Thomaidis¹

¹ Laboratory of Analytical Chemistry, Department of Chemistry, University of Athens, Panepistimiopolis Zografou, 15771 Athens, Greece
² Bruker Daltonik GmbH, Fahrenheitstr. 4, 28359 Bremen, Germany

e-mail: ntho@chem.uoa.gr



- Objectives**
- Development of a Database including analytes compatible with GC-APCI-QToF systems.
 - Development of a methodology for the wide scope screening of multiclass pesticides in orange
 - Analysis of real samples and investigate where the concentration of pesticides is greater
 - Highlight the potential of GC-APCI-HRMS as a promising alternative to GC-MS/MS for the determination of pesticides

Database development

Workflow example

1. Investigation of Ionization: $[M]^+$? $[M+H]^+$? $[M-Cl]^+$? etc

2. EIC $[M+H]^+$

3. MS

4. MS/MS

5. Confirmation of Qualifiers

m/z	RT	formula	name	CAS	Qual1	Qual2	Qual3
330.0804	23.13	C17H14ClFN3O ⁺ 1+	Epoxiconazole	135319-73-2	121.044	192.0282	165.0168
235.0076	21.16	C13H9Cl2 ⁺ 1+	2,4-DDT	789-02-6	318.9424	165.0664	199.0278
180.9373	14.23	C6H4Cl3 ⁺ 1+	a-HCH	319-84-6	172.9632	218.9111	180.9373
366.8541	18.18	C8H9BrCl2O3PS ⁺ 1+	Bromophos Methyl	2104-96-3	142.9902	352.8402	240.8597
271.095	14.12	C10H24O2PS2 ⁺ 1+	Cadusafos	95465-99-9	158.9658	130.938	130.9382
351.9306	17.60	C9H12Cl3NO3PS ⁺ 1+	Chlorpyrifos Ethyl	2921-88-2	199.9207	114.9609	209.9416
220.9532	9.76	C4H8Cl2O4P ⁺ 1+	Dichlorvos	62-73-7	184.9765	144.9818	127.0145
259.8096	9.53	C4Cl6 ⁺ 1+	Hexachlorobutadiene	87-68-3	224.8408	189.8719	140.906
260.8727	13.22	C6HCl4NO2 ⁺ 1+	Tecnazene	117-18-0	202.8759	242.8805	195.9018

>300 Analytes

- 250 Pesticides
- 8 FAMES
- 26 PAHs
- 6 PBDEs
- 7 PCBs
- 8 PCNs

Chlorine Isotopic pattern


3. MS

4. MS/MS

5. Confirmation of Qualifiers

In source Fragments

Compass DataAnalysis



Method

Validation

Linearity - 4 calibration levels:

100 µg/Kg

50 µg/Kg

10 µg/Kg

1 µg/Kg

R² > 0.99 for all the compounds

73 compounds

70 Pesticides

3 IS

2 standard addition levels

2 ppb

10 ppb

Analysis

GC-APCI-QToF

maXis impact

456-GC

APCI II source

Method performance

250/250 pesticides were detected at 100 µg/Kg matrix matched sample

70/70 pesticides were detected at 10 µg/Kg spiked sample

Meet the needs of routine pesticides' screening

Results

Propham 4 µg/Kg

Tetramethrin 4 µg/Kg

Cyhalothrin-lambda 2 µg/Kg

Spirodiclofen 7 µg/Kg

Tebufenpyrad 8 µg/Kg

Phosmet 8 µg/Kg

Chlorpyrifos ethyl 14 µg/Kg

Detected pesticides

Identification example TASQ 1.4

ART: < 0.1 min

Δm/z: < 5 mDa

Isotopic Profile < 100 mSigma

>2 Qualifier ions

Compartment-dependent C_p pesticides

Higher pesticides' concentration in the peel than in other parts of the orange

Peel

Whole orange

Edible part

Juice

28 µg/Kg

14 µg/Kg

< 0.5 µg/Kg

1 µg/Kg

Conclusions

A comprehensive database with Rt, MS and MS/MS information for > 300 compounds was established

GC-APCI-QToF proved a suitable platform for routine pesticide residues' analysis

Significant differences in the concentration of pesticides were observed in the different parts of orange