

Comprehensive 4D workflow for targeted and untargeted screening of PFAS in organisms from different trophic levels utilizing LC-VIP HESI(-)-TIMS-QToF MS



Georgios Gkotsis¹, Dimitrios E. Damalas¹, Maria-Christina Nika¹, Carsten Baessmann², Bob Galvin², Nikolaos S. Thomaidis¹

¹National and Kapodistrian University of Athens, Athens, Greece
²Bruker Daltonics GmbH & Co. KG, Bremen, Germany

Conclusions

- Establishment of an efficient and sensitive combined targeted and untargeted workflow for unveiling of PFAS in complex environmental matrices, such as biota.
- Ion mobility derived CCS values can be used as an additional identification criterion in target analysis, resulting to 4D-identification.
- The ion mobility filtration assists the data processing, leading to increased reliability of results and reduced time for data review.
- Increased sensitivity was observed using VIP-HESI source and mobility filtration.
- An efficient untargeted workflow, combining a suspect list with >4,700 PFAS and novel in-silico tools, was established.

Introduction

Per- and polyfluoroalkyl substances (PFAS)

- easily transportable
- persistent
- non-biodegradable
- bioaccumulated & biomagnified
- "Forever Chemicals"
- adverse effect on human health

> 4,700 registered PFAS plus degradation products

commonly monitored PFAS

PFCAs (n=4-13), PFSAs (n=4-13), ADONA, GenX

*restricted within EU (water DIRECTIVE (EU) 2020/2184)

emerging PFAS

Methods

generic sample preparation protocol

Androulakis et al 2022

Simultaneous extraction of PFAS from different sub-classes (different physicochemical properties)

RPLC - VIP HESI(-) - TIMS - QTOF MS

Target Screener 4D

acquisition modes:

- bbCID
- PASEF

Target analysis

Suspect screening

Non-target screening

Target Analysis

m/z	formula	name	rt	Qual1	Qual2	Qual3	Qual4	ccs	oneOverK0
298.9430	C4HF9O3S	Perfluorobutanesulfonic acid (L-PFBS)	5.68	98.9558	79.9574			133.83	0.6390
448.9334	C7HF15O3S	Perfluorheptanesulfonic acid (L-PFHpS)	8.64	168.9894	79.9574	98.9558		160.00	0.7752
488.9302	C8HF17O3S	Perfluoroctanesulfonic acid (L-PFOS)	9.47	79.9574	98.9558			168.84	0.8204
368.9786	C7F15	Perfluorooctanoic acid (L-PFOA) Fragm 369	8.63	168.9894	368.9786	218.9862	118.9926	157.25	0.7570
328.9743	C6H4F9	4,2 Fluorotetramer sulfonic acid (4,2 FTS)	6.52	306.9881	80.9652			150.81	0.7228
428.9679	C8H4F9	6,2 Fluorotetramer sulfonic acid (6,2 FTS)	8.59	406.9617	80.9652	388.9554		168.81	0.8166
528.9615	C10H4		10.23	506.9553	80.9652			186.00	0.9051

Screening results





Formula	MRSQC	Δm/z [mDa]	ΔRT [min]	mSigma	Ion Score	ΔCCS [%]
C ₄ H ₉ F ₉ O ₃ S	✓	0.27	-0.02	25.6	++	-0.36
C ₄ H ₇ F ₇ O ₂	✓	0.83	-0.03	22.0	++	-0.01
C ₁₂ H ₁₉ F ₁₉ O ₂	✓	0.50	-0.02	3.8	++	-0.02
C ₁₂ H ₁₇ F ₁₇ O ₂	✓	-0.13	-0.01	17.9	++	0.31
C ₁₂ H ₁₅ F ₁₅ O ₂	✓	-0.10	-0.03	14.7	++	0.08
C ₁₂ H ₁₃ F ₁₃ O ₂	✓	0.58	-0.14	1.1	++	-0.12
C ₁₀ H ₁₁ F ₁₁ O ₂	✓	0.45	-0.02	1.8	++	-0.01
C ₁₀ H ₉ F ₉ O ₂	✓	0.33	-0.20	4.3	++	1.61

ASQ

12 PFAS sub-classes
 ✓ 60 PFAS
 ✓ 21 PFAS I.S. representative isotopically labeled I.S.
 ΔCCS typically ≤ 1%

Mass accuracy
 Retention time
 Sigma value
 Qualifier ions
 CCS value agreement
 additional identification criterion

Results



Suspect Screening

MS/MS spectra annotation

Ion Formula	m/z meas.	Δ m/z [mDa]	int.	IntCov. [%]
[C2F5]	118.999	-7.3	18393.838	50.2
[C3F7]	168.989	-0.5	15737.872	42.9
[C6F13]	318.978	0.8	2545.803	6.9

PFAS chemical space

Kendrick Mass Defect (KMD) plot repeating unit: CF₂ ~1,850 features

Blank subtraction ~430 features

procedural blank otter sample

Structure

in-silico tools

predicted MS/MS spectrum

predicted CCS value

accurate mass

isotopic profile

9.9%

5.3%

497.9451

PRORISKPFAS suspect list containing 4,777 PFAS

Norman SusID	Name	Formula	CAS_RN	SMILES	SMILES	Monoisotopic mass	StdInChI	StdInChIKey
4632	NORMAN PRORISKPFAS	Cyclohexanaminium p...	462.96311	9.80	164.5	973.8	11.0	-2.8
4633	NORMAN PRORISKPFAS	Cyclohexanaminium p...	412.96311	9.00	155.6	967.7	8.9	-2.7
4634	NORMAN PRORISKPFAS	Sodium 3,3,4,4,5,5,6,6,7,7...	426.96814	9.06	167.9	924.6	16.8	-1.9
4635	NORMAN PRORISKPFAS	Sodium perfluorodeca...	512.95991	10.51	173.6	957.7	15.5	-2.1
4636	NORMAN PRORISKPFAS	Perfluorundecanoic ac...	562.95648	11.11	182.4	942.1	10.3	-1.5
4637	NORMAN PRORISKPFAS	1H,1H-Perfluoropropyl...	202.01288	1.67	133.6	920.0	81.0	6.9
4638	NORMAN PRORISKPFAS	Perfluorooctanesulfon...	488.93205	9.78	167.8	914.2	2.6	-2.5
4639	NORMAN PRORISKPFAS	Perfluorooctanesulfon...	497.94568	11.05	166.3	906.6	869.1	15.9
4640	NORMAN PRORISKPFAS	Perfluorooctanesulfon...	498.92986	10.06	167.8	906.6	617.1	5.0
4641	NORMAN PRORISKPFAS	Silane, trimethyl(tridec...	391.01917	1.79	170.6	877.0	27.0	24.8
4642	NORMAN PRORISKPFAS	Silane, trimethyl(tridec...	391.01912	1.20	170.8	877.0	28.4	-0.6
4643	NORMAN PRORISKPFAS	(4-Cyclohexylphenyl)di...	344.16081	5.11	183.3	842.0	78.0	47.1
4644	NORMAN PRORISKPFAS	Cyclopentane, 1,3,4,4...	192.28940	5.83	117.5	800.0	100.0	7.2
4645	NORMAN PRORISKPFAS	Perfluorooctanesulfon...	594.96278	10.39	171.2	800.0	100.0	10.0
4646	NORMAN PRORISKPFAS	Perfluorooctanesulfon...	268.98312	7.07	123.8	788.0	100.0	7.8
4647	NORMAN PRORISKPFAS	Perfluorooctanesulfon...	362.96965	8.11	146.8	788.0	953.8	4.9