

# Spatial imaging parallel reaction monitoring with parallel accumulation serial fragmentation (iprm-PASEF) investigating metabolically resilient CD38<sup>-/-</sup> mice

Charles A. Schurman<sup>1\*</sup>, Joanna Bons<sup>1\*</sup>, Prasanna Vadhana Ashok Kumar<sup>1</sup>, Jingji Fang<sup>1</sup>, Andrea Roberts<sup>1</sup>, Genesis Vega Hormazabal<sup>1</sup>, Rebecca R. Riley<sup>1</sup>, Nannan Tao<sup>2</sup>, Eric Verdin<sup>1</sup>, Birgit Schilling<sup>1</sup>

<sup>1</sup> Buck Institute for Research on Aging, Novato, CA, USA,

<sup>2</sup> Bruker, San Jose, CA, USA,




## Overview


Mass spectrometry imaging (MSI) is a rapidly advancing technology that provides mapping of the spatial molecular landscape of tissues for a variety of analytes. One common method employed is matrix-assisted laser desorption/ionization (MALDI)-MSI, however, confident *in situ* identification and accurate quantification of analytes with MALDI-MSI remains challenging without extensive workflows. Here, we present the first study combining trapped ion mobility spectrometry (TIMS)-based parallel accumulation-serial fragmentation (PASEF) with MALDI ionization for targeted **imaging parallel reaction monitoring (iprm-PASEF)**<sup>1,2</sup> to investigate lipids and metabolites in liver tissues from wild-type and CD38 knockout mice (CD38<sup>-/-</sup>). The novel iprm-PASEF workflow enabled confident identification and differentiation of lipid isomers at the MS2 level and revealed greatly increased NAD<sup>+</sup> levels and decreased *adenosine diphosphate ribose (ADPR)*, a by-product of NAD<sup>+</sup> hydrolysis, in CD38<sup>-/-</sup> livers. This approach provides highly confident, specific, and robust MS2-based identification and quantification of fragment ions in spatial MSI experiments. Additionally, the innovative iprm-PASEF opens unprecedented opportunities for spatial metabolomics and lipidomics, offering spatially resolved insights into molecular mechanisms.

## NAD<sup>+</sup>, CD38, and Aging

NAD<sup>+</sup> (nicotinamide adenine dinucleotide) metabolism plays a critical role in maintaining cellular health and is a central metabolite for physiological and cellular processes. Aging and age-related diseases and cellular conditions, such as cellular senescence, are associated with declines in NAD<sup>+</sup>. CD38, is a potent enzymatic consumer of NAD<sup>+</sup> and is highly expressed in Kupffer cells, the resident macrophages of the liver sinusoids, where it drives cellular senescence and inflammation with aging. Notably, mice lacking CD38 (CD38<sup>-/-</sup>) display higher mitochondrial NAD<sup>+</sup> levels in the liver exhibiting a resistance to age-related declines in NAD<sup>+</sup>, suggesting improved metabolic regulation.<sup>3</sup> Here, we utilized iprm-PASEF to investigate the spatial distribution of metabolites and lipids in the liver of CD38<sup>-/-</sup> mice and wild type (WT) controls.




**Figure 1:** CD38 is a potent enzymatic consumer of NAD, catalyzing NAD<sup>+</sup> hydrolysis into its enzymatic by-products ADPR, cADPR and nicotinamide. Mice lacking CD38 (CD38<sup>-/-</sup>) showed improved metabolism and lower cellular senescence in the liver, possibly through the beneficial metabolic effects of NAD. CD38<sup>-/-</sup> liver was expected to exhibit elevated NAD<sup>+</sup> levels with a reciprocal repression of its related by-products.




**Figure 2: iprm-PASEF Workflow.** A) Using the timsTOF fleX mass spectrometer MALDI-trapped ion mobility spectrometry (TIMS)-MS survey acquisitions were collected for MS1 data. Survey imaging scans were processed in Bruker SCILS Lab using the T-ReX3 feature finding algorithm and the MetaboScape software. From identified features, a list of candidate precursor ions, including lipids and metabolites, was generated. B) Targeted iprm-PASEF assays were developed and optimized in SCILS Lab and timsControl software, extracting information on 1/K<sub>0</sub>, m/z and region of interest from the MALDI-TIMS-MS1 scans. Precursor ions were isolated by ion mobility, followed by m/z isolation and *in situ* fragmented in a sequential manner to collect MS2 data in a spatial manner.

## iprm-PASEF enables spatially resolved fragment ion distribution




**Figure 3: iprm-PASEF enables spatial fragment mapping** A) In CD38<sup>-/-</sup> mouse liver, the composite CCS-enabled survey MALDI-MS scan was acquired in negative ionization mode and detected a feature at m/z 885.55 and 1/K<sub>0</sub> 1.47 VS/cm<sup>2</sup> that was highly abundant across the entire tissue section. B) This precursor ion candidate was selected for targeted iprm-PASEF isolation and fragmentation. The generated MS/MS spectrum showed species-specific fragment ions used for both accurate molecular identification and quantification of Phosphatidylinositol PI(18:0/20:4). C) The identified fragment ions were spatially mapped for representative images of wild-type (WT) and CD38<sup>-/-</sup> mouse liver tissue sections based on their m/z. The spatial distribution of the PI(18:0/20:4) fragment ions matched the distribution of the corresponding precursor ion.

## iprm-PASEF resolves lipid isomers




**Figure 4: Detection of characteristic fragment ions for the identification of lipid isomers** A) A zoom into the m/z range 857-861 of the composite survey MALDI-MS scan collected on CD38<sup>-/-</sup> mouse liver tissue showed a profile resembling an isotopic pattern (left). Potential lipid isomer discovery and annotation for the most intense ion at m/z 857.52 was initially given as PI(14:0/22:4) by MetaboScape software. B) The iprm-PASEF fragmentation pattern of the precursor ion [M-H]<sup>-</sup> at m/z 857.52 in the mobility isolation window 1.4308 ± 0.0093 Vs/cm<sup>2</sup> showing the detection of two characteristic fragment ions, m/z 255.23 (pink) and m/z 303.23 (orange), as evidenced for the identification of PI 36:4/PI 16:0/20:4. Lack of characteristic fragments for PI(14:0/22:4) (m/z 227.20 red, & m/z 331.26 blue) confirmed the absence of this species.

## iprm-PASEF identifies NAD and related metabolites



**Figure 5: Identification of NAD<sup>+</sup> and related metabolites in WT and CD38<sup>-/-</sup> mouse livers by iprm-PASEF.** A) The iprm-PASEF mobilogram and MS/MS spectra of synthetic standards: negative ionization mode with N-(1-naphthyl) ethylenediamine dihydrochloride (NEDC) matrix for cADPR, positive ionization mode with 2,5-dihydroxyacetophenone (DHAP) matrix for ADPR and NAD<sup>+</sup>. B) The iprm-PASEF mobilogram and MS/MS spectrum of the synthetic ADPR standard analyzed in negative ionization mode with NEDC matrix. C) Spatial distribution of ADPR precursor ion [M-H]<sup>-</sup> at m/z 558.06 in MS1 survey scans. D) iprm-PASEF analysis in WT and CD38<sup>-/-</sup> mouse liver tissue found a significant decrease of ADPR in CD38<sup>-/-</sup> mouse liver with MS2-based quantification on the fragment at m/z 348.06. \*\* Student's t-test p-value < 0.01

## iprm-PASEF facilitates *in situ* visualization and quantification of endogenous NAD<sup>+</sup>



**Figure 6. Quantification and *in situ* visualization of endogenous NAD<sup>+</sup> fragment ions in WT and CD38<sup>-/-</sup> mouse livers by iprm-PASEF.** A. The iprm-PASEF fragmentation pattern of the endogenous NAD<sup>+</sup> precursor ion [M+H]<sup>+</sup> at m/z 664.12 in the mobility isolation window 1.1224-1.1400 Vs/cm<sup>2</sup> showed the detection of fragment ions at m/z 136.06, m/z 232.08, m/z 348.07, and m/z 428.04. B. Spatial distribution of the four NAD<sup>+</sup> fragment ions detected by iprm-PASEF in one WT and one CD38<sup>-/-</sup> mouse liver replicate. Each fragment ion was quantified *in situ* using the pixel intensity within each tissue section, which all demonstrated up-regulated level in CD38<sup>-/-</sup> vs WT mouse livers. The statistical analysis was performed by comparing the MS2 intensity of each fragment ion in each individual pixel that compose the CD38<sup>-/-</sup> vs WT image. \*\*\*\* Student's t-test p-value < 0.0001 (n = 1885 pixels for the WT image and n = 2300 pixels for the WT image; only one mouse and one image replicate per condition). Cleavage sites of identified NAD<sup>+</sup> fragments are mapped in the NAD<sup>+</sup> structure.

## Conclusions

- iprm-PASEF is a powerful application of MALDI-MSI that allows for MS2 level identification and quantification of metabolites and lipids in place during imaging experiments.
- Iprm-PASEF can differentiate between multiple related lipid species.
- Iprm-PASEF allows for further confidence identification of preliminarily identified analytes through fragment level analysis.
- CD38<sup>-/-</sup> mice were confirmed to have both elevated NAD<sup>+</sup> levels and consequently repressed ADPR levels in liver tissue, demonstrating their metabolic resilience.

## Acknowledgments

**Supporting Grants:** NIAMS R21 AR084303 (PI Schilling), NIA P01 AG066591 (PI Ellerby), NIA T32 AG000266 (PI Ellerby)

## References

- [1] - Li, M. J. *et al. Rapid Commun Mass Spectrom* **39**, e10006 (2025).
- [2] - Sanchez, L., Shepherd, R. A. & Luu, G. T. *bioRxiv*, (2025).
- [3] - Covarrubias, A. J. *et al. Nat Metab* **2**, 1265-1283 (2020).