


# Microanalysis with high spectral resolution: the power of QUANTAX WDS for SEM



Bruker Nano Analytics, Berlin, Germany  
Webinar, September 14, 2017



# Presenters



Dr. Jörg Silbermann

Product Manager WDS,  
Bruker Nano Analytics, Berlin, Germany



Dr. Michael Abratis

Sr. Application Scientist WDS,  
Bruker Nano Analytics, Berlin, Germany

# Microanalysis with high spectral resolution: the power of QUANTAX WDS for SEM



## OUTLINE:

- QUANTAX WDS – an overview
- The XSense™ spectrometer: working principle and spectral resolution
- EDS vs. WDS: the need for high spectral resolution
- Sample measurement data and application examples focusing on high spectral resolution

# QUANTAX WDS

## System Components



QUANTAX WDS: integral part  
of the QUANTAX family




XSense™ WD spectrometer



Signal processing  
unit SVE 6





ESPRIT 2



# XSense WD Spectrometer Setup and Working Principle



- Parallel Beam Optic (PBO) transforms X-rays diverging from the sample into a parallel beam
- Polychromatic beam undergoes Bragg diffraction at flat analyzer crystal
- Angle  $\Theta$  between beam and crystal surface and crystal lattice constant  $2d$  determines the energy that passes through to the detector
- X-ray detection by a flow proportional counter



Bragg equation:  $n\lambda = 2d \sin(\theta)$


# XSense WD Spectrometer

## Spectral resolution




- Spectral resolution usually defined via the *full width at half maximum* (FWHM) of an elemental line peak
- Natural line widths are in the 0.2-3 eV range
- Peak broadening due to spectrometer effects:
  - In a PB-WDS:
    - Imperfections of analyzing crystals/multi-layers (crystal defects, inter-layer diffusion, waviness of bilayer boundaries)
    - Imperfect parallelization of beam due to **(A)** aberrations of the optic, **(B)** optical misalignment.

Natural line width:



Peak broadening due to spectrometer effects:



# XSense WD Spectrometer

## Spectral resolution



- Spectral resolution usually defined via the *full width at half maximum* (FWHM) of an elemental line peak
- Natural line widths are in the 0.2-3 eV range
- Peak broadening due to spectrometer effects:
  - In a PB-WDS:
    - Imperfections of analyzing crystals/multi-layers (crystal defects, inter-layer diffusion, waviness of bilayer boundaries)
    - Imperfect parallelization of beam due to **(A)** aberrations of the optic, **(B)** optical misalignment.


### XSense: Minimization of **A**:

- Use of a grazing incidence mirror optic which (in comparison with a polycapillary-based optic) produces a highly parallel beam of low divergence

Grazing incidence mirror optic:



low beam divergence



# XSense WD Spectrometer

## Spectral resolution



- Spectral resolution usually defined via the *full width at half maximum* (FWHM) of an elemental line peak
- Natural line widths are in the 0.2-3 eV range
- Peak broadening due to spectrometer effects:
  - In a PB-WDS:
    - Imperfections of analyzing crystals/multi-layers (crystal defects, inter-layer diffusion, waviness of bilayer boundaries)
    - Imperfect parallelization of beam due to **(A)** optic figure/slope errors, **(B)** optical misalignment.



### XSense: Minimization of **B**:

- Spectrometer equipped with 3-axis optic positioning unit + powerful optical alignment software algorithm

grazing incidence optic + auto-optic alignment:

**highest possible resolution from a PB-spectrometer**

3-axis optic positioning unit




# QUANTAX WDS

## EDS vs. WDS




- In SDD-based EDS peak broadening is much more pronounced, resolutions are in the 40-120 eV range
- Limited resolution of EDS leads to frequent peak overlaps, mainly in the low energy range

Typical EDS peak overlap:



- XSense features substantially higher spectral resolution (typically 4-20 eV FWHM)

XFlash™  
EDS detector



XSense™  
WD spectrometer



**WDS ideally complements EDS in demanding applications, where resolution is critical**

# QUANTAX WDS

## Resolving common overlaps in EDS microanalysis




Element and line	Interferences with	$\Delta$ eV	Samples or applications where the overlaps are found
Cu-L	Na-K $\alpha$	18	Biological samples (grid)
As-L	Na-K $\alpha$	79	Biological samples (stain or fixative)
Ag-L	Cl-K $\alpha$	10	Biological samples (stain or fixative)
Ru-L	S-K $\alpha$	54	Biological samples (stain or fixative)
Os-M	Al-K $\alpha$	5	Biological samples (stain or fixative)
U-M	K-K $\alpha$	22	Biological samples (stain or fixative)
Sr-L $\alpha$	Si-K $\alpha$	31	Silicates (feldspars in particular)
Y-L $\beta$	P-K $\alpha$	18	Phosphates
Y-L $\beta$	Zr-L $\alpha$	46	Silicates (zircon), oxides (zirconia)
S-K $\alpha,\beta$	Mo-L $\alpha$ ; Pb-M $\alpha$	14; 38	Minerals, lubricants, sulfides, sulfates
Ti-K $\beta$	V-K $\alpha$	20	Steels, Fe-Ti oxides
V-K $\beta$	Cr-K $\alpha$	13	Steels
Cr-K $\beta$	Mn-K $\alpha$	47	Steels
Mn-K $\beta$	Fe-K $\alpha$	87	Steels
Fe-K $\beta$	Co-K $\alpha$	128	Steels, magnetic alloys
Co-K $\beta$	Ni-K $\alpha$	169	Steels, hard surfacing alloys
W-M $\alpha,\beta$	Si-K $\alpha,\beta$	35	Semiconductor processing
Ta-M $\alpha,\beta$	Si-K $\alpha,\beta$	27	Semiconductor processing
Ti-K $\alpha$	Ba-L $\alpha$	45	Optoelectronics, silicates
As-K $\alpha$	Pb-L $\alpha$	8	Pigments

Overlaps known from biological, geological and material sciences and industries

Modified after Goldstein et al. (2007). Scanning Electron Microscopy and X-Ray Microanalysis. Springer

# QUANTAX WDS and EDS


## Geological samples I: Pb sulfide



$\Delta$  S-K $\alpha$  – Pb-M $\alpha$ : 38 eV

# QUANTAX WDS and EDS


## Geological samples II: Mo sulfide



$\Delta S-K\alpha - Mo-La: 14 \text{ eV}$

# QUANTAX WDS and EDS


## Geological samples III: REE phosphates



$\Delta P-K\alpha - Y-L\beta$ : 18 eV;  $\Delta P-K\alpha - Y-L\alpha$ : 92 eV

# QUANTAX WDS and EDS


## Geological & material sciences: Zr oxide



$\Delta$  Zr-L $\alpha$  – Y-L $\alpha$ : 120 eV;  $\Delta$  Zr-L $\alpha$  – Y-L $\beta$ : 46 eV

# QUANTAX WDS and EDS


## Application in material science II: glass



$\Delta \text{Na-K}\alpha - \text{Zn-L}\alpha: 28 \text{ eV}; \Delta \text{Na-K}\alpha - \text{Zn-L}\beta: 5 \text{ eV}$

# QUANTAX WDS and EDS

## Material science III: semiconductors



$\Delta \text{Si-K}\alpha - \text{Ta-M}\alpha$ : 28 eV;  $\Delta \text{Si-K}\alpha - \text{Ta-M}\beta$ : 27 eV

# QUANTAX WDS and EDS


## Comparison of WDS and EDS resolution



Element	Atomic No.	X-ray line	Energy [keV]	FWHM EDS [eV]	WDS diffractor	FWHM WDS [eV]	Resolution improvement
Si	14	K $\alpha$	1.740	75	PET	3.5	21x
P	15	K $\alpha$	2.014	77	PET	5	15x
S	16	K $\alpha$	2.307	85	PET	7	12x
Y	39	L $\alpha$	1.922	82	PET	6.3	13x
Zr	40	L $\alpha$	2.042	83	PET	7.2	12x
Mo	42	L $\alpha$	2.293	87	PET	9.5	9x
Ta	73	M $\alpha$	1.712	71	PET	6	12x
W	74	M $\alpha$	1.775	74	PET	6.4	12x
Hg	80	M $\alpha$	2.195	80	PET	9	9x
Pb	82	M $\alpha$	2.345	91	PET	11.9	8x


# QUANTAX WDS and EDS

## Resolution vs. deconvolution I



Deconvolved EDS spectrum of tungsten silicide (WSi<sub>2</sub>)

# QUANTAX WDS and EDS Resolution vs. deconvolution I



Highly resolving WDS spectrum of tungsten silicide ( $WSi_2$ )

# QUANTAX WDS and EDS

## Resolution vs. deconvolution II




Quantitative results in atomic percentages (5kV)

Sample	Element	Stoichiometry	EDS <sup>1</sup>	EDS <sup>2</sup>	WDS
<b>MoS<sub>2</sub></b>	Mo	33.3	39.0	34.5	33.9
	S	66.7	61.0	65.5	66.1
<b>WSi<sub>2</sub></b>	W	33.3	20.2	32.6	33.4
	Si	66.7	79.8	67.4	66.6


<sup>1</sup>standardless, <sup>2</sup>standard-based

# QUANTAX WDS and EDS

## High spectral resolution for mapping



Combined WDS and EDS mapping



Combined WDS and EDS spectra

Lead accumulator showing surface deposits of  $PbSO_4$


# Summary of today's WDS Webinar



- Bruker QUANTAX PB WDS on SEM
- XSense WDS facilitates high spectral resolution analyses
- Applications include biological, geological and material sciences and industries
- Deconvolution methods cannot replace true spectral resolution
- High spectral resolution is important for qualitative and quantitative analyses as well as mapping
- QUANTAX WDS is a powerful tool for scientific and industrial applications

## Are There Any Questions?

Please type in the questions you might have  
in the Q&A box and press *Send*.



Innovation with Integrity