

# MALDI mass spectrometry imaging guided LC-MS on the same instrument

Frédéric Dewez<sup>1,2</sup>, Janina Oejten<sup>3</sup>, Corinna Henkel<sup>3</sup>, Romano Hebel<sup>3</sup>, Michael Herfs<sup>4</sup>, Ron M.A Heeren<sup>1</sup>, Edwin De Pauw<sup>2</sup> and Benjamin Balluff<sup>1</sup>

<sup>1</sup> Maastricht Multimodal Molecular Imaging Institute (M4i), University of Maastricht, The Netherlands

<sup>2</sup> Mass Spectrometry Laboratory (MSLab), University of Liège, Belgium


<sup>3</sup> Bruker Daltonik GmbH, Bremen, Germany

<sup>4</sup> Laboratory of Experimental Pathology, GIGA-CancerUniversity of Liège, Belgium


## Introduction

Mass spectrometry imaging (MSI) is an analytical technique for the unlabeled and multiplex analysis of molecular spatial distributions in biological tissues. Unfortunately, due to lack of purification steps, the analytical depth of MSI for a comprehensive *in situ* molecular characterization is still limited [1]. Liquid chromatography mass spectrometry (LC-MS), in contrast, gives the possibility to obtain proteomic or metabolomic data for thousands of molecules even from smallest amounts of samples [2]. In the presented project, the goal is to combine MSI and LC-MS by designing a pipeline where MSI can guide a laser microdissection system (LMD) to accurately isolate regions of interest (ROIs) and perform microproteomics on the same instrument.

## I. Lipid MSI of breast cancer tissue section




## II. Segmentation of tumor areas




MSI experiment of lipids was performed on a fresh frozen breast tumor section mounted onto a PEN (PolyEthylene Naphthalate) membrane slide compatible with LMD. The data were acquired in positive ionization mode at 50 μm pixel size on a tims TOF flex (Bruker Daltonik GmbH, Bremen, Germany). The MSI data were then imported into SCiLS Lab MVS 2019c (Bruker Daltonik) for data analysis. K-means segmentation (k=3) was performed on the annotated tumor areas by a pathologist.

## III. Image processing of segmentation




In order to increase the viability of microdissection, the segmentation data of the tumor areas was first exported from SCiLS Lab (.csv) and then imported into Matlab R2018a for image processing. First, a smoothing was performed by *opening* the image with a 2x2 square as structuring element (*imopen*). Then, small objects composed of less than 30 pixels based on a 4-connected neighborhood were removed (*bwareaopen*) and holes were filled in the 8-connected neighborhood (*imfill*). Finally, the boundaries (as polygonal areas) of the MSI segments were determined and up-scaled to the resolution of the optical image for a later coregistration to the LMD system.

## IV. Coregistration of MSI data to LMD



The coordinates of the MSI segments were then transferred to the LMD system (LMD 7000, Leica Microsystems GmbH, Wetzlar, Germany) by recalculating the boundary coordinates using one common reference point that is visible in the LMD as well as in the digital optical image [3].

## V. Microproteomics analysis of microdissected MSI segments



## Conclusions

Here we present a pipeline where MSI can spatially guide the LMD for subsequent microproteomics molecular characterization of regions of interest. The tims TOF flex is an ESI instrument associated with a MALDI source. Combining with this workflow, this will bring together the MSI spatial dimension to the molecular information using X-Omics analyses on the very same instrument for a more comprehensive molecular characterization of *in situ* biological processes.

## References

- McDonnell LA, Heeren RMA. Imaging mass spectrometry. *Mass Spectrom Rev.*
- Alberts D et al. MALDI imaging-guided microproteomic analyses of heterogeneous breast tumors—a pilot study. *Proteomics Clin Appl.*
- Dewez F et al. Precise co-registration of mass spectrometry imaging, histology, and laser microdissection-based omics. *Anal Bioanal Chem.*

## Acknowledgements

This research was financially supported in part through the LINK program of the Province of Limburg (The Netherlands). FD received support from the Joint Imaging Valley program of the University of Liège and the Maastricht University.