

Webinar

Investigation of Concrete by Means of micro-XRF



Dr. Andrew H. Menzies and Dr. Max Buegler
 Bruker Nano Analytics, Berlin, Germany
 Webinar, June 3rd of 2020



Introduction

Presenters / Moderators



Andrew Menzies, PhD

Sr. Application Scientist Geology and Mining,
Bruker Nano Analytics, Berlin, Germany



Max J.L. Bügler, PhD

Applications Specialist,
Bruker Nano Analytics, Berlin, Germany

Construction Materials - Concrete

Outline



- Brief introduction
 - M4 TORNADO
 - Micro-XRF
- Construction Materials - Concrete
 - Alkali-Silicate- (ASR) and Alkali-Carbonate-Reaction (ACR)
 - AMICS (mineralogy)
 - SEM-EDX
- Mapping of relevant elements
- Analytical approaches
- Comparison of anodes materials
- Practical demonstration
- Summary

M4 TORNADO

Micro-XRF Spectrometer



30 W micro-focus Rh tube with polycapillary lens
for excitation spot sizes $< 20 \mu\text{m}$ (for Mo-K α)

Optional 40 W micro-focus W tube with collimator
for excitation of 'heavy' elements, embedded in
lighter matrices

Up to two Silicon drift detectors (SDD)
with 30 or 60 mm² active area each energy
resolution $< 145 \text{ eV}$
(for Mn-K α @ 130 kcps throughput)




Sealed sample chamber with adjustable
pressure between 1 mbar and atmospheric
pressure for detecting elements down to Na

Sample stage with measurable area of 200 mm x 160 mm, maximum
sample height 120 mm, maximum sample weight 5 kg, and sample
stage speed up to 100 mm/s, minimum step size 4 μm



M4 TORNADO Instrument Design

Focusing by Poly-Capillary Lens




Poly-capillary lens collects large solid angle of tube radiation

Concentrates beam at sample into small spot

Small spot size guarantees **high brilliance** for data collection

M4 TORNADO Instrument Design


Design and Safety



HyperMap Data Accumulation Position-Tagged Spectrometry (PTS)



Data mining both on- and offline



What is a HyperMap?

Complete X-ray spectrum for every measured pixel saved.

Generates 4D data cube

2D – Spatial Distribution

3D – Spectral Distribution (E)

4D – Intensity of every spectral channel

Measurement time/pixel down to 1 ms

Introduction

XRF on Different Scales and Fields



X-TRACE
(on SEM)

19 mm



M4 Tornado

25 mm



M6 Jetstream

37.5 cm




Tracer

19 mm



Information from different depths of the sample

38.7 mm



Trace element sensitive

46.5 cm



No sample preparation

Concrete Applications



- Concrete has been used as building material since the Mesopotamian ages and became very popular with Greek and Roman builders
- It withstands immense pressure and has a high durability
- It can be poured into almost any shape
- If reinforced with steel cable or mesh as pre- or post-tensioned it can bridge wide distances



Concrete Weathering



- Concrete can last thousands of years under the right conditions
- Nowadays we frequently find images of cracked and weathered concrete
- As concrete is a mixture of many complex chemical compounds its durability depends on environmental conditions
- Saltwater for example leads to chemical reactions that change the volume and lead to stress and strain in the structure



Concrete

Alkali-Silica Reaction (ASR)



- Aggregates containing certain forms of silica will react with alkali hydroxide in concrete to form a gel that swells as it adsorbs water from the surrounding cement paste or the environment.
- These gels can induce enough expansive pressure to damage concrete.
- The final reaction can be schematically represented by:





- Which is a reaction that will terminate during maturation of the concrete as the acidity goes down

Concrete

Alkali-Silica Reaction (ASR)



- In the presence of alkaline hydroxides (NaOH / KOH) the pH level is significantly increased, and the reaction can continue. In the presence of such alkaline hydroxides a sodium silicate can form:



- The hygroscopic sodium silicate will expand as it absorbs water.
- This expansive pressure in the solid matrix may induce spalling of material, formation of macroscopic cracks, and a consequently a loss of structural integrity.
- However, this ASR can be controlled using certain supplementary cementitious materials.

Concrete

Alkali-Silica Reaction (ASR)



In order of decreasing potential to cause alkali-silica reaction.

Important Mineral Phases:

- Amorphous silica: volcanic glass, opal,
- Unstable silica polymorphs (tridymite and cristobalite),
- Chalcedony,
- Other forms of cryptocrystalline silica,
- Deformed quartz
- Strained quartz
- Quartz with crystallization imperfections,
- Quartz in crystalline forms



Grains of aggregate with different mineral composition and morphology:

- a) initial condition,
- b) with ASR - induced characteristic micro-cracks

Concrete

Alkali–Carbonate Reaction (ACR)




- Dedolomitization, the breaking down of dolomite, is normally associated with expansion.
- The reaction in this case can be schematically represented by:



- In another step the sodium-carbonate in combination with calcium hydroxide (portlandite) can react again into :



- Under the presence of calcium hydroxide, the dolomite reacts as follows:




- This reaction and subsequent crystallization of brucite (**Mg(OH)₂**) may cause considerable expansion.

Advanced Mineral Identification and Characterization System (AMICS)



- AMICS software package applies finger printing to spectra for pixels in order to identify the minerals in a mapping
- The system comes with a preloaded comprehensive set of typical spectra for more than 2000 minerals for quick starting a custom database
- For each material system an individual database is developed as samples are characterized
- Based on the sample the database is extended by standard spectra created from the identified regions of a specimen.
- Statistical tools to calculate mass fractions and areal coverage across the specimen






AMICS


How does AMICS work?



Micro-XRF Analysis: Hypermap
Mineralogy Workflow Schematic



Classification



- Acquires spectrum for each segment

- Identifies minerals based on spectral information


Concrete Hypermap Example

Starting with AMICS



With the preloaded database the concrete sample was analyzed and a according database of minerals with according spectra was defined.

The analysis shows the spatial distribution of each of the identified minerals along the sample surface.



Concrete Sample Mineralogy

Visualizing Results with AMICS




Resulting modal composition of the investigated area.

Relatively large unknown portion due to edge area around grains, that cannot be properly assigned to either one phase.

Higher resolution scan and longer exposure can reduce the amount of unknown material.

Modal Mineralogy



AMICS: Relevant Mineralogical Information



Modal Mineralogy: How much of each mineral is present.

Elemental Assay: How much of each element is present (Whole Rock).

Elemental Distribution: How is the element of interest (EOI) distributed in each mineral? E.g. Cr in Chromite vs Pyrope Garnet.

Particle (Grain) Size Distribution: Classify sizes of particles (grains).
Provide P-80, P-50 and P-20 values

Mineral Association: Identify how the minerals are associated with each other, e.g. Chromite and Garnet and Clinopyroxene etc.

Grain Shape Factor: The shape of the grain, i.e. euhedral, elongated.

Mineral Density Distribution: Classify densities of minerals. Identify how the minerals are distributed among the densities.



Academic Research EDX and X-TRACE



As the spot in micro-XRF is in the range of 15 μm for higher energies and about 30 μm for elements such as Si, to resolve the micrometer and sub-micrometer scale an alternative method is required. The common approach is an e-beam excitation of the fluorescence in a SEM.


To correlate this information with micro-XRF results, X-TRACE provides an X-ray excitation on SEM. This allows for direct correlation of micro-XRF results with the sub-micrometer features in heterogeneous samples.

This adds complimentary information at a superior resolution, its high sensitivity for light elements down to boron. Where micro XRF outperforms an EDX for elements with $Z > 20$ in detection limits, below that threshold the EDX has superior sensitivity and detection limits.




<https://www.understanding-cement.com/alkali-silica.html>

Analyzing Concrete Measurement Conditions




- Map of of 56.4 mm by 10.8 mm with 20 μm pixel size and 10 ms pixel time was measured on a standard M4 TORNADO utilizing 2x30 mm² detectors and a Rhodium anode
- 4 h scan on a region from left to right (thin green frame)
- Gradient in Cl distribution (orange) clearly visible

Map information	
Mapping parameter	
Width:	2820 pixel
	56,4 mm
Height:	540 pixel
	10,8 mm
Pixel Size:	20 μm
Total number of pixel:	1522800 pixel
Acquisition parameter	
Frame count:	1
Pixel time:	10 ms/pixel
Measure time:	4:13 h
Overall time:	4:48 h
Tube parameter	
High voltage:	50 kV
Anode current:	600 μA
Filter:	Empty
Optic:	Lens
Chamber at:	Vacuum 20,1 mbar
Anode:	Rh
Detector parameter	
Selected detectors:	1,2



Analyzing Concrete


Enhanced Visualization and Dimensions



- False color representation of the chlorine concentration to enhance dynamic to the human eye for small changes
- The chemical gradient of Cl is clearly visible by the color changing accordingly to the element's fluorescence intensity
- With the *drawing tool* based on the visualization, a **line** has been drawn (shown above) and is 11.5 mm long, which gives an indication of the depth of Cl diffusion into the concrete

Analyzing Concrete


Cl Penetration Depth from Line Profile




- A line object can be drawn into the Hypermap, creating a line profile along the sample.
- Additionally, the line was widened to average points perpendicular to the line. The Cl signal is higher in the first ~ 11.5 mm (from the left).
- The strong modulation is a result of the large agglomerates in the matrix.

Analyzing Concrete

... a more Quantitative Approach: Objects




Objects with a width of 2 mm and almost the full height are created next to each other. Their center position with respect to the sample surface are determined with the help of the "draw tools". The individual objects are saved as new individual maps, using the "save all objects" function.






Analyzing Concrete

Phase Analysis and Concentrations




- Load each object
- Select only a signature element of the cement (here P)
- Press auto phase
- Select the cement phase (phases and sum them up if necessary)
- Create sum spectrum of the according (sum-)phase
- Quantify the according spectrum with appropriate method
- Move the spectrum to the project (or save to disk)
- Import (drag and drop) all sum spectra into point mode
- Export results table to Excel → plot




Measuring Concrete

Taking a Map and Visualizing Intensity



MAP INFORMATION	
Mapping parameters	
Width:	940 pixel
	70,5 mm
Height:	799 pixel
	59,893 mm
Pixel Size:	75 μ m
Total number of pixel:	751060 pixel
Acquisition parameters	
Frame count:	1
Pixel time:	15 ms/pixel
Measure time:	2:35 h
Overall time:	3:51 h
Stage speed:	5,0 mm/s
Stage position (X,Y,Z):	96.968;89.294;78.717 mm
Tube parameter	
High voltage:	50 kV
Anode current:	600 μ A
Filter:	Empty
Optic:	Lens
Collimator diameter:	0
SpotSize:	20
Chamber at:	Air 2 mbar
Flow rate:	--- l/min
Anode:	Rh
Detector parameters	
Selected detectors:	1,2
Max. pulse throughput:	275000 cps



- Drill core of concrete exposed to NaCl was mapped
- The intensity in a Region-Of-Interest (ROI) for Cl-K-line (as proxy for Na) is displayed
- Intensity is found virtually everywhere

Measuring Concrete





Signal-to-Noise-Ratio (SNR)



- A look at the spectrum shows a strong overlap of the Cl lines with backscattered Rh-L-line intensity
- Fast algorithm can perform deconvolution of 18 elements in 350 000 spectra in 50 sec, in order to separate Cl from Rh intensity
- Rh-anode yields poor SNR due to overlap, suppressing the Rh-L-line by a filter removes the need for deconvolution at the cost of overall signal and particularly for the light elements

Region of Interest (ROI)

Deconvoluted signal




Measuring Concrete

Alternative use of Ag-Anode






- To avoid the need for deconvolution, or requirement for a filter, alternative anode material should be considered
- For the given analytical question silver offers a much more efficient and potent excitation
- The backscattered Ag-L-line intensity not only does not overlap with the Cl, but it provides a much more efficient excitation of Cl by the 2.82 keV absorption edge
- No filter or deconvolution required




ROI of Cl with Rh-anode

ROI of Cl with Ag-anode



Measuring Concrete

Is there an Ideal Source?



- A direct comparison shows that the Ag excitation comes at a cost as well, it overlaps strongly with K
- As the Ag-L-lines are transmitted more efficiently, it provides better signal for elements with $6 < Z < 18$, such as Cl, S, P, etc.
- It is possible to equip the M4 TORNADO with two sources, but only one can be equipped with a poly-capillary lens

From the above, it becomes clear, there is no ideal anode for all questions! In the decision for an excitation source the priority in terms of elements must be considered!

Concrete

Why XRF? From Construction Site to Lab



- XRF allows for fast measurement of a surface on site using a handheld XRF, such as the tracer
- On indication of contaminants a drill core can be taken for more in-depth analysis
- An area of up to 19 cm x 16 cm can be mapped with an M4 TORNADO, or an area of up to 80 cm x 60 cm utilizing an M6 JETSTREAM
- The penetration of catalysts, such as NaCl, into the concrete can be visualized and semi-quantitatively analyzed
- For research purposes it can be analyzed even further utilizing EDX in conjunction with X-TRACE, allowing to correlate micrometer and sub- μ m-scale element distributions

Live Demonstration

Phase Analysis and Concentrations



Concrete Summary




- Micro-XRF can provide compositional and mineralogical information across a range of scales
- Fundamental parameter (FP) based quantification allows for quick assessment of composition in points
- Mapping of mineralogy using AMICS yields information about possible ASR and ACR reactions as well as matrix and aggregate compositions
- Mapping of element intensities on a drill core sample can provide visualization of alkaline hydroxide penetration into the volume
- Smart analysis of mapping allows for semi-quantitative analysis of Cl and K concentration in the matrix
- Phase analysis allows for quick estimation of cement to aggregate ratio as well as for semi-quantitative analysis

Literature

M4 TORNADO in Concrete Research



Since in 2015 the first M4 was purchased for the purpose of research in concrete more than 130 publications can be found, and multiple Masters and PhD thesis utilizing the M4 in the field were published. Some example from various groups:



G. Bonifazi et al., *"Hyperspectral imaging applied to the identification and classification of asbestos fibers,"* 2015 IEEE SENSORS, Busan, 2015, pp. 1-4, doi: 10.1109/ICSENS.2015.7370458.

K. De Weerd et al., *"Comparing chloride ingress from seawater and NaCl solution in Portland cement mortar,"* Cement and Concrete Research, Volume 115, 2019, pp. 80-89, doi: 10.1016/j.cemconres.2018.09.014.

H. Morillas et al., *"Multianalytical approach to evaluate deterioration products on cement used as consolidant on lithic material: The case of Tello Obelisk, Lima (Peru),"* Microchemical Journal, Volume 139, 2018, pp. 42-49, doi: 10.1016/j.microc.2018.02.017.

G. Plusquellec et al., *"Determining the free alkali metal content in concrete – Case study of an ASR-affected dam,"* Cement and Concrete Research, Volume 105, 2018, Pages 111-125, doi: 10.1016/j.cemconres.2018.01.003.

A. Dehghan et al., *"Application of X-ray microfluorescence for the determination of chloride diffusion coefficients in concrete chloride penetration experiments,"* Construction and Building Materials, Volume 148, 2017, pp. 85-95, doi: 10.1016/j.conbuildmat.2017.05.072.

T. Danner et al., *"Long-term Influence of Concrete Surface and Crack Orientation on Self-healing and Ingress in Cracks – Field Observations,"* Nordic Concrete Research, 58(1), 1-16. doi: <https://doi.org/10.2478/ncr-2018-0001>

Questions, Thoughts or Comments?

If you have questions or want to contact us during the Webinar, please **type your questions**, thoughts, or comments in the **Q&A box** and **press Submit**.

We ask for your understanding if we do not have time to discuss all comments and questions within the session.

Any unanswered questions or comments will be answered and discussed by e-mail or in another Webex session.



For more information, please contact us:


Bruker Nano GmbH

info.bna@bruker.com

Or

Max.Buegler@bruker.com

Andrew.Menzies@bruker.com



Innovation with Integrity