



# Bruker's IR Biotyper<sup>®</sup>: Paving the way to real-time outbreak control

Authors: Miriam Cordovana and Markus Kostrzewa, Bruker Daltonics GmbH & Co. KG (Bremen)

## Abstract

*Klebsiella pneumoniae*, a gram-negative bacterium of considerable clinical significance, has emerged as a formidable global health threat over recent decades. It is associated with a wide array of infections, from pneumonia and urinary tract infections to septicemia and hospital-acquired infections, posing a substantial challenge to healthcare systems worldwide. Given its notable propensity to acquire antibiotic resistance and thrive in diverse environments, *K. pneumoniae* requires timely epidemiological surveillance to facilitate infection control. In this study, the performance of the IR Biotyper<sup>®</sup>, a Bruker Fourier transform infrared (FT-IR) spectroscopy system, was evaluated for *K. pneumoniae* isolates typing [1].

## Introduction

*Klebsiella pneumoniae* plays a significant role in the occurrence of infections acquired both in community and hospital settings. The emergence of various multidrug-resistant (MDR) strains within high-risk *K. pneumoniae* clonal lineages, associated with distinct sequence types (ST), frequently leads to outbreaks, posing a substantial threat to public health.

Effectively tracking, comprehending, and mitigating *K. pneumoniae* outbreaks, along with the circulation of specific clones, presents a formidable challenge for infection control measures. Fourier transform infrared (FT-IR) spectroscopy has demonstrated its efficacy in bacterial typing across different intraspecies levels. This technique capitalizes on the distinctive infrared spectrum of each bacterial strain, serving as a unique fingerprint signature that facilitates precise identification.


Bacterial strain typing plays a pivotal role in investigating outbreaks by revealing clonality and possible transmission pathways as well as by establishing relatedness between clinical isolates and environmental sources. Selecting an appropriate typing method involves careful assessment of typeability, ease-of-use, interpretability, speed, reproducibility, cost-effectiveness and the availability of expertise. The IR Biotyper is a spectroscopic microorganism typing system based on FT-IR technology that achieves discriminatory power comparable to many molecular-based methods, delivering results within three hours, while being cost-efficient [2].



## Results

The study found that by applying HCA, the isolates were discriminated into six clusters, with the largest one comprising all ST307 isolates, and five smaller ones, corresponding to the non-ST307 isolates (Figure 1, [1]). IR Biotyper partitioning was in perfect concordance with the K locus (AR = 100%), in very good concordance with ST (AR = 99.1%), and lower concordance with the O locus (AR = 91%). PCA analysis showed a clear clustering of the outbreak isolates, as well as the clustering of the other isolates accordingly to the K locus (Figure 2, [1]).

**Figure 2**  
2D PCA scatterplot [1].  
Each geometric form corresponds to a spectrum,  
and the color corresponds to the locus K.



## Conclusion

In summary, the data demonstrated that the IR Biotyper proved to be a novel and reliable method to reveal clonality among *K. pneumoniae* isolates and could represent a faster, easier and more cost-effective alternative to molecular methods, in the context of outbreak detection.

Despite IR Biotyper clustering showing a slightly weaker correlation with ST than with the K locus, for outbreak isolates this correlation seems strong enough to propose the IR Biotyper as a reliable real-time tool.

The combination of the IR Biotyper's typing results with its low cost and very short turnaround time suggests that it is a promising tool for strain typing that could make real-time outbreak investigation a reality.

Click [here](#) for more information about Bruker Microbiology & Infection Diagnostics.

## References

- [1] Miriam Cordovana et al., 33rd ECCMID conference 2023, P0773: Fourier-transform infrared spectroscopy reliably differentiates outbreak isolates of *Klebsiella pneumoniae*.
- [2] Rakovitsky, N. et al. Fourier transform infrared spectroscopy is a new option for outbreak investigation: a retrospective analysis of an extended spectrum-beta-lactamase producing *Klebsiella pneumoniae* outbreak in a neonatal intensive care unit. *J Clin Microbiol*, 2020, 58(5): <https://doi.org/10.1128/jcm.00098-20>
- [3] Quintelas, C. Ferreira, E C. Lopes, J A. Sousa, C. An overview of the evolution of infrared spectroscopy applied to bacterial typing. *Biotechnol*, 2018, 13(1). <https://doi.org/10.1002/biot.201700449>
- [4] Comparing Partitions. Online tool for quantitative assessment of classification agreement, 2011. <http://www.comparingpartitions.info>

Not for use in clinical diagnostic procedures.  
Please contact your local representative for availability in your country.

IR Biotyper® is a registered trademark of the Bruker group of companies.

**Online information**  
[bruker.com/microbiology](https://bruker.com/microbiology)



**Bruker Daltonics GmbH & Co. KG**

Bremen · Germany  
Phone +49 (0) 421-2205-0

[info.md@bruker.com](mailto:info.md@bruker.com)

**Bruker Scientific LLC**

Billerica, MA · USA  
Phone +1 (978) 663-3660

