


## MagneticTwister™ Option for NanoTracker™ 2 Optical Tweezers

The NanoTracker™2 optical tweezers system comes with a large number of translational degrees of freedom allowing the user to manipulate the stage, sample chamber and optical traps independently with high precision. With the MagneticTwister™ option, an additional rotational degree of freedom becomes available. The three-dimensional manipulation possibilities provided by the optical trap are complemented by the ability to apply torque and rotate magnetically sensitive particles in defined steps.

### Torque Application and Twisting of Magnetic Particles

The custom-designed magnets in the MagneticTwister™ generate a homogeneous field of 700 mT at a distance of 1 mm in the sample chamber. Depending on the particle's size and susceptibility, forces in the Piconewton range can be applied. Both magnet-sample distance (and thus field strength) and rotation of the magnetic field can be adjusted via the NanoTracker™ 2 control software.



**Fig. 1:** Image of the unmounted MagneticTwister™ option. It replaces the detection objective holder and is easily aligned with the trap and camera optics. Controlled by an electric stepper motor, the magnets can be rotated 360° with an accuracy of 4° using the provided control software.

### Undisturbed Bright Field and Fluorescence Microscopy

Due to the tubular design, camera-assisted bright field microscopy remains available for the localization of the particle throughout the experiment. Being fully compatible with the approved modular design of the NanoTracker™ 2, no optical ports are obstructed. Thus, all fluorescence techniques that can be routinely combined with JPK's NanoTracker™ 2 (Epifluorescence, CLSM, TIRF) can also be used in combination with the MagneticTwister™ option.

### Single Molecule Torsion

By the use of functionalized chamber and magnetic bead surfaces, single molecules like DNA or chromatin can be manipulated and held in a fixed configuration by the optical trap(s) provided by the NanoTracker™ 2. Rotation of the magnets will simultaneously apply a torque to the trapped bead which is transferred to the clamped molecule. This will alter its mechanical pre-tension which, in the case of DNA, has profound effects on the molecule's ability to coil [1] or to interact with enzymes and other molecules [2]. In order to optically evaluate the orientation of the magnetic handle, functionalized nanorods with attached superparamagnetic beads can be bound to the molecule of interest [3] and then be exposed to the rotating magnetic field.

### Literature

- [1] W. Li, „Impact of DNA Twist Accumulation on Progressive Helical Wrapping of Torsionally Constrained DNA“, *Phys. Rev. Lett.*, Bd. 109, Nr. 21, 2012.
- [2] F. Kouzine, S. Sanford, Z. Elisha-Feil, und D. Levens, „The functional response of upstream DNA to dynamic supercoiling in vivo“, *Nature Structural & Molecular Biology*, Bd. 15, Nr. 2, S. 146–154, Jan. 2008.
- [3] A. Celedon, I. M. Nodelman, B. Wildt, R. Dewan, P. Searson, D. Wirtz, G. D. Bowman, und S. X. Sun, „Magnetic Tweezers Measurement of Single Molecule Torque“, *Nano Letters*, Bd. 9, Nr. 4, S. 1720–1725, Apr. 2009.