



# **MICRO-XRF**

# Analyzing thin layers under vacuum

# Application Note # XRF 457

# Introduction

The analysis of thin layers or coatings is a common task in micro-XRF spectrometry. Both the non-destructive operation of the method and the ability of X-rays to penetrate into sample and obtain information on the material beneath the surface make this method attractive for the purpose of analyzing single or multiple layers. The special challenge in analyzing the samples discussed here is that both layer (aluminum) and substrate (silicon) are light elements, which requires measurement under vacuum, because otherwise the air in the beam path between sample and detector would absorb the low energy radiation emitted by the sample. Additionally, this lab report compares manual and automated analysis using Auto-Point.

### Instrumentation

The measurements were performed with a Bruker M4 TORNADO. This micro-XRF spectrometer is equipped with a large vacuum sample chamber and uses a focused X-ray beam (spot size < 20  $\mu$ m) to induce fluorescence in the sample. This signal is analyzed with an energy dispersive detector. The M4 TORNADO combines high spatial resolution with fast data processing and a high speed motorized XYZ-stage for sample positioning.



#### Figure 1

Sample 1 video images with top: manual measurement locations, bottom: Auto-Point measurement locations.

### Results

As the numbers in Table 1 and Table 2 show, the thickness of the layers is around 2  $\mu$ m, as it is supposed to be. Although it is difficult to compare the results of the manual and the automated measurements due to the variations in the measured values and measurement locations, there is an indicator that the sample 2 layer is thicker than the sample 1 layer. The comparable standard deviations show that the shorter acquisition time of automated measurements suffices for reliable results.

The results obtained with the M4 layer quantification routine are verified by measurements on a fracture edge in the scanning electron microscope, as shown in Figure 2.

## Conclusion

This example shows that micro-XRF is a fast and accurate means to measure layer thicknesses at arbitrary locations on a sample non-destructively. The vacuum option extends the analytical range to light elements in layers and substrate.



Figure 2 SEM image of a fracture edge of layer (bright) and substrate.

Measurement point manual	Sample 1 Thickness / µm	Sample 2 Thickness / µm
Al-Si-1 1	1.16	2.22
Al-Si-1 2	1.50	2.71
Al-Si-1 3	1.94	2.18
Al-Si-1 4	2.10	2.77
Al-Si-1 5	2.15	2.21
Al-Si-1 6	2.25	2.64
Al-Si-1 7	-	2.56
Al-Si-2 1	1.79	2.20
Al-Si-2 2	1.90	2.17
Al-Si-2 3	2.04	2.22
Al-Si-2 4	2.15	2.18
Al-Si-2 5	2.20	2.26
Al-Si-2 6	2.24	2.16
Al-Si-3 1	1.26	1.85
Al-Si-3 2	0.91	1.64
Al-Si-3 3	1.55	-
Mean ± standard dev.	1.81 ± 0.42	2.26 ± 0.29

Measurement point auto	Sample 1	Sample 2 Thickness / µm
	Thickness / µm	
Line_1_1	0.98	1.52
Line_1_2	1.67	2.21
Line_1_3	1.76	2.73
Line_1_4	1.83	2.58
Line_1_5	1.87	2.24
Line_1_6	1.92	2.22
Line_1_7	1.97	2.60
Line_1_8	2.02	2.76
Line_1_9	2.07	2.52
Line_1_10	2.12	2.22
Line_1_11	2.16	2.25
Line_1_12	2.17	2.74
Line_1_13	2.18	2.78
Line_1_14	2.20	2.25
Line_1_15	2.23	2.21
Line_1_16	2.26	2.54
Line_1_17	2.24	2.83
Line_1_18	2.19	2.51
Line_1_19	2.17	2.09
Line_1_20	1.97	1.70
Mean ± standard dev	2.00 ± 0.92	2.37 ± 0.34

# Table 2

Table 1

Results of the manual

measurements.

Results of the Auto-Point measurements.

#### Author

Dr. Roald Tagle, Senior Application Scientist micro-XRF, Bruker Nano GmbH

Bruker AXS Info.baxs@bruker.com

