

EDULAB FOR STUDENTS: FOURIER 80

Munching on Sweets: The monitoring of catalytic hydrolysis of sucrose via invertase through benchtop NMR

An enzyme kinetic experiment with a sweet taste!

Authors & Affiliation:

Ronald Soong, William Wolff, Katrina Steiner, Kiera Ronda, Katelyn Downey, Peter Costa, Andre Simpson

Environmental NMR Center, University of Toronto, 1265 Military Trail, Toronto, ON, Canada, M1C 1A4

Experiment Hashtag #: #MunchingOnSweets #Educate2Resonate

Keywords:

Enzyme, Catalysis, Kinetics

Target group:

Undergraduate, Analytical Chemistry,

Biochemistry, Food Chemistry

Objectives:

- The determination of the rate of reaction via benchtop NMR.
- The understanding of enzyme catalyzed reactions.
- The processing and interpretation of NMR spectra.
- Introduction to relative quantitation by NMR.

Background of the Experiment:

Invertase is an enzyme that catalyzes the hydrolysis of sucrose into an equimolar mixture of glucose and fructose. This process is commonly used in the food industry to manufacture 'invert syrup' (Figure 1).




Figure 1 The Hydrolysis of sucrose to glucose and fructose by invertase.

Since NMR is both non-selective and time-resolved, it is an ideal platform for non-invasive reaction monitoring, especially in the case where the products or reactants are not detectable by other techniques (e.g., UV-vis spectroscopy). Since benchtop NMR systems do not require cryogenics, this permits them to be sited in a fume hood or directly attached to a chemical reactor to monitor chemical reactions. This has made it a potent tool for monitoring reaction progress and characterizing the processes and mechanisms of a reaction.

NMR is inherently quantitative: each signal is proportional to the number of nuclei in the sample, and thus the concentration. To obtain absolute quantitation, all of the nuclei must be given the time to return to thermal equilibrium which can be a time consuming process. A more convenient approach to reaction monitoring is relative quantitation, in which relative changes in peak areas can be monitored over time. This allows for rate constants to be obtained within shorter time periods. In this experiment, the hydrolysis of sucrose by invertase will be tracked by relative quantitation.

Preparation & Prerequisite

This lab is designed to take ~ 6 hrs as a two-week experiment, and it is assumed that students have basic knowledge of 1D NMR and basics of spectral interpretation. The experiments are designed to be ideally completed in groups of 3-6 students. This investigation aims to demonstrate key NMR concepts, including interpretation of 1D spectra, relative quantitation, and first order kinetics. Prior to carrying out these experiments, it is strongly recommended to be familiar with basic 1D NMR processing.

This information is readily available in the version 001 Fourier EduLab Students Guide, which can be found on the USB stick delivered with the Fourier 80. In addition, having a basic understanding of processing software such as MestreNova is strongly recommended.

To perform this experiment, a properly installed and adjusted Fourier 80 system with TopSpin Software is required. Fourier 80 equipped with pulse field gradient is optional. In addition, a 10 mL volumetric flask & stopper, 10 uL micropipettes, 1000 uL micropipettes and vortex should be available.

Experimental Setup:

- Food grade invertase, purchased online
- Sucrose
- 5 mm disposable NMR tube with cap
- Pulse program: zg, zgpgw5

Glossary

Hydrolysis:

Hydrolysis is any chemical reaction in which a molecule of water breaks one or more chemical bonds

NMR:

Spectroscopic analytical technique based on radio frequency-induced transitions between energy levels that atomic nuclei adopt in an external magnetic field as a result of their own magnetic moment

UV-vis spectroscopy:

UV spectroscopy or UV-visible spectrophotometry refers to absorption spectroscopy or reflectance spectroscopy in part of the ultraviolet and the full, adjacent visible regions of the electromagnetic spectrum

Abbreviations

NMR:

Nuclear Magnetic Resonance

UV-vis Spectroscopy:

Ultraviolet-visible spectroscopy

Results & Discussion:

To complete the exercise, answer the following questions:

1. The protons attached to the alpha carbon were used to track the progress of the reaction. Why weren't hydroxyl groups chosen to monitor the reaction?

2. Aside from the relaxation delay (D1), what changes could be made to the experiment conditions to obtain absolute quantitation of glucose and sucrose?

3. Fit the results to a first-order kinetics decay curve. Is the time resolution sufficient to establish a rate constant? What could be changed about the experiment to enhance the time resolution while obtaining a comparable sensitivity?

4. These experiments were run in 100% D₂O. What would happen if these experiments were to be run in pure H₂O?

5. Advanced Question. Sample results of the invertase experiment are shown in Figure 2. After production the glucose anomeric peak decreases. Can you work out why? A hint the peak plotted and monitored is the peak for α -glucose. The peak for β -glucose (if applicable) would be under the water. Do your experimental results match the sample data? If not, why?




Figure 2. (A-B) A series of NMR spectra as a function of time showing only the anomeric doublet of the (A) glucose and (B) sucrose. (C) A graph of normalized integrals of sucrose and glucose as a function of time. (D) A spectrum showing a mixture of glucose and sucrose during enzymatic hydrolysis reaction.

Key Take Home Messages:

This experiment demonstrates the ability of NMR to monitor an enzymatic reaction. In this case, the enzymatic hydrolysis of sucrose into glucose via invertase. Through this experiment, you learned the following NMR concepts:

- Processing of NMR spectra and integration of NMR resonances.
- Analysis of NMR spectra.
- Using relative quantitation to determine kinetic parameters by NMR.

References:

1. Kehlbeck, J.D.; Slack, C.C.; Turnbull, M.T.; Kohler, S.J.; Exploring the Hydrolysis of Sucrose by Invertase Using Nuclear Magnetic Resonance Spectroscopy: A Flexible Package of Kinetic Experiments. **J. Chem. Educ.** 2014, 91, 5, 734–738.
2. Vang, J.Y.; Her, C.; Krishnan, V.V.; NMR based real-time enzyme kinetics on estimating the inhibitory effect of sucralose in the enzymatic conversion of sucrose. **Biophys. Chem.** 2021, 268, 106495
3. Her, C.; Singh, J.; Krishnan, V.V.; Effect of sucralose on the enzyme kinetics of invertase using real-time NMR spectroscopy and progress curve analysis. **Carbohydr. Res.** 2018, 455, 5-9

