



Software

BioPharma Compass

Accelerate Routine BioPharma
Characterization

Innovation with Integrity

Faster insights in your BioPharma data

BioPharma Compass® increases the productivity of biologics characterization labs with integrated workflows that take full advantage of the benefits of high isotopic fidelity data, complementary ionization modes and trapped ion mobility.

Workflows for multi attribute peptide and protein screening, peptide mapping, oligonucleotide sequence verification as well as top-down sequencing and critical reagent QC are available. This empowers users with the tools to rapidly turn MS data into actionable insights.

Simple interface to get productive quickly


- Single GUI for acquisition and processing
- Entire process contained in a single workflow

Powerful algorithms

- Optimized for high isotopic fidelity data
- Multi attribute screening with collisional cross section (CCS)
- Fast PASEF® peptide maps with CCS-enabled selectivity

Easier implementation of GLP guidelines

- State-of-the-art data security
- Version tracking for all editable items
- User management and audit trails





Algorithms built for instruments with **high isotopic fidelity** deliver sub-ppm mass accuracy.

Support for timsTOF, QTOF and MALDI-TOF offers users a comprehensive range of options to optimize **resolution, speed and selectivity** for their application and budget.



Confident characterization of oligonucleotides and RNA

RNA and oligonucleotides are an increasingly important class of molecules. Whether an RNA molecule is used directly as API or as a tool for gene editing, it is necessary to comprehensively characterize its primary sequence, modifications and impurities.


In particular the MS2 analysis of these compounds raises unique interpretation challenges due to the complexity of the fragmentation pattern and the limited number of nucleic acids building blocks. Embedded in BioPharma Compass, **OligoQuest™** provides BioPharma scientists with the tools to comprehensively verify oligonucleotides and related impurities by accurate mass and CID fragmentation.

Comprehensive tools for RNA and oligonucleotide sequence verification


- Flexible definition of building blocks nomenclature and nucleotides modifications
- Deconvolution and peak picking optimized for high isotopic fidelity data
- Interpretation of termini and internal MS2 fragments enabled by high mass accuracy and fidelity, and resolution of fragment ions.



Oligonucleotide LC-MS analysis



Deconvolution and monoisotopic annotations




MS2 based impurities confirmation




“Isomeric Oligonucleotides can be analyzed with OligoQuest™ allowing to identify base exchanges. Already the first version could be included in our day-to-day workflow.”

Team at Axolabs GmbH Kulmbach, Germany



OligoQuest workflow

Confident annotation of 1 Da C>U variants using high isotopic fidelity data



d	1.6	4.3	-0.7	-3.9	3.2	4.3	5.3						1.9	4	5.5	3.2		2.5	2.6	5.1				
c	3	3.8	3.1	0.9	0.6	4.2	3	4.1	2.8	3.6	2.9	5.9	-0.2	1.8				4.4	5	2.6	4.6			
b	6.2	3.2	-0.6	0.8		4.8	4.1	2.6	2	3.6			2.2			1.2		2.4	5.4	3.3				
a	5.1	2.9	3.7			4.8		2.7	2.8				3.6			1.9					3.2			
a-B	3.2	3.5	2.8	2.8	3.3	2		4	5								3.4	3	0.7	4.7	4.8	4.1	4.6	1.7
internal >																								#1
5' Index	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Sequence	g	c	a	c	g	c	g	u	g	c	u	u	u	u	g	c	a	c	g	c	g	u	g	c
3' Index	24	23	22	21	20	19	18	17	16	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1
internal overlap		#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1	#1
internal <	#1																							
w	5.1	4.3	5	4.9	5	4.6	3.2	4.1	4	1.9						4.1	4.3	1.5	1.3	2.7	2.9		1.6	-2.7
x	4.6						1.8	3	5.9	2.9	3.6	2.8	4.1	5.2	4.2	5.9		6.1		5.4	3			
y				5.3		4.4	5.5	1.2		2.2	3.6	2	3.3	4.1	5.1	4.7	2.6	3.9	4.4	6.2				
z	3.2			4.5		1.9	3.7		3.6	2.8	2.7		4.8	2.9		-1.3		1.8	5.1					
Sequence Validation	16	15	15	15	17	16	18	16	15	12	12	15	15	17	16	15	15	14	14	15	14	14	14	14



Routine characterization of protein therapeutics


In recent years, monoclonal antibody-based therapeutics have been highly represented in both new drug approvals and clinical pipelines. For the related development projects it is essential to characterize the molecule's primary sequence as well as identify and track modifications that are identified as critical quality attributes. The intact analysis of mAbs or related domains requires high mass accuracy in order to confidently confirm the molecule's identity.

BioPharma Compass is optimized to take advantage of high isotopic fidelity data for intact mass analysis.


Heterogeneities of interest can be further characterized by top-down analysis or peptide mapping.

Confident mAb characterization


- Deconvolution and deisotoping optimized for high isotopic fidelity data
- Simulation of theoretical isotopes and butterfly views for visual comparison
- Full automation from sample submission to report




Clone Screening



Intact, reduced and subunit analysis




Peptide mapping






mAb light chain measured with high isotopic fidelity for sub-ppm mass accuracy

Intact protein and drug distribution analysis




Tools optimized to gain the most insights in your intact and subunit analysis for characterization or clone screening. Carry out multiple attribute method (MAM) analysis at the intact level with speed and confidence.


- Algorithms designed to take advantage of data with high isotopic fidelity.
- Flexible peak annotation system
- Support for top-down analysis the investigation of unknown peaks



Automated annotation of clipped fragments









Top-down sequencing by ESI or MALDI



Peptide level analysis

From database search for HCP analysis to in-silico digest for peptide mapping and development of multiple attribute monitoring methods, fully characterize and quantify your molecule's impurities and heterogeneities.

- Interactive sequence coverage maps simplify data review
- Annotated chromatograms and MS2 spectra to verify primary sequence and identify modifications
- Easy conversion of peptide map data into MAM screening methods for EIC based quantitation
- Easy conversion of peptide map data into MAM screening methods



Fast sample screening


Row	Result	Position	Base Peak Mr [Da]	Δ Base Peak Mr [Da]	Sample Name	Matrix	Operator
1		B3	555.2642	-0.0051	Leu-Enkephalin	CCA	AR
2		C3	756.3875	-0.0043	Bradykinin fragment	CCA	AR
3		D3	898.4560	-0.0101	Angiotensin III	CCA	AR
4		E3	1045.5256	-0.0090	Angiotensin II	CCA	AR
5		F4	1295.6731	-0.0044	Angiotensin I	CCA	AR
6		G3	1362.7093	15.9811	Substance P	CCA	AR
7		H3	1551.6370	-18.0325	Glu-Fib B	CCA	AR
8		I3	1618.7885	-0.0266	Bombesin	CCA	AR
9		J3	1757.8914	-0.0339	Renin Substrate	CCA	AR
10		K3	2092.0380	-0.0409	ACTH clip 1-17	CCA	AR
11		L3	2464.1436	-0.0475	ACTH clip 18-39	CCA	AR

Confident detection of analytes based on accurate mass, collisional cross section and retention time.


Flexible workflow can be applied to the QC of oligonucleotides, synthetic peptides, glycan mixtures or raw materials.

Rapidly analyze large sample batches with automated analysis and data processing. Results are easily visualized with traffic light based attributes reporting.

CCS aware screening of leucine and isoleucine containing peptides



Released glycan profiling with MALDI



Peptides and oligonucleotides screening (ESI/MALDI)



Accelerate your routine BioPharma assays

BioPharma Compass simplifies data acquisition and processing for users of Bruker QTOF, timsTOF and MALDI instruments with an interface optimized for the routine characterization of therapeutics proteins and oligonucleotides.

- Automation of routine characterization tasks from acquisition to report from a single interface
- Integrated tools assist with data integrity requirements
- Workflows supporting major analytical tasks for proteins and oligonucleotides analysis

Bruker Daltonics is continually improving its products and reserves the right to change specifications without notice. © BDAL 05-2023, 1901123

For Research Use Only. Not for use in clinical diagnostic procedures.

Bruker Switzerland AG

Fällanden · Switzerland
Phone +41 44 825 91 11

ms.sales.bdal@bruker.com - www.bruker.com

Bruker Scientific LLC

Billerica, MA · USA
Phone +1 (978) 663-3660

Online information
bruker.com

