MALDI HiPLEX-IHC: Highly Multiplex & Multiomic Tissue Imaging with Photocleavable Mass-Tagged Probes Gargey B. Yagnik; Ziying Liu; Zhi Wan; Kenneth J. Rothschild, Mark J. Lim AmberGen AmberGen Incorporated, 313 Pleasant Street, Watertown, MA, 02472

Introduction

Multiomics approaches stand to transform our ability to understand human disease by providing a comprehensive picture of the immense complexity, heterogeneity and interplay within biological systems. However, "omic" data must not forfeit spatial information through homogenization of the biospecimen, but rather, must be put into a spatial context at the multi-cellular and tissue levels. This is best achieved by tissue imaging approaches, yet current methods lack the multiplexity and flexibility to obtain multiomic scale information. For example, immunohistochemistry (IHC) provides an important and widely used tool for researchers and pathologists to image multiple protein biomarkers in tissue specimens. However, standard fluorescence based IHC is generally limited to 3-5 different biomarkers (hyperspectral/multispectral methods <10). IHC also cannot image small molecules in tissue specimens such as metabolites and drugs. While mass spectrometry (MS) is a proven proteomic/metabolomic tool and the advent of mass spectrometric imaging (MSI) has extended MS to the spatial dimension, it is generally limited to untargeted analysis of small molecules and peptides. To overcome this barrier, MALDI HiPLEX-IHC (MALDI-IHC) uses probes such as antibodies, lectins and oligonucleotides conjugated to photocleavable mass-tags (PC-MTs) for MSI of targeted proteins and other macromolecules in cells and tissues. MALDI HiPLEX-IHC significantly exceeds the multiplexity of both fluorescence and previous cleavable mass-tag based methods, without the need for iterative cycling procedures. Moreover, we have combined on the same tissue section untargeted MSI of endogenous small molecules with this targeted protein MSI approach for a truly unique multiomic capability. Novel dual-labeled fluorescent MALDI HiPLEX-IHC probes extend the utility of this new approach to allow co-registration of low-plex but high spatial resolution fluorescence images with MALDI HiPLEX-IHC images of the same tissue section, such as for cell segmentation purposes. Overall, this approach stands to transform the fields of "omics" based research and discovery, tissue pathology, tissue diagnostics, drug development, therapeutics and precision medicine.

Instrumentation Bruker rapifleX & timsTOF





Example Immunofluorescence Validations



*Dual-Labeled Probes: Same Tissue Section & Same ROI Shown •For Fluorescent-Only Probes: Adjacent Sections & Different ROIs Shown •Fluorophores used for Probes:

- DyLight 650: Myelin, NeuN, Syn, NRGN, PanCK, CD20, αSMA, HER2 • DyLight 594: Vimentin
- DyLight 550: CD3, Aß42
- Fluorescein: GLUT1



Human Tonsil

Dual-Labeled Fluorescent Probes for Multimodal Imaging

1.	Pan-CK
2.	CD3
3.	CD4
4.	CD8
5.	CD20
6.	CD45RO
7.	FoxP3
8.	ER
9.	PR
10.	Her2
11.	dsDNA
12.	CD68
13.	Ki67
14.	ACTA2
15.	VIM
16.	FAP
17.	PDPN
18.	PDGFRB
19.	PTEN

23-Plex Imaging with Antibody & Lectin Probes



19 Antibody Probes



[Yagnik, Liu et al. (2021) J Am Soc Mass Spectrom 32(4): 977-988 pubs.acs.org/doi/10.1021/jasms.0c00473] **Funding:** Work funded in part by SBIR grant R44 CA236097 from the NIH (NCI) to AmberGen For more information, please contact us at info@ambergen.com or go to our website at www.ambergen.com





Breast Cancer Whole Section

- Positive Ion, Reflector
- 10 kHz Laser Rate
- 20 µm Resolution

Conclusions

• Multiplex: High-Plex (up to 200) Probe-Based MALDI Imaging of Intact Proteins • **Multiomic:** MALDI Imaging of label free small molecules and intact proteins on same tissue sample • Multimodal: Fluoresence and MS images on same tissue sample with Dual-Labeled Probes