

Maximizing information content in ion mobility-enhanced DIA using overlapping, ion mobility-encoded quadrupole windows

Ute Distler¹; Mateusz Krzysztof Łacki¹; Michał Piotr Startek^{1,2}; David Teschner³; Sven Brehmer⁴; Jens Decker⁴; Oliver Raether⁴; Andreas Hildebrandt³; Stefan Tenzer^{1,5}

¹University Medical Center Mainz, Mainz, Germany; ²University of Warsaw, Warszawa, Poland; ³Institute of Computer Science, Johannes Gutenberg University, Mainz, Germany; ⁴Bruker Daltonics GmbH & Co.KG, Bremen, Germany; ⁵Helmholtz Institute for Translational Oncology (HI-TRON), German Cancer Research Center (DKFZ), Mainz, Germany



INTRODUCTION


Cycling deterministically through segments of a predefined precursor m/z range, data-independent acquisition approaches provide a comprehensive record of all detectable precursor and fragment ions by isolating and fragmenting populations of different precursor ions. The recently introduced diaPASEF method utilizes the correlation of molecular weight and ion mobility in a trapped ion mobility device (timsTOF Pro 2) to extend sensitivity and specificity of established targeted data extraction workflows by the additional ion mobility dimension. Here, we discuss a novel DIA scan mode, which uses mobility-specific micro-encoding of overlapping quadrupole windows to maximize information content in DIA acquisitions.

METHODS

Total protein extracts from HeLa cells were digested and analyzed via nanoLC coupled to a timsTOF Pro 2. Spectral libraries were generated from triplicate PASEF-DDA acquisitions of high pH reversed-phase fractions of each sample. Overlapping ion-mobility-dependent quadrupole windows were defined based on spectral library data to provide >99% library coverage over the full mass range. Datasets were acquired on a modified timsTOF Pro 2 platform in both diaPASEF and several variations of the novel acquisition scheme (varying cycle time, quadrupole window size and overlap) with and without collision energy to evaluate and optimize acquisition workflow performance. Precursor-fragment deconvolution tools were trained on unfragmented HeLa datasets.

DATA ACQUISITION CONCEPT

diaPASEF mobility-scan DIA



RESULTS




Figure 1: Mobility-scan DIA produces unique data patterns and increases sensitivity by overlapping diagonals.




Figure 2: Mobility-scan DIA data complexity on fragment ion level. DIA analysis of 200 ng HeLa digest on a modified TimsTOF Pro2 platform. Extracted ion chromatograms of the indicated mass ranges were generated in Data Analysis viewer. Discernible features are highlighted.




Figure 3: Low-intensity mobility-scan DIA data pattern enable precursor position modelling. Algorithmic prediction of the most likely precursor ion position based on a data occupancy matrix enables to pinpoint precursor ion positions with an accuracy of a few Da.




Figure 4: Database search results of deconvoluted mobility-scan DIA data. Deconvoluted mobility-scan DIA data from a HeLa digest were exported as .mgf files and submitted to Mascot Search. Exemplary peptide identifications are shown.

	DDA1	DDA3	DDA5	DDA6	DDA8	DDA9	DDA10	MSDIA	count	DDA	
b3	2402.9	7136.4	459.76	448.28	255.36	170.52	50.001	482.28	17728	8	
b4	802.82	4029.7	488.96	316.39	29.195	28.042	73.934	184.24	5090	8	
b6-NH3	392.58	1400.8	137.57	218.76	129.03	66.787	104.03	118.07	2378	8	
y10	39485	138130	5478.0	15052.6	4388.1	3257.0	106.68	3531.3	121624	8	
y8	2857.4	11812	345.2	984.8	221.89	255.36	231.04	148.66	16779	8	
y6	2014.8	5093.6	130.26	349.87	132.03	68.579	50.714	37.041	5202	8	
y6	3936.4	10116	399.31	1518.63	175.72	300.07	131.93	48.535	15108	8	
y7	2533.6	10035	228.47	638.9	224.72	183.46	48.246	35.355	11959	8	
y8	4452.4	15756	479.13	1066.8	532.38	405.54	93.317	205.98	21864	8	
b2	2093.4	7265.4	438.35	133.21	49.797	35.962	207.79			7	
y5	1634.9	7101.5	317.82	522.32	176.15	52.81	79.635	11836		7	
y6	7213.9	26118	622.92	2620.7	504.96	1010.3	1298.5	39181		7	
b9	578.22	3414.4	254.247	38.824			19.307	58.056	3541	6	
b7	882.74	3798.1	219.591	283.45	38.36			154.66	2556	6	
y10-NH3	5574.7	16047		372.06	755.11	419.4	118.61		11893	6	
b6-NH3	498.55	1702.4		57.324	20.401	26.124	0.266		2238	6	
b5	915.36	3124.3	297.7	263.26			33.599		2915	5	
b6-NH3	669.57	2508.9	68.134	73.981			72.791	2458		5	
b5-NH3			106.6	111.64			73.777	42.717		4	
b8-H2O	168.01	1250.3		10.573	53.153			2812		4	
b9	956.08	2221		92.37	250.89			3119		4	
b9-H2O	238.98	2098.8		128.63	96.742			2732		4	
b9-NH3	487.51	712.93		4.5817			39.956		2888	4	
y10-H2O	1185.8	135.71		100.85				38.452		4	
y11	5961.5	18546		1534.65			580.11		6400	4	
y11[2+]				17.314			117.68		136.37	13446	4

Figure 5: Exemplary fragment ion coverage provided by mobility-scan DIA compared to eight reference DDA runs.

CONCLUSIONS

- Novel quadrupole scanmode provides the basis for comprehensive, fast and highly sensitive DIA
- Using overlapping quadrupole selection windows, the method can theoretically provide a 2.5-fold increase in fragment ion sensitivity compared to diaPASEF
- To enable in-depth data processing of resulting high-complexity datasets, we developed algorithms for multidimensional peak detection and to classify whether the intensity-distribution of a given fragment matches the one of a target precursor to generate highly specific pseudo-MSMS spectra on the basis of multidimensional deconvolution (RT, IMS, quad window), which can be converted into common MS/MS data formats such as mgf and searched directly with established tools like PEAKS and Mascot

ACKNOWLEDGEMENTS

This work was supported by grants from the the BMBF (MSCORESYS, DIASyM FKZ 161L0217A), as well as the Forschungszentrum Immuntherapie (FZI).