

Enhanced 4D workflows using TIMS for advancing small molecule research

Aiko Barsch¹; Niels Goedecke¹; Matthias Szesny¹; Sven W. Meyer¹; Viola Jeck¹; Mohamed Elsadig¹; Michael Krause¹; Birgit Schneider¹; Carsten Baessmann¹

Bruker Daltonics GmbH, & Co KG
Fahrenheitstraße 4, 28359 Bremen,
Germany

Introduction

Small molecule analyses in metabolomics, pharma and applied research often face similar challenges requiring the analysis of complex samples for detection, identification and (semi-) quantification of known and unknown analytes across a broad mass and analyte concentration range. Recent advances in high resolution LC-MS in combination with trapped ion mobility separation (TIMS) address the need for combining rapid analyses with confident identification. The ability of TIMS to separate co-eluting isobaric and isomeric analytes and measure collisional cross sections (CCS) can increase confidence in compound identification. The TIMS-enabled PASEF data acquisition mode provides extensive MS/MS precursor coverage in every injection. Here we present an optimized workflow for 4D metabolomics using LC-MS and PASEF data acquisition capable of capturing analytes across a wide mass range in a single analysis, with dedicated data processing for broad profiling/screening applications.

Methods

Samples: Human urine; centrifuged and filtered (0.22µm). Methanolic human plasma (SRM 1950) extract. Custom standard mix including Restek Pharmaceutical Mix 1, Agilent Forensic Toxicology Comprehensive Mix, Imidazole, Tylosine A (both Sigma Aldrich).

LC: Elute UHPLC, Intensity Solo C18 column (Bruker).

- LC gradient according to T-ReX LC-QTOF Solution (Bruker), allows matching of retention times for >600 metabolomics relevant compounds [1].
- LC gradient according to TargetScreener 15 min gradient allows for annotation of > 2800 drugs, pesticides or mycotoxins [2].

MS: timsTOF Pro 2 (Bruker) equipped with VIP-HESI source

Acquisition: Optimized broad range PASEF and bbCID acquisitions

Software: TASQ and MetaboScape 2022b, preliminary version (Bruker).

Target Lists for annotation in MetaboScape:

- Target List 1: Compounds contained in Unified CCS Compendium [3] containing CCS, appended with retention times (RT) from Bruker HMDB Personal Library 2.0 and MS/MS references from Bruker HMDB Personal Library 2.0, Bruker MetaboBase Personal Library 3.0 and NIST 2020, respectively.
- Target List 2: Metabolites present in serum derived from HMDB 5.0 [4] including name, molecular formula and structure but no assigned MSMS or CCS reference values.


References

- <https://www.bruker.com/en/applications/academia-life-science/metabolomics/metabolomics-solution.html>
- <https://www.bruker.com/en/products-and-solutions/mass-spectrometry/ms-solutions/targetscreener.html>
- <https://doi.org/10.1039/C8SC04396E>
- <https://hmdb.ca/>
- <https://doi.org/10.1186/1471-2105-11-148>
- <https://doi.org/10.1186/s13321-016-0115-9>

Note: HMDB and CCS Compendium are not Bruker products.

Results

A) LC-TIMS-MS/MS for broad profiling in metabolomics, pharma and applied research



The novel data acquisition utilizes TIMS stepping and provides broad profiling and screening/quantitation in a single acquisition for target samples.


4D-Metabolomics™ data acquisition of a human urine sample spiked with a custom reference standard mix highlights wide mass and mobility transfer: Spiked compounds detected from 69 m/z (Imidazole) to 916 m/z (Tylosin A).

This optimized transfer methods were combined with two complementary MS/MS acquisition modes. Parallel accumulation serial fragmentation (PASEF) (see B) and broad band CID (bbCID) (see C) which are routinely used for profiling and screening/quantitation workflows, respectively.

B) Profiling and annotation: DDA workflows with PASEF & MetaboScape


For routine profiling of complex samples the optimized transfer method was combined with parallel accumulation serial fragmentation (PASEF) data dependent MS/MS acquisition. PASEF provides fast MS/MS acquisition, increasing the depth of coverage for all small molecules. For methanolic extracts of human plasma this provided an MS/MS coverage of ~60-70% of extracted features. TIMS intrinsically generates ion mobility information which is automatically recalibrated and transformed into CCS values in MetaboScape® by the T-ReX® 4D feature extraction algorithm. CCS values in combination with MS/MS, retention time, accurate mass, and isotopic pattern fit allow for annotation of known target compounds with high confidence.

1. Annotation of knowns:



The automatic annotation of known metabolites using Target List 1 including reference CCS values from CCS Compendium, retention times (from Bruker HMDB 2.0 Library) and reference MSMS spectra, allowed for confident assignment of e.g. Indole-3-acetate. Mass error was below 1 ppm, retention time deviation below 0.1 min vs. HMDB 2.0 reference, MS/MS score >990 and CCS below 1% vs. CCS Compendium reference [3].

2. Annotation of unknowns:




Annotating the human plasma data set with Target List 2 of compounds reported to be present in serum (derived from HMDB 5.0 [4]) enabled tentative annotation of further metabolites. Note: This Target List could be appended with RTs, MS/MS spectra and/or CCS values from reference repositories as highlighted in B1. In case these are not readily available, MetaboScape can perform automatic CCS prediction and MS/MS matching based on InChI encoded structures. This is based on the CCS-Predict Pro model and MetFrag [5, 6], respectively. Here it permitted tentative annotation of > 50 extracted features (not shown) with high scoring (<1 ppm precursor mass, <20 mSigma, >800 MetFrag MSMS score, < 2% CCS vs. prediction). The example highlighted shows the tentative annotation of gamma-Glutamyl-tyrosine MetFrag MS/MS score 992 (max. 1000); delta CCS vs prediction 0.3%. Next step: Validate by comparison to reference standard.

C) Screening and quantitation: DIA workflows with bbCID and TASQ

Data independent broadband CID (bbCID) acquisition continuously cycles between low collision energy and elevated collision energy (bbCID) MS/MS. This provides seamless collection of accurate mass precursor, true isotopic pattern intensities and fragment qualifier ion data in a single analysis. In combination with the optimized method for broad transfer, TIMS-bbCID provides enhanced quantitation capabilities by mobility separating target compounds from co-eluting background noise.

Evaluating this data and considering high resolution diagnostic qualifier ion information and retention time information from the TargetScreener HR solution [2] means false positive detections are minimized, delivering maximum confidence in screening result.

A dilution series acquired from reference standards using TIMS-bbCID data was investigated using the TASQ software solution and linear dynamic ranges of quantitation were determined. The example shown for Diazepam shows >3.0 orders of magnitude linear dynamic range (0.1 to 250 ppb). The residuals plot shows deviations of <10% for all concentration levels, proving the capability of the established method to enable quantitation across wide dynamic ranges.



Summary

We developed a CCS-enabled, single shot data acquisition for routine small molecule applications. This new workflow can be used for both targeted and untargeted samples and combines the confidence of high MS/MS coverage with CCS, resulting in a powerful workflow for broad screening, profiling, quantitation and ID capabilities for small molecule research.

Conclusion

- The CCS-enabled, single shot data acquisition routines provide deeper insights for metabolomics, pharma and applied research.
- The CCS-enabled profiling workflow provides high confidence for automated compound annotation.
- The CCS-enabled screening/quantitation provides more than three orders of magnitude in linear dynamic range, separating target compounds from co-eluting background noise.
- Both workflows (profiling using PASEF MS/MS and screening/quantitation using TIMS-bbCID MS/MS) provide up to five confidence criteria for evaluation in MetaboScape and TASQ, respectively:

- Accurate mass
- Isotopic pattern fit
- Retention time
- MS/MS information
- CCS values

4D-Metabolomics