

Mercury Analysis in Food, Environmental and Medical Samples with TXRF



Bruker Nano GmbH, Berlin, Germany
Webinar 09.12.2020

The background is a blue-toned graphic. On the left, a portion of the periodic table is shown with elements like Ca, Sc, Ti, Sr, Y, Zr, Ba, La, Hf, Ra, and Ac. In the center, the text "XFlash® Technology" is displayed above a spectral plot with peaks labeled with elements like V, Cr, Mn, Fe, Ni, Cu, Zn, Ga, Ge, As, Se, Br. Below the plot are labels for X-ray lines: $K\alpha$, $K\beta$, $L\alpha$, $L\beta$, $M\alpha$, and $M\beta$. On the right, another portion of the periodic table is shown with elements like V, Cr, Mn, Fe, Co, Ni, Cu, Zn, Ga, Ge, As, Se, Br, Ta, W, Re, Os, Ir, Pt, Au, Hg, Tl, Pb, Bi, Po, At, Rn, Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu, Th, Pa, U, Np, Pu, Am, Cm, Bk, Cf, Es, Fm, Md, No, and Lr. The large text "TXRF" is positioned on the right side. At the bottom right, there is a 3D rendering of a TXRF instrument with a sample tray and a detector.

Welcome



Today's Topics

- TXRF – how does it work?
- Mercury – Where does it come from and where does it go?
- TXRF analysis of mercury
- Comparison with Atomic Spectroscopy methods
- Interactive Q & A

Speakers

Dr. Hagen Stosnach
Applications Scientist TXRF
Berlin, Germany



Dr. Armin Gross
Global Product Manager TXRF
Berlin, Germany




TXRF – How does it work?

Principles of total reflection X-ray fluorescence (TXRF) spectroscopy



Total reflection X-ray fluorescence spectroscopy



Beam angle: $0^\circ / 90^\circ$

- Samples must be prepared on a reflective media
- Polished quartz glass or polyacrylic glass disc
- Dried to a thin layer, or as a thin film or microparticle

Principles of total reflection X-ray fluorescence spectroscopy




Samples for TXRF

- Powders: Direct preparation or as suspension
- Liquids: Direct preparation
- Always as a thin film, micro fragment or suspension of a powder
- Necessary sample amount: Low μg respectively μl range

Simple quantification

- ➔ Matrix effects are negligible due to thin layer
- ➔ Quantification is possible by internal standardization




Principles of total reflection X-ray fluorescence spectroscopy



In TXRF the samples are prepared as thin films or layers

- Matrix effects are negligible
- Quantification is possible





Principles of total reflection X-ray fluorescence spectroscopy



In TXRF the samples are prepared as thin films or layers

- Matrix effects are negligible
- Quantification is possible


- TXRF detects elements from Na(11) to U(92)
- The element sensitivities depend on the atomic number
- The sensitivity factors are calibrated ex works
- Quantification requires the addition of one standard element



Mercury – Where does it come from and where does it go?

Mercury

Background information



Mercury

Background information




Sources for mercury uptake

- Fish & other seafood, non-alcoholic beverages and composite food are the most important contributors to inorganic mercury dietary exposure

Toxicity of mercury

- The critical target for toxicity is the kidney
- Other targets include the liver, nervous system, immune system, reproductive and developmental systems
- Methylmercury is able to enter the hair follicle, and to cross the placenta as well as the blood-brain and blood-cerebrospinal fluid barriers, allowing accumulation in hair, the foetus and the brain



Mercury

Background information



Regulations

Maximum levels for mercury in certain foods
(Commission Regulation (EC) No 1881/2006)

3.3	Mercury	Maximum level (mg/kg wet weight)
3.3.1	Fishery products ⁽²⁶⁾ and muscle meat of fish ⁽²⁴⁾ ⁽²⁵⁾ , excluding species listed in 3.3.2. The maximum level applies to crustaceans, excluding the brown meat of crab and excluding head and thorax meat of lobster and similar large crustaceans (<i>Nephropidae</i> and <i>Palinuridae</i>)	0,50
3.3.2	Muscle meat of the following fish ⁽²⁴⁾ ⁽²⁵⁾ : anglerfish (<i>Lophius species</i>) atlantic catfish (<i>Anarhichas lupus</i>) bonito (<i>Sarda sarda</i>) eel (<i>Anguilla species</i>) emperor, orange roughy, rosy soldierfish (<i>Hoplostethus species</i>) grenadier (<i>Coryphaenoides rupestris</i>) halibut (<i>Hippoglossus hippoglossus</i>) marlin (<i>Makaira species</i>) megrin (<i>Lepidorhombus species</i>) mullet (<i>Mullus species</i>) pike (<i>Esox lucius</i>) plain bonito (<i>Orcynopsis unicolor</i>) poor cod (<i>Tricopterus minutes</i>) portuguese dogfish (<i>Centrosymnus coelolepis</i>) rays (<i>Raja species</i>) redfish (<i>Sebastes marinus</i> , <i>S. mentella</i> , <i>S. viviparus</i>) sail fish (<i>Istiophorus platypterus</i>) scabbard fish (<i>Lepidopus caudatus</i> , <i>Aphanopus carbo</i>) seabream, pandora (<i>Pagellus species</i>) shark (all species) snake mackerel or butterfish (<i>Lepidocybium flavobrunneum</i> , <i>Ruvettus pretiosus</i> , <i>Gempylus serpens</i>) sturgeon (<i>Acipenser species</i>) swordfish (<i>Xiphias gladius</i>) tuna (<i>Thunnus species</i> , <i>Euthynnus species</i> , <i>Katsuwonus pelamis</i>)	1,0

Mercury

Background information



Regulations

- US EPA Safe Drinking Water Act:
Maximum level 2 $\mu\text{g/l}$
- German drinking water regulation:
Maximum level 1 $\mu\text{g/l}$
- International norms for beverages:
Maximum levels 5 – 50 $\mu\text{g/l}$





Mercury analysis by TXRF - Instrumentation

Mercury analysis Instrumentation



S2 PICOFOX

- Mo tube, 50 kV/1000 μ A
- 60 mm² XFlash SDD
- 25 position sample changer



1 H Hydrogen																	2 He Helium				
3 Li Lithium	4 Be Beryllium															10 B Boron	11 C Carbon	12 N Nitrogen	13 O Oxygen	14 F Fluorine	16 Ne Neon
11 Na Sodium	12 Mg Magnesium															13 Al Aluminium	14 Si Silicon	15 P Phosphorus	16 S Sulphur	17 Cl Chlorine	18 Ar Argon
19 K Potassium	20 Ca Calcium	21 Sc Scandium	22 Ti Titanium	23 V Vanadium	24 Cr Chromium	25 Mn Manganese	26 Fe Iron	27 Co Cobalt	28 Ni Nickel	29 Cu Copper	30 Zn Zinc	31 Ga Gallium	32 Ge Germanium	33 As Arsenic	34 Se Selenium	35 Br Bromine	36 Kr Krypton				
37 Rb Rubidium	38 Sr Strontium	39 Y Yttrium	40 Zr Zirconium	41 Nb Niobium	42 Mo Molybdenum	43 Tc Technetium	44 Ru Ruthenium	45 Rh Rhodium	46 Pd Palladium	47 Ag Silver	48 Cd Cadmium	49 In Indium	50 Sb Antimony	51 Sn Tin	52 Te Tellurium	53 I Iodine	54 Xe Xenon				
55 Cs Cesium	56 Ba Barium	57 La Lanthanum	72 Hf Hafnium	73 Ta Tantalum	74 W Tungsten	75 Re Rhenium	76 Os Osmium	77 Ir Iridium	78 Pt Platinum	79 Au Gold	80 Hg Mercury	81 Tl Thallium	82 Pb Lead	83 Bi Bismuth	84 Po Polonium	85 At Astatine	86 Rn Radon				
87 Fr Francium	88 Ra Radium	89 Ac Actinium																			
			L Lanthanides	58 Ce Cerium	59 Pr Praseodymium	60 Nd Neodymium	61 Pm Promethium	62 Sm Samarium	63 Eu Europium	64 Gd Gadolinium	65 Tb Terbium	66 Dy Dysprosium	67 Ho Holmium	68 Er Erbium	69 Tm Thulium	70 Yb Ytterbium	71 Lu Lutetium				
			Ac Actinides	90 Th Thorium	91 Pa Protactinium	92 U Uranium	93 Np Neptunium	94 Pu Plutonium	95 Am Americium	96 Cm Curium	97 Bk Berkelium	98 Cf Californium	99 Es Einsteinium	100 Fm Fermium	101 Md Mendelevium	102 No Nobelium	103 Lr Lawrencium				

Mercury analysis Instrumentation



S4 T-STAR

- Mo- and W-tube, 50 kV/1000 μ A
- Mo-K, W-L and W-Brems excitation
- 90 position sample changer



1 H Hydrogen																	2 He Helium
3 Li Lithium	4 Be Beryllium															10 Ne Neon	
11 Na Sodium	12 Mg Magnesium															18 Ar Argon	
19 K Potassium	20 Ca Calcium	21 Sc Scandium	22 Ti Titanium	23 V Vanadium	24 Cr Chromium	25 Mn Manganese	26 Fe Iron	27 Co Cobalt	28 Ni Nickel	29 Cu Copper	30 Zn Zinc	31 Ga Gallium	32 Ge Germanium	33 As Arsenic	34 Se Selenium	35 Br Bromine	36 Kr Krypton
37 Rb Rubidium	38 Sr Strontium	39 Y Yttrium	40 Zr Zirconium	41 Nb Niobium	42 Mo Molybdenum	43 Tc Technetium	44 Ru Ruthenium	45 Rh Rhodium	46 Pd Palladium	47 Ag Silver	48 Cd Cadmium	49 In Indium	50 Sb Antimony	51 Sn Tin	52 Te Tellurium	53 I Iodine	54 Xe Xenon
55 Cs Cesium	56 Ba Barium	57 La Lanthanum	72 Hf Hafnium	73 Ta Tantalum	74 W Tungsten	75 Re Rhenium	76 Os Osmium	77 Ir Iridium	78 Pt Platinum	79 Au Gold	80 Hg Mercury	81 Tl Thallium	82 Pb Lead	83 Bi Bismuth	84 Po Polonium	85 At Astatine	86 Rn Radon
87 Fr Francium	88 Ra Radium	89 Ac Actinium															
			L Lanthanides	58 Ce Cerium	59 Pr Praseodymium	60 Nd Neodymium	61 Pm Promethium	62 Sm Samarium	63 Eu Europium	64 Gd Gadolinium	65 Tb Terbium	66 Dy Dysprosium	67 Ho Holmium	68 Er Erbium	69 Tm Thulium	70 Yb Ytterbium	71 Lu Lutetium
			Ac Actinides	90 Th Thorium	91 Pa Protactinium	92 U Uranium	93 Np Neptunium	94 Pu Plutonium	95 Am Americium	96 Cm Curium	97 Bk Berkelium	98 Cf Californium	99 Es Einsteinium	100 Fm Fermium	101 Md Mendelevium	102 No Nobelium	103 Lr Lawrencium



Mercury analysis by TXRF – Water and environment

Mercury analysis

Water and food samples



Background

Drinking water:

- Strict public limit values for various metals

Element	Limit value (µg/l)
Cu	2000
Mn	50
Fe	200
Ni	200
As	10
Cd	3
Hg	1
Pb	10
U	10

Example:
German drinking water regulation

Mercury analysis

Water and food samples




Background

High-end laboratories for drinking water

- Very high sample throughput
- Only analysis of drinking water
- In most cases application of ICP-MS

Water analysis in environmental research

- Normal sample throughput
- Water with different matrices and various content of solid matter
- Demand for the analysis of biomonitors
- Sometimes demand for mobile on-/near-site analysis
- TXRF more suitable than other methods like ICP-MS



Mercury analysis

Water



Sample preparation

- Aliquotation of 1 ml sample
- Addition of EDTA- or thiourea solution for Hg-fixation
- Addition of gallium solution for internal standardization
- Homogenization, transfer of 10 μ l onto sample carrier
- Drying
- Measurement



Mercury analysis

Water



Results

	LLD Water (µg/l)		LLD Sewage (µg/l)	
	S2 PICOFOX	S4 T-STAR	S2 PICOFOX	S4 T-STAR
Manganese	0,62	0,35	1,4	0,88
Iron	0,49	0,27	1,1	0,72
Nickel	0,33	0,18	0,73	0,49
Copper	0,29	0,17	0,63	0,44
Arsenic	0,15	0,10	0,27	0,21
Cadmium	n. d.	51	n. d.	64
Mercury	0,22	0,11	0,42	0,26
Lead	0,20	0,11	0,37	0,23
Uranium	0,34	0,28	0,34	0,30


Mercury analysis

Soil



Background

- Heavy metal distribution in soils is often very inhomogeneous
- Typical TXRF analysis of soils only analysed minute sample amounts
→ only screening analysis is possible
- For soil samples the application of handheld XRF is more suitable.
These also fulfill several international norms (e.g. US EPA Method 6200-05, EN 15309—07, ISO 13196:2013)





Mercury analysis by TXRF – Food and beverages

Mercury analysis

Non-alcoholic beverages



Iced tea - Motivation

Iced tea is a widely consumed non-alcoholic beverage and challenging for trace element analysis because of the high sugar matrix

Sample preparation

- Spiking of an iced tea samples with As, Cd, Hg and Pb
- Internal standardization with Ga and Pd
- Addition of polyvinyl alcohol solution
- Homogenization
- Preparation of 10 μ l on quartz glass sample carriers







Mercury analysis

Non-alcoholic beverages



Results – S4 T-STAR






Mercury analysis

Non-alcoholic beverages



Results – S2 PICOFOX



Mercury analysis

Non-alcoholic beverages



Results

($\mu\text{g/L}$)	LLD S2 PICOFOX	LLD S4 T-STAR
Arsenic	3,4	2,9
Cadmium	n.d.	132
Mercury	4,3	3,1
Lead	4,1	2,8

Mercury analysis

Solid food



Background

Dried food

- On a world-wide scale fish and rice are two dominant food sources
- For test purposes the certified reference materials „DORM-3“ (fish protein) and “NIST 1568a” (rice flour) were tested. Both are available as freeze-dried powders.



Mercury analysis

Solid food



Sample preparation solid food

- Grinding ($< 25 \mu\text{m}$)
- Weighing of about 100 mg into plastic tubes
- Suspension in 5 ml aqueous Triton X100 solution
- Addition of internal standard
- Homogenization, treatment in ultrasonic bath
- Preparation of 10 μl suspension on sample carrier
- Drying



Mercury analysis

Solid food



Results

	LLD S2 PICOFOX (mg/kg)		LLD S4 T-STAR (mg/kg)	
	NIST 1568a	DORM-3	NIST1568a	DORM-3
Arsenic	0,075	0,074	0,074	0,067
Cadmium	n. d.	n. d.	0,269	0,3
Mercury	0,084	0,098	0,080	0,072
Lead	0,086	0,099	0,072	0,065



Mercury analysis by TXRF – Medical analysis

Mercury analysis

Medical samples



Urine - Motivation

- Toxic metals like As and Hg are normally present in blood and serum at ppt-levels. This is too low for TXRF and is normally done by ICP-MS (TXRF is widely used in medical research for these samples for other elements like Mn, Cu, Zn and Se)
- As the kidney is the main secretion part, analysis of toxic metals in urine is of high interest
- Especially when dealing with sampling in rural regions, when intoxications must be monitored over longer time periods or samples must be stored and transported the analysis of urine samples becomes important



Mercury analysis

Medical samples



Urine - Samples

- Seronorm Urine CRM

Urine – Sample preparation

Same as for all high-matrix liquid samples

- Aliquotation
- Internal standardization
- Addition of polyvinyl alcohol solution
- Transfer of 10 μ l to quartz glass sample carrier
- Drying (here on air at room temperature)


Mercury analysis

Medical samples



Urine - Results

	LLD Seronorm Urine Level 2 (µg/l)	
	S2 PICOFOX	S4 T-STAR
Arsenic	2,9	2,8
Cadmium	n. d.	420
Mercury	3,6	3,1
Thallium	3,5	2,9
Lead	3,4	2,7
Bismuth	3,4	2,6



Mercury analysis

Medical samples



Hair - Motivation

- Hair is easier to collect, transport and store than other medical samples like blood or serum
- As elements are fixed in the hair, element control over longer time period is possible

Hair – Samples

- Certified reference sample NIES 13 (National Institute for Environmental Studies)
- 5 human hair samples (cleaned with acetone and ultrapure water)



Hair – Sample preparation by suspension

- Grinding
- Weighing of about 100 mg
- Suspension in aqueous Triton X100 solution
- Internal standardization
- Homogenization
- Transfer of 10 µl suspension to quartz glass carrier
- Drying

Hair – Sample preparation by acid digestion on hot plate


- Grinding
- Weighing of about 100 mg
- Digestion at 120 ° C for 2h in 1 ml HNO₃/H₂O₂ mixture (3:1), filled up to 5 ml with pure water
- Internal standardization
- Homogenization
- Transfer of 10 µl suspension to quartz glass carrier
- Drying

Mercury analysis Medical samples



Hair - Results

LLD	S2 PICOFOX (mg/kg)		S4 T-STAR (mg/kg)	
	Digest.	Susp.	Digest.	Susp.
Manganese	0,40	0,61	0,15	0,31
Nickel	0,20	0,27	0,059	0,12
Arsenic	0,11	0,19	0,036	0,076
Cadmium	n. d.	n. d.	4,7	12,3
Mercury	0,13	0,23	0,040	0,084
Lead	0,13	0,22	0,040	0,085





Mercury analysis by TXRF – Analysis of cosmetics

Mercury in cosmetics - Motivation

- Mercury is a common ingredient found in skin lightening soaps and creams
- It is also found in other cosmetics, such as eye makeup cleansing products and mascara
- Mercury can be found in cosmetics exists in two forms:
 - inorganic (e.g. ammoniated mercury), which is used in skin lightening soaps and creams
 - organic mercury compounds (thiomersal ethyl mercury and phenyl mercuric salts), which are used as preservatives



Mercury in cosmetics - Motivation

The allowed mercury content in cosmetics is strictly regulated in many countries

Regulatory body	Hg limits for cosmetics other than eye area products
European Union	Banned
Many African nations	Banned
United States FDA	< 1 mg/kg
Health Canada	< 3 mg/kg
Phillipines FDA	< 1 mg/kg
Regulatory body	Hg limits for eye area products
European Union	< 70 mg/kg
United States FDA	< 65 mg/kg

Sample preparation

- The complex and variable matrix of cosmetics (oil, fats, minerals and numerous organic and inorganic detergents) make the complete digestions for techniques like ICP-MS extremely difficult
- For TXRF a rather simple suspension preparation can be applied for all types of cosmetics

Measurements

- Measurements were performed with a S2 PICOFOX, applying Mo-K excitation and measurement times of 1000 s

Mercury analysis Cosmetics



Results

Table 3 Comparison of determination results (mg/kg) of cosmetics samples by different methods

样品 Samples	Pb		As		Hg	
	本方法 TXRF	其它方法 Other method	本方法 TXRF	其它方法 Other method	本方法 TXRF	其它方法 Other method
GBW09305	39.4	37.2 (Reference value)	9.1	9.9 (Reference value)	-	-
GBW09306	-	-	-	-	0.64	0.67 (Reference value)
GBW09307	-	-	-	-	1.33	1.35 (Reference value)
洁面膏 Face wash gel	1.25	1.20 (ICP-MS)	-	-	1.20	1.16 (ICP-MS)
护臀膏 Bottom Butte	-	-	1.30	1.24 (ICP-MS)	-	-
日霜 Day cream *	0.62	0.60 (Added value)	0.58	0.60 (Added value)	0.57	0.60 (Added value)
美白霜 Whitening cream	6.71	6.84 (ICP-MS)	-	-	-	-

* : 空白样品加标 (Spiked blank sample)。

Wang, Q. et al. (2018): „Determination of Lead, Arsenic and Mercury in Cream Cosmetics by Total Reflection X-Ray Fluorescence Spectrometry Using Suspension Sampling“, Chinese Journal of Analytical Chemistry, 4, 517-523

Mercury analysis Cosmetics



Summary

- TXRF is suitable for the analysis and control of Hg (and other toxic metals) in all types of cosmetics
- The applied sample preparation is fast, easy and does not need any toxic or dangerous chemicals

Outlook

- Detailed work on the analysis of toxic metals in cosmetics (optimisation of sample preparation procedures, development of SOPs, publications) in cooperation with external partners are planned for Q1 2021



Mercury analysis by TXRF - Summary


Mercury analysis

Summary



Environmental analysis

- TXRF is suitable for the analysis of Hg and other metals in water samples from fresh water to highly contaminated sewages
- All different kinds of biomonitors (plants, algae, biofilms etc.) can be analysed for environmental research and control



Food and beverages

- Beverages can easily be analysed into the sub mg/l range for Hg and other toxic elements like As and Pb
- Fast screening for toxic elements in solid foodstuff is possible after an easy sample preparation, avoiding complex digestion procedures



Mercury analysis

Summary




Medical research

- Hg and other metals can easily be analysed in urine and hair samples applying fast and easy preparation techniques

Cosmetics

- Hg and other metals like As and Pb can be analysed in various types of cosmetics for control of regulatory limits avoiding time consuming and expensive preparation techniques








Comparison with Atomic Spectroscopy

Feature comparison TXRF versus ICP




Topic	Feature	Benefit	
Suspensions, solid samples	direct analysis without sample digestion possible	saves time and avoids use of hazardous chemicals	
Solid, valuable samples	non-destructive analysis possible (semi-quantitative)	samples are not destroyed and can be used for other purposes	
Linear range	high linear quantify-cation range from sub-ppb to low %	avoids multiple calibrations or measurements; full results after one run	

Feature comparison TXRF versus ICP






Topic	Feature	Benefit
Media, disposables	no consumption of media (e.g. carrier gases, chemicals, standards), no need for periodic replacement of parts (lamps, nebulizer)	cost saving; no special infrastructure (gas lines) required
Quantification	by internal standardization, instrument calibrated ex works	simple; time saving, no daily calibration required
Memory effects	No memory effects	ICP: intense purging of the system, risk of carryover from previous samples



Feature comparison TXRF versus ICP



Topic	Feature	Benefit	
System maintenance and cleaning	QC system check automated, no cleaning of specific parts required	ICP requires daily/weekly system checks	 A close-up photograph of an ICP (Inductively Coupled Plasma) sample introduction system. It shows a glass nebulizer chamber with a yellow band around its neck, connected to a glass tube. The chamber contains a white, foamy sample. Handwritten text on the chamber reads "Cleared KIC 15-07-09".
Ease of operation	simple, suitable for academic training	technicians are able to use TXRF; ICP typically used by (scientific) experts	 A photograph of two scientists in a laboratory. A woman in the foreground is wearing safety glasses and a white lab coat, looking intently at a small plant in a container. A man in the background is also in a lab coat and is looking at a computer monitor. There are various lab equipment like test tubes and beakers visible.
Footprint	compact benchtop instrument, foot print 693 x 528 mm	saves lab space	 A photograph of a lab technician in a white lab coat and safety glasses, wearing blue gloves. She is working at a lab bench, handling a small vial. The bench is cluttered with various lab supplies, including bottles, pipettes, and containers. Shelves with more lab equipment are visible in the background.

Cost comparison TXRF versus ICP-MS



Total (5 years)	S4 T-STAR	ICP-MS	Remarks
Installation	136.900 €	188.100 €	incl. peripheral devices
Operation costs	64.600 €	179.000 €	disposables, gas, media, standards etc.
Man hours	82.000 €	134.000 €	100 k€/a, 220 working days, 20 samples/d
Total	283.500 €	501.100 €	
Costs / sample	12,89 €	22,78 €	

Sources: Automotive study 2015
 EPA study
 ICP-OES cost calculator Spectro
 Discussion forums
 Bruker data

More detailed PPT and Excel sheet available

Technology comparison

Day-to-day checks

Technology comparison

Daily checks



ICP	TXRF
Setup before analysis	
Inspect torch for injector blockage	Check system status
Check nebulizer for blockage	Select excitation
Inspect peristaltic pump tubing	
Check exhaust system operating	
Rinse between samples (ICP-MS)	
After analysis / end of day	
Aspirate rinse solution	
Release pressure and detach pump	
Empty waste vessel	
Wipe down exterior surfaces	Wipe off exterior surfaces
Leave system in stand-by mode	Automatic system stand-by

Technology comparison

Weekly/monthly checks



ICP	TXRF
Clean torch and test plasma ignites after reconnection	Sensitivity test (automatic!)
Clean nebulizer frequently (ICP-MS)	Start auto beam adjustment, when sensitivity is decreasing
Check sample introduction tubing, O-rings	
Inspect cone or snout and clean, if needed; check vacuum after reconnection	
Check/replace graphite gasket (ICP-MS)	
Inspect torch bonnet (radial ICP, -MS)	
Clean spray chamber	Check auto beam adjustment
Check/clean nebulizer (ICP-OES)	
Inspect induction coil	
Check/clean air filter for cooling air	
Check/clean water level and air filter on water chiller	


Any Questions?

Please **type in** the questions you may have for our speakers in the **Questions form** and click **Submit**



Thank you for your attention!

www.bruker.com
www.s2picofox.com
www.s4tstar.com



www.bruker.com