

## Nanomechanical Testing Journal Club

*Keep apprised of recent research*

### BEoL Stack Failure Mode Evaluation Utilizing Micromechanical Testing and FEM Modelling

by Jendrik Silomon, André Clausner, and Ehrenfried Zschech

*2022 IEEE International Symposium on the Physical and Failure Analysis of Integrated Circuits (IPFA), Singapore, Singapore, 2022, pp. 1-8*

This article is paywalled: <https://ieeexplore.ieee.org/document/9915739>

Hello researcher,

Thank you for reading this month's edition of Bruker's [Nanomechanical Testing Journal Club](#).

This month I'd like to highlight a [paper](#) by André Clausner and co-authors that focuses on the **methodology of inducing stresses (and eventually failure) in microelectronic devices to understand failure modes, damage initiation, and damage propagation**.

This work summarizes and puts into perspective some of the authors' earlier, experimental work on both shearing and tensile loading of copper pillars that are bumped onto aluminum pads. The device stack also includes a polyimide underfill and multiple metallization layers and low-*k* material. I would encourage any readers interested in microelectronic device failure to look back at some of these earlier works.

In this article, finite element modeling (FEM) is employed to understand stresses within the stack. These FEM results are analyzed in the context of micro-mechanical loading results from Bruker's [Hysitron TI 950 TriboIndenter](#) equipped with a [3D OmniProbe](#) high-bandwidth transducer.

For shear tests, a flat-sided probe is brought into contact with the side of a Cu bump. If the probe is lifted, a longer lever arm is created at the Cu-Al pad interface. This long lever arm results in crack initiation at this interface, as determined by FIB slices for both failed and sub-critical loading. If the lever arm is shortened – both the Cu-Al interface crack and low-*k* – metallization delamination occurs. This failure mode, described in this work as cratering, results in these two cracks growing to meet and causing failure well into the 4<sup>th</sup> and 5<sup>th</sup> metallization layers. Layers M1-M3 delaminate and are ejected from the damage region.

Since both long and short lever arm modes tend to cause damage initiation at the Cu-Al interface, the authors use the FEM to look at the possibility of changing the Al pad size. They do find that the Al pad should be larger than the Cu pillar to increase the

load at failure, but further work is required to optimize the structure to prevent damage to the device.

I encourage you to spend some time with this paper, as this is a great example of **using experimental and modeling results to make engineering decisions on real products**.



### **Douglas Stauffer, Ph.D.**

Sr. Manager of NI Applications Development  
Bruker Nano Surfaces and Metrology

[Douglas.Stauffer@bruker.com](mailto:Douglas.Stauffer@bruker.com)


## **Bruker Technology Featured in this Article**



### **Hysitron TI 980 TribolIndenter**

The world's most powerful nanomechanical and nanotribological test system.

[Learn More](#)



### **3D OmniProbe**

Bringing nanoindentation and microindentation testing techniques into the world of 3D characterization.

[Learn More](#)

## **Engage with the Nanomechanical Testing Community**



### **Revisit Last Month's Journal Club**

[Pop-In Identification in Nanoindentation Curves with Deep Learning Algorithms](#)



### **Attend Nanobrücken 2023**

Register for Bruker's annual nanomechanical testing conference. May 23-25 in Saarbrücken, Germany.

[Learn More and Register](#)



### **Contact a Bruker Nanomechanical Testing Expert**

Get in touch with an expert to discuss your challenges and schedule a personalized virtual demo.

[Contact Expert](#)

**Bruker Nano Surfaces and Metrology**  
Hysitron Nanomechanical Test Instrumentation  
9625 West 76th Street  
Minneapolis, MN 55344 USA  
+1 952-835-6366  
[productinfo@bruker.com](mailto:productinfo@bruker.com)

© 2023 Bruker Corporation