

CCS-aware wide-scope target screening utilizing LC-TIMS-HRMS and a new heated ESI source – The answer to environmental and human biomonitoring challenges

Konstantina S. Diamanti¹; Reza Aalizadeh¹; Dimitrios E. Damalas¹; Georgios O. Gkotsis¹; Bob Galvin²; Artem Filipenko³; Carsten Baessmann⁴; Nikolaos S. Thomaidis¹

¹National and Kapodistrian University of Athens, Athens, Greece; ²Bruker UK Ltd., Coventry, United Kingdom; ³Bruker Daltonics, Billerica, MA; ⁴Bruker Daltonics GmbH & Co. KG, Bremen, Germany

Introduction

Thousands of chemicals are dispersed throughout the environment and may threaten ecosystems and health. Through environmental and human biomonitoring, the exposure to persistent, bioaccumulative and toxic substances can be assessed. Nowadays, various chemical classes can be reliably identified and quantified, thanks to target screening using high-resolution mass spectrometry coupled to chromatographic techniques. However, the sample complexity often results in high matrix effects and hampers the detection of compounds with low concentration (suppressing ionization efficiency of compounds below signal-to-noise ratio). Reliable detections are also prevented by the complexity of MS/MS data produced by data-independent acquisition (DIA) modes. Therefore, apart from the established target screening workflow (retention time, MS1 and MS2), a new technology is required to identify xenobiotics with high confidence.

Methods

Reference standard curves & standard addition curves at a wide concentration range

Spiked samples: biota (bird eggs) & human biospecimens (urine)



Sample preparation: Generic protocols for simultaneous extraction of semi-polar to polar xenobiotics (e.g. pesticides and drugs)

Instrumentation: LC-TIMS-QTOFMS with optimized broad mass and mobility bbCID mode (DIA)



Data treatment: Targeted analysis using TargetScreener 4D platform for the screening of xenobiotics (precursor ion formula, RT, MS1 and MS2 qualifier ions, ion mobility-derived CCS values)



Results




Figure 1: Distribution of 198 detected spiked compounds into groups based on the ratio of the peak area using VIP-HESI source to the peak area using ESI source at $C_{\text{instrumental}} = 5 \mu\text{g/L}$

Ionization efficiency: ESI vs VIP-HESI




Figure 2: Presentation of a representative set of 25 compounds spiked at $C_{\text{instrumental}} = 25 \mu\text{g/L}$ that cover a wide range of the following chemical properties: (1) m/z , (2) logKow, (3) RT, (4) logIE (developed for ESI), (5) CCS




Based on the results presented in Figure 1:

- 50% of the spiked compounds have >2 times higher response using VIP-HESI source in comparison with ESI.
- 35% of the spiked compounds have no significant differences in response with both sources (ratio 1-2).
- 15% of the spiked compounds show decreased response or are not detected using VIP-HESI source in comparison with ESI.

By utilizing VIP-HESI source, lower method LODs can be achieved for half of the compounds.

Incorporation of TIMS in LC-VIP-HESI-QTOF MS


Example of ethopabate (veterinary drug) spiked in an egg sample at $C_{\text{instrumental}} = 10 \mu\text{g/L}$




Data processing with ion mobility filtering
The background signal caused by the presence of co-eluting isobaric analytes and matrix components is deconvoluted.
Full-scan and bbCID EICs of even higher quality

Data processing with timsOFF vs timsON


Example on the identification of co-eluting isomeric compounds in a spiked urine sample



Example on the effect of ion mobility dimension in target screening using a spiked egg sample



Overall results of the detected compounds spiked at a urine sample at $C_{\text{instrumental}} = 5 \mu\text{g/L}$



Conclusions

- VIP-HESI source offers enhanced sensitivity for a great number of compounds, achieving lower LODs in many cases.
- The addition of TIMS dimension in LC-HRMS improves signal-to-noise ratio.
- Ion mobility filtering during data processing provides higher quality chromatograms and full-scan MS and bbCID MS/MS spectra.
- The overall identification confidence is increased by adding CCS values.

More reliable identification of xenobiotics and lower method LODs are achieved in environmental and human biomonitoring studies.

Ion Mobility: Applications I