


# Do it yourself - Metal analysis in biological and medical samples



Bruker Nano Analytics, Berlin  
Webinar, February 27<sup>th</sup>, 2019



# Welcome



## Speakers

Dr. Hagen Stosnach  
Applications Scientist TXRF  
Berlin, Germany



Dr. Po-Wah So  
Senior lecturer in biomedical  
imaging and spectroscopy  
King's College London



Dr. Armin Gross  
Global Product Manager TXRF  
Berlin, Germany



# Itinerary



## Part I: Introduction and background

- Motivation
- TXRF spectroscopy

## Part II: Po-Wah So

- TXRF-based element profiling in Alzheimer's Disease

## Part III: Protein analysis

## Part IV: Summary/Q & A



# Part I: Introduction and background

# Introduction

## Motivation



### Metals play a crucial role for cellular and subcellular functions

- The divalent cations  $Zn^{2+}$  ,  $Ca^{2+}$  and  $Mg^{2+}$  prevent cytotoxicity and in vivo antagonize Cd-induced carcinogenesis.
- Lack of body iron is common in cancer patients and it is associated with complications in surgery and in animal experiments. The transport of iron and other metal ions by the blood plasma is achieved through the formation of protein complexes.
- Copper is recognized as an essential element and is primarily associated with copper-dependent cellular enzymes triggering redox reactions or for the stabilization of protein particulates.
- In addition it is well known that a number of toxic metals can inhibit enzyme functions, which leads to serious health issues.

# Introduction

## Motivation



### Common methods for metal analysis


- Atomic absorption spectroscopy (AAS)
  - Inductively-coupled plasma optical emission spectroscopy (ICP-OES)
  - Inductively-coupled plasma mass spectrometry (ICP-MS)
    - High demands on laboratory environment and staff
    - Necessary sample amounts are often too high
- ⇒ Samples are sent to central or service labs for analysis
- ⇒ restricted numbers of analyzed samples
  - ⇒ high costs
  - ⇒ long waiting time for samples

# Introduction

## TXRF spectroscopy



### Total reflection X-ray fluorescence spectroscopy



Beam angle:  $0^\circ$  /  $90^\circ$

- Samples must be prepared on a reflective media
- Polished quartz glass or polyacrylic glass disc
- Dried to a thin layer, or as a thin film or microparticle


# Introduction

## TXRF spectroscopy



### Sample preparation of minute liquid samples

- Aliquotation of sample (10  $\mu$ l – 10 ml)
- Addition of internal standard (mono-element solution)
- Homogenization (mixing)
- Transfer of 10  $\mu$ l sample to sample carrier
- Drying






# Introduction

## TXRF spectroscopy



### Sample preparation of minute liquid samples ( $< 10 \mu\text{l}$ )

- Preparation of sample droplet on sample carrier
- Direct addition of internal standard onto sample droplet
- Drying
- ⇒ Results are not as reproducible as for external mixing but in most cases still sufficient
- ⇒ Also possible: Weighing of microsamples directly on carrier and subsequent internal standardization and ashing or digestion




# Introduction

## TXRF spectroscopy



### Sample preparation of solid or very viscous samples

- Weighing of sample (amount depending on balance working range)
- Suspension in water or alcohol (ultrasonic bath optional)
- Addition of internal standard
- Homogenization
- Transfer of 10  $\mu\text{l}$  sample to sample carrier
- Drying



# Introduction

## TXRF spectroscopy




### S2 PICOFOX

- Mo-tube 50 kV/1000  $\mu$ A
- 60 mm<sup>2</sup> XFlash SDD
- 25 position sample changer



### S4 T-STAR

- Mo tube, 50 kV/1000  $\mu$ A
- W-tube, 50 kV/1000  $\mu$ A
- Monochromator system for Mo-K, W-L and W-Brems monochromatisation
- 60 mm<sup>2</sup> XFlash SDD
- 90 position sample changer



# Introduction TXRF spectroscopy



## Element range S2 PICOFOX

1 <b>H</b> Hydrogen																	2 <b>He</b> Helium						
3 <b>Li</b> Lithium	4 <b>Be</b> Beryllium																	5 <b>B</b> Boron	6 <b>C</b> Carbon	7 <b>N</b> Nitrogen	8 <b>O</b> Oxygen	9 <b>F</b> Fluorine	10 <b>Ne</b> Neon
11 <b>Na</b> Sodium	12 <b>Mg</b> Magnesium																	13 <b>Al</b> Aluminium	14 <b>Si</b> Silicon	15 <b>P</b> Phosphorus	16 <b>S</b> Sulphur	17 <b>Cl</b> Chlorine	18 <b>Ar</b> Argon
19 <b>K</b> Potassium	20 <b>Ca</b> Calcium	21 <b>Sc</b> Scandium	22 <b>Ti</b> Titanium	23 <b>V</b> Vanadium	24 <b>Cr</b> Chromium	25 <b>Mn</b> Manganese	26 <b>Fe</b> Iron	27 <b>Co</b> Cobalt	28 <b>Ni</b> Nickel	29 <b>Cu</b> Copper	30 <b>Zn</b> Zinc	31 <b>Ga</b> Gallium	32 <b>Ge</b> Germanium	33 <b>As</b> Arsenic	34 <b>Se</b> Selenium	35 <b>Br</b> Bromine	36 <b>Kr</b> Krypton						
37 <b>Rb</b> Rubidium	38 <b>Sr</b> Strontium	39 <b>Y</b> Yttrium	40 <b>Zr</b> Zirconium	41 <b>Nb</b> Niobium	42 <b>Mo</b> Molybdenum	43 <b>Tc</b> Technetium	44 <b>Ru</b> Ruthenium	45 <b>Rh</b> Rhodium	46 <b>Pd</b> Palladium	47 <b>Ag</b> Silver	48 <b>Cd</b> Cadmium	49 <b>In</b> Indium	50 <b>Sb</b> Antimony	51 <b>Sn</b> Tin	52 <b>Te</b> Tellurium	53 <b>I</b> Iodine	54 <b>Xe</b> Xenon						
55 <b>Cs</b> Cesium	56 <b>Ba</b> Barium	57 <b>La</b> Lanthanum	72 <b>Hf</b> Hafnium	73 <b>Ta</b> Tantalum	74 <b>W</b> Tungsten	75 <b>Re</b> Rhenium	76 <b>Os</b> Osmium	77 <b>Ir</b> Iridium	78 <b>Pt</b> Platinum	79 <b>Au</b> Gold	80 <b>Hg</b> Mercury	81 <b>Tl</b> Thallium	82 <b>Pb</b> Lead	83 <b>Bi</b> Bismuth	84 <b>Po</b> Polonium	85 <b>At</b> Astatine	86 <b>Rn</b> Radon						
87 <b>Fr</b> Francium	88 <b>Ra</b> Radium	89 <b>Ac</b> Actinium																					
			L Lanthanides	58 <b>Ce</b> Cerium	59 <b>Pr</b> Praseodymium	60 <b>Nd</b> Neodymium	61 <b>Pm</b> Promethium	62 <b>Sm</b> Samarium	63 <b>Eu</b> Europium	64 <b>Gd</b> Gadolinium	65 <b>Tb</b> Terbium	66 <b>Dy</b> Dysprosium	67 <b>Ho</b> Holmium	68 <b>Er</b> Erbium	69 <b>Th</b> Thulium	70 <b>Yb</b> Ytterbium	71 <b>Lu</b> Lutetium						
			Ac Actinides	90 <b>Th</b> Thorium	91 <b>Pa</b> Protactinium	92 <b>U</b> Uranium	93 <b>Np</b> Neptunium	94 <b>Pu</b> Plutonium	95 <b>Am</b> Americium	96 <b>Cm</b> Curium	97 <b>Bk</b> Berkelium	98 <b>Cf</b> Californium	99 <b>Es</b> Einsteinium	100 <b>Fm</b> Fermium	101 <b>Md</b> Mendelevium	102 <b>No</b> Nobelium	103 <b>Lr</b> Lawrencium						

Impossible to analyse

Analysed using K-lines

Difficult to analyse

Analysed using L-lines

# Introduction TXRF spectroscopy



## Element range S4 T-STAR

1 <b>H</b> Hydrogen																	2 <b>He</b> Helium
3 <b>Li</b> Lithium	4 <b>Be</b> Beryllium											5 <b>B</b> Boron	6 <b>C</b> Carbon	7 <b>N</b> Nitrogen	8 <b>O</b> Oxygen	9 <b>F</b> Fluorine	10 <b>Ne</b> Neon
11 <b>Na</b> Sodium	12 <b>Mg</b> Magnesium											13 <b>Al</b> Aluminium	14 <b>Si</b> Silicon	15 <b>P</b> Phosphorus	16 <b>S</b> Sulphur	17 <b>Cl</b> Chlorine	18 <b>Ar</b> Argon
19 <b>K</b> Potassium	20 <b>Ca</b> Calcium	21 <b>Sc</b> Scandium	22 <b>Ti</b> Titanium	23 <b>V</b> Vanadium	24 <b>Cr</b> Chromium	25 <b>Mn</b> Manganese	26 <b>Fe</b> Iron	27 <b>Co</b> Cobalt	28 <b>Ni</b> Nickel	29 <b>Cu</b> Copper	30 <b>Zn</b> Zinc	31 <b>Ga</b> Gallium	32 <b>Ge</b> Germanium	33 <b>As</b> Arsenic	34 <b>Se</b> Selenium	35 <b>Br</b> Bromine	36 <b>Kr</b> Krypton
37 <b>Rb</b> Rubidium	38 <b>Sr</b> Strontium	39 <b>Y</b> Yttrium	40 <b>Zr</b> Zirconium	41 <b>Nb</b> Niobium	42 <b>Mo</b> Molybdenum	43 <b>Tc</b> Technetium	44 <b>Ru</b> Ruthenium	45 <b>Rh</b> Rhodium	46 <b>Pd</b> Palladium	47 <b>Ag</b> Silver	48 <b>Cd</b> Cadmium	49 <b>In</b> Indium	50 <b>Sb</b> Antimony	51 <b>Sn</b> Tin	52 <b>Te</b> Tellurium	53 <b>I</b> Iodine	54 <b>Xe</b> Xenon
55 <b>Cs</b> Cesium	56 <b>Ba</b> Barium	57 <b>La</b> Lanthanum	72 <b>Hf</b> Hafnium	73 <b>Ta</b> Tantalum	74 <b>W</b> Tungsten	75 <b>Re</b> Rhenium	76 <b>Os</b> Osmium	77 <b>Ir</b> Iridium	78 <b>Pt</b> Platinum	79 <b>Au</b> Gold	80 <b>Hg</b> Mercury	81 <b>Tl</b> Thallium	82 <b>Pb</b> Lead	83 <b>Bi</b> Bismuth	84 <b>Po</b> Polonium	85 <b>At</b> Astatine	86 <b>Rn</b> Radon
87 <b>Fr</b> Francium	88 <b>Ra</b> Radium	89 <b>Ac</b> Actinium															
			L Lanthanides	58 <b>Ce</b> Cerium	59 <b>Pr</b> Praseodymium	60 <b>Nd</b> Neodymium	61 <b>Pm</b> Promethium	62 <b>Sm</b> Samarium	63 <b>Eu</b> Europium	64 <b>Gd</b> Gadolinium	65 <b>Tb</b> Terbium	66 <b>Dy</b> Dysprosium	67 <b>Ho</b> Holmium	68 <b>Er</b> Erbium	69 <b>Th</b> Thulium	70 <b>Yb</b> Ytterbium	71 <b>Lu</b> Lutetium
			Ac Actinides	90 <b>Th</b> Thorium	91 <b>Pa</b> Protactinium	92 <b>U</b> Uranium	93 <b>Np</b> Neptunium	94 <b>Pu</b> Plutonium	95 <b>Am</b> Americium	96 <b>Cm</b> Curium	97 <b>Bk</b> Berkelium	98 <b>Cf</b> Californium	99 <b>Es</b> Einsteinium	100 <b>Fm</b> Fermium	101 <b>Md</b> Mendelevium	102 <b>No</b> Nobelium	103 <b>Lr</b> Lawrencium

Impossible to analyse

Analysed using K-lines (Mo-K excitation)

Analysed using K-lines (W-Brems excitation)

Difficult to analyse


Analysed using L-lines (Mo-K excitation)

# Introduction

## TXRF spectroscopy



### Sensitivity S2 PICOFOX vs. S4 T-STAR





Part II: Dr. Po-Wah So

TXRF elemental profiling in Alzheimer's disease



# TXRF-based Elemental Profiling in Alzheimer's Disease

**Dr Po-Wah So**

*Senior Lecturer in Biomedical Imaging and Spectroscopy*

**Head, Phenomics Lab**

Po-wah.so@kcl.ac.uk  @powahso

*Department of Neuroimaging,*

*Maurice Wohl Clinical Neuroscience Institute,*


*Institute of Psychiatry, Psychology and Neuroscience*






# Causes of Alzheimer's Disease

Alzheimer's Brain Cells



Amyloid and Tau


Metal Dyshomeostasis



Braidy et al., 2014. Front. Aging Neurosci., <https://doi.org/10.3389/fnagi.2014.00138>

## ***Aim***

**Discriminate  
between control and  
Alzheimer's Disease  
by plasma TXRF  
elemental profiling**



**Blood Metals/Elements in Alzheimer's Disease**


**No consensus**



**Traditional elemental/metal analysis methods**


matrix effects


limited elements measured concomitantly

# Typical TXRF Spectrum of Human Plasma












# Percentages of Elements in Plasma




Control


Alzheimer's Disease



Significant difference between control and Alzheimer's Disease




**Percentage  
s of  
Elements in  
Plasma**




**Control**

**Alzheimer's  
Disease**



Significant difference between control and Alzheimer's Disease





Absolute levels


No difference



Alzheimer's Disease




Alzheimer's Disease



# Iron Deficiency In Alzheimer's Disease

**Fe**



**Calcium**

**Iron**

**Zinc**

**Phosphorous**

Elemental profiling discriminates between control and AD

Crucial to measure panel of elements

**Conclusions**

*Ashraf, Stosnach, Parkes, Hye, Lovestone and So, to be published in Scientific Reports*

# Acknowledgements

Hagen Stosnach


Azhaar Ashraf

Harry G. Parkes,

Abdul Hye,

John Powell

Simon Lovestone



*Samples from AddNeuromed Consortium*



## Part III: Protein analysis

# Protein analysis


## Introduction



### Metal ions and enzymes

- Metal ions are important for the biological function of enzymes
- Various modes of metal-protein interaction: metal-, ligand-, enzyme-bridge complexes
- Metals serve as electron donors or acceptors, Lewis acids or structural regulators

(Riordan JF.: "The role of metals in enzyme activity.", Ann Clin Lab Sci. 1977 Mar-Apr;7(2):119-29



C. Arnaud: CEN, January 7, 2008  
Volume 86, Number 1, p. 8

# Protein analysis

## Introduction



### Examples of metal ions in enzymes

Role	Metals	Protein
Oxygen transport and storage	Fe, Cu	Haemoglobin
Electron transport	Fe, Cu	Cytochromes
Nitrogen fixation	Fe, Mo, V	Nitrogenase
Oxygen atom transfer	Mo, W	Oxidases, Reductases
Alkyl group transfer	Co	Vitamin B <sub>12</sub>
Hydrolysis	Zn, Cu, Mn	Hydroxylases, Peptidases
Storage and transport	Fe, Cu, Zn	Ferritin, Metallothioneins

# Protein analysis

## Feasibility tests



### Samples

Certified reference materials

- "BCR 273 – Single cell protein"
- "BCR 274 – Single cell protein"

500 mg powderous sample were suspended in 25 ml pure water





# Protein analysis

## Feasibility tests




### Sample preparation

First approach:  
mixing of microsamples in **vials**

- 100  $\mu$ l sample + internal standard  
+ 10  $\mu$ l Sc (10 mg/l)  
+ 10  $\mu$ l Ga (10 mg/l)



- Preparation of
  - 10  $\mu$ l
  - 5  $\mu$ l
  - 3  $\mu$ l
  - 1  $\mu$ l
  - 0,5  $\mu$ l



# Protein analysis


## Feasibility tests



### Sample preparation

Second approach:  
mixing of microsamples on **discs**

- 5  $\mu$ l sample + 5  $\mu$ l Sc/Ga solution (20 mg/l)
- 3  $\mu$ l sample + 3  $\mu$ l Sc/Ga solution (20 mg/l)
- 1  $\mu$ l sample + 1  $\mu$ l Sc/Ga solution (20 mg/l)



# Protein analysis

## Feasibility tests



### S4 T•STAR

- Mo tube, 50 kV/1000  $\mu$ A
- W-tube, 50 kV/1000  $\mu$ A
- 60 mm<sup>2</sup> XFlash SDD
- 90 position sample changer
- Mo-K excitation, 1000 s
- W-L excitation, 1000 s
- W-Brems, 1000 s





# Protein analysis

## Feasibility tests



- Accuracy and reproducibility for the light elements P (S and Cl) fluctuates with deposited sample amount  $\Rightarrow$  effect of sample height and size
- Other major and trace elements (K, Ca, Mn, Fe, Cu, Zn and Se) are not affected
- Reproducibility of samples, mixed in vessels is better than for those, directly mixed on the carrier
- The detection limits are directly related to the deposited amount




# Protein analysis

## Application example



### Background

- Ribonucleotide reductase (RNS) catalyses the only known de novo pathway for the production of all four deoxyribonucleotides that are required for DNA synthesis
- Fe is required for function in the aerobic, class I RNR in all eukaryotes and many bacteria
- A dinuclear metal site has been viewed as necessary to generate and stabilize the catalytic radical, essential for RNR activity
- This study<sup>\*)</sup> describes a group of RNR proteins that possess a metal independent stable radical



<sup>\*)</sup>: Srinivas, V. et al. (2018): „Metal-free ribonucleotide reduction powered by a DOPA radical in Mycoplasma pathogens“, Nature, Vol. 563, 416 - 420

# Protein analysis

## Application example



### Samples

- Proteins (*MfR2* active, *MfNrdI*, *MfR1*, *EcR2a*) and buffer solutions
- Control sample: metal-free *E. coli* class 1a R2 protein, Fe-reconstituted by incubation and subsequent removal of all unbound Fe

### Preparation

- Addition of a gallium standard solution (2 mg/l) with a volume ratio of 1:1 with a few  $\mu$ l sample
- Preparation on quartz glass carriers

### Measurements

- Threefold preparation and measurement with a S2 PICOFOX system for 1000 s


# Protein analysis

## Application example



### Results

- None of the *MfRNR* proteins contains a substantial amount of metal
- The *EcR2a* contains in the order of two metal ions per monomer, also after the desalting step




# Protein analysis

## Application example



### Conclusions

- The cumulative amount of transition metals is less than 4 % per protein. It is therefore not possible that a metal ion is required to stabilize the observed radical species
- Comprehensive structural, EPR, UV-vis absorption, TXRF and mass spectrometric data support the hypothesis that the novel radical species is metal-independent





## TXRF offers an ideal analytical solution for elemental analysis in medical and biological research



- Analysis of small sample amounts in the low  $\mu\text{l}$ -range
- Simultaneous analysis of main- and trace elements
- Simultaneous analysis of other important samples types like buffer solutions
- Instruments can be operated in normal laboratory environments (small footprint, no external gases or cooling water necessary)
- Moderate analytical demands on laboratory staff
- Low analytical and lifetime costs


# Q & A



## Any Questions?

Please **type in** the questions you may have for our speakers in the **Questions Box** and click **Submit**





Innovation with Integrity

Thank you for your attention!

[www.bruker.com](http://www.bruker.com) / [www.s4tstar.com](http://www.s4tstar.com)