

Compositional analysis of archaeological ceramics by Handheld XRF



Examples of pottery sherds showing the complexity of target ceramic materials



Archaeological ceramics in situ




BUFFALO STATE
The State University of New York

Prof. Aaron Shugar
Art Conservation Department,
SUNY Buffalo State College



Dr. Nigel Kelly
Senior Market Application Scientist,
Bruker Nano Analytics


Overview of the Talk



- Basics of X-ray Fluorescence (a review)
- Analysis of Ceramics
 - Complexity for analysis
 - Challenges for analysis – complexity and variability
 - Approaches to achieving high quality data
 - Examples of the issues
- Quantification
 - Types of quantification
 - Basics on developing a calibration and choosing standards
 - Examples of some attempts
- Suggestions for obtaining quality data


THE ELECTROMAGNETIC SPECTRUM

How does light affect molecules and atoms?




The X ray Reactions in the Sample

- **Compton:** inelastic scattering
- **Raleigh:** elastic scattering
- **Bragg:** coherent scattering
- **Electron ejection:** X ray fluorescence
- **Inter element effects:**
 - Secondary fluorescence
 - Sample absorption
 - Density effects
- Reaction probability as a function of incoming x ray energy
- **Attenuation** as function of x ray energy and material type




- Since the *electronic energy levels for each element are different*, the energy of X-ray fluorescence peak can be correlated to a specific element


Electrons, the Photoelectric Effect and Fluorescence




Electrons, the Photoelectric Effect and Characteristic Fluorescence



Excitation / Emission of Characteristic X-ray Radiation



The Photoelectric Effect and Absorption Edges



- Sudden jumps in absorption (edges) occur where the energy of the incident X-rays exceeds the binding energies of a particular electron state.
- The energy of a particular transition is always less than the energy of edge corresponding edge.
- Absorption is additive –photoelectric absorption from earlier edges is superimposed and hence the fluorescence spectrum will contain emission from all edges up to the energy of the incident radiation.

Electronic Transitions – Notation (Example: Barium)

Line	Siegbahn	Electron Binding Energies (keV)	Energy of emitted radiation (keV)	Wavelength (Å)
K-L ₂	Kα ₂	37.441, 5.634	31.807	0.390
K-L ₃	Kα ₁	37.441, 5.247	32.194	0.385
K-M ₂	Kβ ₃	37.441, 1.137	36.304	0.342
K-M ₃	Kβ ₁	37.441, 1.062	36.379	0.341
K-N ₂	Kβ' ₅	37.441, 0.192	37.249	0.333
K-N ₃	Kβ'' ₂	37.441, 0.180	37.261	0.333
L ₃ -M ₁	Lα ₁	5.247, 1.293	3.954	3.136
L ₃ -M ₄	Lα ₂	5.247, 0.796	4.451	2.786


$$E = \frac{h \times c}{\lambda} = E_{\text{keV}} = \frac{12.4}{\lambda(\text{Å})} = \frac{1.24}{\lambda(\text{nm})}$$

or

$$\lambda(\text{Å}) = \frac{12.4}{E_{\text{keV}}} \quad \lambda(\text{nm}) = \frac{1.24}{E_{\text{keV}}}$$

$$\text{Energy (keV)} = 12.4 \div \text{wavelength (Å)}$$

XRF Spectrum of Pure (98+%) Barium



X-RAY DATA BOOKLET

Center for X-ray Optics and Advanced Light Source

Lawrence Berkeley National Laboratory


- [Introduction](#)
- [X-Ray Properties of Elements](#)
- [Electron Binding Energies](#)
- [X-Ray Energy Emission Energies](#)
- [Fluorescence Yields for K and L Shells](#)
- [Principal Auger Electron Energies](#)
- [Subshell Photoionization Cross-Sections](#)
- [Mass Absorption Coefficients](#)
- [Atomic Scattering Factors](#)
- [Energy Levels of Few Electron Ions](#)
- [Periodic Table of X-Ray Properties](#)
- [Synchrotron Radiation](#)
- [Characteristics of Synchrotron Radiation](#)
- [History of X-rays and Synchrotron Radiation](#)
- [Synchrotron Facilities](#)
- [Scattering Processes](#)

UPDATE October 2009


A new version of the booklet is now available!
Please visit the [CXRO](#) web site to obtain your copy.

[X-Ray Data Booklet \(October 2009\)\[pdf\]](#)


<http://xldb.lbl.gov/>



Voltage setting and region of interest



Excitation / Emission of Characteristic X-ray Radiation



Theory: X-ray identification (open architecture)

- **“organic”** elements are not identified in spectrum (absorbed by air)
- **Low Z** elements only show K lines (L lines too low to be detected)
- **Medium Z** elements show both K and L lines
- **High Z** elements show L and some M lines (K can't be excited)

Bruker Tracer 5g
Can see F with He flush

1 H																	2 He
3 Li	4 Be											5 B	6 C	7 N	8 O	9 F	10 Ne
11 Na	12 Mg											13 Al	14 Si	15 P	16 S	17 Cl	18 Ar
19 K	20 Ca	21 Sc	22 Ti	23 V	24 Cr	25 Mn	26 Fe	27 Co	28 Ni	29 Cu	30 Zn	31 Ga	32 Ge	33 As	34 Se	35 Br	36 Kr
37 Rb	38 Sr	39 Y	40 Zr	41 Nb	42 Mo	43 Tc	44 Ru	45 Rh	46 Pd	47 Ag	48 Cd	49 In	50 Sn	51 Sb	52 Te	53 I	54 Xe
55 Cs	56 Ba	57 La	72 Hf	73 Ta	74 W	75 Re	76 Os	77 Ir	78 Pt	79 Au	80 Hg	81 Tl	82 Pb	83 Bi	84 Po	85 At	86 Rn
87 Fr	88 Ra	89 Ac	104 Rf	105 Db	106 Sg	107 Bh	108 Hs	109 Mt	110	111	112	113	114				

Na can be identified with special instrumentation setup or with Tracer 5g

58 Ce	59 Pr	60 Nd	61 Pm	62 Sm	63 Eu	64 Gd	65 Tb	66 Dy	67 Ho	68 Er	69 Tm	70 Yb	71 Lu
90 Th	91 Pa	92 U	93 Np	94 Pu	95 Am	96 Cm	97 Bk	98 Cf	99 Es	100 Fm	101 Md	102 No	103 Lr

Ceramics – complexity for analysis and approaches to achieving high quality data



Ceramic types

Earthenware



<http://www.mfa.org/artemis/zoom.asp?oid=315436&pic=SC37245.fpx&tem=1>



China



<http://netra.glendale.cc.ca.us/ceramics/bowkakiemonplate.html>

Ceramic Bodies	Porosity	Firing Temp	Comment
• Terra-cotta	High ~30%	well below 900	Unglazed, coarse fabric
• Earthenware	~10-25%	900-1200	Glazed or unglazed, non vitrified
• Stoneware	~0.5-2%	1200-1350	Glazed or unglazed, vitrified body
• China	~<1%	1100-1200	White fabric, vitrified
• Porcelain	Near 0%	1300-1450	Hard body, translucent, fine fabric

• Ceramic composed mainly of aluminosilicates (kaolinites smectites, montmorillonite and vermiculites)



Typical elements of interest in ceramics analysis

(Yellow = major, Red = trace)

Lower Z

Periodic Table of the Elements

Atomic Number

Symbol

Name


Atomic Mass

Atomic mass values reflect the IUPAC accepted values as of 09/2013.
 Masses expressed in [x,y] format show the lower and upper limit of atomic mass depending on the physical and chemical history of the element.
 Masses expressed in <x> format are the mass numbers of the longest-lived isotope for elements with no stable nucleus.

1 1A 1A 1 H Hydrogen [1.00784;1.00811]	2 IIA 2A 4 Be Beryllium 9.0121831(5)											13 IIIA 3A 5 B Boron [10.806;10.821]	14 IVA 4A 6 C Carbon [12.0096;12.0116]	15 VA 5A 7 N Nitrogen [14.00643;14.00728]	16 VIA 6A 8 O Oxygen [15.99903;15.99977]	17 VIIA 7A 9 F Fluorine [18.998403163(6)]	18 VIIIA 8A 2 He Helium 4.002602(2)		
3 Li Lithium [6.938;6.997]	11 Na Sodium [22.98976928(2)]	12 Mg Magnesium [24.304;24.307]	21 Sc Scandium 44.955908(5)	22 Ti Titanium 47.867(1)	23 V Vanadium 50.9415(1)	24 Cr Chromium 51.9961(6)]	25 Mn Manganese [54.938043(5)]	26 Fe Iron [55.845(2)]	27 Co Cobalt [58.933194(4)]	28 Ni Nickel [58.9334(4)]	29 Cu Copper [63.546(3)]	30 Zn Zinc [65.38(2)]	31 Ga Gallium [69.723(1)]	32 Ge Germanium [72.630(6)]	33 As Arsenic [74.921595(5)]	34 Se Selenium [78.971(8)]	35 Br Bromine [79.90179(9)]	36 Kr Krypton [83.798(2)]	
19 K Potassium [39.0983(1)]	20 Ca Calcium [40.078(4)]	37 Rb Rubidium [85.4678(3)]	38 Sr Strontium [87.62(1)]	39 Y Yttrium [88.90584(2)]	40 Zr Zirconium [91.224(2)]	41 Nb Niobium [92.90637(2)]	42 Mo Molybdenum [95.95(1)]	43 Tc Technetium <98>	44 Ru Ruthenium [101.07(2)]	45 Rh Rhodium [102.90550(2)]	46 Pd Palladium [106.42(1)]	47 Ag Silver [107.8682(2)]	48 Cd Cadmium [112.4114(4)]	49 In Indium [114.818(1)]	50 Sn Tin [118.710(2)]	51 Sb Antimony [121.760(1)]	52 Te Tellurium [127.60(3)]	53 I Iodine [126.90447(3)]	54 Xe Xenon [131.293(6)]
55 Cs Cesium [132.90545196(9)]	56 Ba Barium [137.327(7)]	57-71 Lanthanide Series	72 Hf Hafnium [178.49(2)]	73 Ta Tantalum [180.94789(2)]	74 W Tungsten [183.84(1)]	75 Re Rhenium [186.207(1)]	76 Os Osmium [190.23(3)]	77 Ir Iridium [192.217(3)]	78 Pt Platinum [195.084(9)]	79 Au Gold [196.966569(5)]	80 Hg Mercury [200.592(3)]	81 Tl Thallium [204.382,204.385]	82 Pb Lead [207.2(1)]	83 Bi Bismuth [208.98040(1)]	84 Po Polonium <209>	85 At Astatine <210>	86 Rn Radon <222>		
87 Fr Francium <223>	88 Ra Radium <226>	89-103 Actinide Series	104 Rf Rutherfordium <261>	105 Db Dubnium <268>	106 Sg Seaborgium <271>	107 Bh Bohrium <272>	108 Hs Hassium <278>	109 Mt Meitnerium <276>	110 Ds Darmstadtium <281>	111 Rg Roentgenium <289>	112 Cn Copernicium <285>	113 Uut Ununtrium unknown	114 Fl Flerovium <289>	115 Uup Ununpentium unknown	116 Lv Livermorium <293>	117 Uus Ununseptium unknown	118 Uuo Ununoctium unknown		
			57 La Lanthanum [138.90547(7)]	58 Ce Cerium [140.116(1)]	59 Pr Praseodymium [140.90766(2)]	60 Nd Neodymium [144.242(3)]	61 Pm Promethium <145>	62 Sm Samarium [150.36(2)]	63 Eu Europium [151.964(1)]	64 Gd Gadolinium [157.25(3)]	65 Tb Terbium [158.92535(2)]	66 Dy Dysprosium [162.500(1)]	67 Ho Holmium [164.93033(2)]	68 Er Erbium [167.256(3)]	69 Tm Thulium [168.93422(2)]	70 Yb Ytterbium [173.054(5)]	71 Lu Lutetium [174.9668(1)]		
			89 Ac Actinium <227>	90 Th Thorium [232.0377(4)]	91 Pa Protactinium [231.03588(2)]	92 U Uranium [238.02891(3)]	93 Np Neptunium <237>	94 Pu Plutonium <244>	95 Am Americium <243>	96 Cm Curium <247>	97 Bk Berkelium <247>	98 Cf Californium <251>	99 Es Einsteinium <252>	100 Fm Fermium <257>	101 Md Mendelevium <258>	102 No Nobelium <259>	103 Lr Lawrencium <262>		

Heterogeneity of Ceramics

- Variable particle size based on –
 - Fabric type
 - Temper
 - Mineral inclusion
 - Firing temperature – vitrification
 - Porosity
 - Etc..
- Potential for variable surface chemistries
- Increased Matrix effect possible



Importance of homogenization


Specimen Heterogeneity
(Variable particle size and distribution)




i.e.

- Porosity
- Coarse ceramic body


Homogenize



Homogenous Sample
(Powdered or Fused)




Specimen Heterogeneity
(Layered Structure)




i.e.

- Corroded Material
- **Surface alterations**
- Layered structures

Homogenize



Homogenous Sample
(Powdered or Fused)




Heterogeneous materials allow for additional attenuation altering resulting peak height – thus changing the reported chemistry

X-ray physics - Quantification

Physics calculations based on homogenous materials that are infinitely thick.

"In practice, however, there are two major constraints that may prevent this ideal circumstance (*quantification*) from being achieved, these being **sample size** and **samples heterogeneity**" (Jenkins 1999: 141).

Samples prepared from heterogeneous materials **require significant** pre-treatments prior to analysis



- T = Sample Thickness
- x = Path Length
- d = Penetration Depth
- ψ_1 = Incident Angle
- ψ_2 = Take-Off Angle
- ψ = Scattering Angle
- μ = Mass Absorption Coefficient
- ρ = Density

- Transmitted Beam Intensity: $I(\lambda_0) = I_0(\lambda_0) \exp(-\mu\rho x)$
- Incoherent Scatter Wavelength: $\lambda_c - \lambda_0 = 0.0243 [1 - \cos\psi]$
- Fluorescence Penetration Depth: $d = x \sin\psi_2$

Figure 1.3. Interaction of the primary X-ray beam with the sample.
From Jenkins 1999

Change in surface chemistry due to deposition

- First pointed out by Sayre, Dodson and Burr Thompson in 1957.
- Alkali and alkali-earth metals are most reactive with water and typically the most mobile— Leaching of cesium and strontium (among others) noted in low-fired ceramic.
- Reduction of Ca, Cs, Rb, K, and Na at surface.
- Vitreous phases are prone to breakdown and leaching as glass would.
- Dependent on many factors - pH, redox conditions, organic activity etc...
- Mobile soluble salts.



DESALINATING ARCHAEOLOGICAL CERAMICS

December 20, 2013 · by Sarah Klopf ·


<https://rammleventisproject.wordpress.com/2013/12/20/desalinating-ceramics/>

Schwedt, A., Mommsen, H. and Zacharias, N., 2004 *Archaeometry*, 46(1), pp.85-101


Buxeda iGarrigóset al. 2001, 360–2.

Golitko, M., Dudgeon, J.V., Neff, H. and Terrell, J.E., 2012 *Archaeometry*, 54(1), pp.80-100.

Sample 1	Al	Si	P	K	Ca	Fe
Crushed cupped	1.4	8.3	0.13	1.35	15.78	2.65
Cleaned	3.46	12.85	0.29	2.07	21.05	2.74
Red Slip	2.91	13.43	0.41	3.71	8.1	6.09
Crushed bagged	-0.37	3.51	0.02	1.18	14.18	2.44



Effect of Moisture on calibration




- Dry = Green: Wet = Red normalized to Rh La
- Note the increased area under peak for all light elements with dry sample. Trace elements of Mg, Al and P **highly attenuated** by the addition of water.
- The resulting quantification would be a **52%** decrease for **Al** and a **47 %** decrease for **P**.

Overcoming heterogeneity and surface alterations

- Create standards (reference materials) that
 - a) **Match** the composition (dynamic range and all elements present),
 - b) have **similar matrix** and **particle size** to the samples.
- Only when they are prepared in the same way can you hope to create a more linear relationship used for calibration.
- When standards (reference materials) and samples are prepared the same way, **ONLY then**, can you achieve accuracy and precision in calibration.
- There is **no substitute** for good sample preparation. If the sample is poorly prepared **NO AMOUNT OF CORRECTION** can improve the analytical results. Time and resources should be devoted to sample preparation to benefit fully from an XRF spectrometer.

Quantification – background, issues and application to ceramics



XRF quantification

Assumed factors (infinitely thick samples)

- Samples are **homogenous**
- Samples have similar – **small particle size**
- Samples have **flat smooth surfaces**
- Samples have similar **density** to each other **AND** the reference materials used for calibration.
- Samples are '**infinitely thick**' (for bulk samples)
- Samples are **large** enough to cover X-ray beam
- Samples are **dry**
- Samples are **not porous**

Dependent on several factors:


- Analyte **concentration**
- The **matrix**
- Accompanying elements (**matrix effects**)
- Sample type (solid, liquid, powder, fused etc...)
- Method of **sample preparation**
- Shape, thickness
- Instrumentation setup (geometry, spot size etc..)




Typical XRF metal reference standards: thick, flat, smooth and homogenous.

Typical analytical range

- Measured readily
- Measured with difficulty




Quantification: Peak height directly related to concentration



	Mn	Fe	Co	Ni	Cu	Zn	As	Pb	Ag	Cd	Sn	Sb
Red	0.72	1.20	0.14	0.50	78.96	0.49	0.06	0.26	0.10	0.15	16.05	0.66
Green				0.00	97.08		2.90	0.00			0.01	

Methods offered to achieve quantification



Newer approaches include:
Machine learning and
Artificial intelligence algorithms

- Fundamental parameters
 - Assumes:
 - Homogenous matrix
 - Values total 100% (usually)


Empirical calibrations

- Compton normalization
 - Assumes:
 - Lighter element composition
 - Values provided in ppm or volume/area (usually)


Theoretical calibrations

Accuracy, Precision and Bias


- **Accuracy** : Accuracy is a measure of how **close** the measured value is to the true value.
- **Precision**: Precision is a measure of the **agreement** among a group of individual measurements. (How close repeat measurements are to one another.)
- **Bias**: Bias is due to **systematic errors** such as a change in voltage since calibration or wrong calibration constants would introduce a constant error into each measurement. Bias can be reduced by calibration and carefully following established measurement procedures.




1
Unbiased
Imprecise
Inaccurate



2
Biased
Imprecise
Inaccurate



3
Biased
Precise
Inaccurate





4
Unbiased
Precise
Accurate

Matrix Corrections – Compton Normalization

- Based on the Compton scattering where Compton scatter is used as an internal standard and a matrix correction.

$$I / I_{\text{Compton}} \sim C (+ \text{ corrections})$$

- Matrix related corrections are calculated based on known concentrations and are determined based on calibration using standards (i.e., Fe, As_2O_3).
- The intensity of the Compton scatter can be used to obtain an estimate of the absorption coefficient of the specimen.
- Best for $Z > 23$ and for samples that are relatively low density and samples containing relatively low concentrations of elements (i.e., soil, **potentially ceramics**).



Matrix Corrections – Fundamental Parameters

- Fundamental parameters (FP) calibrations are typically very accurate if the intensities of all matrix elements can be determined, or if the matrix of a sample is known.
- Calibration is **KEY** to success – either through pure-element standards, or through theoretical calculations (standardless analysis).
- Can greatly improve results by using (one or many) standard reference material(s) during calibration.


Matrix Corrections – Fundamental Parameters

Advantages


- Method can be applied to any sample type
- Do not need to use standards (expensive to purchase)

Disadvantages

- Do not consider all physical interactions (i.e. tertiary fluorescence, scatter of low Z elements etc.)
- Strongly dependent on **consistent** mass attenuation coefficients, fluorescence yield, instrument geometry tube energy etc...
- Can only be applied if **ALL** the elements in the sample are identified (low Z elements like C, O, N should not be present).
- Often assumes normalization to 100%
- Not reliably consistent between instrumentation (even with same calibration).



FP misidentification/peak overlap



Reported Results (in ppm)



Ti	Cr	Fe	Zn	As	Ag	Cd	Sn	Pb
18992	2360	2472	436	50637	2776	764	1654	476703

Matrix Corrections – Influence coefficient algorithms

- Uses the relationship between counts and concentration (simple linear calibration curve). Matrix effects can affect a good linear calibration curve but can be compensated for with 'correction' models called empirical or theoretical coefficients.

$$C \sim I (+ \text{ corrections})$$

- This method can be very accurate – only limited by experimental design.
- Can provides **more accurate results** than FP and Compton Normalization.
- Preparation of reference materials manufactured of the element(s) of interest in a matrix that **closely approximates** that of the unknown samples is **necessary**.
- The matrix correction and calibration are clearly seen so interpretation of the method is easy.




(a) Lachance-Trail algorithm and (b) Claisse-Quintin algorithm - Sitko and Zawisza, (2012)

Spectral Interference Effects on calibration


Spectral interference – peak overlap – can lead to increased error in the calibration curve as well.

Think of $\text{CuK}\beta$ and $\text{ZnK}\alpha$, or $\text{PbL}\alpha$ and $\text{AsK}\alpha$.

This error is improved with increased resolution – SDD detector ($\sim 145\text{eV}$) vs SiPin detector ($\sim 205\text{eV}$).



Increased resolution in new detectors to help with peak overlap



Difference in resolution between a standard Pin diode detector (Green) and an SDD detector (Red). Typical pin diode resolution at **205 eV** while SDD at **145 eV** at MnKa

Issues related to achieving quality quantification

- Matrix Issues
- Surface Roughness
- Particle Size Effects
- Heterogeneity
- Mineralogical Effects
- Curvature
- etc...


Poor Analytical Accuracy



Section of a Chalcolithic crucible showing high porosity, organic and inorganic temper and potential alteration to surface chemistry


Matrix Effects on calibration

- Matrix effects are caused by the absorption and enhancement of characteristic x-ray photons within the sample.
- Can cause alteration to linear relationship established with pure element calibrations.




Relationship between radiation intensity of Fe and weight fraction of Fe:

- Curve A – negligible matrix effects
- Curve B – FeCr,
- Curve C – FeNi,
- Curve D – FeMn



Remember - XRF is a surface sensitive analytical technique



- No excitation within the upper part of the samples(white)
- Central portion of the sample will be excited by the emitted radiation will be absorbed within the sample (light yellow)
- The measured radiation comes from the surface layers (yellow)
- This 'escape potential' can be measured based on mass attenuation coefficient

$$\frac{\ln(I_0/I)}{u/p \times \text{density}}$$

Escape Potential in Different Matrices (dependent on density)

Line	Energy	Graphite	Glass	Iron	Lead
Cd KA1	23.17 keV	14.46 cm	8.20 mm	0.70 mm	77.30 μm
Mo KA1	17.48	6.06	3.60	0.31	36.70
Cu KA1	8.05	5.51 mm	0.38	36.40 μm	20.00
Ni KA1	7.48	4.39	0.31	29.80	16.60
Fe KA1	6.40	2.72	0.20	* 164.00	11.10
Cr KA1	5.41	1.62	0.12	104.00	7.23
S KA1	2.31	116.0 μm	14.80 μm	10.10	4.83
Mg KA1	1.25	20.00	7.08	1.92	1.13
F KA1	0.68	3.70	1.71	0.36	0.26
N KA1	0.39	0.83	1.11	0.08	0.07
C KA1	0.28	* 13.60	0.42	0.03	0.03
B KA1	0.18	4.19	0.13	0.01	0.01

Approx. Density:

~2.09–2.23 g/cm³

~2.5 g/cm³

~7.874 g/cm³

~11.34 g/cm³

Photon Escape Potential for Ceramics

ELEMENT	PHOTON ESCAPE DEPTH	UNIT
Si	24	μm
Ca	70	μm
Fe	271	μm
Rb	0.19	cm
Sr	0.27	cm
Y	0.3	cm
Zr	0.38	cm


Table 3: Approximate photon escape potential for several elements of interest (based on a density of 2.3 g/cm³). Note that the higher atomic number elements have the potential to escape from a much larger volume of ceramic.

Attenuation of Mylar and Polyethylene based on polymer thickness for low Z elements

Al Si P


K Ca

Mn Fe




Density (g/cm³): Mylar = 1.4, Polyethylene = 0.93

Sample size and curvature



Effect on spectrum with change in sample size and curvature.




Success in analysis – suggestions

Sample preparation
Collecting data

Options for Ideal Sample and Reference Material Preparation

- Sample should be ground to uniform particle size
- Ideally the particle size should be much smaller than the analyzed layer depth



- Powdering, Pressing and Fusing are all options.

Rocks, minerals, ores, ceramics, . . .



Crushing,
pulverizing



Pouring

Pressing

Fusing



Sample preparation to improve homogeneity

- Each has benefits and disadvantages.
- Maintaining consistency in your methodology is key to success.

Powdered samples











Pressed pellets



Borax glass beads / Fusion



Comparison of **WD-XRF**, **ED-XRF**, **pXRF-author built** **Pressed Pellet Cal**, **pXRF- Bruker mudrock general cal**




Suggestions for Success – Collecting Data

- Take advantage of X-ray Physics – Change your instrument setting (kV, current, filters, time, and atmosphere) to enhance your data collection.
- Low Z elements – long count times, low kV, relatively high current, He flow, no filter.
- Mid – high Z elements – long count times, high kV, relatively high current, no He, optimal filter.

Table 4


Analytical protocols used to evaluate performance of clay/sediment and 'mudrock' PXRF calibrations.

	CAIS Clay/Sediment calibration	Bruker mudrock calibration
Low Z elements	200 s count time	180 s count time
	15 kV/20 μ A	15 kV/20 μ A
	He flow	Vacuum
	No filter	No filter
Mid Z elements	200 s count time	60 s count time
	40 kV/30 μ A	40 kV/30 μ A
	'green' filter:	'yellow' filter:
	300 μ m Al/20 μ m Ti/150 μ m Cu	300 μ m Al/20 μ m Ti



Creating a Calibration – Cross Analysis

- In addition to proper sample preparation – obtain quality **comparable data**
 - NAA and/or ICP-MS can provide reliable compositional data to build a cross reference set.
- Retain some samples for a **Validation Set** to test your calibration
 - Comparing your cal to samples used to build your calibration is NOT recommended.
- Different XRF manufacturers have different options for calibration programs (i.e. EasyCal – Bruker), or use a self standing software (i.e. Cloudcal - v3.0 Drake, B.L. 2018. CloudCal v3.0. GitHub. <https://github.com/leedrake5/CloudCal>. doi: 10.5281/zenodo.2596154)



New investigations in ceramics analysis

Scientific investigation of ceramics materials from Honduras



(Shugar, Schortman and Urban, NSF award # 1733890)

- Aim – investigate creating two different calibration models

- Sanded flat and cleaned fragments.
- Thin films

Thin samples


- No matrix effects
- Intensity is linear function of mass per unit area of analyte




Scientific investigation of ceramics materials from Honduras

(Shugar, Schortman and Urban, NSF award # 1733890)

- 430 samples from archaeological context in Honduras
- Preparation
 - Sanded flat (both sides) and cleaned fragments.
 - Thin films of 0.02 grams from drilled cores.
- Analysis
 - Using empirical calibrations and new theoretical ones




Sanded and
cleaned




Drilled
and
thin
filmed

Solid 45 kV, 60 uA Cu/Ti/Al filter 60 sec.

Thin Film 45 kV, 60 uA Cu/Ti filter 60 sec.





Thin Film and Solid 10 kV, 120 uA no filter No He flush 60 sec.





Preliminary Results

- Solid samples trace element calibration curves for Rb and Sr
 - Lucas-Tooth Calibration with total counts normalization



Preliminary Results


- Thin Film samples Major elements calibration curves for Si and Ca
 - Lucas-Tooth Calibration with time normalization



Newer Machine learning algorithms

- Forest algorithm
- Neural Networks
- XGBoost


- Rahman, Ashfaque, et al. "A machine learning approach to find association between imaging features and XRF signatures of rocks in underground mines." *2015 IEEE SENSORS*. IEEE, 2015.
- Radtke, Martin. "Machine learning for direct quantification of XRF measurements." (2019).
- Heginbotham, Arlen, Robert Erdmann, and Lee-Ann C. Hayek. "The dating of French gilt bronzes with ED-XRF analysis and machine learning." *Journal of the American Institute for Conservation* 57.4 (2018): 149-168.



XGBoost calibration curve for ceramics solid trace scan for Na and Al – courtesy of Lee Drake

Potentially accept qualitative data comparison

- Faster results – reduced cost – can always apply calibration in the future
- Look at the area under peak for elements of interest.
- In many cases, grouping can be extracted.



Area under peak
elemental ratios to
define ceramic grouping

Summary

- Come prepared –
 - Know what your samples is composed of
 - If you need quantification –
 - you must create reference materials that fully resemble your unknowns
 - Cross analyze your unknowns using a different technique (i.e. ICP)
 - Be ready to properly prepare your samples
 - Reduce the variables
 - Be prepared to accept qualitative data
 - Ideal for survey of large assemblages
- The Future is bright!
 - The development of XRF is continuing
 - Newer technologies are pushing the boundaries of what is possible
 - Development of new algorithms is making calibration of complex systems easier



New Tracer 5g – Graphene window on detector allows for the identification of F