


Non-targeted metabolomics CCS-enabled annotation workflows for ID and detailed annotation confidence reporting






Metabolomics 2021 – Poster 244

Overview

Aiko Barsch; Matthias Szesny; Ulrike Schweiger-Hufnagel; Sofie Weinkouff; Nikolas Kessler
Bruker Daltonics GmbH & Co. KG,
28359 Bremen, Germany

Lack of exhaustive repositories providing retention time, MS/MS spectra and CCS values hamper the annotation and identification (ID) of target compounds in Metabolomics research. Researchers require a solution that automatically and transparently annotates features. Here we present a solution for automatic annotation of targets with up to 5 confidence criteria. Additional, tentative annotation of knowns and predicted known compounds is supported by automatic *in-silico* fragmentation and CCS prediction based on the novel CCS-Predict Pro model. Paired with customizable annotation quality scoring and visualization the presented annotation workflows enable researchers to assess and report ID level confidence suitable for the study, as recommended by Schymanski *et al.* [1].

A) 4D-Metabolomics™: RP-LC-TIMS-MS/MS featuring PASEF provides five indicators of data quality



LC: Elute UHPLC
MS: timsTOF Pro 2

- Retention time
- Accurate mass
- True isotopic pattern

4. TIMS (Trapped Ion Mobility Separation) = reproducible **CCS values**

5. PASEF® data acquisition provides clean MS/MS spectra by **MOMA** = Mobility Offset Mass Aligned

B) High confidence annotation of knowns

Enabled by Target List containing CCS values for metabolites derived from CCS Compendium [2], appended by RT and MS/MS references

ID criteria	m/z	Isotopic pattern	RT	CCS	MS/MS
Target List	✓	✓	Reference	Reference	Reference

Unified CCS Compendium [2] 1438 compounds
HMDB Metabolite Library BRUKER Daltonics >600 RTs
Unified CCS Compendium [2] CCS values 1438 compounds
MS/MS spectra for > 800 compounds
MetaboBASE Personal Library BRUKER Daltonics
NIST Library BRUKER Daltonics > 27,000 compounds

Name	Molecular Formula	RT [min]	CCS [M+Na] ⁺	CCS [M+H] ⁺	CCS [M-H] ⁻	Library Compound
N-Acetylsertotonin	C ₁₂ H ₁₄ N ₂ O ₂	6.44	153.8	149.8	150.6	N-Acetylsertotonin
5-Hydroxyindoleacetate	C ₁₀ H ₉ NO ₃	6.29	148.6	141.1	137.3	5-Hydroxyindoleacetic acid
D-Glucuronolactone	C ₁₂ H ₁₆ O ₆	6.13	146.0			D-Glucurono-6,3-lactone
3,4-Dihydroxyphenylacetate	C ₈ H ₈ O ₄	5.73				3,4-Dihydroxybenzoic acid

Name	Molecular Formula	AQ	AQ Details	Δm/z [ppm]	mSigma	ΔRT	ΔCCS [%]	MS/MS score
Acetaminophen	C ₉ H ₉ NO ₂	5	131.4	132.7	Acetaminophen			
Acetaminophen	C ₉ H ₉ NO ₂	5	131.4	132.7	Acetaminophen			
Phenylacetate-L-glutamine	C ₁₃ H ₁₆ N ₂ O ₄	4	1.899	9.5	0.16	0.8	969.4	
Xanthurate	C ₁₀ H ₇ NO ₄	4	1.437	9.2	0.16	0.6	747.0	
5'-Methylthioadenosine	C ₁₁ H ₁₃ N ₅ O ₃ S	4	-0.083	8.8	0.19	1.3	947.9	
3-Methyl-L-histidine	C ₇ H ₁₁ N ₃ O ₂	4	-0.779	9.0	0.16	1.2	918.9	

Structure of annotated compound, here Acetaminophen, enables triggering of metabolite prediction for this drug compound (see D).

Annotation Quality Scoring (AQ): Access the confidence in compound annotation with a glimpse!

- 1 Mass accuracy
- 2 RT deviation
- 3 mSigma value
- 4 MS/MS score
- 5 CCS deviation

Higher confidence (green icons)
Lower confidence (red icons)

C) Tentative annotation of expected knowns

Target List based annotation for tentative ID based on HMDB 4.0 [3] derived urine metabolites

ID criteria	m/z	Isotopic pattern	RT	CCS	MS/MS
Target List	✓	✓		Prediction - CCS-Predict Pro	<i>in-silico</i> fragmentation - MetFrag

hmdb 4306 compounds
InChI encoded structures enable:
• CCS prediction
• *in-silico* fragmentation

Name	Molecular Formula	InChI	CCS	MS/MS
1-Methylhistidine	C ₇ H ₁₁ N ₃ O ₂	InChI=1S/C7H11N3O2/c1-10-3-5(9-11)-2		
1,3-Diaminopropane	C ₃ H ₁₀ N ₂	InChI=1S/C3H10N2/c4-2-1-3-5/h1-3		
2-Ketobutyric acid	C ₄ H ₆ O ₃	InChI=1S/C4H6O3/c1-2-3(5)4(6)7/h		
2-Hydroxybutyric acid	C ₄ H ₈ O ₃	InChI=1S/C4H8O3/c1-2-3(5)4(6)7/h		
Glycylleucine	C ₈ H ₁₆ N ₂ O ₃	InChI=1S/C8H16N2O3/c1-5(2)3-6(8)		

D) Tentative annotation of BioTransformer [4] predicted knowns

MetaboScape enables triggering of local BioTransformer based prediction of metabolites based on target structures. Visualization allows for the review of the generated metabolism tree and automatically matched compounds in the Feature Table:

BioTransformer: Predicted Metabolite Structures

Annotations were found for metabolites highlighted in boxes

Name	Molecular Formula	AQ	AQ Details	Δm/z [ppm]	mSigma	ΔRT	ΔCCS [%]	MS/MS score
Acetaminophen	C ₉ H ₉ NO ₂	5	131.4	132.7	Acetaminophen			
Acetaminophen + CH ₂ O	C ₉ H ₉ NO ₃	4	-1.088	6.1	0.14	0.0	902.2	
Acetaminophen + C ₂ H ₄ O ₄	C ₁₁ H ₁₃ NO ₆	4	0.590	2.2		-0.3	1000.0	
Acetaminophen + C ₂ H ₃ O ₂	C ₁₁ H ₁₁ NO ₃	4	-0.068	18.8		0.3	968.0	
Acetaminophen + O ₂ (1)	C ₉ H ₉ NO ₃	4	0.485	15.9		-2.7	987.3	

Icon indicates: Matching to reference MS/MS spectrum
Icon indicates: CCS matching to predicted value by CCS-Predict Pro

Combined, AQ and AQ Details enable to easily assess confidence in ID and subsequent custom reporting of ID confidence levels

Visual inspection of MetFrag generated *in-silico* fragment structures enables to readily investigate tentative structure annotations.

Note: HMDB and CCS Compendium are no Bruker products.

Methods

- Sample:** Human Urine; centrifuged and filtered (0.22µm), 2µl injected, 3 replicates
- LC:** Elute UHPLC, Intensity Solo C18 column (Bruker).
 - Gradient: Acetonitrile / Water based LC gradient according to T-ReX LC-QTOF Solution (Bruker), allows matching of retention times for >600 compounds
- MS:** timsTOF Pro 2 (Bruker)
- Acquisition:** PASEF positive mode
- Software:** MetaboScape 2022a, preliminary Version (Bruker).
- Libraries:**
 - Bruker HMDB Metabolite Library 2.0
 - Bruker MetaboBASE Personal Library 3.0
 - Bruker NIST 2020 Mass Spectral Library
- CCS reference values:** CCS Compendium [2]

Summary

The described workflows provide automated annotation routines across different confidence levels for known and predicted metabolites, building on trapped ion mobility separation and MS/MS spectra.

- References:
- <https://doi.org/10.1021/es5002105>
 - <https://doi.org/10.1039/C8SC04396E>
 - <https://hmdb.ca/>
 - <https://doi.org/10.1186/1471-2105-11-148>
 - <https://doi.org/10.1186/s13321-016-0115-9>
 - <https://doi.org/10.1186/s13321-018-0324-5>

Conclusions

- MetaboScape provides users highest confidence in ID by matching up to **5 confidence criteria** including **reproducible CCS values**
- MetaboScape's novel **CCS-Predict Pro** model enables small molecule CCS prediction based on the molecular structures of target compounds
- Automatic *in-silico* fragmentation and CCS prediction** (in case of lacking reference MS/MS spectra and / or CCS values) provides users **higher confidence in annotations**
- Reporting of **Annotation Quality Scores** and **Annotation Quality Details** enables researchers to readily **assess and report their ID level confidence**

4D-Metabolomics